

# °MEQU °M WARMER SYSTEM - STARTER KIT PLUS

## A Portable Blood and IV Fluid Warming System

°MEQU has developed an IV-blood warmer for emergency transfusions. It's a portable fluid warmer which warms up blood from cold to body temperature within seconds.

### Fast set-up time and intuitive use

The set-up time is less than 30 seconds, and the user-friendly design ensures that the °M Warmer System cannot be assembled incorrectly.

### High flow rates and large warming capacity

Blood and other infusion fluids can be heated from 5°C to 37°C at flowrates up to 150 mL/min. One charged battery can heat up to 2 liters of cold fluids 5°C or 4 liters of fluids at room temperature ( 21°C) to body temperature.

### Low weight and compact design

The °M Warmer (5pc) has a compact design and weighs from 1.8 lbs. (820 g) making it portable and ideal for prehospital use.

### Adhesive for secure attachment close to the infusion site

The single-use warmer's small size (3.94 in. x 2 in. x 0.79 in. / 10 cm. x 5cm. x 2 cm.) allows it to be attached to the patient close to the infusion site, using its integrated adhesive pad. This secures infusion site and reduces heat loss in the tubing.

## FEATURES/SPECIFICATIONS

- **Easy Setup:** User-friendly design; Standard Luer Locks
- **Lightweight:** Battery with warmer only 1.8 lbs.
- **Instant Warming:** Warms fluids in less than 10 seconds
- **Meets TCCC Guidelines:** Warms fluids to body temperature at flowrates up to 150 mL/min.
- **Large Warming Capacity:** Warms up to 2L of cold fluids on one fully charged battery
- **No Maintenance:** Wipe clean and recharge after use

- **Field-proven:** In use since 2018 in both Military, HEMS, EMS and Hospitals
- **Dust and Splashproof:** IP54 battery and warmer
- **For the entire continuum of care:** °M Warmer can stay with the patient for up to 72 hours
- **Adhesive Patch:** For secure attachment close to the infusion site
- **Interoperable:** Seamlessly supporting continuity of care across the entire patient journey



**°MEQU products are restricted to DoW and Federal Purchases only. Please contact Customer Service for details.**

*°Caution: Federal (USA) law restricts this device to sale by or on the order of a practitioner licensed to use or order use of the device or Federal (USA) law prohibits dispensing without a prescription.*

# °MEQU

## DIMENSIONS

- °M Warmer: W 1.97 in. x L 3.94 in. x H 0.79 in. (5 cm x 10 cm x 2 cm)
- Power Pack\*: W 3.54 in. x L 7.09 in. x H 1.38 in. (9 cm x 18 cm x 3.5 cm)  
» Cable length: 47.24 in. (120 cm)
- AC Charger: W 2.95 in. x L 4.61 in. x H 1.73 in. (7.5 cm x 11.7 cm x 4.4 cm)  
» Cable length: 66.14 in. (168 cm)

## WEIGHT

- °M Warmer: 0.24 lbs. (110g.)
- Power Pack\*: 1.57 lbs. (710g.)
- AC Charger: 0.55 lbs. (250g.)

## °MEQU °M WARMER SYSTEM - STARTER KIT PLUS

ITEM #  
35-0048

*\* Explore other °MEQU products at [www.NARescue.com](http://www.NARescue.com)*



**Learn more about this and other Products with a Mission® at [www.NARescue.com](http://www.NARescue.com)**



**NORTH AMERICAN RESCUE®**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

# °MEQU °M WARMER SYSTEM - STARTER KIT PLUS

## all KITS INCLUDE:



## Operating Parameters

- Output temperature range: 39°C (+/- 3°C at flowrates from 2-150ml/min)
- Max fluid warming volume on one fully charged Power Pack Plus (output temperature of min. 36°C:
  - » Input temperature 5°C: Power Pack\*: 2L
  - » Input temperature 21°C: Power Pack\*: 4L
- Priming volume: 3.5 ml (including tubing and luer locks)
- Maximum inlet pressure: 300 mmHg

## Environmental

- Continuous Operation:
  - » Relative Humidity RH: 15-95%
  - » Pressure: 60-106kPa
  - » Temperature: Charging the Power Pack Plus: 10°C - + 40°C  
Fluid warming mode (any power source connected): 0°C - + 40°C
- Storage and transportation:
  - » Relative Humidity RH: 15-95%
  - » Pressure: 60-106kPa
  - » Temperature: -20°C - + 50°C
- Transient Operating Conditions:
  - » Relative Humidity RH: 15-95%
  - » Pressure: 60-106kPa
  - » Temperature: -20°C - + 40°C

- Shock and vibration: Power Pack: ISO 1789 + A2 Sec. 6.4
- Basic safety and essential performance: ISO 60601-1-12
- Electrical compliance: ISO 60601-1



Step-by-step guide for the MEQU °M Warmer System

Learn more about this and other Products with a Mission® at [www.NARescue.com](http://www.NARescue.com)



## °MEQU °M WARMER SYSTEM - STARTER KIT PLUS

ITEM #

35-0048

\* Explore other °MEQU products at [www.NARescue.com](http://www.NARescue.com)

**\*MEQU products are restricted to DoW and Federal Purchases only. Please contact Customer Service for details.**

\*Caution: Federal (USA) law restricts this device to sale by or on the order of a practitioner licensed to use or order use of the device or Federal (USA) law prohibits dispensing without a prescription.



**NORTH AMERICAN RESCUE®**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277