

# CYCLONE® 1000 POCKET BVM®

1000 mL

**"... In the case of impaired ventilation and uncorrectable hypoxia with an oxygen saturation less than 90%, the casualty should be ventilated with a 1,000 mL resuscitator bag valve mask."**

- Airway Management in TCCC, CoTCCC guideline update, 01/2024

The compact, fully functional adult BVM that reduces storage space by up to 75% compared to traditional manual resuscitators.

The Cyclone® 1000 Pocket BVM® is today's most compact, easy-to-use BVM available to provide ventilatory support to patients who are not breathing or who are in respiratory distress. Its unique collapsible design is quick-to-deploy, easy-to-operate and does not need to be manually inflated prior to use.

To help reduce rescuer fatigue, this BVM requires minimum pressure for compression of the bag to ventilate and quickly recovers its shape during use. The bag also has a textured grip to reduce slippage from a gloved hand. Packaged in a green, rugged, hard plastic, twist-top case, the Cyclone® Pocket BVM® comes complete with oxygen tubing for delivering supplemental oxygen when possible.

The Cyclone® 1000 Pocket BVM is a single-use, disposable device.

*The compact, fully functional adult BVM that reduces storage space by up to 75% compared to traditional manual resuscitators.*



## CYCLONE® 1000 POCKET BVM®

ITEM #	NSN	Rx ONLY
10-0081	6515-01-737-9112	

Learn more about this and other **Products with a Mission®** at [www.NARescue.com](http://www.NARescue.com)



# CYCLONE® 1000 POCKET BVM®

1000 mL


## How to deploy the improved Cyclone Pocket BVM:



▲ Separate the Face Mask from the folded Resuscitator Bag



◀ Pull both end valves (Front patient valve housing (air/O<sub>2</sub> exit port) and rear air/O<sub>2</sub> inlet valve) outwards in opposite directions while twisting the Patient Valve



◀ Attach the Face Mask to the Patient Valve Housing

## PRODUCT FEATURES

- Disposable, full Adult BVM that delivers up to 400-600 cc of air
- Compact adult, single-use, disposable BVM that reduces storage space by up to 75% of traditional BVM's
- Quick to deploy and does not need to be manually inflated prior to use
- Comes with 6 ft 6 in. of oxygen tubing for delivering supplemental oxygen when available
- Compact and collapsible BVM in a green rugged, hard plastic, twist-top case
- Top quality materials for top quality performance
- Includes a textured grip to reduce slippage from a gloved hand

## SPECIFICATIONS

- Tidal (Balloon) Volume: 1000 mL
- Suitable for use with the following body mass range: > 88 lbs.
- Dead space: < 6 mL
- Delivered oxygen concentrations under various test conditions: 55% (at 2 L/min) to 100% (at 8 L/min)
- Recommended operating temperature: 0°F to +122°F (-18°C to +50°C) to EN ISO 10651-4:2002
- Recommended storage temperature: < Tested at -40°F and +140°F (-40°C and + 60°C) according to EN ISO 10651-4:2002

## DIMENSIONS

- Packaged: Diameter 5 in x H 2.75 in.
- Weight: 16.18 oz.

Learn more about this and other **Products with a Mission®** at [www.NARescue.com](http://www.NARescue.com)

