

Intended Use

The BOA® Constricting IV Band provides circumferential pressure on either of the upper extremities for the purpose of obtaining peripheral venous access.

Warning

The BOA® should be used only by trained medical professionals. Improper use could result in injury to casualty. Reuse of this device will degrade efficacy, resulting in adverse casualty reaction, including potential death. Continually monitor casualty to ensure BOA® is functioning per medical protocols. In the event of a BOA® malfunction, follow local protocols and report any serious incident to North American Rescue or authorized representative and the competent authority of the Member State.

Contraindications

Within the European Union (EU), the BOA® Constricting IV Band is not indicated for pediatrics and pregnant women due to EU regulatory requirements.



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

35 Tedwall Court • Greer, South Carolina 29650

Phone: (864) 675-9800 • Fax: (864) 675-9880

BOA®-Constricting IV Band

REF 30-0009 NSN: 6515-01-537-2611


REF 30-0071 BOA XL



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

35 Tedwall Ct. • Greer SC 29650 • USA



MIDI Europa GmbH
30855 Hannover-Langenhagen
Germany




Instructions for Use



1. To find the correct size, place the BOA[®] (in its relaxed state) gently around the arm. It should measure 1/2 to 2/3 around the extremity.




2. Stretch BOA[®] straight out and place as high as possible around upper arm.



3. Keep fingers underneath the connectors as you secure in place to prevent pinching of the skin.



4. Press the "Quick Release" button and pull to remove the BOA[®].



5. Monitor casualties with compromising clinical conditions. Evacuate casualty to secondary treatment facility, advising follow on care to monitor casualty for allergic reactions or infections.

6. In the event of device issues refer to medical protocols.

Harmonized Standard Symbols:



Lot number



Single use



Medical Device



Consult instructions for use



Manufacturer



Authorized Representative



Device part number



Not for use in MRI



Date of manufacture



Expiration Date



Keep Away from Sunlight



Keep Dry