

TOURNIQUETS | FIRST AID KITS | TRAUMA KITS | RESCUE TOOLS | LITTERS

LE

LAW ENFORCEMENT
Catalog



Since 1996, North American Rescue® has focused on **increasing survivability** on both the battlefield and the homefront. With our innovative medical/rescue equipment, NAR® continues to lead the way by delivering the **tailored, mission-critical healthcare solutions** needed by today's first responders to counter the rapidly evolving threats of the 21st century. See what a lifesaving difference our combat proven make in the prehospital environment!



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS WITH A MISSION™



BE THE DIFFERENCE™

www.NARescue.com


CONTENTS

TOURNIQUETS		3-4
FIRST AID KITS		4 -16
TRAUMA KITS		17-20
MASS CASUALTY KITS		21-22
EXTRAS		23

DISCOUNT PRICING

Contact Customer Service to have your account upgraded to discount pricing. Shop your pricing on our full catalog of products, available online 24-7 at www.NARescue.com

888.689.6277





CUSTOMIZABLE SOLUTIONS

BUILT TO MEET YOUR REQUIREMENTS

If you haven't found exactly what you need, our team of *subject matter experts* at North American Rescue® are standing by to guide you to the *ideal lifesaving solution* tailored to meet all your specifications.

Whether it's the essential gear needed for individuals or customized kits for multiple casualties, a *North American Rescue® Specialty Kit* is your answer to the precise aid kit you and your organization need to *increase survivability* should injury occur.




Call or Email us TODAY to discover all your options for the tailored lifesaving solution YOU can live with!



Explore all the possibilities at
NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277



Massive Hemorrhage | Tourniquets

C-A-T® Combat Application Tourniquet Gen 7

- **Single route buckle** decreases application time to help decrease blood loss, allows for slack removal for proper TQ application
- **Applied same way** to both upper and lower extremities
- True **one-handed tourniquet**
- **Scientifically proven** 100% effective in occluding blood flow in both upper and lower extremities by the USAISR and recommended by CoTCCC



The FASTEST, SAFEST, Most EFFECTIVE Prehospital Field Tourniquet!



COMBAT APPLICATION TOURNIQUET

SPECIFICATIONS

Packaged: H 6.5 in. x W 2.4 in. x D 1.5 in.
Device: L 37.5 in.
Weight: 2.7 oz.

COLOR	STYLE	ITEM #	NSN #
Black	Tactical	30-0001	6515-01-521-7976
Orange	Rescue	30-0023	6515-01-612-7129
Blue	Trainer	30-0033	6910-01-560-2972



online video available

< scan to view



RIGID C-A-T® TQ CASE with Shirt Shield

- Allows carriage of a Gen 7 C-A-T® Tourniquet on your duty/patrol/pistol belt (Tourniquets sold separately)
- Shirt Shield protects uniform shirts from undue friction wear
- Blade-Tech TekLok Belt Attachment fits up to a 2.25 in. belt
- Can be mounted either vertically or horizontally, or to a Safariland Drop-leg Shroud
- Injection Molded Nylon Case stands up to rugged military usage
- Cross medical identifier patch included

RIGID C-A-T TQ CASE

COLOR	ITEM #
Black	80-0944
Black, Basket Weave	80-0962
Black, Polished	80-0963
Coyote	80-0945
Ranger Green	80-0946
Navy	80-1161

SPECIFICATIONS

H 5.88 in. x W 2.5 in. x D 2 in.
Weight: 5.36 oz

Explore other TQ carrier options available at NARescue.com

RIGID C-A-T® INSIDE-THE-WAISTBAND TOURNIQUET CASE



- Keeps your C-A-T® Tourniquet close at hand for rapid deployment
- Holds the C-A-T® Tourniquet (Gen 7) utilizing friction (Tourniquet sold separately)
- Designed with built-in Shirt Shield that creates a barrier between your body, shirt, and the exposed portion of your tourniquet
- Injection Molded Case to stand up to rugged LE/MIL usage
- Drainage hole at bottom of case
- Offers a range of adjustment
- Available Mounting Options – 1.5" and 1.75" IWB Overhook (included) designed, manufactured, and assembled in the USA.

RIGID CAT TQ CASE

COLOR	ITEM #
Black	80-1488
Gray	80-1489

SPECIFICATIONS

Packaged: L 5 in. x W 3.25 in. x D 1.625 in.
Weight: 3.0 oz

LEATHER C-A-T® HOLDERS



- Designed to meet the specific duty gear needs of law enforcement officers - for use on daily duty belt to fit within organizational uniform standards
- Securely fastens to your duty belt with easy one-handed access to allow you to take control of extreme medical situations in seconds
- Constructed of black full grain leather in different finishes, attachment options, and snaps colors

- High quality snaps ensure secure placement on your duty belt and easy accessibility to your tourniquet
- Ideal for the C-A-T® Gen 7 but also accommodates other popular windlass TQs
- Attachment options: Clip, Loop or Hook & Loop Strap attachment on select finishes and snap colors



◀ **with Strap Attachment**
with hook and loop allows direct attachment to a radio strap.



Rigid Case MOLLE ATTACHMENT

ITEM

ZZ-1696

SPECIFICATIONS

H 4 in. x W 1 in. x D 0.5 in.
Weight: 0.4 oz.



SPECIFICATIONS

L 5.75 in. x W 2 in. x D 2 in.
Weight: 3.2 oz.

UNIFORM DUTY BELT LEATHER TQ HOLDER

FINISH	ATTACHMENT	SNAP	ITEM #
High Gloss	Clip	Silver	30-0194
Full Grain	Clip	Silver	30-0193
Basket Weave	Loop	Silver	30-0195
Full Grain	Loop	Silver	30-0197
Basket Weave	Strap Attach	Silver	30-0228
High Gloss	Strap Attach	Silver	30-0229

Nylon C-A-T® HOLDER

- Tailored specifically for your C-A-T®, to protect it from exposure to the elements
- PALS/MOLLE-style connectors for versatile attachment to the duty vest, gear or belt (mounted horizontally or vertically) for rapid & easy access
- Easy-open pull elastic tab with tourniquet ID patch - "TQ"



NYLON TQ HOLDER

COLOR	ITEM #
Black	30-0057
Coyote	30-0044
Multi-Cam	30-0059
Red	30-0128
Ranger Green	30-0208
RealTree®	30-0207
OD Green	30-0053

SPECIFICATIONS

L 6.25 in. x W 1.75 in. x D 1.5 in.
Weight: 2.17 oz

LASER CUT C-A-T® HOLDER

- Constructed with laser cut laminated Cordura® nylon
- Easy Pull Elastic Cover
- PALS/MOLLE-style connections
- Accommodates GEN 6 and GEN 7 C-A-T® Tourniquets
- Weighs less than traditional TQ holders

SPECIFICATIONS

L 7.25 in. x W 2.75 in. x D 1.5 in.
Weight: 1.29 oz

LASER CUT TQ HOLDER

COLOR	ITEM #
Black	30-0297
Coyote	30-0298
Ranger Green	30-0299
OCP	30-0306
Red	30-0551



IFAK Individual First Aid Kits



Individual Patrol Officer Kit **IPOK**

- Individual hemorrhage control kit
- Compact, vacuum-sealed & packaged to fit into a pocket, vest pouch or individual bag
- Rugged, durable pouch features NAR's signature Red Tip Technology® ensuring quality & performance

SPECIFICATIONS

H 4 in. x W 6.5 in. x D 2.75 in.
Weight: 7.4 oz



IPOK:



IPOK - INDIVIDUAL PATROL OFFICER KIT

HEMOSTATIC/GAUZE	ITEM#
Wound Packing Gauze	80-0167
Combat Gauze, LE	80-0168

ETAK

Enhanced Trauma Aid Kit

Designed to pair with a separate C-A-T® Tourniquet & holder (not included) for a great option to add more capability to previously owned kit

- Compact IFAK platform accessible from a belt, duty belt, assault vest/body armor or any equipment with PALS/MOLLE attachment points.
- Featuring multiple interior elastic loops keep your kit organized and secure
- Exterior tear away panel is easily removable for added mobility and access to your kit
- MED patch included

Every ETAK includes these core items:



Bag only option also available at NARescue.com


PLUS+

BASIC also includes:



with **COMBAT GAUZE** also includes:

ADVANCED also includes:



ETAK ENHANCED TRAUMA AID KIT

COLOR	LEVEL	ITEM #
● Black	Basic	80-1168
● Coyote	Basic	80-1169
● Ranger Green	Basic	80-1171
● MultiCam	Basic	80-1170
● Black	Combat Gauze LE	80-1172
● Coyote	Combat Gauze LE	80-1173
● Ranger Green	Combat Gauze LE	80-1175
● MultiCam	Combat Gauze LE	80-1174
● Black	Advanced	80-1176
● Coyote	Advanced	80-1177
● Ranger Green	Advanced	80-1179
● MultiCam	Advanced	80-1178

SPECIFICATIONS

H 5.75 in. x W 4 in. x D 1.75 in.
Weight: 0.65 lb

IPOK HOLDER

- Protects your medical kit from the elements
- Enclosure can be secured by Fastex snap buckle, hook-and-loop or both



IPOK HOLDER

COLOR	ITEM#
● Black	80-0205T
● Coyote	80-0207
● OD Green	80-0254T
● Red	80-0358T

SPECIFICATIONS

H 6.5 in. x W 6 in. x D 3.5 in.
Weight: 4.1 oz

VISOR MED POUCH

- Securely mounts to a vehicle visor
- Orange webbing handle to grab for easy deployment
- Pouch with easy-open buckle



VISOR MED POUCH

COLOR	ITEM#
● Black	80-0269
● Coyote	80-0270
● Red	80-0360

SPECIFICATIONS

H 6.5 in. x W 5.5 in. x D 3 in.
Weight: 5.4 oz

HEADREST IFAK

- Designed to provide personnel with an easily accessible, quick to deploy medical pouch that conveniently mounts to a vehicle's headrest
- Fully loaded with emergency medical equipment for point-of-wounding treatment until definitive care is available
- Durable nylon construction features multiple interior elastic loops allowing for the contents to be organized and secured
- Easy to grab side handles allow for easy deployment
- IFAK pouch only is also available
- The pouch design allows for the contents of an individual first aid kit or vacuum sealed kit

SPECIFICATIONS

H 5.25 in. x
W 7.25 in. x D 4.25 in.
Weight: 1.63 lb



HEADREST IFAK

BASIC QTY 2 Gauze, Wound Packing	-OR- BASIC with Combat Gauze QTY 1 Gauze, Wound Packing	QTY 2 Combat Gauze® LE	QTY 2 C-A-T® Tourniquet	QTY 2 ETD® 4 in. Mini Responder
---	---	---------------------------	-------------------------------	------------------------------------

QTY 2 HyFin® Vent Compact, Twin	QTY 1 MicroShield CPR Face Mask	QTY 1 Survival Blanket	QTY 20 1 in. x 3 in. Bandages	QTY 1 Triangular Bandage	PAIR 2 Responder Gloves Black, Lg
---------------------------------------	---------------------------------------	------------------------------	-------------------------------------	--------------------------------	---

HEADREST IFAK


COLOR	OPTIONS	ITEM #
● Black	Basic	80-1091
● Coyote	Basic	80-1150
● Black	Combat Gauze LE	80-1151
● Coyote	Combat Gauze LE	80-1152
● Black	Bag Only	80-1090T
● Coyote	Bag Only	80-1149T

M-FAK LCL *Mini First Aid Kit - Laser Cut Laminate*

- **NEW ITEM:** NAR's most compact, versatile Mini First Aid Kit - featuring Laser Cut Laminate MOLLE
- First Responder's most requested critical point-of-wounding medical equipment for treating penetrating, blast or other traumatic injuries in the line of duty
- The LCL M-FAK is double-spaced for the MOLLE attachment points on the rear of the bag for maximum adjustability and super secure mounting to any platform
- Half the size of our popular T.O.R.K platform



▲ NEW!
LCL M-FAK Bag features laser cut MOLLE on front, side & rear

SPECIFICATIONS 
H 6 in. x W 3 in. x D 3 in.
Weight: 13 oz

M-FAK LCL

COLOR	LEVEL	ITEM #
● Black	Basic	80-1667
● Coyote	Basic	80-1668
● OCP	Basic	80-1669
● Red	Basic	80-1671
● RGR Green	Basic	80-1672
● Black	Combat Gauze LE	80-1673
● Coyote	Combat Gauze LE	80-1674
● OCP	Combat Gauze LE	80-1675
● Red	Combat Gauze LE	80-1677
● RGR Green	Combat Gauze LE	80-1678
● Black	Advanced	80-1379
● Coyote	Advanced	80-1680
● OCP	Advanced	80-1681 <small>Rx ONLY</small>
● Red	Advanced	80-1683
● RGR Green	Advanced	80-1684



Bag only option also available at NARescue.com


M-FAK LCL BASIC:



M-FAK LCL with COMBAT GAUZE:




M-FAK LCL ADVANCED:



ROO® M-FAK LCL *Mini First Aid Kit - Laser Cut Laminate*

- **NEW ITEM:** Super compact, rugged nylon platform featuring **Laser Cut Laminate MOLLE**
- Essential point-of-wounding equipment needed for treating penetrating, blast or other traumatic injuries
- Unique front pouch that easily opens with the pull of a tab for **rapid one-handed access to your tourniquet**
- The LCL ROO® M-FAK is double-spaced for the MOLLE attachment points on the rear of the bag.



SPECIFICATIONS 
 H 7 in. x W 3.75 in. x D 4 in.
 Weight: 13 oz

ROO® M-FAK LCL

COLOR	LEVEL	ITEM #
● Black	Basic	80-1703
● Coyote	Basic	80-1704
● OCP	Basic	80-1705
● Red	Basic	80-1707
● Ranger Green	Basic	80-1708
● Black	Combat Gauze LE	80-1709
● Coyote	Combat Gauze LE	80-1710
● OCP	Combat Gauze LE	80-1711
● Red	Combat Gauze LE	80-1713
● RGR Green	Combat Gauze LE	80-1714
● Black	Advanced	80-1715
● Coyote	Advanced	80-1716
● OCP	Advanced	80-1717
● Red	Advanced	80-1719
● RGR Green	Advanced	80-1720

Rx ONLY


ROO® M-FAK LCL BASIC:



ROO® M-FAK LCL with COMBAT GAUZE:



ROO® M-FAK LCL ADVANCED:



Bag only option also available at NARescue.com

▲ **NEW!** ROO® bag includes laser cut MED ID Cross and double-spaced LCL MOLLE

IFAK Individual First Aid Kits

M-FAK Mini First Aid Kit - LE

- Contains the essential point-of-wounding equipment needed for treating penetrating, blast and other traumatic injuries
- Super compact, rugged nylon platform that allows attachment both vertically (MOLLE backing) and horizontally (3 in. belt loop)
- Clam shell configuration utilizes (2) two main sleeves that open on both ends for easy access

Bag only option also available at NARescue.com

SPECS


H 6 in. x W 3 in. x D 3 in.
Weight: 13 oz



LE M-FAK | ROO® M-FAK BASIC or BASIC with Combat Gauze:



LE M-FAK | ROO® M-FAK ADVANCED:



M-FAK MINI FIRST AID KIT - LE

COLOR	LEVEL	ITEM #
Black	Basic	80-0494
Coyote	Basic	80-0495
Gray	Basic	80-0872
Multi-Cam	Basic	80-0496
OD Green	Basic	80-0497
Ranger Green	Basic	80-0996
Black	Combat Gauze LE	80-0498
Coyote	Combat Gauze LE	80-0499
Gray	Combat Gauze LE	80-0873
Multi-Cam	Combat Gauze LE	80-0500
OD Green	Combat Gauze LE	80-0501
Ranger Green	Combat Gauze LE	80-0998
Black	Advanced	80-0502
Coyote	Advanced	80-0503
Gray	Advanced	80-0874
Multi-Cam	Advanced	80-0504
OD Green	Advanced	80-0505
Ranger Green	Advanced	80-0997



Compact Officer Response Emergency Kit CORE KIT

Designed to go where you go, this compact medical kit will fit in a BDU pocket, vest pouch or Go-Bag and provides personnel with essential life-saving bleeding control equipment that can be kept on their person. The CORE Kit is packaged in a compact, rugged, clear vacuum-sealed bag with easy-to-open tear notches.

- Compact design with the essential life-saving equipment needed to address moderate and severe traumatic bleeding
- Contains the #1 CoTCCC Recommended Prehospital Tourniquet and Hemostatic Gauze
- The Best Tourniquet paired with the Best Hemostatic Gauze and an effective pressure bandage at the Best Price
- Meets the Tactical Emergency Casualty Care (TECC) Guidelines for hemorrhage control

CORE KIT:



SPECS

H 7 in. x W 3.5 in. x D 2.25 in. | Weight: 6.4 oz

CORE KIT

ITEM #
85-0834

ROO® M-FAK Mini First Aid Kit

Rapid one handed access to your tourniquet without ever having to open the full kit!

Contains the same packout options as the M-FAK:
Basic, Basic with Combat Gauze® LE and Advanced

ROO M-FAK MINI FIRST AID KIT

COLOR	LEVEL	ITEM #
● Black	Basic	80-1049
● Coyote	Basic	80-1050
● OD Green	Basic	80-1051
● Ranger Green	Basic	80-1052
● Black	Combat Gauze LE	80-1053
● Coyote	Combat Gauze LE	80-1054
● OD Green	Combat Gauze LE	80-1055
● Ranger Green	Combat Gauze LE	80-1056
● Black	Advanced	80-1057
● Coyote	Advanced	80-1058
● OD Green	Advanced	80-1059
● Ranger Green	Advanced	80-1060



NEW!
ROO® bag includes Mini 1 in. PVC ID Patch

Bag only option also available at NARescue.com

L-FAK Lumbar First Aid Kit

- Durable 500-denier exterior carrier
- Compact design to reduce space needed for life-saving equipment
- Designed for wear in the small of the back
- Multiple interior elastic loops for securing gear
- Rapid access to medical gear
- C-A-T® tourniquet sold separately



Every L-FAK includes these core items:




SPECIFICATIONS

Pkgd - H 4.75 in. x W 5.75 in.
Weight: 0.88 lbs
Adv. Weight: 0.94 lbs



PLUS+ BASIC
also includes:



COMBAT GAUZE®
also includes:



ADVANCED also includes:



LUMBAR FIRST AID KIT

COLOR	LEVEL	ITEM #
● Black	Basic	85-3627
● Coyote	Basic	85-3633
● Ranger	Basic	85-3630
● Black	w/ Combat Gauze™	85-3628
● Coyote	w/ Combat Gauze™	85-3634
● Ranger	w/ Combat Gauze™	85-3631
● Black	Advanced	85-3629
● Coyote	Advanced	85-3635
● Ranger	Advanced	85-3632

Bag only option also available at NARescue.com

IFAK Individual First Aid Kits



EDC ANKLE TRAUMA KIT

OPTIONS	ITEM #
Wound Packing	80-1157
BCD	80-1158
Carrier Only	30-0294T

EDC ANKLE TRAUMA KIT

- Lightweight, low profile, versatile carrier for essential trauma and EDC equipment with quick, single-handed, ambidextrous access
- Separate concealed compartment for essential EDC items or ID
- Designed to accommodate a C-A-T® Tourniquet and gauze/bandage
- Advanced materials that provide user comfort and long-term wearability
- Individual hook and loop pouches for each carried component
- One size fits most with convenient hook and loop adjustability
- Universal magazine pocket

SPECIFICATIONS

L 17.75 in. x H 4 in. x D 0.25 in.
Weight: 0.657 lb

EVERY DAY CARRY ANKLE TRAUMA KIT



REFLEX™ SYSTEM ROLL & CARRIER

- Primary construction of heavy-duty nylon laminate allows rapid deployment of medical supplies
- Features High Abrasion-Resistant Neoprene Kevlar® (HANK™) composite grab handles for long-term durability
- Can mount vertically or horizontally on belts 1.5 in.- 2.25 in., and horizontally on MOLLE
- Carry handle and shock cord loops allow the roll to be attached to a patient, or hung from an I.V. rack

- "Mini" MOLLE allows for precise fit on belts as narrow as 1.5 in. or two-row MOLLE panels
- Features loop panels on med roll and carrier for label patches (Medical Cross Patch included)
- Two-piece system designed for ambidextrous access
- Medical contents sold separately

REFLEX SYSTEM ROLL & CARRIER

COLOR	ITEM #
Black	80-1117
LE Blue	80-1118
Coyote	80-1119
Gray	80-1120
Multi-Cam Blk	80-1121
Multi-Cam	80-1122
OD Green	80-1123



SPECIFICATIONS

L 8 in. x W 3 in. x H 5.5 in.
Weight: .58 lb



Carrier

Roll

REFLEX SYSTEM IFAK

COLOR	LEVEL	WGT	ITEM #
● Black	Basic	1.58 lb	80-1127
● LE Blue	Basic	1.58 lb	80-1128
● Coyote	Basic	1.58 lb	80-1129
● Gray	Basic	1.58 lb	80-1130
● Multi-Cam Black	Basic	1.58 lb	80-1131
● Multi-Cam	Basic	1.58 lb	80-1132
● OD Green	Basic	1.58 lb	80-1133
● Black	Combat Gauze	1.6 lb	80-1134
● LE Blue	Combat Gauze	1.6 lb	80-1135
● Coyote	Combat Gauze	1.6 lb	80-1136
● Gray	Combat Gauze	1.6 lb	80-1137
● Multi-Cam Black	Combat Gauze	1.6 lb	80-1138
● Multi-Cam	Combat Gauze	1.6 lb	80-1139
● OD Green	Combat Gauze	1.6 lb	80-1140
● Black	Advanced	1.61 lb	80-1141
● LE Blue	Advanced	1.61 lb	80-1142
● Coyote	Advanced	1.61 lb	80-1143
● Gray	Advanced	1.61 lb	80-1144
● Multi-Cam Black	Advanced	1.61 lb	80-1145
● Multi-Cam	Advanced	1.61 lb	80-1146
● OD Green	Advanced	1.61 lb	80-1147

REFLEX™ IFAK SYSTEM

- Developed with direct input from active-duty medical personnel, the ReFlex™ System IFAK comes fully loaded with essential medical equipment to treat life threatening injuries in the field of operation
- Organized to put all medical supplies within easy reach
- Heavy-duty nylon laminate allows for rapid deployment
- HANK™ composite grab handles for long-term durability
- Carry handle and shock cord loops allow the roll to be attached to a patient, or hung from an I.V. rack
- "Mini" MOLLE allows for precise vertical/horizontal mount on belts 1.5 in. - 2.25 in.
- Mounts horizontally on MOLLE



SPECIFICATIONS

L 8 in. x W 3 in. x H 5.5 in. | Weight: 1.6 lb

Every REFLEX SYSTEM IFAK includes these core items:

 QTY 1 C-A-T® Tourniquet	 QTY 1 ETD® 4 in. Flat	 QTY 1 Wound Packing Gauze	 QTY 1 HyFin® Vent, Twin Pk	 QTY 1 Survival Blanket
 QTY 1 NAR® Flat Tape, 2"	 QTY 1 Trauma Shears, Lg	 PAIR 2 Responder Black Gloves, LG	 QTY 2 Combat Casualty Card	 QTY 1 Permanent Marker, Large

PLUS+ BASIC

also includes:

 QTY 1 Wound Packing Gauze

COMBAT GAUZE®

also includes:

 QTY 1 Combat Gauze® Z-fold

ADVANCED also includes:

 QTY 1 Combat Gauze® Z-fold	 QTY 1 Pre-lubed NAR® Naso, 28F	 QTY 1 ARS® Needle, 10G
-----------------------------------	---------------------------------------	-------------------------------

ReFlex™ System VEHICLE MOUNT

- Attaches ReFlex™ Med Roll to the back of a vehicle headrest
- All-nylon construction with built-in elastic tourniquet loop



REFLEX MOUNT

COLOR	ITEM #
● Coyote	ZZ-1369
● Gray	ZZ-1370

SPECIFICATIONS

L 5.5 in. x W 0.125 in. x H 6 in.
Weight: 0.21 lb



IFAK Individual First Aid Kits

TORK® LCL Tactical Operator Response Kit (Laser Cut Laminate)

The NEW Tactical Operator Response Kit (TORK® LCL) is designed to meet the need for a compact individual first aid kit to address leading causes of potentially preventable death in a tactical environment. Now featuring Laser Cut Laminate (LCL) MOLLE.

- Essential life-saving tools to address injuries caused by penetrating or blast trauma
- Compact Individual First Aid Kit (IFAK) providing Basic Life Support for immediate Self-Aid/Buddy-Aid
- Includes the added advantage of a hemostatic dressing for extra capability in the control of life threatening bleeding
- The LCL TORK® is double-spaced for the MOLLE attachment points on the rear of the bag for maximum adjustability and super secure mounting to any platform
- Also attaches to NAR's Quick Detach Panel (sold separately)
- Signature Red-Tip Technology® rapid opener for fast and easy, one-pull tab opening of the kit.

NEW!
TORK® LCL bag features double-spaced Laser Cut Laminate MOLLE

TORK LCL

COLOR	LEVEL/HEMOSTATIC	ITEM #
Black	Basic	80-1685
Coyote	Basic	80-1686
OCP	Basic	80-1687
Red	Basic	80-1689
RGR Green	Basic	80-1690
Black	Basic w/Combat Gauze	80-1691
Coyote	Basic w/Combat Gauze	80-1692
OCP	Basic w/Combat Gauze	80-1693
Red	Basic w/Combat Gauze	80-1695
RGR Green	Basic w/Combat Gauze	80-1696
Black	Advanced w/Combat Gauze	80-1697
Coyote	Advanced w/Combat Gauze	80-1698
OCP	Advanced w/Combat Gauze	80-1699
Red	Advanced w/Combat Gauze	80-1701
RGR Green	Advanced w/Combat Gauze	80-1702


SPECIFICATIONS

H 7 in. x W 4.5 in. x D 2.75 in.
Weight: 1 lb. 5 oz.

TORK LCL BASIC / BASIC with Combat Gauze



TORK LCL ADVANCED



Individual First Aid Kits IFAK

TORK TACTICAL OPERATOR RESPONSE KIT

COLOR	LEVEL/HEMOSTATIC	ITEM #
Black	Basic	85-0008
Coyote	Basic	85-0009
OD Green	Basic	85-0054
Multi-Cam	Basic	85-0012
Ranger Green	Basic	85-2301
Red	Basic	85-1805
MTC Alpine	Basic	80-1633
Black	Basic w/Combat Gauze	80-1097
Coyote	Basic w/Combat Gauze	80-1098
OD Green	Basic w/Combat Gauze	80-1101
Multi-Cam	Basic w/Combat Gauze	80-1100
Ranger Green	Basic w/Combat Gauze	80-1102
Red	Basic w/Combat Gauze	80-1108
MTC Alpine	Basic w/Combat Gauze	80-1635
Black	Advanced w/Combat Gauze	85-0048
Coyote	Advanced w/Combat Gauze	85-0049
OD Green	Advanced w/Combat Gauze	85-0053
Multi-Cam	Advanced w/Combat Gauze	85-0052
Ranger Green	Advanced w/Combat Gauze	85-2302
Red	Advanced w/Combat Gauze	85-1806
MTC Alpine	Advanced w/Combat Gauze	80-1636
7 options online	Advanced w/ChitoGauze	NARescue.com

TORK® Tactical Operator Response Kit

- Streamlined, compact IFAK designed to MOLLE to the assault vest or body armor
- Includes essential life-saving equipment to be used at the point of wounding
- For self-aid/buddy-aid injuries caused by penetrating or blast trauma

TORK ReSUPPLY KITS*

CONTENT OPTIONS	ITEM #
Basic	80-1037
Basic Combat Gauze	80-1103
Advanced Combat Gauze	80-1039
Advanced ChitoGauze	80-1038



*ReSupply Kits include contents only. Nylon bag is sold separately.

TORK BASIC/BASIC with Combat Gauze®:



SPECIFICATIONS

H 7 in. x W 4.5 in. x D 2.75 in.
Weight: 1 lb 5 oz

TORK ADVANCED:



Bag only options also available at NARescue.com

IFAK Individual First Aid Kits

RIG SERIES EAGLE IFAK

- Tourniquet attached for quick deployment by banding it to inside of front opening flap
- Protective hood with quick pull latch release helps protect IFAK components and includes sewn in embroidered medical cross identifier
- Detachable (hook & loop backed) shock-cord inner panel adjusts to secure internal components and is easily removable for rapid access



RIG SERIES EAGLE IFAK

COLOR	LEVEL	ITEM #
● Black	Basic	80-0509
● Coyote	Basic	80-0510
● Multi-Cam	Basic	80-0511
● Ranger Green	Basic	80-0991
● Black	Combat Gauze	80-0517
● Coyote	Combat Gauze	80-0518
● Multi-Cam	Combat Gauze	80-0519
● Ranger Green	Combat Gauze	80-0993
● Black	Advanced	80-0513
● Coyote	Advanced	80-0514
● Multi-Cam	Advanced	80-0515
● Ranger Green	Advanced	80-0992

RX ONLY

Bag only option also available at NARescue.com

RIG SERIES EAGLE BASIC:



RIG SERIES EAGLE with COMBAT GAUZE LE:



SPECIFICATIONS

H 8 in. x W 4 in. x D 3 in.
 Kit Packed Weight: 1.18 lb
 Empty Bag Weight: 7 oz

RIG SERIES EAGLE ADVANCED:



Individual First Aid Kits IFAK

Provides several options for carrying essential life-saving equipment



RAID KIT V2

COLOR	ITEM#	OPTION	NSN#
BLK	80-1545	Basic	----
COY	80-1546	Basic	----
MTC	80-1547	Basic	----
RGR	80-1548	Basic	----
BLK	80-1549	Basic w/ CG	----
COY	80-1550	Basic w/ CG	----
MTC	80-1551	Basic w/ CG	----
RGR	80-1552	Basic w/ CG	----
BLK	80-1625	Basic w/ CELOX™	----
COY	80-1626	Basic w/ CELOX™	----
MTC	80-1627	Basic w/ CELOX™	----
RGR	80-1628	Basic w/ CELOX™	----
BLK	80-1553	Advanced	----
COY	80-1554	Advanced	----
MTC	80-1555	Advanced	----
RGR	80-1556	Advanced	----
BLK	80-1629	Adv. w/ CELOX™	----
COY	80-1630	Adv. w/ CELOX™	----
MTC	80-1631	Adv. w/ CELOX™	----
RGR	80-1632	Adv. w/ CELOX™	----

Rx only

RAID KIT V2

Releasable Abdominally Integrated Drop kit

Every RAID Kit V2 includes these core items:

 QTY 1 C-A-T® Tourniquet	 QTY 1 NAR® Gauze S-Rolled	 QTY 1 ETD® 4 in. Dressing	 QTY 1 HyFin® Vent Compact, Twin Pk	 PAIR 1 Nitrile Gloves, Lg	 QTY 1 NAR Trauma Shears, Lg
-----------------------------------	-------------------------------------	-------------------------------------	--	-------------------------------------	---------------------------------------

- RAID V2 pouch's unique take on the industry standard hanger combined with the utility of a fanny pack allows user to decide how to wear it based on their needs during self-aid/buddy-aid
- Designed to meet the need for an individual first aid kit with multiple carry options
- First Spear Tube makes reattaching the pouch as easy as removal

PLUS+
BASIC also includes: OR BASIC with COMBAT GAUZE® also includes: OR BASIC with CELOX™ GAUZE™ also includes: OR ADV. with COMBAT GAUZE® also includes: OR ADV. with CELOX™ GAUZE™ AND

 QTY 1 NAR® Gauze S-Rolled	 QTY 1 Combat Gauze® Z-Fold	 QTY 1 CELOX™ Gauze Z-fold	 QTY 1 Combat Gauze® Z-Fold	 QTY 1 CELOX™ Gauze Z-fold	 QTY 1 Pre-lubed NAR® Naso, 28F	 QTY 1 10G ARS® Needle Decomp.
-------------------------------------	--------------------------------------	-------------------------------------	--------------------------------------	-------------------------------------	--	---

SPECIFICATIONS

H 8 in. x W 7 in. x D 2 in.
Weight: 3 lbs 3.2 oz.



MINI MEDIC KIT

- Compact, multi-purpose trauma aid bag with vital medical equipment for point of wounding treatment of multiple injuries
- Two external pockets with two external tourniquet holders
- One external pouch with additional drop-down pouch for use with breaching tools or storage for a hydration pack
- Expansion zipper to increase load capacity
- External MOLLE on sides and bottom of bag / Internal MOLLE rows on inside back of bag
- Three clear tear away pouches on the inside main compartment
- Two internal mesh pouches and one larger internal pouch on inner lid
- Internal tape roll strap



Inside Back Compartment for storing extra equipment

Every MINI MEDIC KIT includes these core items:

QTY 2 C-A-T® TQ Holder	QTY 4 C-A-T® Tourniquet	QTY 4 Gauze S-Rolled	QTY 4 ETD™ 4 in. Flat
QTY 2 HyFin® Vent. Twin Pk	QTY 2 BurnTec 4"x4" Dressing	QTY 2 SAM® Splint II	QTY 1 NAR QuikLitter LITE
QTY 2 Surgical Tape 2" Flatfold	QTY 1 Mini Tape Multipurpose	QTY 1 Trauma Shears, Lg	QTY 1 Scissor Leash™
QTY 4 Combat Casualty Card	QTY 2 Permanent Marker, Large	QTY 1 Casualty Reference Card	QTY 1 Medical Cross Patch

Bag only option also available at NARescue.com

Breaching Equipment Sold Separately

SPECIFICATIONS

H 14 in. x
W 11 in. x D 5 in.
(Depth expandable to 7.25 in.)
Weight: 12.5 lb
13 lb (adv.)

MINI MEDIC KIT

COLOR	LEVEL	ITEM #
● Black	Basic	80-0979
● Coyote	Basic	80-0980
● OD Green	Basic	80-0981
● Ranger Green	Basic	80-1003
● MultiCam	Basic	80-1231
● Black	Combat Gauze	80-0985
● Coyote	Combat Gauze	80-0986
● OD Green	Combat Gauze	80-0987
● Ranger Green	Combat Gauze	80-1005
● MultiCam	Combat Gauze	80-1232
● Black	Advanced	80-0982
● Coyote	Advanced	80-0983
● OD Green	Advanced	80-0984
● Ranger Green	Advanced	80-1004
● MultiCam	Advanced	80-1233

PLUS+ BASIC also includes:

QTY 4 Gauze S-Rolled

COMBAT GAUZE also includes:

QTY 4 Combat Gauze Z-fold

QTY 2 Pre-lubed NAR® Naso, 28F

QTY 1 BAC-Pack Bougie Aided Cric

ADVANCED also includes:

QTY 2 ARS® Needle, 14g

QTY 4 Combat Gauze Z-fold

QTY 1 Cyclone Pocket BVM

Rx ONLY

Patrol **VEHICLE TRAUMA KIT**

- Provide personnel with an easily accessible, quick to deploy medical kit from your Patrol Vehicle
- Modular storage system comes organized with an array of medical products ranging from basic first aid to more intermediate trauma care
- One side covers all your agency needs to administer basic first aid and community assistance
- Second side offers a first time combination of Basic First Aid and Basic Trauma medical items
- Exterior includes MOLLE-style system to accommodate extra equipment as needed



PATROL VEHICLE TRAUMA KIT:

 QTY 1 C-A-T® Tourniquet	 QTY 1 ETD® 6in. Dressing	 QTY 1 ETD® Abd./Stump	 QTY 2 Responder Compressed Gauze	 QTY 1 Responder Compressed Gauze	 QTY 1 Combat Gauze® LE	
 QTY 1 CPR Pocket Resuscitator	 QTY 1 HyFin® Vent Chest Seal, Twin Pk	 QTY 1 Survival Blanket	 QTY 1 BurnTec® 4"x4" Dressing	 QTY 1 SAM® Splint II	 QTY 2 Triangular Bandage	 QTY 4 4" x 4" Gauze Pad 2 PK
 QTY 2 4" Elastic Wrap Bandage	 QTY 1 Eye Wash Solution, 4oz.	 QTY 2 PES Eye Shield	 QTY 1 Surgical Tape, 2 in. Mini	 QTY 1 NAR Trauma Shears, Lg	 PAIR 3 Nitrile Gloves, Lg	 QTY 25 Flexible Fabric Bandages
 QTY 1 Splinter Forceps, 4.5"	 QTY 1 Pocket Penlight	 QTY 1 Permanent Marker, Large				

PATROL VEHICLE TRAUMA

LEVEL	ITEM #
Basic	80-0542
Combat Gauze LE	80-0859

SPECIFICATIONS

H 7 in. x W 12 in. x D 5 in.
Weight: 6 lb



Trauma Kits

RANGE TRAUMA AID KIT *Hard Case*

- Water-tight case maintains the integrity of your medical products when you need them most
- Kit contents arranged in a modular system featuring rugged lid storage organizer and removable nylon pouches
- Labeled internal pouches for ease of use



RANGE TRAUMA AID

ITEM

85-0889

SPECIFICATIONS

H 9.5 in. x W 13.5 in. x D 6 in.
Weight: 8 lb 12 oz

RANGE TRAUMA AID KIT HARD CASE:



School Resource Officer Kits

SRO CRISIS RESPONSE KIT

- Equip School Resource Officers with the essential life-saving tools to begin treating injured victims until medical personnel arrive
- Empowers an SRO (School Resource Officer) to administer self-aid in the event he or she is injured or provide immediate care to injured casualties during a school incident
- Multiple carry options

SRO CRISIS RESPONSE KIT:



Bag only option
also available at
NARescue.com

SPECIFICATIONS

H 7.75 in. x W 13 in. x D 9 in.
Weight: 5 lb 16 oz

SRO CRISIS RESPONSE KIT

ITEM

85-0417

I-ROK *Individual Rugged Outdoor Kit*



- Contains essential first-aid supplies to treat basic injuries as well as control major bleeding
- Storage case is crushproof, dustproof and waterproof up to a depth of 2 meters
- Flexible lid organizer helps keep smaller medical supplies neatly stored
- Built-in pressure valve keeps water and dust out while balancing air pressure
- Available in 3 hardcase colors and 2 sizes with the option of Basic or with Bleeding Control Dressing (hemostatic bandage) versions

Every I-ROK contains the following items:



SPECIFICATIONS

H 7 in. x W 4.5 in.
x D 2.5 in.
Weight: 2.48 lb

I-ROK

COLOR	LEVEL	ITEM #
Black	Basic	82-0031
OD Green	Basic	82-0033
Orange	Basic	82-0034
Black	w/BCD	82-0035
OD Green	w/BCD	82-0037
Orange	w/BCD	82-0038
Black	w/CELOX	82-0066
OD Green	w/CELOX	82-0068
Orange	w/CELOX	82-0069



I-ROK XL

SPECIFICATIONS

H 8.9 in. x W 5.4 in. x D 3 in.
Weight: 3 lb

I-ROK XL includes these extra items:



I-ROK XL

COLOR	LEVEL	ITEM #
Black	Basic	82-0049
OD Green	Basic	82-0051
Orange	Basic	82-0052
Black	w/BCD	82-0053
OD Green	w/BCD	82-0055
Orange	w/BCD	82-0056
Black	w/CELOX	82-0070
OD Green	w/CELOX	82-0072
Orange	w/CELOX	82-0073

Multiple Injury Trauma Kits

MCI-TKD *Mass Casualty Incident Throw Kit Dispenser*

- Contains an inner main compartment specifically designed for quick deployment of up to eight (8) point of injury kits
- One-handed quick release, silent zippers
- G-Hook backpack straps can switch from 2 shoulder straps to bandolier style.
- Easy backpack carry for mass casualty events
- Quickly deploy bleeding control kits in emergency situations
- Side pockets have C-A-T® Tourniquets and Emergency Trauma Dressings (ETD®) for maximum preparedness for multiple injuries
- Front pockets for Litter, Breathing, and wound packing scenarios



Bag only option also available at NARescue.com

MCI-TKD

COLOR	LEVEL	ITEM #
Black	Basic	80-1749
Red	Basic	80-1750
Black	w/ BCD	80-1751
Red	w/ BCD	80-1752
Black	w/ CELOX™	80-1753
Red	w/ CELOX™	80-1754

Every MCI-TKD includes these core items:

 QTY 8 ORG or BLK C-A-T® Tourniquet	 QTY 8 NAR® Gauze, Wound Packing	 QTY 8 ETD® 4 in. Responder, Flat
 QTY 4 HyFin® Vent Compact, Twin Pk	 QTY 1 NAR® QuikLitter LITE	 QTY 2 Responder Shears, Lg

SPECIFICATIONS

H 21 in. x W 13 in. x D 8 in.
Weight: 15 lbs.



PLUS+

BASIC
also includes:



w/BCD
also includes:



w/CELOX™
also includes:



Multiple Injury Trauma Kits



MCI-TKD Bleeding Control Throw Kit includes:

 QTY 1 C-A-T® Tourniquet	 QTY 1 ETD® 4 in. Responder, Flat	 QTY 1 HyFin® Vent Compact, Twin Pk	 QTY 1 Survival Blanket	 PAIR 1 Responder Blue Gloves, LG	 QTY 1 Instructions Bleeding Control
--	--	--	--	--	--

PLUS+

BASIC
also includes:



QTY 1
NAR® Gauze,
Wound Packing

w/BCD
also includes:



QTY 1
QuikClot®
BCD

w/CELOX™
also includes:



QTY 1
5' CELOX™
Hemostatic



Rescue Tools & Litters



Leatherman® RAPTOR® MULTI-TOOL

High-performance multi-tool equipped with a variety of mission critical tools essential to emergency medical practitioners in the field

1. 420HC Stainless Steel Folding Medical Shears
2. Strap Cutter
3. Oxygen Tank Wrench
4. Ring Cutter
5. Carbide Glass Breaker
6. 5 cm Ruler
7. Holster with Clip

RAPTOR MULTI-TOOL

HANDLE	ITEM #	NSN #
● Black	70-0014	5110-01-626-8766
● Coyote	70-0265	----
● Blu/Blk	70-0264	----
● Red/Blk	70-0266	----
● Blu/Ylw	70-0326	----

SPECS


Shear Length: 1.9 in.
Closed Length: 5 in.
Weight: 5.8 oz

Leatherman® RAPTOR® RESPONSE SHEARS

1. Stainless Steel folding shears
2. Integrated ring cutter
3. 5 cm ruler
4. Oxygen tank wrench
5. Lanyard connection in handle
6. Removable pocket clip

RAPTOR RESPONSE SHEARS

HANDLE	ITEM #
● Cement	70-0328
● Crimson	70-0329
● Navy	70-0330



SPECS

Length: Closed 4.3 in.
Open 7.25 in.
Shear Length: 1.9 in.
Width: 1.7 in | 4.3 cm
Overall Thickness: .85 in.
Weight: 5.5 oz

Benchmade® 7 HOOK PERSONAL SAFETY CUTTER

- Saves time – more than 1.5 minutes faster than traditional trauma shears when cutting 2 articles of clothing.
- Extensively tested and widely used by the US military for medical and egress purposes
- Maintains sharpness for hundreds of cuts (up to 500+ cuts through seatbelt test material)
- 440C cutting surface and opening with handle to provide a more compact carry
- Includes a MOLLE-compatible carrying sheath



BENCHMADE

HOOK SAFETY CUTTER

ITEM #	STYLE
ZZ-0652	7 Hook

SPECS

Length: 4.3 in.
Weight: 1.5 oz

NAR® QUIKLITTER LITE

- Perfect balance of size, weight, cost, and capability
- Features ten large carry handles with a load capacity tested to 750 lb
- Made of lightweight, heavy-duty nylon
- Vacuum-packed to conserve space and fit in a cargo pocket
- Latex-free & water-resistant



QUIKLITTER

SIZE	ITEM #
LITE*	60-0071

*Nylon Construction

SPECIFICATIONS

Payload 750 lb.
Litter Bed: L 74 in. x W 27 in.
Weight: 12 oz

UNO EMERGENCY SPLINT

- Modular Design
- Stepless Adjustment
- Small and Lightweight
- Compatible with TCCC Equipment
- Pack of 3 includes additional ankle strap for stabilizing the ankle and enables the device to provide traction



SPECIFICATIONS

Packaged - L 21.50 in. x W 10.55 in. x D 1.50 in.
Unpackaged - L 15.75 in. x W 2.17 in. x D 2.36 in.
Weight: 0.91 lbs (sgl)
2.27 lbs (3 pk)

UNO EMERGENCY SPLINT

ITEM #	OPTION
50-0045	Single
50-0046	Pack of 3



NORTH AMERICAN RESCUE®

Education and Training

North American Rescue Education and Training (NAR® Training)

has roots that run deep and wide in the world of pre-hospital medicine and technical rescue. Our team was built from the ground up on the vision and experience of seasoned military medics seeking to fill the training gaps commonly encountered in the operational environment.


Since then, we have continued to broaden our vision to include other invaluable first responder perspectives such as pre-hospital healthcare providers in civilian, law enforcement, emergency medical services and community public safety initiatives.

CORE Competencies

- Continuing Medical Education (CMEs)
- Prolonged Casualty Care
- Austere Medicine
- Technical Rescue
- Emergency Medical Services
- Critical Care Flight Medicine



NAR Training delivers a broad range of curricula which frequently include **Hands-on Skills Stations** and **Field Scenarios** designed to develop both competence and confidence. Using **Professional Role Players, Moulage Artists, Subject Matter Experts, High Fidelity Manikins, Re-animated Cadavers** and **State of the Art Facilities**, our team delivers unparalleled medical instruction to today's heroes.



Training Locations:

- NAR Training
- Training Partners & Affiliates

ACCREDITATIONS:



American Heart Association®


Now SAME DAY Cards Issued At End Of Course!



Explore training opportunities online

NAReducation.com ➤ **Phone: 888.689.6277** ➤ **Email: NARTraining@NARescue.com**

CAT[®] COMBAT APPLICATION TOURNIQUET



**THE FASTEST,
SAFEST, MOST
EFFECTIVE
PREHOSPITAL FIELD
TOURNIQUET**

TOURNIQUETS | FIRST AID KITS | TRAUMA KITS | RESCUE TOOLS | LITTERS



NORTH AMERICAN RESCUE[®]
www.NARescue.com • 888.689.6277

35 Tedwall Court | Greer, SC 29650 | USA

PRODUCTS WITH A MISSION[®]



BE THE DIFFERENCE[®]

www.NARescue.com