

NAR BP/Stethoscope Combo Kit

PN: 20-0057

Special Features:

- Ultra-compact blood pressure cuffs in a set of 4 (adult, large adult, child, infant)
- Ambidextrous palm gauge with high-impact casing and trigger valve
- One-tube bladder with quick-release luer slip-connector
- ADC stethoscope with extra-large bell and ultra-sensitive diaphragm

Dimensions:

- Case Closed - L 6.75 in. x W 9 in. x H 5 in.
- Weight: 2 lb

The NAR BP/Stethoscope Combo Kit is an ultra-compact, multi-pocketed unit including a stethoscope and set of 4 blood pressure cuffs which are instantly accessible and easily stored in most trauma kits. This combo kit features a Stethoscope from ADC with extra large bell and ultra-sensitive diaphragm coupled with a set of four rugged nylon cuffs. The cuff set features adult, large adult, child and infant size blood-pressure cuffs and one-tube bladder with quick-release luer slip-connector, as well as an ambidextrous palm gauge with high impact casing and trigger valve.



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277