

Black Talon®

Ultimate Nitrile Gloves

ITEM# 70-0001 (SM.)
NSN# 6515-01-515-0195

ITEM# 70-0002 (MED.)
NSN# 6515-01-515-0200

ITEM# 70-0003 (LRG.)
NSN# 6515-01-515-0197

ITEM# 70-0004 (XL.)
NSN# 6515-01-515-0190

Special Features:

- Textured for Improved Wet and Dry Gripping
- Longer and Thicker than Standard Examination Gloves
- Extended, Beaded Cuff for Quick and Easy Donning
- Form Fitting
- Non-sterile



Specifications:

- Length: 11.8"
- Cuff Thickness (mils): 4.7
- Palm Thickness (mils): 6.0
- Finger Thickness (mils): 8.0

Black Talon® Ultimate Nitrile gloves are made of 100% nitrile material which is free of the allergens found in standard latex gloves. They are engineered to give our users high performance by offering both maximum protection and unsurpassed tactile sensitivity. This careful attention to detail includes textured fingertips to help with fine motor skills, beaded cuffs for quick and easy donning and a thickened palm for durability and reliability.

These ambidextrous gloves are packaged 50 pairs per box (Extra large, 45 pairs per box); 10 boxes per case. Designed for single-use only.



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277