

Hasty Harness

PN: 90-0001

Multi-Purpose Applications:

- Hoisting and Lowering
- Dragging
- Assisted Carry

Designed to Allow the Rescuer to:

- Maintain Engagement with Threat
- Utilize Improved Body Mechanics (Maintain an Anatomically Correct Posture) to Rapidly Extract Larger Weight
- Maintain Natural Postures & Positions that Preclude Shifting of Personal Protective Equipment (PPE) that either Obviates that Equipment or Exposes Additional Unprotected Body Regions

Decrease Manpower Needs in Zeroed in Location

Note:

- Designed to be Utilized by Trained Professionals Only
- Not Meant to be Shock Loaded
- All items must be visually inspected for any deterioration of stitching and webbing both before and after each use.



MULTI-PURPOSE HARNESS OPTIONS FOR HIGH THREAT EXTRACTIONS...

This resourceful tool is ideal for critical adaptability during high threat rescues. Constructed of a circumference of 22 feet of 1" Tech Tape whose tensile strength of 4200 lbs. outperforms its military spec counterpart in strength, flexibility, knotability and durability.



NORTH AMERICAN RESCUE®

www.NARescue.com 888.689.6277