

High Risk Warrant Casualty Kit

ITEM# 80-0121 (BLK)

ITEM# 80-0122 (ODG)

Contents:

- 1 x Dragon® Handle System Extraction Device
- 10 x Black Talon® Nitrile Trauma Gloves (5 pr.)
- 2 x Nasopharyngeal Airway 28F with Lubricant
- 2 x Petrolatum Gauze (3 in. x 18 in.)
- 2 x ARS® Needle Decompression Kit (14 G x 3.25 in.)
- 1 x HyFin® Chest Seal
- 2 x C-A-T® (Combat Application Tourniquet®)
- 4 x ETD™ 6 in. Emergency Trauma Dressing
- 4 x S-Rolled Gauze™ (4.5 in. x 4.1 yd)
- 1 x Tactical Suction Device
- 1 x Cyclone® Bag Valve Mask
- 1 x Trauma Shears (7.25 in.)
- 2 x Triage Card (NAR T2 Tag®)

Special Features:

- Can Be Worn as a Belly Bag or Slung over the Shoulder
- MOLLE Style Attachments
- Trauma Shears Storage on Front for Quick Access
- Extra Storage Sleeve on Back contains the Dragon® Handle System to facilitate Immediate High Threat Extraction

Dimensions:

- Cased Closed: H 7.75 in. X W 13 in. X D 9 in.

Weight:

- 4 lb 13 oz



High Risk Warrant Service operations account for the majority of all SWAT call-outs nationwide. Due to the rapid dynamic nature of these missions Tactical Health Care Professionals must focus on mobility, adaptation and efficiency.

The High Risk Warrant Casualty Kit™ from North American Rescue allows the Tactical Healthcare Professional the ability to rapidly extract casualties and initiate care for the top three causes of death from penetrating trauma in the tactical environment (60% extremity bleeding / 33% tension pneumothorax / 6% airway obstruction). Offering a large central compartment and multiple attachment points, this kit maximizes storage space while maintaining a compact profile. Each component is designed to be effective under survival stress and potentially reduce some of the physiological pitfalls associated with it.



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277