

RaptorTM

IV Securing Device

PN: 30-0006 NSN: 6515-01-523-3317

Special Features:

- Minimizes equipment for Tactical Healthcare Professionals
- Works well in conjunction with Saline Locks
- Secures IV tubing on wet skin
- Sold 6 per box



Engineered specifically to quickly secure an IV line in the most difficult conditions. Once vascular access is achieved, the RaptorTM secures the IV tubing to the patient without tape.

Dimensions:

- Packaged (Individual): L 3.75 in. x W 1 in.
- Device Length (Fully Extended): L 19 in.
- Weight: 1.9 oz each
- Packaged (Box of 6): L 6 in. x W 4.8 in. x D 0.5 in.



NORTH AMERICAN RESCUE[®]

www.NARescue.com • 888.689.6277