




NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277
PRODUCTS WITH A MISSION®



*Increasing
Survivability
on the Battlefield
and Beyond...*

PRODUCTS WITH A MISSION®

ICONS



MADE IN USA: Compliant with the Buy American Act mandating "substantially all" materials, manufacturing & labor are of U.S. content



CE: Meets the requirements of the European Medical Device Directive and appropriate Quality System standards



LATEX FREE: Product is considered latex free and therefore safe for latex sensitive individuals



STERILE: Product is sterile - sterility guaranteed unless packaging is opened or damaged



RED TIP TECHNOLOGY®: Specifically branded with NAR's signature Red Tip Technology® packaging. Product packaging has been engineered for less-than-optimal, often hostile conditions that exist when our equipment is needed the most



NEW PRODUCTS: Marks the newest products throughout the catalog, making them easily recognizable as our newest innovations



Table of Contents

volume 18 • issue 1

M.A.R.C.H. Products

Products designed to address Massive Hemorrhage, Airway Management, Respiratory Support, Circulation/Vascular Access and Hypothermia Prevention

4

Combat Casualty Response Kits • CCRK®

Kits designed to meet the needs/challenges of our customers and to address the leading causes of potentially preventable death

22

DoD Bleeding Control Response Kits • D-BCRK

Individual kits, packs and stations designed to equip bystanders and first responders with the essential medical gear needed to help save lives

48

TAC-EVAC

The most effective and highest quality evacuation platforms and accessories for use in the tactical environment

52

Medical/Surgical Products

Medical and surgical equipment for use in austere environments, during prolonged field care situations and in forward operating field hospitals

60

Rescue Equipment

Includes an array of rescue products ranging from lightweight, versatile Pegasus® Ropes to tactical lights and accessories

74

Treatment Support

From personal protection products to splinting and immobilization products to help from the point of injury throughout the entire continuum of care

86

Training Products & Services

NAR Training Course Overview • Products designed to train users on TCCC protocols to reduce preventable death at the point of wounding

97

Product Index

Product offerings indexed by name, NSN number, and Item number for quick reference and information

111



M•A•R•C•H



pg 6-9



pg 10-12



pg 13-15



pg 16-20



pg 21



MASSIVE HEMORRHAGE | Tourniquets

C-A-T® Combat Application Tourniquet - Gen 7

- Scientifically proven 100% effective in occluding blood flow in both upper and lower extremities by the USAISR and recommended by CoTCCC
- Affordable, lightweight, rugged design for true circumferential pressure featuring patented free-moving internal band, proprietary Red Tip Technology® elliptical tip, and new single routing buckle for simplified training with single protocol for all applications
- Reinforced windlass and optimized windlass securing system with writeable time stamp to document application time

SPECIFICATIONS

Packaged: H 6.5 in. x W 2.4 in. x D 1.5 in.

Device: L 37.5 in.

Weight: 2.7 oz



online video
available

< scan to view



ORDERING INFORMATION

COLOR	STYLE	ITEM #	NSN #
Black	Tactical	30-0001	6515-01-521-7976
Orange	Rescue	30-0023	6515-01-612-7129
Blue	Trainer	30-0033	6910-01-560-2972



SPECIFICATIONS

Packaged: H 5.5 in. x W 7 in. x D 3.75 in.

Overall Device: L 56.75 in. | Weight: 1 lb 9.6 oz

Buttocks Pad: L 19 in. x W 6.5 in.

Compression Pad: L 4.85 in. x W 2.5 in. x D 11 in.

T-Handle: H 4.62 in. x W 4.04 in. x D 0.63 in.

ORDERING INFORMATION

ITEM #	NSN #
30-0088	6515-01-616-5841

online video
available

scan to view >



JETT® Junctional Emergency Treatment Tool

- Immediately treat high extremity wounds in the thigh/groin region at or near the point of injury
- Superior performance during patient transport - not affected by changing atmospheric pressure and temperature like a pneumatic device
- Simultaneous occlusion of blood flow to both lower limbs that will not impede respiration
- Dual compression pad assemblies for rapid bilateral and/or unilateral application based on the injury
- Mesh sleeves help pre-position the pads
- Ideal for the prehospital environment
- Fits waists up to 50 inches

Also available for
separate purchase:

JETT Pack

#80-0356 Coyote

#80-0355 Black

Nylon pouch
designed to carry
the JETT® on the outside of an aid bag





RIGID C-A-T® TQ CASE *with Shirt Shield*

- Allows carriage of a Gen 7 C-A-T® Tourniquet on your belt while protecting your shirt from undue wear caused by friction from the windlass clip
- Blade-Tech TekLok Belt Attachment fits up to a 2.25 in. belt
- Can be mounted either vertically or horizontally
- Injection Molded Nylon Case stands up to rugged usage
- Cross medical identifier patch included



ORDERING INFORMATION

COLOR	ITEM #
Black	80-0944
Coyote	80-0945
Ranger Green	80-0946

SPECIFICATIONS

H 5.88 in. x W 2.5 in. x D 2 in.
Weight: 5.36 oz

Explore other TQ carrier options available at NARescue.com



Nylon C-A-T® HOLDER



- PALS/MOLLE-style connectors provide versatile attachment capabilities (mounts horizontally or vertically)
- Includes easy-open pull elastic tab and tourniquet "TQ" ID patch
- Compatible with both Gen 6 & Gen 7 C-A-T® Tourniquets (sold separately)

SPECIFICATIONS

L 6.25 in. x W 1.75 in. x D 1.5 in.
Weight: 1.26 oz

ORDERING INFORMATION

NATO NSN*

COLOR	ITEM #	NSN #
Black	30-0057	----
Coyote	30-0044	----
MultiCam	30-0059	----
◆ Ranger Gr.	30-0208	----
OD Green	30-0053	6530-01-605-9432*

Cummerbund C-A-T® HOLDER

- Allow attachment to tactical armor such as assault and plate carriers
- Utilizes 4 in. elastic coupled with hook and loop attachment to external cummerbund flap
- Rapid, easy access with ambidextrous one handed TQ removal
- Compatible with both Gen 6 & Gen 7 C-A-T® Tourniquets (sold separately)



SPECIFICATIONS

L 5 in. x W 4 in. x D 0.25 in.
Weight: 0.6 oz

ORDERING INFORMATION

COLOR	ITEM #
◆ Black	30-0212
◆ Coyote	30-0213



Elastic C-A-T® HOLDER

- PALS/MOLLE-style strap with tuck finish eliminating the need for a button snap
- Versatile attachment to body armor, gear or belt for rapid & easy access
- Easy-open quick release hood & extended elastic body for TQ release outwards under pressure
- Compatible with both Gen 6 & Gen 7 C-A-T® Tourniquets (sold separately)

SPECIFICATIONS

L 7 in. x W 1.75 in. x D 1.5 in.
Weight: 2 oz

ORDERING INFORMATION

COLOR	ITEM #
◆ Black	30-0184
◆ Coyote	30-0185
◆ MultiCam	30-0186





MASSIVE HEMORRHAGE | Hemostatics

COMBAT GAUZE® QuikClot Z-Fold Hemostatic

- Sterile, non-woven Z-fold gauze impregnated with kaolin for temporary external control of bleeding from penetrating trauma
- Effective, simple and safe to use on wounds not amenable to tourniquet application
- Flexible & pliable, contours to all wounds
- Does not generate heat
- Individually wrapped in easy-rip, military-grade pouch
- Features an embedded x-ray detectable strip

SPECIFICATIONS



Standard Gauze: L 12 ft x W 3 in.

Packaged: H 5 in. x W 5.5 in. x D 0.625 in. Weight: 0.8 oz

XL Gauze: L 4 yd x W 4 in. (2 ply)

Packaged: H 6.5 in. x W 6.5 in. x D 0.75 in. Weight: 1.75 oz

ORDERING INFORMATION

SIZE	ITEM #	NSN #
Standard	30-0039	6510-01-562-3325
XL	30-0070	6510-01-600-4705



See page 105 for information on the Combat Gauze® Trainer Version

CHITOGAUZE®

SPECIFICATIONS



Gauze: L 12 ft x W 3 in.

Packaged: H 6 in. x W 5 in. x D 0.65 in.

Weight: 0.8 oz

ORDERING INFORMATION

ITEM #	NSN #
30-0049	6510-01-591-7740



NUSTAT®

- Patented dual-fiber technology controls bleeding through absorption at the injury site, which results in a physical aggregation of blood cells and clotting factors
- Compact packaging allows for convenient storage and rapid application

SPECIFICATIONS



Pads: 4 in. x 4 in. (Pk of 2)

Tactical: W 2 in. x L 36 in.

ORDERING INFORMATION

OPTION	ITEM #
◆ 4 in. x 4 in.	30-0236
◆ Tactical	30-0170



CELOX™

SPECIFICATIONS



Gauze: L 10 ft x W 3 in.

Packaged:

H 6.75 in. x W 4.25 in. x D 0.75 in.

Weight: 2.1 oz

ORDERING INFORMATION

ITEM #	NSN #
30-0065	6510-01-619-8636





Emergency Trauma Dressing **ETD™**

- Premium cotton construction for max quality and absorption
- High quality elastic wrap with a sterile non-adherent pad
- Durable C-clasp securing device



SPECIFICATIONS

4 in. ETD: L 70 in. x W 4 in.

Pad: H 6 in. x W 4 in

Weight: 2.6 oz

Std Pkg: L 4.5 in. x W 2.5 in. x D 175 in.

Flat Pkg: L 4.5 in. x W 3.5 in. x D 113 in.

Abd/Stump ETD:

L 70 in. x W 6 in.

Weight: 5oz

Pad: H 12 in. x W 12 in.

Pkg: L 5.25 in. x W 6 in. x D 151 in.

6 in. ETD: L 70 in. x W 6 in.

Pad: H 6 in. x W 6 in.

Weight: 4 oz

Std Pkg: L 6 in. x W 2.5 in. x D 21 in.

Flat Pkg: L 6 in. x W 3.5 in. x D 1 in.

8x10 inch ETD:

L 70 in. x W 6 in.

Weight: 4oz

Pad: H 8 in. x W 10 in.

Pkg: L 4.13 in. x W 7 in. x D 1 in.

ORDERING INFORMATION

SIZE	STYLE	ITEM #	NSN #
4 inch	Standard	30-0031	6510-01-558-4108
4 inch	Flat	30-0093	----
6 inch	Standard	30-0032	6510-01-558-4114
6 inch	Flat	30-0094	----
Abd/Stump	----	30-0012	6510-01-541-8121
8x10 inch	----	30-0165	6510-01-654-6801



**online
video
available**

< scan to view

COMPRESSED GAUZE

- Low-cube, ultra-compact sterile gauze
- Premium rolled 100% cotton gauze for maximum absorption

SPECIFICATIONS



Gauze:

W 4.5 in. x L 41 yd

Packaged:

L 3 in. x W 2.5 in. x D 1 in.

Weight: 12 oz



ORDERING INFORMATION

ITEM #	NSN #
30-0052	6510-01-503-2117

S-ROLLED GAUZE

- Unique dispenser, allows controlled linear feed application
- Premium 100% cotton

SPECIFICATIONS



Gauze:

W 4.5 in. x L 41 yd

Packaged:

L 3 in. x W 2.5 in. x D 1 in.

Weight: 12 oz



ORDERING INFORMATION

ITEM #	NSN #
30-0003	6510-01-605-2234

CRAVAT (Triangular Bandage)

Bandage: 37 in. x 37 in. x 52 in.

Pkg: 3.5 in. x 2 in. x 0.63 in. Wgt: 1.6 oz



ORDERING INFORMATION

ITEM #	NSN #
30-0089	6510-00-201-1755



TACTICAL COMPRESSION WRAP



Wrap: L 38 in. x W 5 in. Wgt: 6oz

Pkg: 5.4 in. x 2.7 in. x 2 in.

ORDERING INFORMATION

ITEM #	NSN #
30-0016	6510-01-614-0938



AIRWAY MANAGEMENT



BAC-PACK™ *Bougie-Aided Cricothyroidotomy Pack*

- All-in-one cricothyroidotomy kit for use in any environment
- Pre-inserted bougie in tube allows for a more rapid, easy procedure
- Tube meets TCCC recommendations for bougie-aided open surgical technique
- PenBlade #10 Safety Scalpel included for safety and ease of use
- Small, durable pouch keeps the items organized and ready to use
- Lightweight and ruggedized
- Components configured by order of use
- PenBlade safety scalpel, tube, syringe all provided sterile
- Patent pending

INCLUDED IN THE BAC-PACK:






SPECIFICATIONS

Packaged: H 8 in. x W 4 in. x D 1.75 in.
Weight: 5.2 oz
U.S. Patents: Pending


ORDERING INFO

ITEM #

◆ 10-0052



NAR Tracheostomy Kit with BOUGIE INTRODUCER



SPECIFICATIONS

Tube diameter:
Inner: 6.0mm | Outer: 8.0mm
Tube length: 67.0mm
Packaged:
L 5.5 in. x H 3.625 in. x D 0.5 in.
Weight: 0.84 oz

- Convenient device that gives the user a pre-inserted bougie for introduction of the tube into the trachea and indication of correct placement prior to insertion
- Adjustable for maximum user preference
- Latex free tube features 6.0mm inner diameter and inflatable cuff
- Contains the same preferred flexible tube as NAR's Tactical CricKit®

ORDERING INFO

ITEM #

◆ 10-0048

TACTICAL CRICKIT[®] *Tactical Surgical Cricothyroidotomy Kit*

- Contains all features required by CoTCCC for a surgical airway kit in the tactical setting
- Ideal length & size tracheal tube with securing device
- Configured in the sequence of use, ideal for high-stress, low-light, austere environments
- Lightweight, ruggedized and minimal in cube size
- Included rolled pouch provides a clean working surface



INCLUDED IN THE TACTICAL CRICKIT:



SPECIFICATIONS

Packaged: H 6.5 in. x W 2.5 in. x D 1.25 in.
Unrolled: H 6.5 in. x W 7.5 in. x D 1.25 in.
Weight: 3.7 oz

ORDERING INFORMATION

ITEM #	NSN #
10-0017	6515-01-540-7568



NAR TRACHEAL HOOK




ORDERING INFORMATION

ITEM #	NSN #
10-0020	6515-01-562-9734

SPECIFICATIONS

L 4.2 in. x W 0.5 in. x D 0.8 in.
Weight: 0.5 oz

NAR TACTICAL SUCTION DEVICE



SPECIFICATIONS

L 9.75 in. x W 3.75 in. x D 3 in. | Wgt: 7.3 oz
Suction Capacity: 1000 mL
Vacuum Force: 100 mmHg

ORDERING INFORMATION

ITEM #	NSN #
10-0018	6515-01-540-7206

NASOPHARYNGEAL AIRWAY

with lubricant

SPECIFICATIONS

L 6 in. x W 16 in. x D 0.5 in.
Weight: 0.5 oz

ORDERING INFORMATION

SIZE	ITEM #	NSN #
28F	ZZ-0034	6515-01-529-1187



22F-36F available
at NARescue.com

KING LTS-D[™] Supraglottic Airway Kit

- Unique dual channel design, anatomical shaped distal tip/cuff & open second lumen
- Tapered ramp design provides additional ventilation outlet

SPECIFICATIONS

Folded: H 8 in. x W 11.43 in. x D 1.92 in.
Unfolded: H 8 in. x W 13.36 in. x D 1.22 in.
Weight: 5.7 oz

ORDERING INFORMATION

CLR	SIZE	ITEM #	NSN #
YEL	3	10-0025	----
RED	4	10-0026	6515-01-591-1578
PUR	5	10-0027	6515-01-591-0397

INCLUDED WITH KING LTS-D[™]





AIRWAY MANAGEMENT

AMBU® RDIC *Mark III Resuscitator*

- Genuine Ambu Mark III reusable resuscitator equipped with impermeable elastomeric components on the compression bag and patient valve (butyl rubber)
- Specialized inlet filtration system, including C2 canister, for filtering contaminated air

SPECIFICATIONS

L 14 in. x W 7.5 in. x H 7 in.
Weight: 5.4 lb

ORDERING INFO

CE

ITEM #	NSN #
10-0047	6515-01-338-6602



CYCLONE® *Improved Pocket BVM*

- Unique collapsible design packaged 25% smaller than conventional BVMs
- Requires minimum pressure for compression of the bag to ventilate



SPECS

Diameter:
5.3 in. x H 2.9 in.
Weight: 17.4 oz

ORDERING INFO

CE

ITEM #	NSN #
10-0019	6515-01-568-0193

AMBU® SPUR II BVM

SPECIFICATIONS

L 11.63 in. x W 5 in. x H 5 in.
Weight: 11.2 oz

ORDERING INFORMATION

CE

ITEM #	NSN #
10-0043	6515-01-519-3612



Special Operations

LARYNGOSCOPE SET

- Perform direct visual laryngoscopy even in the most severe environments
- Lightweight, compact kit to equip medical care providers in the field with tools to secure an airway
- Contains 2 single-use fiber optic blades, textured handle & removable light source all rolled up in a rugged nylon pouch

SPECIFICATIONS

Open: L 8 in. x W 5.5 in. x D 2 in.
Rolled: L 8 in. x W 2 in. x D 2 in.
Weight: 10 oz

ORDERING INFORMATION

ITEM #	NSN #
10-0001	6515-01-450-9790



INCLUDED IN THE LARYNGOSCOPE SET:



QTY 1 Pouch, Nylon



QTY 1 Light, Pocket Exam



QTY 1 #2 Miller Blade



QTY 1 #3 Mac Blade



QTY 1 Handle Layrngoscope



QTY 1 Spare Bulb




QTY 2 Size AA Batteries

NAR Fiber Optic LARYNGOSCOPE BLADES

ORDERING INFORMATION

CE

AVAILABLE SIZES	ITEM #	NSN# available on all sizes offered
Miller #1, #2, #3 & #4	See NARescue.com	See NARescue.com for details
Mac #2, #3 & #4	See NARescue.com	See NARescue.com for details






HYFIN® Chest Seals

Twin Pack

- Two individually packaged, vented/valve-type chest seals in one unique package designed for treatment of both entry/exit wounds or multiple penetrating injuries.
- Patented 3 channel pressure relief vents
- Ruggedized dual perforated pouch designed to allow the provider to open only one dressing at a time, as needed
- Advanced adhesive technology designed to provide a superior seal



ORDERING INFORMATION

ITEM #	NSN #
10-0037	6510-01-642-6212



SPECIFICATIONS

Device: 6 in. x 6 in. (excludes tab) Wgt: 2.5 oz
 Packaged: H 7.5 in. x W 4.5 in. x D 0.25 in.
 Unfolded: H 7.5 in. x W 9 in. x D 0.13 in.




online video
available

< scan to view

HyFin® Vent Individual

ORDERING INFORMATION

ITEM #	NSN #
10-0029	6510-01-624-0840



SPECIFICATIONS

Device: 6 in. x 6 in. (excludes tab) Wgt: 1.31 oz
 Packaged: H 4.63 in. x W 7.5 in. x D 0.1 in.



HyFin® Occlusive Individual

ORDERING INFORMATION

ITEM #	NSN #
10-0015	6510-01-532-8019



SPECIFICATIONS

Device: 6 in. x 6 in. (excludes tab) Wgt: 1.31 oz
 Packaged: H 4.63 in. x W 7.5 in. x D 0.1 in.



USMC HyFin® Chest Seal COMBO PACK

- Versatile combo pack designed for the treatment of both entry/exit or multiple penetrating injuries to the chest
- Includes one HyFin® Occlusive & one HyFin® Vented Chest Seal packaged in a clear resealable pouch

SPECIFICATIONS

Packaged: L 9 in. x W 5 in. x D 0.2 in.
 Each Chest Seal (when deployed):
 H 6 in. x W 6 in. x D 0.1 in.
 Total Weight: 2.9 oz
 US Patent 7,504,549 &
 Patent(s) Pending



USMC COMBO PACK:



ORDERING INFORMATION

ITEM #	NSN #
85-1045	6510-01-654-6787





RESPIRATORY SUPPORT

*Tension
Pneumothorax*

10 ga x 3.75 inch

SPEAR™ Simplified Pneumothorax Emergency Air Release

- Latest innovation in Needle Decompression technology
- Designed based on known gaps in the management of tension pneumothorax
- Enhanced needle tip for improved insertion control
- User friendly features like spin-lock connection, depth markings, and detachable low-pressure one-way valve
- Fenestrations along catheter for additional air flow
- Familiar, rugged, pen-style case for ease of storage and deployment
- Allows for either lateral decompression or anterior needle thoracostomy
- Catheter made of flexible kink-resistant material with a soft tip
- Single use only. Supplied sterile



**online video
available**

< scan to view

SPECIFICATIONS


Packaged:
L 7.625 in. x Diameter 0.75 in.
Catheter: 10Ga x 3.75 in.
Weight: 1 oz
US Patents: Pending

ORDERING INFORMATION

GAUGE	ITEM #	NSN #
◆ 10 ga	10-0051	---



NOTE: Not for sale in Canada



*Versatility that allows for
both Lateral and Anterior
Decompression Procedures...*

**10 gauge 3.75 in.
FLEXIBLE CATHETER with
DEPTH GAUGE**

DISTAL FENESTRATIONS

SOFT SAFETY TIP

**NEEDLE with
MODIFIED TIP**

Lateral Approach

Anterior Approach

3.25 inch **ARS® NEEDLE DECOMPRESSION KIT**

- Penetrate the pleural space in 99% of casualties presenting the signs and symptoms of a tension pneumothorax
- Quality, reliable devices with 3.25 in. length needle/catheter
- Capless flash chamber provides immediate confirmation of correct placement
- Patented design packaged in an easy-to-open, crush-proof protective tube featuring easy-to-identify textured twist top with clip

SPECIFICATIONS

Pkg: H 5.7 in. x W 0.7 in. x D 0.8 in.
Weight: 0.6 oz

ORDERING INFORMATION



GAUGE	ITEM #	NSN #
14 ga	ZZ-0056	6515-01-541-0635
10 ga	ZZ-0298	6515-01-673-1701

U.S. Design Patents: D584,409 S; D595,847S; D584,410S; U.S. Patent 7,874,426B2; © North American Rescue,® LLC., 127820, 127822, 127823, 001002372-0001, 001013940-0001, 001013940-0002

NOTE: Not for sale in Canada



online video available

< scan to view



INTRO

M.A.R.C.H

KITS

TAC-VIAC

MED-SURG

RESCUE

TREATMENT SUPPORT

TRAINING

INDEX



Meets all requirements for US Army Task: 081-831-1011

NAR SALINE LOCK KIT



Essential components for establishing a saline lock for safer, quicker and more effective IV cannulation in the most demanding environments

ORDERING INFORMATION

ITEM #	NSN #
30-0013	6515-01-537-4094

SPECIFICATIONS

Pkg: H 5.5 in. x W 3.25 in. x D 1 in.
Weight: 13 oz



SALINE LOCK KIT CONTENTS:



NAR NEEDLELESS SALINE LOCK KIT

IV safety catheter & needle-free valve port reduce or eliminate accidental needle-stick exposures



ORDERING INFORMATION

ITEM #	NSN #
30-0046	6515-01-607-1943

SPECIFICATIONS

Pkg: H 5.5 in. x W 3.25 in. x D 1 in.
Weight: 13 oz



NEEDLELESS SALINE LOCK KIT CONTENTS:



BOA® IV CONSTRICTING BAND

- Direct peripheral IV access utilizing a revolutionary rolling action for transfusing venous blood to the distal veins while still allowing arterial flow
- Quick Click connector & two-finger Quick Release improves "time on task"
- XL version stretches approximately 50% longer than the standard BOA®

SPECIFICATIONS

Standard: L 7.5 in. to 18 in. XL: L 10 in. to 27 in.
Packaged: L 4 in. x W 3 in. x D 0.9 in. Weight: 1.5 oz
US Patents: 5,304,202 & 5,607,448



online video available

< scan to view

ORDERING INFORMATION

SIZE	ITEM #	NSN #
Standard	30-0009	6515-01-537-2611
XL	30-0071	----



SHARPS SHUTTLE

SPECIFICATIONS

L 6.38 in. x Diameter 1.11 in.
Weight: 1.5 oz

ORDERING INFORMATION

UNIT	ITEM #	NSN #
Single	30-0014	6515-01-518-8551
6 pack	30-00146	6515-01-619-5671



RAPTOR® IV Securing Device

SPECIFICATIONS

Extended: L 19 in. x W 1 in.
Weight: 19 oz

ORDERING INFORMATION

UNIT	ITEM #	NSN #
6 pack	30-0006	6515-01-523-3317



G3 POWER DRIVER by EZ-IO®

- Battery-powered solution for providing rapid vascular access to humerus & tibia in less than 10 seconds (Not for sternal use)
- Trigger guard prevents accidental activation
- Molded Hard Side Carrying Case (sold separately)

ORDERING INFORMATION

UNIT	ITEM #
Driver	30-0075
Case	30-0077

SPECIFICATIONS


H 4.5 in. x W 6.5 in. x D 2.5 in.
Weight: 15 oz



NOTE: EZ-IO Trainers also available - see page 108

MANUAL EXTREMITY NEEDLE SET

Manual I/O insertion at 6 extremity sites by EZ-IO®



SPECIFICATIONS

Device: L 6 in. x W 3.25 in.
Weight 1.4 oz

ORDERING INFO

ITEM #
30-0078

MANUAL NEEDLE SET:



POWER DRIVER NEEDLE SET

for use with the EZ-IO® G3 Power Driver


ORDERING INFORMATION

COLOR	SIZE	ITEM #
Pink	15mm/PD	30-0081
Blue	25mm/AD	30-0083
Yellow	45mm/LD	30-0085

SPECS

Packaged:
H 4.75 in.
W 4.25 in.
D 0.7 in.
Weight: 1 oz

POWER DRIVER NEEDLE SET:



T.A.L.O.N.™ MANUAL NEEDLE SET

Manual I/O insertion at 7 extremity sites - including the sternum

SPECIFICATIONS

Device: L 6 in. x W 3.25 in.
Weight 1.5 oz

ORDERING INFORMATION

SIZE	ITEM #
38.5mm	30-0132

T.A.L.O.N. MANUAL NEEDLE SET:





QUANTUM™ Blood & Fluid Warming System

- **Lightweight** The Quantum™ has the lightest operational weight of any portable blood & fluid warming system (Total Operational Weight: 22 oz)
- **Powerful** Warms two units of 39°F refrigerated whole blood at 100mL/min+ on a single charge
- **Compact** Fewest components to carry (three)
- **Fast** Warm-up time of only 24 seconds
- **Intuitive** The thermal transfusion/infusion set looks, handles, and sets up/primes like any conventional IV set
- **Safe** Tested safe for blood, blood products and IV solutions
- **Effective** The system is plug-and-play for ease of setup under stressful conditions such as tactical field care (only two connection points, with no plug orientation or set-up order requirements – system status constantly visible at a glance)
- **Complete** Does not require additional standard transfusion or infusion set – all other warmers do
- **Tactical** Only system FDA cleared allowing easy disable of audible alerts to protect operators in combat

Quantum Charger

Quantum Battery


Quantum Controller

TTS-B Thermal Transfusion Set - Blood (Blood, blood products & IV solutions)

TIS Thermal Infusion Set (IV solutions)

PRODUCT WEIGHTS

Total Operational Weight: 22 oz
(Quantum™ Controller, Battery & TTS-B)
Battery: 16 oz
Controller: 3 oz
Thermal Transfusion Set: 2.7 oz
Thermal Infusion Set: 2.3 oz
Charger: 22 oz



OPERATING PARAMETERS

IV FLUIDS

- Input Temp: 68°F
- Output Temp: 100.4°F (38°+/- 2°C)
- Flow-Rate: 2 - 200mL/min (+)
- Warm-up Time: 24 seconds
- Warming Capacity with a Single Battery: 1700 mL

BLOOD or BLOOD PRODUCTS

- Input Temp: 39°F
- Output Temp: 100.4°F (38°+/- 2°C)
- Flow-Rate: 2 - 100mL/min (+)
- Warm-up Time: 24 seconds
- Warming Capacity with Single Battery: 2 units (+)

*NOTE:
Complete System includes
1 x Controller
1 x Battery
1 x Charger
2 x TTS-B

ORDERING INFORMATION

COMPONENTS	ITEM #
♦ Complete System*	35-0006
♦ TTS-B (Pack of 2)	35-0005
♦ TIS (Pack of 2)	35-0004
♦ Controller	35-0003
♦ Charger	35-0002
♦ Battery	35-0001

Experience the **Quantum** **DIFFERENCE** The CAPACITY vs. CAPABILITY solution...



*Blood & Fluid
Warming System*

QuantumTM

Close the **Fluid Resuscitation Capability Gap** that requires a **Normothermic Whole Blood** delivery system for far forward, near point-of-injury use

FAST | EFFECTIVE | SAFE | INTUITIVE | LIGHTWEIGHT | POWERFUL

REQUIREMENT:

- ▶ Delivery of 2 units of 39°F (4°C) Whole Blood at 100ml/min with an output temperature of 100.4°F (38°± 2°C) on a single battery charge
- ▶ Warm-up time - 24 seconds ▶ Lightweight & Portable - 22 oz total

We aren't playing the game. We're changing it.

Explore all NEW Online Educational Materials/Videos available at
NARescue.com/Quantum



NORTH AMERICAN RESCUE[®]

www.NARescue.com • 888.689.6277

REV101819

CIRCULATION & Vascular Access



FBTK *Field Blood Transfusion Kits*

- Comprehensive kit that provides all the needed supplies to allow the user to rapidly perform a field blood transfusion
- Equipped to allow for verifying the blood types of the donor and recipient and collecting blood samples from the donor for additional testing

ORDERING INFORMATION

OPTIONS	ITEM #
FBTK Basic	83-0022
Blood Surveillance	83-0025
FBTK Complete*	83-0000

FBTK *Basic:*

Donor Module:

- 5 x Alcohol Prep Pad
- 2 x 2-pack Gauze Pad - 4 in. x 4 in.
- 1 x IV Constricting Band - 1 in. x 18 in.
- 1 x Clear Adhesive Dressing - 4 in. x 4.75 in.
- 1 x Blood Bag Label

Recipient Module:

- 5 x Alcohol Prep Pad
- 1 x IV Catheter - 18 ga x 1.25 in.
- 2 x 2-pack Gauze Pad - 4 in. x 4 in.
- 1 x IV Constricting Band - 1 in. x 18 in.
- 1 x Clear Adhesive Dressing - 4 in. x 4.75 in.
- 1 x IV Catheter - 16 ga x 1.25 in.

Additional Supplies:

- 1 x roll 1 in. Surgical Tape
- 2 x pair Black Nitrile Gloves, Large
- 1 x Black Permanent Marker
- 1 x Kelly Curved Hemostat - 5.5 in.
- 1 x Single Blood Bag (450ml)
- 1 x 550 Cord - 9.5 in.
- 1 x Beaded Security Tie - 11 in.
- 1 x IV Set Y-Type with Filter 170-260 UM
- 1 x Questionnaire Card - DD572 Modified

FBTK *Blood*

Surveillance Module:

Donor Module:

- 1 x IV Hypodermic Needle - 18 ga x 1.5 in.
- 1 x IV Syringe, 10cc
- 1 x Eldon Card
- 1 x Interlink Injection Site Needle
- 1 x Blood Transfer Device
- 2 x Vacutainer, Gold Top (5ml)
- 1 x Vacutainer, Purple Top (5ml)
- 2 x Vacutainer, Pearl Top (5ml)

Recipient Module:

- 1 x 3 Way Stopcock
- 1 x Eldon Card
- 1 x Interlink Injection Site Needle

*FBTK *Complete:*

Contains all components of **FBTK Basic** & **Blood Surveillance Module**



HYPOTHERMIA Prevention & Management

HPMK® Hypothermia Prevention & Management Kit

- Includes reinforced Heat Reflective Shell (HRS™) of 4-ply, composite fabric with a protected non-conductive reflective layer
- Built-in hood, integrated fluid absorption pad & top-to-bottom taper providing excellent thermal insulation
- Continuous 1.5 in. Velcro® closures provide rapid 360° access to the casualty with minimal exposure to the elements
- Self-heating, 4-cell liner (Ready-Heat™ Blanket) provides 10 hours of continuous dry heat up to 125°F (52°C)
- Vacuum-packaged for a low-cube footprint

SPECIFICATIONS

Shell Pkg: H 6.75 in. x W 10.5 in. x D 5.5 in.
 Shell Open: L 78 in. x W 43 in.
 US Patent: 7,766,950 B2
 Weight: 3 lb 8 oz



HPMK® CONTENTS:



ORDERING INFORMATION

ITEM #	NSN #
80-0027	6515-01-532-8056

NOTE: Litter is sold separately

NOTE: Not for sale in Canada

HRS® Heat Reflective Shell



As an individual unit, reinforced 4-ply, composite shell with built-in hood, weight tested to support 250+ lb.

ORDERING INFORMATION

ITEM #	NSN #
ZZ-0051	6532-01-627-7888

SPECIFICATIONS

Shell Packaged: L 10 in. x W 6.75 in. x D 2.5 in.
 Open: L 78 in. x W 43 in. US Patent: 7,766,950 B2
 Weight: 1 lb 13 oz



READY HEAT™ BLANKET



As an individual unit, the self-heating, 4-cell liner provides 10 hours of continuous dry heat after reaching peak temperature approximately 15-20 minutes after opening

SPECIFICATIONS

Open: L 48 in. x W 36 in.
 Weight: 1 lb 11 oz



ORDERING INFORMATION

ITEM #
ZZ-0052

EMERGENCY SURVIVAL WRAP

Designed to reflect heat back to the body

ORDERING INFORMATION

ITEM #
ZZ-0834

SPECIFICATIONS

Unfolded: 60 in. x 90 in.
 Weight: 3.7 oz







KITS

COMBAT CASUALTY RESPONSE (CCRK®)

- Warfighter/Operators • Combat Lifesavers
- 1st Responders • Tactical Healthcare Professionals
- DoD Bleeding Control Response Kits **pages 24-52**



SIRK *Supplemental IFAK ReSupply Kits*

KIT CAPABILITIES INCLUDE:



SIRK *Gen 2*


SIRK GEN 2:



QTY
2 C-A-T®
Tourniquet



QTY
1 ETD™
4 in. Flat



QTY
1 Combat
Gauze Z-fold



QTY
1 Gauze,
Compressed



QTY
1 28F Naso
Airway & Lube



QTY
1 HyFin®
Vent. Twin Pk



QTY
1 Fox Eye
Shield w/garter



PAIR
2 Nitrile
Gloves, Lg



QTY
1 Mini Duct
Tape, 2 in.



QTY
1 Model 7
Rescue Hook



QTY
1 Permanent
Marker, Large



QTY
1 Combat
Casualty Card



85-1899 ONLY
QTY
1 ARS®
Needle, 14g

SPECIFICATIONS

H 7 in. x W 5 in. x D 2 in.
Weight: 14 oz

SIRK GEN 2 ORDERING INFO

ARS OPTION	ITEM#
w/o ARS Needle	85-1738
with ARS Needle	85-1899

- Vacuum-sealed first aid kits ideal for re-stocking an existing IFAK that contains the essential items to treat leading causes of preventable combat death
- Designed for small cube space in rugged, clear package allowing contents to be easily visualized

SIRK *Gen 1*

SIRK GEN 1:



QTY
1 C-A-T®
Tourniquet



QTY
1 ETD™ 6 in.
Dressing



QTY
1 Gauze,
Compressed



QTY
1 Choice of
Hemostatic



QTY
1 28F Naso
Airway & Lube



QTY
1 HyFin® Chest
Seal, Individual



PAIR
2 Nitrile
Gloves, Lg



QTY
1 Surgical
Tape, 1 in.



QTY
1 PES
Eye Shield



QTY
1 Medevac
9 Line Card



SPECIFICATIONS

H 7 in. x W 5 in. x D 2 in.
Weight: 14 oz

SIRK GEN 1 ORDERING INFO

HEMOSTATIC	ITEM#
Combat Gauze	80-0291
ChitoGauze	80-0225**



IPOK *Individual Patrol Officer Kit*

KIT CAPABILITIES INCLUDE:



- Individual hemorrhage control kit small enough to be carried on your person
- Compact, vacuum-sealed & packaged to fit into a BDU pocket, vest pouch or individual bag
- Rugged, durable pouch features NAR's signature Red Tip Technology® ensuring quality & performance

IPOK:



QTY
1 C-A-T®
Tourniquet



QTY
1 ETD™
4 in. Flat



QTY
1 Choice of S-Rolled Gauze
or Hemostatic Dressing



PAIR
1 Nitrile
Gloves, Lg


IPOK ORDERING INFO

HEMOSTATIC/GAUZE	ITEM#
S-Rolled Gauze	80-0167**
Combat Gauze, LE	80-0168
ChitoGauze	80-0223**



SPECIFICATIONS

H 4 in. x W 6.5 in. x D 2.75 in.
Weight: 7.4 oz





M-FAK *Mini First Aid Kit*

KIT CAPABILITIES INCLUDE:



- Contains the essential point-of-wounding equipment needed for treating penetrating, blast and other traumatic injuries
- Super compact, rugged nylon platform that allows attachment both vertically (MOLLE backing) and horizontally (3 in. belt loop)
- Clam shell configuration utilizes (2) two main sleeves that open on both ends for easy access

ORDERING INFORMATION

COLOR	LEVEL	ITEM #
Black	Basic	80-0494
Coyote	Basic	80-0495
Multi-Cam	Basic	80-0496
OD Green	Basic	80-0497
Ranger Green	Basic	80-0996
Black	Combat Gauze LE	80-0498
Coyote	Combat Gauze LE	80-0499
Multi-Cam	Combat Gauze LE	80-0500
OD Green	Combat Gauze LE	80-0501
Ranger Green	Combat Gauze LE	80-0998
Black	Advanced	80-0502
Coyote	Advanced	80-0503
Multi-Cam	Advanced	80-0504
OD Green	Advanced	80-0505
Ranger Green	Advanced	80-0997




SPECIFICATIONS

H 6 in. x W 3 in. x D 3 in.
Weight: 13 oz

M-FAK BASIC:



M-FAK with COMBAT GAUZE:



M-FAK ADVANCED:



INTRO

M-A-R-C-H

KITS

TAC-EVAC

MED-SURG

RESCUE

TREATMENT SUPPORT

TRAINING

INDEX



EAGLE IFAK

Compact and versatile Individual First Aid Kit

KIT CAPABILITIES INCLUDE:



- Can connect to the assault vest/body armor or to a duty/pistol belt
- Contains the essential equipment needed for self-aid/buddy-aid to address the leading causes of potentially preventable death in a tactical environment
- Quick Pull Tab easily opens the kit with just one hand
- Unique, detachable shock-cord inner panel that can be adjusted to your components
- Inner panel is removable for easy access to your kit

SPECIFICATIONS

H 8 in. x W 4 in. x D 3 in.

Weight: 1 lb 2 oz

Bag only option
also available at
NARescue.com



ORDERING INFORMATION

COLOR	LEVEL	HEMOSTATIC	ITEM #
Black	Basic	None	85-0177
Coyote	Basic	None	85-0178
Multi-Cam	Basic	None	85-0179
OD Green	Basic	None	85-0180
Black	Advanced	Combat Gauze	85-0189
Coyote	Advanced	Combat Gauze	85-0190
Multi-Cam	Advanced	Combat Gauze	85-0191
OD Green	Advanced	Combat Gauze	85-0192
4 options online	Advanced	ChitoGauze	See NARescue.com

EAGLE IFAK BASIC:



EAGLE IFAK ADVANCED:





RIG Series **EAGLE IFAK**

New, low profile IFAK with rapid access to your C-A-T® and other equipment



- Tourniquet attached for quick deployment by banding it to inside of front opening flap
- Protective hood with quick pull latch release helps protect IFAK components and includes sewn in embroidered medical cross identifier
- Detachable (hook & loop backed) shock-cord inner panel adjusts to secure internal components and is easily removable for rapid access

ORDERING INFORMATION

COLOR	LEVEL	ITEM #
Black	Basic	80-0509
Coyote	Basic	80-0510
Multi-Cam	Basic	80-0511
OD Green	Basic	80-0512
Ranger Green	Basic	80-0991
Black	Combat Gauze LE	80-0517
Coyote	Combat Gauze LE	80-0518
Multi-Cam	Combat Gauze LE	80-0519
OD Green	Combat Gauze LE	80-0520
Ranger Green	Combat Gauze LE	80-0993
Black	Advanced	80-0513
Coyote	Advanced	80-0514
Multi-Cam	Advanced	80-0515
OD Green	Advanced	80-0516
Ranger Green	Advanced	80-0992

Bag only option also available at NARescue.com

SPECIFICATIONS

H 8 in. x W 4 in. x D 3 in.
Kit Packed Weight: 118 lb
Empty Bag Weight: 7 oz


RIG SERIES EAGLE BASIC:



RIG SERIES EAGLE with COMBAT GAUZE LE:



RIG SERIES EAGLE ADVANCED:





OCHO IFAK

Compact and Low Visibility IFAK

KIT CAPABILITIES INCLUDE:



- Compact, low visibility, point-of-wounding IFAK designed for your duty/pistol belt
- Purposely designed to look like a non-tactical, accessory bag
- Equipped with a 3 in. belt loop attachment to accommodate a duty/pistol belt
- Constructed with rip stop nylon, elastic loop organizational storage and easy open zipper

SPECIFICATIONS

L 7 in. x H 2.25 in. x D 3.25 in.

Weight: 13 oz

ORDERING INFORMATION

COLOR	LEVEL	ITEM #
Black/Gray	Basic	80-0374
Blue/Gray	Basic	80-0375
Red/Gray	Basic	80-0915
Black/Gray	Combat Gauze	80-0921
Blue/Gray	Combat Gauze	80-0922
Red/Gray	Combat Gauze	80-0917
Black/Gray	Advanced	80-0372
Blue/Gray	Advanced	80-0373
Red/Gray	Advanced	80-0916

See OCHO Accessory Bag only option on page 50

OCHO IFAK BASIC:



OCHO IFAK with COMBAT GAUZE:



OCHO IFAK ADVANCED:





Tactical Operator Response Kit **T.O.R.K.**



- Streamlined, compact IFAK designed to MOLLE to the assault vest or body armor
- Includes essential life-saving equipment to be used at the point of wounding
- For self-aid/buddy-aid injuries caused by penetrating or blast trauma

ORDERING INFORMATION

COLOR	LEVEL	HEMOSTATIC	ITEM #
Black	Basic	None	85-0008
Coyote	Basic	None	85-0009
Digital Camo	Basic	None	85-2300
Multi-Cam	Basic	None	85-0012
OD Green	Basic	None	85-0054
Red	Basic	None	85-1805
Ranger Green	Basic	None	85-2301
Black	Advanced	Combat Gauze	85-0048
Coyote	Advanced	Combat Gauze	85-0049
Digital Camo	Advanced	Combat Gauze	85-0050
Multi-Cam	Advanced	Combat Gauze	85-0052
OD Green	Advanced	Combat Gauze	85-0053
Ranger Green	Advanced	Combat Gauze	85-2302
Red	Advanced	Combat Gauze	85-1806
7 options online	Advanced	ChitoGauze	See NARescue.com

T.O.R.K. BASIC: CE




SPECIFICATIONS

H 7 in. x W 4.5 in. x D 2.75 in.
Weight: 1 lb 5 oz



T.O.R.K. ADVANCED:



Bag only options
also available at
NARescue.com

INTRO

M.A.R.C.H.

KITS

TAC/EVAC

MED-SURG

RESCUE

TREATMENT SUPPORT

TRAINING

INDEX



USCG IFAK

Meets specific requirements of the US Coast Guard

Versatile modular system can be worn by the operator as a drop-down leg rig, duty-belt add-on or assault vest/body armor attachment (via MOLLE-style connectors)

ORDERING INFO

COLOR	ITEM #	NSN #
Coyote	80-0158	6545-01-587-8513
Black	80-0149	6545-01-572-7248
Resupply*	85-0104	----

USCG IFAK:



QTY
1 C-A-T®
Tourniquet



QTY
1 ETD™ 6 in.
Dressing



QTY
1 Gauze
S-Rolled



QTY
1 Combat
Gauze Z-fold



QTY
2 28F Naso
Airway & Lube



QTY
2 Gauze,
Petrolatum



QTY
1 HyFin®
Occlusive, Sgl.



PAIR
2 Nitrile
Gloves, Lg



QTY
1 PES
Eye Shield



QTY
1 Tape, 2 in.
Gecko-Grip

SPECIFICATIONS

H 8.25 in. x W 7.1 in. x D 5.5 in
Weight: 2 lb 10 oz

Bag only option also available at NARescue.com

*Resupply does not include the Nylon IFAK Bag

KIT CAPABILITIES INCLUDE:



MARITIME ASSAULT KIT

Meets specific requirements of the Navy Expeditionary Combat Command

Designed to meet the unique & specific requirements of the Navy Expeditionary Combat Command for a small, efficient IFAK

SPECIFICATIONS

H 7.4 in. x W 6.9 in. x D 5 in Weight: 1 lb 13 oz

MARITIME ASSAULT KIT:



QTY
1 C-A-T®
Tourniquet



QTY
1 ETD™ 6 in.
Dressing



QTY
1 Gauze
S-Rolled



QTY
1 Combat
Gauze Z-fold



QTY
1 28F Naso
Airway & Lube



QTY
1 HyFin®
Vent. Twin Pk



QTY
1 ARS®
Needle, 14g



PAIR
2 Nitrile
Gloves, Lg



QTY
1 PES
Eye Shield



QTY
1 T2® Tag
Triage Card



QTY
1 Casualty
Reference Card



QTY
1 Model 7
Rescue Hook

KIT CAPABILITIES INCLUDE:



ORDERING INFORMATION

COLOR	ITEM #	NSN #
Multi-Cam	80-0245	----
Coyote	80-0015	6545-01-529-4187
OD Green	80-0014	6545-01-531-6626
Black	80-0013	6545-01-529-4185
Digital Camo	80-0058	----

Bag only option also available at NARescue.com

INDIVIDUAL KIT

Designed for the Warfighter/Operator

Designed with all the tools for the Warfighter/Operator to mitigate extremity bleeding, tension pneumothorax & airway obstruction

SPECIFICATIONS

H 8.5 in. x W 8.75 in. x D 4.5 in. Weight: 1 lb 12 oz

ORDERING INFORMATION

COLOR	ITEM#	NSN#
Coyote	80-0005	----
OD Green	80-0004	----
Black	80-0003	6530-01-529-1134
Digital Camo	80-0091	----

Bag only option also available at NARescue.com



KIT CAPABILITIES INCLUDE:



INDIVIDUAL KIT:



MEDIC LEG RIG

Rapid Access and Organization of Critical Equipment

Designed for Tactical Healthcare Professionals configurable based on individual skill set & mission requirement

SPECIFICATIONS

H 9.5 in. x W 9 in. x D 7.75 in. Weight: 3 lb 9 oz

ORDERING INFORMATION

COLOR	ITEM#	NSN#
Multi-Cam	80-0249	----
Coyote	80-0048	6545-01-543-3170
OD Green	80-0047	6545-01-543-3174
Black	80-0046	6545-01-543-3185
Digital Camo	80-0087	----


Bag only option also available at NARescue.com



KIT CAPABILITIES INCLUDE:



MEDIC LEG RIG:





SQUAD KIT

Enhance the survivability of the Assault Squad

KIT CAPABILITIES INCLUDE:



- Contoured profile, large central compartment, multiple attachment points & front trauma-shear storage
- Multiple-load carriage capabilities (fanny pack, over the shoulder, attachments for 3-day assault packs & MOLLE-style webbing)



Bag only option also
available at NARescue.com

SPECIFICATIONS



H 7.75 in. x W 13 in. x D 9 in.

Weight: 4 lb 6.5 oz

ORDERING INFORMATION

COLOR	ITEM #	NSN #
Multi-Cam	80-0251	----
Coyote	80-0039	6545-01-543-2956
OD Green	80-0038	6545-01-543-3159
Black	80-0037	6545-01-543-5721
Digital Camo	80-0093	6545-01-543-5733

SQUAD KIT:





MARINER KIT

Meets requirements of EMAT



- Designed to meet the rigorous requirements of the Emergency Med Assist Team (EMAT) Response Bag
- Standard casualty response bag for the Submarine Forces
- Rugged rescue orange nylon bag for high visibility, easy transport & convenient storage
- Compact, contoured profile with large central compartment, multiple attachment points & front trauma-shear storage
- Designed to be worn as a belly bag or slung over the shoulder



SPECIFICATIONS


H 7.75 in. x W 13 in. x D 9 in.
Weight: 4 lb 6.5 oz

ORDERING INFORMATION

ITEM #	NSN #
80-0107	6545-01-543-5770

Bag only option
also available at
NARescue.com

MARINER KIT:





NAR-4™ CHEST POUCH

Rapid access to critical medical equipment from your chest

KIT CAPABILITIES INCLUDE:



- Auxiliary NAR-4 companion kit streamlined to provide rapid access to critical basic & advanced equipment from the chest
- Designed to be secured to the front of the body armor using a MOLLE-style attachment system

- Ambidextrous trauma shear accessibility

SPECIFICATIONS

H 7.75 in. x W 9.75 in. x D 4.5 in

Weight: 2 lb 4 oz

ORDERING INFORMATION

COLOR	HEMOSTATIC	ITEM #
Multi-Cam	Combat Gauze	80-0247
Coyote	Combat Gauze	80-0175
OD Green	Combat Gauze	80-0289
Black	Combat Gauze	80-0173
Digital Camo	Combat Gauze	80-0176
5 options online	Chito Gauze	See NARescue.com

NAR-4 CHEST POUCH:



QTY
2 C-A-T®
Tourniquet



QTY
2 ETD™ 6 in.
Dressing




QTY
2 Gauze
S-Rolled



QTY
2 Choice of
Hemostatic



QTY
2 ARS®
Needle, 14g



QTY
2 28F Naso
Airway & Lube



QTY
1 Tactical
Crickit®



QTY
1 HyFin®
Vent. Twin Pk



PAIR
3 Nitrile
Gloves, Lg



QTY
1 Trauma
Shears, Lg



QTY
1 Scissor
Leash™



QTY
2 PES
Eye Shield



QTY
2 Combat
Casualty Card



QTY
1 Permanent
Marker, Large

Bag only option also
available at NARescue.com





NAR-4™ AID KIT

Configurable for individual skill set and mission requirements



Bag only option
also available at
NARescue.com



- Low-profile, load carriage for Tactical Healthcare Professionals
- Modular system with 5 segregated removable equipment pouches to carry equipment consistent with TCCC Guidelines

- Designed to attach directly to the body armor with four rapid attachment straps or carried using shoulder straps

SPECIFICATIONS



H 19 in. x W 15 in. x D 7 in.

Weight: 12 lb 8.3 oz

ORDERING INFORMATION

COLOR	HEMOSTATIC	ITEM #
Black	Combat Gauze	80-0181
Coyote	Combat Gauze	80-0183
Digital Camo	Combat Gauze	80-0184
Multi-Cam	Combat Gauze	80-0244
OD Green	Combat Gauze	80-0308
5 options online	ChitoGauze	See NARescue.com

NAR-4 AID KIT:





Medic TRAUMA PACK

Comprehensive trauma management platform

KIT CAPABILITIES INCLUDE:



- 9 removable, color-coded pull bags and 2 removable side compartments provide modular organization
- Over 1,000 cubic inches of gear capacity and MOLLE-style attachment system compatible with removable side pouches

ORDERING INFORMATION

**NATO NSN

COLOR	ITEM #	NSN #
Black	80-0070	6545-01-572-7251
Coyote	80-0072	6545-01-588-1862
Digital Camo	80-0096	----
OD Green	80-0071	6545-01-594-5617**

SPECIFICATIONS

H 19 in. x W 19 in. x D 10 in

Weight: 18 lb 1 oz

MEDIC TRAUMA PACK:

 QTY 2 C-A-T® Tourniquet	 QTY 6 ETD™ 6 in. Dressing	 QTY 2 ETD™ Abdominal	 QTY 4 Gauze S-Rolled	 QTY 2 28F Naso Airway & Lube
 QTY 1 King LT-D Size 4	 QTY 1 Cyclone Pocket BVM	 QTY 1 Tactical Suction Device	 QTY 2 Gauze Petrolatum	 QTY 2 ARS® Needle, 14g
 QTY 1 BOA® IV Band	 QTY 4 Saline Lock Kit	 QTY 1 Sharps Shuttle	 QTY 1 HRS, Heat Reflective Shell	 QTY 1 Tactical Traction Splint
 QTY 4 SAM® Splint II	 QTY 1 Trauma Shears, Lg	 QTY 2 Surgical Tape, 2 in.	 QTY 1 Armadillo Med Case	 QTY 1 Head Lamp
 PAIR 5 Nitrile Gloves, Lg	 QTY 2 PPE Kit	 QTY 6 PES Eye Shield	 QTY 4 T2® Tag Triage Card	 QTY 1 Casualty Reference Card

Bag only option
also available at
NARescue.com

9 removable, color-coded pull bags for modular organization

Warrior Aid & Litter Kit **WALK®**

- Official vehicle life-saving kit of the U.S. Army
- For treatment of multiple trauma casualties and the immediate evacuation of one patient
- Includes full-size Talon® 90C quad-fold litter
- Pocket system holds self or buddy-aid medical supplies
- Back pack straps provide easy transport from the vehicle to the point of need
- Includes Hypothermia Prevention & Management Kit™ to help stabilize patients during transport
- Interfaces with any standard or non-standard TACEVAC asset


ORDERING INFORMATION

COLOR	ITEM #	NSN #
Black	80-0097	----
Coyote	80-0089	6545-01-587-1199
Digital Camo	80-0059	6545-01-532-4962
OD Green	80-0060	----

SPECIFICATIONS

H 22.75 in. x W 16.5 in. x D 11.5 in.

Weight: 29 lb 15 oz

**WARRIOR AID & LITTER KIT:**

Bag only option
also available at
NARescue.com

WALK ReSupply Kit (Item# 85-0020):

Contains WALK components for resupply - excluding nylon bag & litter



NAVAL FIRST AID BOX

Response Kit Designed to stow away perfectly in the Naval First Aid Box

KIT CAPABILITIES INCLUDE:



- First responder rapid access and organization of critical life saving equipment for rendering on scene emergency care at sea
- Designed to meet the rigorous requirements of the U. S. Naval Ships, Military Sealift Command Ships and afloat commercial vessels
- Lightweight, durable slimline pouch is securable using anti-pilferage seals and features a clear window to readily view contents
- Pouch easily fits into the standard Naval First Aid Box (*not included*) and is ideal for restocking the Naval Gun Response Aid Kit as needed

SPECIFICATIONS



H 15 in. x W 6 in. x D 4 in.
Weight: 1 lb 15 oz

ORDERING INFORMATION

OPTION	ITEM #	NSN #
Kit	80-0419	6545-01-648-2705
Pouch only	80-0357	8465-01-666-7328



NAVAL FIRST AID BOX RESPONSE KIT:



Easily fits into the standard Naval First Aid Box (*not included*) and is ideal for restocking the Naval Gun Response Aid Kit (see below) as needed

Bag only option also available at NARescue.com



NAVAL GUN Response Aid Kit

Response Aid Kit Designed for the rigorous requirements of the US Navy

KIT CAPABILITIES INCLUDE:



- Designed for the rigorous requirements of the U.S. Navy for first responder emergency care equipment on board U.S. Naval Ships, Military Sealift Command Ships & afloat commercial vessels
- Ergo dynamic, low profile casualty response system features 2 removable equipment pouches that attach via hook-&-loop to a lightweight, durable outer bag
- Two internal pouches are identical & provide rapid access & organization of critical life saving equipment needed to render on-scene care at sea
- Both outer bag & 2 inner pouches are securable using anti-pilferage seals
- Slim-line, lightweight kit stores easily on a wall and can be carried over the shoulder or bandolier-style



SPECIFICATIONS

H 15 in. x W 6 in. x D 4 in.
Weight: 1 lb 15 oz

See Above for Contents of the Naval First Aid Box Response Kit



ORDERING INFORMATION

OPTION	ITEM #	NSN #
Kit	80-0418	6545-01-648-2690
Bags only	80-0417	8465-01-666-7169

Bag only option also available at NARescue.com



NAR-5 SAR AID KIT

Specifically designed for U.S. Navy Search and Rescue Swimmers



Bag only option also available at NARescue.com

- The SAR Level A Medical Bag designed for U.S. Navy Search and Rescue (SAR) Swimmers
- Low profile & durable response bag featuring 6 internal pouches that attach via hook-and-loop allowing access to the specific pouch needed to provide on-scene medical care to the patient
- Various sized modules can be self labeled for the specific mission set; 4 smaller pouches for hemorrhage control products and 2 larger ones to accommodate airway & intravenous products
- Kit can be carried as a backpack using shoulder straps that stow away until needed

ORDERING INFORMATION NATO NSN*

OPTION	ITEM #	NSN #
Kit	80-0416	-----
Bags only	80-0399	8465-01-669-5796*

SPECIFICATIONS

H 24 in. x W 15 in. x D 7 in.
Weight: 14 lb 10 oz

NAR-5 AID KIT:

 QTY 4 C-A-T® Tourniquet	 QTY 4 ETD™ 4in. Dressing	 QTY 1 ETD™ Abdominal	 QTY 4 Gauze, Compressed	 QTY 1 CPR Mask
 QTY 4 PPE Kit	 QTY 2 1000 cc IV Solution	 QTY 1 Cyclone Pocket BVM	 QTY 4 HyFin® Vent. Twin Pk	 QTY 1 Survival Blanket
 QTY 2 SAM® Splint II	 QTY 4 Triangular Bandage	 QTY 1 Mini Maglight	 QTY 2 ACE Cervical Collar	 QTY 4 Elastic Bandage, 3 in.
 QTY 1 Trauma Shears, Lg	 QTY 1 Surgical Tape, 1 in.	 QTY 12 Flexible Fabric Bandages	 QTY 1 Sterile Burn Sheet	

Naval | USCG BOAT RESPONSE AID KITS

Provides Naval & USCG first responders with critical life saving equipment while at sea



For additional information about these products designed for emergency medical care on military & commercial afloat vessels, contact North American Rescue or visit www.NARescue.com

ORDERING INFORMATION

OPTION	ITEM #	NSN #
Naval Boat Response Aid Kit	80-0420	6545-01-648-3228
Naval Boat Response Aid Box Only	80-0354	8115-01-656-8316
USCG Boat Response Aid Kit	80-0353	6545-01-646-2623

NAVAL BOAT RESPONSE AID KIT:

 QTY 2 C-A-T® Tourniquet	 QTY 1 ETD™ 6 in. Dressing	 QTY 1 ETD™ Abdominal	 QTY 1 CPR Mask	 QTY 1 28F Naso Airway
 QTY 1 Surgical Lubricant	 QTY 1 HyFin® Vent. Twin Pk	 QTY 1 SAM® Splint II	 QTY 1 Dry Sterile Burn Dressing	 QTY 1 Elastic Bandage, 3 in.
 QTY 1 ACE Cervical Collar	 QTY 1 Blizzard Survival Blanket	 QTY 1 Fox Eye Shield w/garter	 QTY 1 Gauze Pads 4"x4" (2pk)	 QTY 1 Trauma Shears, Lg
 QTY 1 Biological PPE Kit				



USMC CLS KIT U.S. Marine Corps Combat Life Saver Kit

KIT CAPABILITIES INCLUDE:



- The Official Combat Lifesaver kit of the U.S. Marine Corps
- Constructed of durable 1000 denier Cordura® solution-dyed nylon
- Sand and water resistant zippers with silent pulls & ITW Nexus Hardware

SPECIFICATIONS

Packaged: L 12 in. x W 10 in. x D 8 in.
Open Length: 37.5 in.
Weight: 6 lb 14 oz



USMC CLS KIT:



Bag only **#ZZ-0252**
also available at NARescue.com

ORDERING INFORMATION

ITEM #	NSN #
80-0159	6545-01-571-4470

USMC CLS ReSupply (#80-0424)

Contains items for resupply in a clear resealable plastic bag - *excludes nylon bag*



KIT CAPABILITIES INCLUDE:



SQUAD RESPONDER KIT

Medical equipment Kit for the US Army Combat Life Saver (CLS)



- Designed for Army Combat Lifesavers
- Medical treatment for 2 to 4 casualties

SPECIFICATIONS

H 14 in. x W 14 in. x D 5 in.
Weight: 4 lb 14 oz

SQUAD RESPONDER KIT:



Bag only **#80-0932**
also available at NARescue.com



Squad Responder ReSupply (#80-2495)

Contains items for resupply in a clear resealable plastic bag - *excludes nylon bag*



ORDERING INFORMATION

ITEM #	NSN #
♦ 85-1764	----



Also available at NARescue.com: **U.S. Army CLS ReSupply - #85-0197** (Combat Gauze) | **#85-0196** (Chito Gauze)



- Constructed of 500D Cordura® solution dyed nylon
- Sand- and water-resistant, self repairing zippers with silent 550 pull cords
- Separate lid compartment for the R.E.E.L. Non-Rigid Litter
- Six segregated removable pouches secured by hook and loop:
 - 4 Internal Bleed Pouches with Red pull handle
 - 1 Airway Pouch with Blue pull handle
 - 1 Accessories Pouch with Coyote pull handle
- 2 internal lid pouches to secure additional medical products
- 2 external side pouches on either side of the kit for additional medical products
- 2 external straps with dual adjust MILSPEC IR reduced buckles to secure the kit in a vehicle or boat

SPECIFICATIONS

Case closed:

H 17 in. x W 16.5 in. x D 10 in.




Weight: 18 lb

ORDERING INFORMATION

ITEM #	NSN #
♦ 80-0889	----

Bag only #80-0888also available at NARescue.com

Lightweight platform to rapidly access and organize all your essential trauma management equipment

EXPEDITIONARY CASUALTY RESPONSE KIT:



RANGE TRAUMA KIT

Addresses injuries common to shooting and training ranges

KIT CAPABILITIES INCLUDE:



- Essential equipment for medical and non-medical personnel to address injuries common to shooting and training ranges
- Treatment capabilities include penetrating injuries, chest wounds, burns, eye injuries, minor lacerations, sprains and fractures

ORDERING INFORMATION

*NATO NSN#

HEMOSTATIC	COLOR	ITEM #	NSN #
Combat Gauze	Orange	80-0299	----
Chito Gauze	Orange	80-0213	6545-01-610-5411*

SPECIFICATIONS

H 7.75 in. x W 13 in. x D 9 in
Weight: 5 lb 10 oz



RANGE TRAUMA KIT:



Bag only option
also available at
NARescue.com

RANGE TRAUMA AID KIT

in Hard Case



- Water-tight case maintains the integrity of your medical products when you need them most
- Kit contents arranged in a modular system featuring rugged lid storage organizer and removable nylon pouches, in the main compartment, with easy-to-read treatment labels

ORDERING INFO

COLOR	ITEM #
Orange	85-0889

SPECIFICATIONS

H 9.5 in. x W 13.5 in. x D 6 in.
Weight: 8 lb 12 oz

RANGE TRAUMA AID KIT HARD CASE:



KIT CAPABILITIES INCLUDE:



K-9 HANDLER IFAK

Designed especially for K-9 and Handler



K-9 HANDLER IFAK:



Bag only option also available at NARescue.com



- Designed for immediate lifesaving treatment to a wounded handler, operator or K-9
- Treatment capabilities range from major trauma typical to operating in tactical or patrol environments to K-9 emergencies

SPECIFICATIONS

H 7 in. x W 4 in. x D 2.75 in

Weight: 1 lb 5 oz

ORDERING INFORMATION

COLOR	HEMOSTATIC	ITEM #
Multi-Cam	----	80-0441
Coyote	----	80-0440
OD Green	----	80-0442
Black	----	80-0439
Multi-Cam	Combat Gauze	80-0445
Coyote	Combat Gauze	80-0444
OD Green	Combat Gauze	80-0446
Black	Combat Gauze	80-0443

KIT CAPABILITIES INCLUDE:



K-9 TACTICAL FIELD KIT

Kit is designed specifically to treat wounded K-9s



K-9 TACTICAL FIELD KIT:

- Designed for the handler or medic to deliver lifesaving treatment to a wounded K-9 at the point of injury
- K-9 treatment capabilities from major trauma typical to operating in tactical or patrol environments to more common K-9 emergencies
- Provides the essential equipment to render care to the injured K9 when veterinary care is not immediately available

SPECIFICATIONS

H 7.75 in. x W 13 in. x D 9 in.

Weight: 5 lb 2 oz



ORDERING INFORMATION

COLOR	HEMOSTATIC	ITEM #
Black	Combat Gauze	80-0300
Coyote	Combat Gauze	80-0301
OD Green	Combat Gauze	80-0303
3 options online	ChitoGauze	See NARescue.com

Also available: **K-9 Litter (#60-0076)**, See p.58

INTRO

M-A-R-C-H

KITS

TAC-EVAC

MED-SURG

RESCUE

TREATMENT SUPPORT

TRAINING

INDEX



MULTI-MISSION TRAUMA PACK

Customizable modular solution ideal for personalized gear organization spanning multiple Mission profiles

- A 3200-cubic inch backpack with (4) compartments
- Multi-mission belt attached for different configurations
- Internal Molle straps covered in hook and loop
- Detachable compartments
- Converts to trauma sheet for vehicles, boats, or aircrafts
- Internal channeling for electrical communications wiring
- Dedicated loops, providing rappelling or lowering of the bag

SPECIFICATIONS

H 25.5 in. x W 11 in. x D 9 in.

Weight: 7 lb 8 oz

ORDERING INFORMATION

COLOR	ITEM #
Black	80-0950
Coyote	80-0949



Tiered Trauma Medical Bag **TTMB**



INTRO

M.A.R.C.H

KITS

TAC-EVAC

MED-SURG

RESCUE

TREATMENT SUPPORT

TRAINING

INDEX



A. ReSupply Bag
Item# 80-0271
(MTC only)




B. NAR-4 Aid Bag
Item# 80-0241 (MTC)

Other colors available
see NARescue.com




C. NAR-4 Chest Pouch
Item# 80-0248 (MTC)
Other colors available



D. Leg Rig Conversion Panel
Item# 80-0272 (MTC only)



E. MOLLE Directional Conversion Panel
Item# 80-0273
(MTC only)



- 5-piece modular ensemble targeting the military medic's "Soldier-as-a-System" requisite
- Includes low-profile backpack-style aid bag, ample re-supply bag, medical chest pouch, leg rig conversion panel & directional MOLLE converter
- Interchangeable internal pouches with removable labels compatible with the NAR-4 system or any hook and loop attachment
- Leg rig converter & directional change MOLLE panel allow the chest pouch to be mounted & carried using a drop-down, leg-attachment style platform

SPECIFICATIONS

- (A) ReSupply Bag: H 20 in. x W 15 in. x D 5 in.
(B) NAR-4 Aid Bag: H 20 in. x W 15 in. x D 5 in.
(C) NAR-4 Chest Pouch: H 7.75 in. x W 9.75 in. x D 4.5 in.
(D) Leg Rig Conversion Panel: L 14 in. x W 6 in.
(E) MOLLE Conversion Panel: L 7.5 in. x W 6 in.

ORDERING INFORMATION

BAG	COLOR	ITEM #
TTMB Complete System	Multi-Cam	80-0274
ReSupply Bag (A)	Multi-Cam	80-0271
NAR-4 Aid Bag (B)	Multi-Cam	80-0241
NAR-4 Chest Pouch (C)	Multi-Cam	80-0248
Leg Rig Panel (D)	Multi-Cam	80-0272
MOLLE Conv. Panel (E)	Multi-Cam	80-0273

NOTE: Bags only all contents sold separately



OCHO BAGS

Low Visibility | Customizable Modular Solutions

Featuring 2 interchangeable mesh aid pouches*, hook and loop fields and hidden compartments – all within purposely designed exteriors that resemble standard off-the-shelf bags to minimize the threat of theft

OCHO BRIEFCASE

SPECIFICATIONS

H 12 in. x W 17 in. x D 4.5 in.
Weight: 2 lb 10 oz

COLOR	ITEM #
Brown	80-0258



OCHO BACKPACK

SPECIFICATIONS

H 19 in. x W 14 in. x D 9 in.
Weight: 3 lb 14 oz

COLOR	ITEM #
Blue/Gray	80-0256
Red/Gray	80-0338
Green/Gray	80-0339
Gray/Black	80-0340



OCHO IFAK/ACCESSORY BAG



COLOR	ITEM #
Blue/Gray	80-0371
Gray/Black	80-0370

*NOTE: OCHO IFAK/
Accessory Bag does
NOT include Mesh
Aid Pouches

SPECIFICATIONS

H 7 in. x W 2.25 in. x D 3.25 in. Weight: 2.75 oz



OCHO SHOULDER BAG

SPECIFICATIONS

H 12 in. x W 14 in. x D 6 in.
Weight: 1 lb 7 oz

COLOR	ITEM #
Blue/Gray	80-0257

OCHO QUICK DETACH GEAR ORGANIZERS



OCHO Shock Cord Panel

SPECIFICATIONS

H 6 in. x W 3 in. x D 0.25 in.
Weight: 116 oz

COLOR	ITEM #
Black	80-0260



OCHO Small Loop Panel

SPECIFICATIONS

H 7 in. x W 5.5 in. x D 0.25 in.
Weight: 183 oz

COLOR	ITEM #
Black	80-0261

Offers additional storage - compatible with any storage solution that has a loop panel attached. Rapid detach capability via red quick-pull tab.



OCHO Mesh Aid Pouch

SPECIFICATIONS

H 7 in. x W 4.5 in. x D 3 in.
Weight: 187 oz

COLOR	ITEM #
Black/Foliage	80-0259



MEDIC TRAUMA SHEET

Visible and easy access to your medical equipment

- Versatile, rugged, lightweight equipment access platform for managing casualties during TACEVAC operations
- Oversized pockets & large grommet attachment points facilitate rapid set-up and equipment retrieval
- Quick attach and release straps allow for secure attachment to the skin of an aircraft or vehicle

SPECIFICATIONS

Open: H 41.5 in. x W 215 in. x D 2.5 in.
Folded: H 20.25 in. x W 215 in. x D 5 in.
Weight: 4 lb 12 oz

Bag only. Contents are sold separately

ORDERING INFORMATION

COLOR	ITEM #	NSN #
Coyote	80-0069	6515-01-540-7603
OD Green	80-0068	6515-01-540-7585
Black	80-0067	6515-01-540-7612
Digital Camo	80-0094	6515-01-540-7579





Horizontal & Vertical Mounting capability available with the 4.5 x 6.5 in. Panel



- Allows the operator the fastest access to their kit during the critical moments after wounding
- Allows any Molle/PALS kit to be detachable
- Gives user faster access to Molle/PALS equipped kits
- Mil-Spec hook & loop and heavy duty snaps to ensure kit security
- Designed to be compatible with all of NAR IFAK's, including a **NEW horizontal mount option** in the medium version (W 4.5 x H 6.5 in.)

Small QUICK DETACH PANEL (3 in. x 6 in.)

- Utilizes only 3 rows by 2 columns on user vest/bag, vertical mount configuration
- Designed specifically for NAR's most compact M-FAK kits

SPECIFICATIONS

W 3 in. x H 6 in. x D 0.5 in.
Weight: 2.8 oz

ORDERING INFORMATION

SIZE	COLOR	ITEM #
3 x 6 in.	Black	80-0893
3 x 6 in.	Coyote	80-0894



Medium DUAL DIRECTION QUICK DETACH PANEL (4.5 in. x 6.5 in.)

- Utilizes 3 rows by 3 columns vertical mount or 2 rows by 4 columns horizontal mount
- Designed for compatibility with all NAR IFAKs such as horizontal mount for M-FAK and vertical or horizontal mount of the TORK

SPECIFICATIONS

W 4.5 in. x H 6.5 in. x D 0.5 in.
Weight: 4.5 oz

ORDERING INFORMATION

SIZE	COLOR	ITEM #
4.5 x 6.5 in.	Black	80-0895
4.5 x 6.5 in.	Coyote	80-0896



Large QUICK DETACH PANEL (9 in. x 6.5 in.)

- Utilizes 3 rows by 7 columns on user vest/bag, vertical mount configuration
- Designed specifically for use with NAR 4 Chest Pouch & Emergency Evacuation Litter (EEL)

SPECIFICATIONS

W 9 in. x H 6.5 in. x D 0.5 in..
Weight: 5.6 oz

ORDERING INFORMATION

SIZE	COLOR	ITEM #
9 x 6.5 in.	Black	80-0897
9 x 6.5 in.	Coyote	80-0898





BLEEDING CONTROL RESPONSE

Individual Kits • 8-Packs • Wall Stations

D-BCRK *Department of Defense* Individual BLEEDING CONTROL RESPONSE KITS

- Individual Bleeding Control Response Kit designed to provide essential equipment that empowers DOD Forces and Federal Employee responders in public areas to take immediate action in stopping life threatening bleeding.
- Contained in a ruggedized, nylon carry bag or vacuum-sealed pouch
- provides intuitive and easy-to-use tools that are proven to help to save lives
- Available configurations: Basic, Intermediate, & Advanced

SPECIFICATIONS

Nylon Bag: H 7 in. x W 5 in. x D 2.75 in.

Vacuum Sealed: H 7 in. x W 5 in. x D 2.75 in.

ORDERING INFORMATION

LEVEL	BAG	WEIGHT	ITEM #
◆ Basic	Red Nylon	15 oz	85-2057
◆ Intermediate	Red Nylon	16.5 oz	85-2058
◆ Advanced	Red Nylon	17.4 oz	85-2059
◆ Basic	Vacuum Sealed	11.7 oz	85-1749
◆ Intermediate	Vacuum Sealed	13.25 oz	85-2066
◆ Advanced	Vacuum Sealed	14 oz	85-2069

Nylon Bag Option:

- Easy to restock in cases of use
- Convenient to carry and easy to open in stressful situations
- Opens flat for ease of use and immediate identification of enclosed devices
- Easily recognizable red nylon bag with attached PVC Bleeding Control Patch

Vacuum Sealed Option:

- Cost-effective alternative to Nylon
- Compact and rugged Vacuum-Sealed pouch
- Large Identifying Card Vacuum-Sealed inside the pouch for ease of identification



Every D-BCRK INDIVIDUAL BLEEDING CONTROL RESPONSE KIT:



PLUS these Individual Components based upon configuration of choice:



D-BCRK 8-Pack Department of Defense Bleeding Control Response Kits (Nylon or Vacuum-Sealed)



- Eight (8) Nylon or Vacuum-sealed Individual Bleeding Control Response Kits inside an easy-to-open nylon carrying case
- Available in Basic, Intermediate, and Advanced configurations

SPECIFICATIONS

Each 8-Pack: H 14 in. x W 11 in. x D 8 in.
Individual Bleeding Control Kit:
H 7 in. x W 5 in. x D 2.5 in.



ORDERING INFORMATION

LEVEL	BAG	WEIGHT	ITEM #
◆ Basic	Red Nylon	8.5 lb	85-2072
◆ Intermediate	Red Nylon	9.3 lb	85-2073
◆ Advanced	Red Nylon	9.7 lb	85-2074
◆ Basic	Vacuum Sealed	8.5 lb	85-1750
◆ Intermediate	Vacuum Sealed	9.3 lb	85-2067
◆ Advanced	Vacuum Sealed	9.7 lb	85-2070

D-BCRK 8-PACK - Vac



D-BCRK 8-PACK - Nylon



See page left for contents of the Individual Kits

D-BCRK Wall Station Department of Defense Bleeding Control Response Kits (Nylon or Vacuum-Sealed)

- Wall Station contains One D-BCRK 8-Pack each containing Eight Nylon or Vacuum-sealed Individual Bleeding Control Response Kits
- Available in Basic, Intermediate, and Advanced configurations
- Includes a durable, clear polycarbonate cabinet with rounded corners provides both physical and environmental protection for medical equipment
- Includes a wire shelf to hold the 8-pack & an external eyelet that allows for a secondary safety seal to easily identify if the station has been opened or tampered with
- External mounting tabs and hardware enable easy installation on any wall
- Clear design allows visibility at any angle


SPECIFICATIONS

Clear Wall Case: H 18 in. x W 15.5 in. x D 10 in.
Each 8-Pack: H 14 in. x W 11 in. x D 8 in.
Individual Bleeding Control Kit:
H 7 in. x W 5 in. x D 2.5 in.

ORDERING INFORMATION

LEVEL	BAG	WEIGHT	ITEM #
◆ Basic	Red Nylon	23.5 lb	85-2075
◆ Intermediate	Red Nylon	24.3 lb	85-2076
◆ Advanced	Red Nylon	24.7 lb	85-2077
◆ Basic	Vacuum Sealed	23.5 lb	85-1751
◆ Intermediate	Vacuum Sealed	24.3 lb	85-2068
◆ Advanced	Vacuum Sealed	24.7 lb	85-2071

D-BCRK CLEAR WALL STATION:



Also available at **NARescue.com:**
D-BCRK Wall Stations in Semi-Recessed Metal Cabinets



CUSTOMIZABLE SOLUTIONS

550 Cord Silent Zipper Pulls

Velcro® Brand Hook & Loop

Sand & Water Resistant, Self-Repairing Zippers

500 Denier CORDURA® Nylon

Mil Spec Nylon Webbing

PALS-MOLLE-style attachment points

Nitrile Trauma Gloves

Elastic Storage Loops

C-A-T® Tourniquet

Hemostatic Dressing

ARS® Needle Decompression Kit

28F Naso Airway with lubricant

6 inch ETD Bandage

HyFin® Vent Chest Seal

...fulfilling YOUR complete Medical & Material logistics requirements!

TRAINING

INDEX



North American Rescue® is **THE market leader** for designing, producing, integrating and assembling tactical field medical sets and kits **supporting the services of the Department of Defense, federal agencies** as well as **state and local law enforcement**.

Over the past several years, our ability to provide rapid response products and services has aided in the delivery of critical battlefield medical care when and where it's needed most with an **end result of saving more lives**.



A comprehensive list of supply chain services affords a unique opportunity to North American Rescue's customers in the procurement, assembly, storage and management of their **complete medical and material logistics** requirements. NAR maintains a 75,000 square foot warehouse facility which consists of 50,000 square feet for receiving and distribution staging and an additional 14,000 square feet solely dedicated to the assembly room where sets, kits and **custom medical solutions are built to order**.

Highly qualified logistics and warehouse staff ensure that orders, inventory and distribution processes **maintain the superior quality standards** associated with the North American Rescue brand. Our integral partnership with the South Carolina Vocational Rehabilitation **NISH Partner** is a pillar in the successful execution of these processes. ■

IAC-EVAC

MED-SURG

RESCUE

TREATMENT SUPPORT

TRAINING

INDEX



TACEVAC

CASUALTY EVACUATION

- FRAMED LITTERS • ACCESSORIES
- NON-RIGID LITTERS
- SPINE BOARDS

see pages 54-59



Rigid Litters

TACEVAC



TALON® II 90C

Collapsible Handle Litter

Advanced, innovative combat evacuation system features portable quad-fold design, retractable handles, six IV-attachment points and two patient securing straps

- Standard 90 in. NATO-compatible evacuation platform
- Developed to U.S. Army requirements
- Certified Safe-to-Fly U.S. Army Rotary Wing

- Retractable, ergonomic handles and automatic locking hinges
- Slip-resistant, chemical-resistant rip-stop fabric tested & approved by Federal & private labs using live warfare agents

SPECIFICATIONS



Folded w/handles collapsed: H 20.5 in. x W 8.5 in. x D 9.5 in.
 Open w/handles collapsed: L 78 in. x W 22.5 in. x D 6 in.
 Open w/handles extended: L 90 in. x W 22.5 in. x D 6 in.
 Weight: 15 lb 6 oz (6.97 kg)
 Static Payload: 1,800 lb (816 kg)
 Ground Clearance: 1.5 in.
 U.S. Patents: 6,842,923 & 5,598,592, Patent(s) Pending

ORDERING INFORMATION

ITEM#	NSN#
60-0002	6530-01-504-9051

RAVEN® 90C

Bi-fold Litter

- The official litter of the U.S. Navy and certified safe-to-fly by the U.S. Air Force
- Bifold design for reliable long-term storage
- Ergonomically designed handles
- Slip-resistant, chemical-resistant, rip-stop fabric tested & approved by Federal & private labs using live warfare agents
- Includes 6 IV-attachment points and 2 patient-securing straps

SPECIFICATIONS



Folded: H 45.5 in. x W 7.5 in. x D 8 in.
 Open: 90 in. x W 22.5 in. x D 6 in.
 Weight: 16 lb 8 oz (7.5 kg)
 Static Payload: 2,000 lb (909 kg)
 Ground Clearance: 1.5 in.
 U.S. Patents: 6,842,923 & 5,598,592, Patent(s) Pending

ORDERING INFORMATION

ITEM#	NSN#
60-0001	6530-01-432-5114

T3™ LITTER

Lightweight, Assault

- Ultra light, ultra strong litter weighing less than 12 lb
- Specifically designed for Special Operations Forces
- Provides maximum maneuverability during TACEVAC in urban terrains
- Overall length of 81 in. enables floor loading on a UH-60 helicopter
- 2 straps for securing patients
- Ergonomically designed handles
- Quad-fold design for maximum cube space storage

SPECIFICATIONS



Folded: H 21 in. x W 8.75 in. x D 7.25 in.
 Open: L 81 in. x W 22.5 in. x D 6 in.
 Weight: 11 lb 13 oz (5.34 kg)
 Static Payload: 1,200 lb (544.3 kg)
 Ground Clearance: 1.5 in.

ORDERING INFORMATION

ITEM#
60-0022



Wheeled LITTER CARRIER

- Lightest wheeled litter cart on the market for single person transport to and from MedEVAC vehicles
- Compatible with the Talon® II 90C, Raven® 90C, MedEvac & T3, in addition to all NATO litters
- Quick and easy to deploy or disassemble for storage

- Locking clamp system keeps the litter firmly in place and allows for different litter widths
- Performs as litter stand in the field due to its unique wheel locking system
- Also available with optional Carrying Case

SPECIFICATIONS

Deployed:
L 56.29 in. x W 27.5 in. x H 36.2 in.
Folded:
L 42.1 in. x W 27.5 in. x H 9.84 in.
Case: L 48 in. x W 29.5 in. x H 12.5 in.
Weight:
Carrier only: 38 lb 8 oz (17.5 kg)
Carrier with Case: 88 lb (39.9 kg)
Payload: 550 lb (250 kg)



ORDERING INFORMATION

*NATO NSN

OPTION	ITEM #	NSN #
Carrier with Case	60-0063	-----
Carrier only	60-0041	6530-01-663-2308*

LITTER STANDS

33 in. Collapsible


- Set of 2 lightweight, compact litter stanchions designed to collapse for convenient storage in a standard issue military carrying case
- Quickly locks open
- Allows the capacity to store a quad-folded Talon®II 90C (sold separately)

SPECIFICATIONS

H 17.95 in. x W 25.5 in. x D 6.92 in.
Extended: 33 in.
Base Width: 25.5 in.
Weight: 10 lb 12 oz

ORDERING INFORMATION

Unit	ITEM#
Set of 2	60-0064



LITTER STANDS

33 in. Standard

- Set of 2 lightweight, foldable stanchions for standard sized litters
- Constructed of a lightweight aluminum frame that conveniently folds flat for ease in transport & storage
- Quickly locks open to safely secure the litter in place

SPECIFICATIONS

H 34 in. x W 25.5 in. x D 6.92 in.
Extended: 33 in.
Base Width: 23.5 in.
Weight: 8 lb 8 oz

ORDERING INFORMATION

UNIT	ITEM#	NSN#
Set of 2	60-0014	6530-01-523-3314



BACKREST TALON® or RAVEN®

- Allows patients on a Talon® 90C or Raven® 90C litter to be placed at 90° or 45° angles
- Constructed of a lightweight aluminum frame with Nomex/Kevlar-blend fabric covering

SPECIFICATIONS

H 31.5 in. x W 29.25 in. x D 1.6 in.
Weight: 2 lb 13 oz

ORDERING INFORMATION

ITEM#	NSN#
60-0027	6530-01-594-7444



IV POLE

Adjustable

Exclusively designed to fit the pre-molded holes of our framed litters

SPECIFICATIONS

H 19-33 in. x W 8.75 in. x D 0.6 in.
Weight: 4.8 oz

ORDERING INFORMATION

ITEM#	NSN#
60-0008	6530-01-504-9054





RIG SERIES CARRIER *Assault Litter Carrier*

- Enhanced carrier with additional strap over the top of the litter to hold it in place during drastic movement and FRIES (Fast Rope Insertion Extraction System) operations
- Upgraded shoulder straps provide comfort for long wear
- Includes pocket in the top flap for storage and three additional handles facilitate packing and unloading of the litter
- Accommodates all NAR quad-fold litter platforms

- Equipment attaching webbing (MOLLE) along the outside of the carrier allows attachment of additional equipment

SPECIFICATIONS

H 22 in. x W 9 in. x D 9 in.

Weight: 3 lb (1.36 kg)

ORDERING INFORMATION ^{*NATO NSN}

COLOR	ITEM #	NSN #
Multi-Cam	60-0057	6530-01-647-8571*
Coyote	60-0056	6530-01-643-9138*



SPECIFICATIONS

Collapsed: H 20 in. x W 9 in. x D 2 in.

Expanded: H 20 in. x W 9 in. x D 10 in.

Weight: 1 lb 7 oz

TALON® CARRIER *Assault Litter Carrier*

- Transport any NAR quad-fold litter (Talon, MedEvac & T3)
- MOLLE-style attachment points
- Includes thin back plates, drain hole, non-slip, adjustable shoulder straps & top-carrying handle

ORDERING INFORMATION ^{*NATO NSN}

COLOR	ITEM #	NSN #
Coyote	60-0011	6530-01-531-6854
OD Green	60-0010	6530-01-504-9055
Black	60-0009	6530-01-531-6907
Digital Camo	60-0025	6530-01-643-9134*

RETENTION STRAP *Litter Support Bar*

Designed to prevent the litter support bar from collapsing

ORDERING INFO

ITEM#

60-0026

SPECS

L 104 in. x W 1 in.

Weight: 2 oz



LITTER BUCKLE *Replacement Set*

SPECIFICATIONS

L 3.75 in. x W 2.38 in. x D 0.5 in.

Weight: 1 oz

ORDERING INFO

ITEM#	NSN#
ZZ-0032	6530-01-435-7337



EQUIPMENT BAG *Store a casualty's personal items*

- Roomy nylon bag with reinforced shoulder straps and info tag to bag & tag casualty's personal items during evacuation
- Lightweight & vacuum-packed for low cube, easy transport

SPECIFICATIONS

H 5 in. x W 5 in. x D 2.5 in.

Weight: 6.25 oz

ORDERING INFO

ITEM#	NSN#
80-0132	8465-01-575-6644





MEDEVAC LITTER



- The ultimate, cost efficient casualty evacuation stretcher for First Responders operating in Active Shooter, Mass Casualty and Disaster incidents
- Compact, collapsible, quad-fold litter is easy to deploy and carry
- Two patient securing buckles, six IV-attachment points, ergonomically designed handles & patented joint configuration with locking leg crossbar

SPECIFICATIONS



Folded: H 21.5 in. x W 8.5 in. x D 7.5 in.
Open: H 84.5 in. x W 22.5 in. x D 4 in.
Weight: 13 lb 8 oz
Static Payload: 1,200 lb (544.3 kg)
Ground Clearance: 1.5 in.

ORDERING INFORMATION

BED COLOR	FRAME COLOR	ITEM #
Black	OD Green	60-0051
Orange	Black	60-0052

OSL™ LITTER

OverSized Litter

- Safely evacuate larger patients with 30% larger surface area compared to standard-sized evacuation platforms
- Bi-fold configuration enables lower cube storage capability
- Includes 6 litter stirrups, 6 padded carrying loops & 3 patient securing straps for firmer patient support as well as better weight distribution on the floor
- Flame-retardant cover & attached storage pouch made of a durable Nomex/Kevlar-blend
- Certified safe-to-fly by the U.S. Air Force

SPECIFICATIONS



Folded: H 45 in. x W 9 in. x D 8 in.
Open: L 90 in. x W 31 in. x D 6 in.
Weight: 22 lb 6 oz (10.2 kg)
Static Payload: 2,000 lb (909 kg)
Ground Clearance: 1.5 in.
U.S. Patents: 5,598,592 & Patent(s) Pending



ORDERING INFORMATION

ITEM#	NSN#
60-0023	6530-01-552-7708

OSL™ LITTER STANDS

for the OverSized Litter™

- Set of 3 lightweight aluminum stanchions accommodate our oversized OSL™ Litter
- Quickly opens and securely locks the litter in place
- Easily folds for transport & storage

SPECIFICATIONS



Folded: H 34 in. x W 25 in. x D 2.5 in.
Extended Height: 32 in.
Base Width: 23.5 in.
Weight: 28 lb 8 oz (12.93 kg)

ORDERING INFO

UNIT	ITEM#
Set of 3	60-0029



OSL™ BACKREST

for the OverSized Litter™

- Position patient at 45° & 90° angles for maximized treatment options & increased comfort
- Constructed of a lightweight aluminum frame with durable Nomex/Kevlar-blend fabric cover
- Compatible with our OSL™ litter

SPECIFICATIONS



H 31.5 in. x W 30.25 in. x D 1.6 in.
Weight: 2 lb 15 oz (1.33 kg)

ORDERING INFO

ITEM#
60-0024





TACEVAC | Non-Rigid Litters

TED Tactical Extrication Device

- Designed to allow two people to carry a casualty to safety
- Lighter, more compact non-rigid evacuation platform designed specially for use in high-threat environments & limited spaces
- Constructed of heavy duty vinyl with 4 sturdy carrying handles
- Packaged in a sturdy, transport-friendly bag with adjustable shoulder strap

SPECIFICATIONS



Folded: H 5 in. x W 10.25 in. x D 3.5 in.
Open: L 43.5 in. x W 19 in.
Weight: 1 lb 8 oz (0.68 kg)
Payload: 600 lb (272 kg)

ORDERING INFORMATION

COLOR	ITEM #	NSN #
OD Green	60-0006	6530-01-492-4114
Black	60-0005	6530-01-531-6808



R.E.E.L. Rapid Expeditionary Evacuation Litter

- Cost-effective, non-rigid casualty evacuation platform designed for quick and effective casualty evacuation
- Rugged, lightweight, nylon construction
- 8 carry handles accommodate mission teams ranging from 2 to 8 rescuers
- Packaged in a compact, small cube, vacuum sealed pouch

SPECIFICATIONS

Packaged - L 13 in. x W 10 in. x D 0.75 in.
Litter (Open) - L 78 in. x W 30 in. | Weight: 1 lb 3 oz
Carrying Capacity: 700 lb

ORDERING INFORMATION

ITEM#

60-0075

NAR K-9 LITTER & Carrying Strap Accessories

- Lightweight, compact litter for the movement of an injured canine
- Litter can be carried in cargo pocket or pouch
- Leg cut outs to safely carry the canine in the litter
- Six large and four small carry handles that allow for one to four person carries
- Optional K-9 Carry Strap & 2 Carabiners for over the shoulder carry of injured K-9 also available



SPECIFICATIONS

Litter: W 35 in. x H 48 in.
Adjustable Carry Strap:
L 9 to 54 in. x W 2 in.
Carabiner open gate
clearance: 27 mm
Packaged:
W 5.5 in. x H 10 in. x D 2.5 in.
Weight: 1.38 lb

ORDERING INFORMATION

OPTIONS	ITEM#
♦ K-9 Litter Only	60-0074W
♦ Carry Strap & Carabiners	ZZ-1004
♦ K-9 Litter with Carry Strap & Carabiners	60-0076

STINGRAY® Poleless Litter

- Ideal for decontamination roller systems & moving wet patients
- Non-skid, flame-retardant monofilament polypropylene fabric that is resistant to mildew, acids, alkalis and penetrating liquids
- Tested by Federal & private labs using live warfare agents
- Includes six heavy-duty carrying handles & 2 casualty securing straps
- Easy to store and clean

SPECIFICATIONS

Folded: H 20 in. x W 12 in. x D 2 in.
Open: L 78 in. x W 24 in.
Weight: 13 oz (0.37 kg)
Payload: 1,000 lb (453.6 kg)



ORDERING INFO

ITEM #	NSN #
60-0004	6530-01-492-4079



RESCUE SLEEVE II

- The industry's ultimate patient securing system refined for enhanced stability & ease of use
- Ideal for rescues requiring vertical and/or horizontal hoisting capabilities

SPECIFICATIONS

Folded: L 25.5 in. x W 13 in. x D 8 in.
Open: L 73 in. x W 18 in. x D 2 in.
Weight: 16 lb (7.26 kg)
Payload: 1,000 lb (453.6 kg)



ORDERING INFORMATION

COLOR	ITEM #	NSN #
Black	50-0011	6530-01-531-6649
Orange	50-0012	6530-01-477-8525

NOTE: Spine board sold separately. Hoisting straps & carabiners not included

SPINE BOARD

- Thin & lightweight with tapered design
- Surface texture to prevent slippage & 14 raised, oversized hand-holds
- Black spine board includes ten rapid attachment pins

SPECIFICATIONS

L 72 in. x W 16 in. x D 1.75 in.
Weight: 14 lb (6.35 kg) Payload: 1,100 lb (498.96 kg)
Float Capacity: 250 lb (113.4 kg)

ORDERING INFORMATION

COLOR	ITEM #	NSN #
OD Green	50-0013	6530-01-490-2487
Orange	50-0014	6530-01-531-6771
Black w/Pins	50-0021	----



NAR CASUALTY IMMOBILIZATION SYSTEM

- Constructed of black polypropylene webbing with hook & loop closures
- Featuring 6 attachment points: 2 "suspender" type, over-the-shoulder straps & 4 adjustable torso/limb straps with adjustable center strap
- No choking "V" configuration by the patient's neck



ORDERING INFO

ITEM #	NSN #
50-0031	6530-01-586-2411

SPECIFICATIONS

Pkg: L 16.5 in. x W 4 in. x D 2.75 in.
Weight: 1 lb 1 oz



NAR HEAD IMMOBILIZER

- Quick, effective head block constructed of soft, closed cell foam with hook & loop closures and padded straps
- Ears and head are accessible without removing the device



SPECIFICATIONS

H 10 in. x W 16 in. x D 5.75 in.
Weight: 2 lb 3 oz



ORDERING INFO

ITEM #	NSN #
50-0020	6530-01-619-1777

NAR CASUALTY RESTRAINT STRAPS

Heavy-duty nylon webbing & convenient seat belt buckle utility

- A. 2 piece Plastic Buckle**
L 5 ft x W 2 in. • Weight: 4 oz
Black: 50-0023 • Orange: 50-0024
- B. 2 piece Plastic Buckle/Speed Clip**
L 5 ft x W 2 in. • Weight: 7.8 oz
Black: 50-0026 • Orange: 50-0025
- C. 2 piece Metal Buckle**
L 5 ft x W 2 in. • Weight: 8 oz
Black: 50-0027 • Orange: 50-0028
- D. 2 piece Metal Buckle/Speed Clip**
L 5 ft x W 2 in. • Weight: 11.7 oz
Black: 50-0029 • Orange: 50-0030



CERVICAL COLLAR ACE®

- Adjustable neck support with foam lining
- Extended back to ensure complete
- Adjusts to fit 16 adult neck sizes

SPECIFICATIONS

H 6.75 in. x W 20.5 in. x D 0.75 in.
Weight: 5.2 oz

ORDERING INFORMATION

COLOR	ITEM #	NSN #
OD Green	50-0010	6515-01-541-8147
White	50-1007	6515-01-650-0036





MED•SURG

MEDICAL | SURGICAL PRODUCTS

- DIAGNOSTICS
- FIELD HOSPITAL
- SURGICAL SUPPLIES

see pages 62-73





HANDHELD WIRELESS ULTRASOUND SCANNERS

from Clarius

• Wireless Freedom

Utilizes secure point-to-point wireless network to connect to the app on the device or devices you choose for imaging. No Internet access required to operate. No cables to manage...

• Point-and-Shoot Ultrasound

High resolution ultrasound images with automated gain and frequency settings designed to be almost as easy as using the camera on your smart phone. No adjustments necessary unless you want to.

• Full Metal Jacket. Waterproof

Designed to withstand unpredictable environments, Sealed and built with a durable magnesium shell that can be completely immersed in liquid

• Rechargeable Batteries

Rechargeable batteries that will last for more than 45 minutes in active scanning & up to 7 days idle

• Images on Virtually Any Smart Device

Operates with most iOS or Android smart devices. Simply download the free Clarius App and you can start scanning with the touch of a button. Clarius can image on multiple devices simultaneously

Every Clarius Scanner includes:

- 2x Rechargeable batteries
- 1x Battery Charger
- 1x Power Cable with Global AC Adapter
- 1x Quickstart Guide



C3 Premium Ultrasound Scanner

- Convex | Virtual Phased Array for scanning all parts of the patient's torso, including the heart

ITEM

20-0089

SPECIFICATIONS

L 6.6 in. x W 3.9 in. x D 16 in.
Weight: 1.21 lb



L7 Premium Ultrasound Scanner

- Linear Array ideal for guiding procedures and to see fine detail

ITEM

20-0090

SPECIFICATIONS

L 6.6 in. x W 3.9 in. x D 16 in.
Weight: 1.21 lb



C7 VET Premium Ultrasound Scanner

- Micro-Convex probe designed for frequencies of 3-10MHz and a depth of 2-15cm for scanning an animal's torso

ITEM

20-0091

SPECIFICATIONS

L 6.6 in. x W 3.9 in. x D 16 in.
Weight: 1.21 lb



For complete specifications and more information on Clarius Scanners please visit www.NARescue.com



VSCAN™ with DUAL PROBE

Pocket-sized Ultrasound Tool from GE

- Directional Color Flow, Voice Recorder, Auto Image Optimization
- Intuitive thumb-wheel navigation with system presets helps speed diagnostic decisions and enables efficient patient triage
- Access to Vscan™ Web Portal / Vscan™ Gateway Software for PC
- Storage capacity: 4GB upgradable to 32GB on microSD or microSDHC cards
- Data storage formats: .jpg for still frames | .mpg for loops | .wav for voice
- **Clinical applications:** Cardiac, Abdominal, Renal, Obstetrics and gynecology, Urology, Fetal, Evaluation of presence of fluid, Peripheral vascular, Procedure guidance for arterial and venous vessels, Imaging guidance for needle/catheter placement, Thoracic/lung, Thyroid and other small organs, Long bone, Hip and knee joints, Pediatric

Kit Contents:

- 1 x GE Vscan™ with Dual Probe Pocket-sized Ultrasound Tool
- 1 x Docking Station with Global AC Adaptor & USB cable
- 1 x Rechargeable Battery
- 1 x 4 GB microSD card
- 1 x User's Manual
- 1 x Soft Case

SPECIFICATIONS

Display unit :
H 5.32 in. x W 2.87 in. x D 1.1 in.
(135 x 73 x 28 mm)

Probe: H 5.08 in. x W 1.5 in. x D 1.1 in.
(129 x 39 x 28 mm)

Weight (unit and probe):
0.96 lb (436 g)

Display: 3.5 in.

Image Resolution: 240 x 320 pixels

ORDERING INFO

ITEM

20-0092



Masimo MightySat™ Rx

Fingertip Pulse Oximeter

- High Resolution Display - Rotating and adjustable color screen for clear real-time display of all measurements
- Plethysmographic Waveform - Provides depiction of the patient's pleth
- Touchpad - Allows for easy access of settings and device features
- Durable, Lightweight Design - Small and portable for use in clinical environments
- Flexible Pad - Easily conforms to finger shape for enhanced comfort & secure application to finger
- Signal IQ (SIQ™) - An indicator of the SpO2 measurement confidence
- Bluetooth Models are compatible with free downloadable Masimo Professional Health App

Kit Contents:

- 1 x MightySat™ Rx Device
- 1 x Lanyard
- 2 x AAA Batteries
- 1 x Carrying Case

SPECIFICATIONS

Device Dimensions :
H 2.9 in. x W 1.6 in. x D 1.2 in.
(7.4 cm x 4.1 cm x 3.0 cm)

Weight (with batteries): 2.56 oz (73 g)

Display: 3.3 cm

ORDERING INFORMATION

OPTION	ITEM #
◆ Standard	20-0093
◆ Bluetooth LE	20-0094
◆ Bluetooth LE/RRp/PVi	20-0095



◆ New Product



MedicECG[®] from Resolute Medical

- Lightest and smallest ECG designed for use in the field - ruggedized, portable, and easy to use
- Automated interpretation of ECG rhythms to alert the user of patient risks
- Results can be transferred through field communication channels
- Useful for early diagnosis of abnormal heart conditions before sending personnel to the field
- Versatile, inexpensive, and highly portable solution for base hospital ECG testing
- Exceeds industry standards for ECG accuracy and quality
- Supplied in muted mil spec colors with crushproof, waterproof, and dustproof case for easy and reliable storage

For complete specifications
and more info, please visit
www.NARescue.com



ORDERING INFORMATION

OPTION	TOUGHPAD	ITEM #
MedicECG Lite	Not Included	20-0083
MedicECG Compact	7 in.	20-0084
MedicECG Elite	10 in.	20-0085

SPECIFICATIONS

Medic ECG[®] Unit:
W 4.06 in. x H 6.26 in. x D 0.64 in.
Weight: 1 lb

EVERY MedicECG[®] KIT Contains:

- 1 x MedicECG[®] 12-lead ECG data acquisition unit
- 12 x Electrode Clips compatible with both snap and tab electrodes
- 12 x Electrodes
- 12 x Electrode Leads
- 1 x Charger
- 1 x USB Cable
- 1 x crushproof, waterproof case
- 1 x Nylon Bag



MedicECG[®] LITE
(no tablet)

SPECIFICATIONS

MedicECG[®] Lite Kit:
W 13.5 in. x H 11.25 in. x D 5.5 in.
Weight with hardcase: 8 lb
Weight in Rucksack: 1.5 lb



MedicECG[®] COMPACT
(with 7 in. FZ-M1
Toughpad[®] tablet)

SPECIFICATIONS

FZ-M1 Toughpad[®] tablet
W 7.9 in. x H 5.2 in. x D 0.7 in.
Weight: 1.2 lb

MedicECG[®] Compact Kit
W 13.5 in. x H 11.25 in. x D 5.5 in.
Weight with hardcase: 8 lb
Weight in Rucksack: 1.5 lb



MedicECG[®] ELITE
(with 10 in. FZ-G1
Toughpad[®] tablet)

SPECIFICATIONS

FZ-G1 Toughpad[®] tablet
W 10.6 in. x H 7.4 in. x D 0.8 in.
Weight: 2.5 lb

MedicECG[®] Elite Kit
W 13.5 in. x H 11.25 in. x D 5.5 in.
Weight with hardcase: 10.5 lb
Weight in Rucksack: 3.8 lb



CardioPAK™ *from Resolute Medical*

- Lightest and smallest ECG available
- Designed for use in the field - ruggedized, portable, and easy to use
- Automated interpretation of ECG rhythms to alert the user of patient risks
- Results can be transferred through field communication channels
- Includes Panasonic Toughpad® for viewing, interpreting, or sending data
- Includes Philips HeartStart FR3 AED with 300 shock capacity and Bluetooth compatibility
- Useful for early diagnosis of abnormal heart conditions before sending personnel to the field
- Versatile, inexpensive, and highly portable solution for base hospital ECG testing

ORDERING INFORMATION

OPTION	TOUGH PAD	ITEM #
CardioPak Compact	7 in.	20-0086
CardioPak Elite	10 in.	20-0087



EVERY CardioPAK™ Kit Contains:

- 1 x MedicECG® 12-lead ECG
- 1 x Panasonic Toughpad® Rugged Tablet
- 1 x Philips HeartStart FR3 AED with pads
- 12 x Electrode Leads
- 12 x Electrode Clips compatible with both snap and tab electrodes
- 12 x Electrodes
- 1 x crushproof, waterproof case
- 1 x Nylon Bag
- 1 x Charger
- 1 x USB Cable



CardioPAK® ELITE (with 10 in. FZ-G1 Toughpad® tablet)

SPECIFICATIONS

Philips HeartStart FR3 AED:
W 5.3 in. x H 8.7 in. x D 2.7 in.
Weight: 3.5 lb

FZ-G1 Toughpad® tablet
W 10.6 in. x H 7.4 in. x D 0.8 in.
Weight: 2.5 lb

CardioPAK™ Elite Kit
W 13.5 in. x H 11.25 in. x D 8.2 in.
Weight with hardcase: 14 lb
Weight in Rucksack: 7.3 lb



CardioPAK® COMPACT (with 7 in. FZ-M1 Toughpad® tablet)

SPECIFICATIONS

Philips HeartStart FR3 AED:
W 5.3 in. x H 8.7 in. x D 2.7 in.
Weight: 3.5 lb

FZ-M1 Toughpad® tablet
W 7.9 in. x H 5.2 in. x D 0.7 in.
Weight: 1.2 lb

CardioPAK™ Compact Kit
W 13.5 in. x H 11.25 in. x D 8.2 in.
Weight with hardcase: 12.7 lb
Weight in Rucksack: 6.0 lb



FIELD CORPSMAN KITS

Portable ENT diagnostics

- Lightweight, durable ENT diagnostic sets for portable field care of the eyes, ears & nose
- A wide selection of high intensity lights & magnifiers for enhancing both day and night examinations

FIELD CORPSMAN KITS:



Additional items *in Deluxe Kit only* ▼



SPECIFICATIONS

H 8 in. x W 4.25 in. x D 2.25 in.
Weight: Basic Kit: 6.2 oz
Deluxe Kit: 11.6 oz
Ophthalmoscope only: 3.7 oz

ORDERING INFORMATION

KIT OPTIONS	ITEM #	NSN #
Basic	20-0001	6515-01-458-6057
Deluxe (with Heine)	20-0002	6515-01-458-6178
Ophthalmoscope only	20-0007	6540-01-491-6651

BP/STETHOSCOPE COMBO KIT *from NAR*

- Ultra-compact, multi-pocketed BP kit that can be stored in most trauma kits
- Blood-pressure cuffs in 4 sizes (adult, large adult, child, infant).
- ADC stethoscope with extra-large bell & ultra-sensitive diaphragm
- Ambidextrous palm gauge

ORDERING INFORMATION

ITEM #	NSN #
20-0057	6515-01-542-7723

SPECIFICATIONS

H 6.75 in. x W 9 in. x D 5 in.
Weight: 2 lb



HEINE MINI 3000 *Ophthalmoscope*

- High quality optics for examining the interior of the eye - in a pocket-size
- Enhanced XHL Xenon Halogen Technology; concentrated white light for perfect illumination & crisp retinal detail

ORDERING INFORMATION

ITEM #	NSN #
20-0007	6540-01-491-6651

SPECIFICATIONS

H 6.75 in. x W 1.25 in. x D 0.9 in.
Weight: 4 oz



STETHOSCOPE *from ADC*


- Best full-range stethoscope for auscultation & comfort
- Superior acoustics for accurate auscultation of heart, lung & blood-pressure sounds

SPECIFICATIONS

Device Length: 30.5 in.
Weight: 13.5 oz

ORDERING INFO

ITEM #
20-0060





Clear Lid
ZZ-0259 only
includes Loop Liner
on the bottom

ARMADILLO CASE

Medication Storage

- Compact solution for keeping critical medications safe, secure and easily accessible
- Features Automatic Pressure Equalization Valve
- Hardware constructed of Marine grade 316 SS
- Case is water-resistant and crush-proof
- Hook & loop attachment system allows storage compartment to be tailored to hundreds of custom configurations

SPECIFICATIONS

Ext: L 7.95 in. x W 5.1 in. x D 2.14 in.
Int: L 6.45 in. x W 3.87 in. x D 1.75 in.
Weight: 10.5 oz

ORDERING INFORMATION

* NATO NSN

LID STYLE	ITEM #	NSN #
Black	ZZ-0081	8145-01-573-2533
Clear	ZZ-0259	6545-01-594-9399*

also available: Half Inch Hook Discs with pressure sensitive adhesive
Set of 250 • Item# ZZ-0258 • NSN# 8145-01-617-9075*

M.U.L.E.

Multi-Use Luggage & Equipment Carrier

- Heavy-duty, water-resistant nylon platform with large compartment and mesh pockets on the inside lid for maximum storage of critical items
- Outside end pockets feature open channels for expanded storage capability
- Suitcase-style carrying handles on top and side; Telescoping handle and wheels for easy transport
- Hidden rucksack shoulder pads

ORDERING INFO

ITEM

ZZ-0212

SPECIFICATIONS

H 30.5 in. x W 18.5 in. x D 11 in.
Weight: 11 lb 10 oz Payload: 70 lb



I&A CASE

Injury & Ailment Case



- Superbly engineered medical equipment organization system
- Lid's clear pockets and hook & loop fasteners allow secure storage as well as superior inventory management
- Polymer O-ring provides dust and water-proof seal
- Automatic Gore-TEX® pressure equalization valve stops moisture from entering the case – also prevents vacuum-lock
- Double-throw latches engineered to tolerate up to 400 lb of pressure
- All-weather, sure-grip handle and reinforced padlock protectors

ORDERING INFO

ITEM

ZZ-0222

SPECIFICATIONS

H 8.5 in. x W 20.5 in. x D 17 in.
Weight: 13 lb 12 oz



PAXLIGHT® LED Field Medical Light

- Compact, lightweight precision lighting solution with the rugged engineering required for use in field hospitals
- LED Technology with intense light fifteen times brighter than traditional halogen field lights
- Patented 3-pod head eliminates shadow interference common with single point light sources
- Energy efficient, cool white color with an exceptionally long life
- Stays cool to the touch, projecting no heat on the patient
- Locking joints make it fully adaptable and secure with adjustable mounting clamp for attachment to hospital bed, litter or litter stand
- Includes a rugged storage case with wheels for easy transport

SPECIFICATIONS

Light in case:
L 47.2 in. x W 16.5 in. x D 6.7 in.

Light Head:
12.8 in. Diameter

Light Arm/Post:
L 78 in. x 2.8 in. Diameter

Weight: Case & Light:
37 lb 8 oz

Weight: Light with Arm/Post:
11 lb 8 oz

Light Output: 50,000 lux
(4590 foot candles)

ORDERING INFO

ITEM #	NSN #
60-0042	6530-01-549-1089



Wheeled Base

ITEM #
60-0045



EZ Light Handles

UNIT	ITEM #
Pack of 10	60-0046

MARK IV Field Hospital Bed

- Lightweight solution designed for patient comfort in the combat support or field hospital/treatment tent
- Durable & easy-to-open, supporting up to 850 lb of static weight
- Designed to be rapidly deployed with minimal training in less than 60 seconds by a single person

SPECIFICATIONS

Folded: L 42 in. x W 17 in. x D 32 in.
Opened: L 79.25 in. x W 33 in. x H 35.25 in.

Weight:
54 lb (Mattress and Frame)
68 lb (Mattress, Frame & Accessories)
U.S. Patents Pending

- Constructed of a bi-fold aluminum frame that is adjustable to 6 positions
- 3 in. polyurethane foam pad with water-resistant decontaminable cover easily cleans with 10% bleach solution or other standard hospital cleaning products
- Includes IV Pole and a removable table platform that rests on its adjustable bed rails

ORDERING INFORMATION

ITEM #	NSN #
66-0019	6530-01-624-7626




microMend™

Skin Closure Devices

- Safe, secure, and rapid skin closure
- Equivalent holding strength to sutures
- Patent pending microstaples anchor the device to the skin
- Three times faster to apply than sutures
- Available in two sizes

SPECIFICATIONS

Standard Device: L 43 mm x W 12 mm
 Small Device: L 28 mm x W 12 mm
 Box: L 5.375 in. x H 2.25 in. x D 2.375 in.
 U.S. Patents Pending



Available in small and standard strips

ORDERING INFORMATION

DEVICE SIZE	UNIT (per box)	ITEM #
◆ Standard	24	30-0203
◆ Standard	48	30-0204
◆ Small	24	30-0230
◆ Standard	8	30-0231
◆ Small	8	30-0232



R-SICK *Rapid Surgical Intervention Convenience Kits*

- Flexible austere and prolonged field care surgical instrument storage and recovery system tray
- For both major and minor surgical procedures in less than ideal or highly mobile environments
- Sequentially set up for rapid "incision, spreading, damage control, and closure"
- Ideal for complex, emergent procedures such as surgical airways, hemorrhage control, field amputation, and REBOA



Instruments are available in Floor & OR Grade

SPECIFICATIONS

Packaged: H 8 in. x W 7 in. x D 3.5 in.
Fully deployed (flat): L 24 in. x W 43 in. x D 0.1 in.
Weight: Kits 3 lb • Bag Only 11 oz

R-SICK CONTENTS:

2 x Weitlaner Retractor, 3 x 4 Blunt (5.5 in.)
1 x Trousseau Tracheal Dilator (5.25 in.)
2 x Adson Forceps, 1 x 2 (4.75 in.)
1 x Ewald Forceps 1 x 2
1 x Plain Splinter Forceps (4.5 in.) Serr. Handle
1 x Kelly Forcep Heavy
1 x Converse Needle Holder (4.24 in.) Long
1 x Olsen-Hegar Needle Holder (7.25 in.) Serr
1 x Stevens Tenotomy Scissors Trbl/Bl
1 x Stevens Tenotomy Scissors Cvdbl/Bl
1 x Scissors, Mayo (5.5 in.) Str
1 x Mayo Dissecting Scissors, Crv (5.5 in.)
1 x Knife Handle (#3)
2 x Gigli Saw Blade (20 in.)
2 x Gigli Saw Handle (3.5 in.) Loop
1 x Hiupp Tracheal Hook, (6.5 in.)
2 x Forcep, Rochester-Pean (7.25 in.) Curv
6 x Scalpel #10 w/Handles

6 x Scalpel #11
6 x Scalpel #15
6 x Blade, #10 Pluss/S(Tefl)
10 x Blade, #15 Pluss/S(Tefl)
2 x Suture Ethilon 4-0 Blk Mono Fs-218
2 x Suture 3-0 Silk
2 x Black Braided Silk Suture, 2-0/2-60 in.
2 x 2-0 Suture Ct-2
4 x 4-0 Silk Suture Blk (30 in.)
1 x Sharps Shuttle
1 x Farabeuf Retractor

ORDERING INFORMATION

INSTRUMENTS	ITEM #
Floor Grade	83-0016
OR Grade	83-0017
Bag Only	80-0871



PENBLADE

Safety Scalpels

- Pen-like activation that uses existing muscle memory
- Blade retraction without repositioning of the hand
- Suture trimming capability when blade is retracted
- High-quality materials and construction

SPECIFICATIONS

Packaged: L 7.5 in. x W 1.7 in. x D 0.45 in.
Unpackaged: L 6 in. x W 0.7 in. x D 0.4 in.
Weight: 0.4 oz

ORDERING INFORMATION

BLADE	COLOR	ITEM #
#10	Blue	ZZ-0877
#11	Orange	ZZ-0878
#15	Green	ZZ-0879



MINOR SURGERY Sets

Minor Surgery Contents:

1 x Tray, Shallow
2 x Kelly Forceps, Straight 5.5 in.
2 x Kelly Forcep, Curved 5.5 in.
1 x Mayo Hegar Needle Holder, 6 in.
1 x Adson Forceps, Serr 4.75 in.
1 x Forceps Tissue, 5 in.
1 x Collier Needle Holder, Straight 5 in.
1 x Handle Surgical Knife #3
1 x Probe W. Eye 5 in.
1 x O.R. Scissors, Sb St 5.5 in.
1 x Wrap, Kimguard 24 x 24 in.

SPECIFICATIONS

Packaged: L 9 in. x W 5 in. x D 1 in.
Weight: 5 oz



ORDERING INFORMATION

ITEM #	UNIT	NSN#
ZZ-0960	Single	-----
ZZ-0959	Pack of 20	-----

SURGICAL SUTURE Kits

Surgical Suture Contents:

1 x Tray, Shallow
1 x Adson Forceps w/Teeth
1 x Webster Needle Holder
1 x Iris Scissors, Straight
1 x Mosquito Hemostat, Curved
1 x Poly Bag w/Twist Tie
1 x Poly-backed Towel
1 x Fenestrated Drape
1 x Alcohol Prep Swab
1 x PVP Prep Swab
1 x Cotton Filled Sponge.

SPECIFICATIONS

Packaged: L 9 in. x W 5 in. x D 1 in.
Weight: 6 oz



Sutures
not included



ORDERING INFORMATION

ITEM #	UNIT	NSN#
ZZ-0958	Single	-----
ZZ-0961	Pack of 20	6515-01-311-0360

TOENAIL REMOVAL Trays

Toenail Removal Contents:

1 x Tray, Half Deep
1 x CSR Wrap, Kimguard 24 in. x 24 in.
1 x Kelly Forcep, Straight 5.5 in.
1 x Kelly Forcep, Curved 5.5 in.
1 x Adson Forceps, 1 x 2 Teeth, 4.75 in.
1 x Iris Scissors, Straight 4.5 in.
1 x Cuticle Nipper 4 in.
1 x Probe w/Eye 5.5 in.
10 x 12-Ply Gauze, 4 in. x 4 in.
1 x Kwiktear Header Bag, 10 in. x 15 in.

SPECIFICATIONS

Packaged: L 9.72 in. x W 5.4 in. x D 2 in.
Weight: 11 oz



ORDERING INFORMATION

ITEM #	UNIT	NSN#
ZZ-0956	Single	-----
ZZ-0957	Pack of 10	6515-01-630-8050



SUTURE REMOVAL KIT #1



QTY
1 KC 100
Kimguard



QTY
1 Littauer
Scissors



QTY
1 Adson
Forceps



QTY
1 Gauze
Pad, 3 x 3

ITEM #

ZZ-0664

SPECIFICATIONS

STERILE

Pkg: L 6 in. x W 3.5 in. x D 0.38 in.
Weight: 2 oz

SUTURE REMOVAL KIT #2



QTY
1 KC 100
Kimguard



QTY
1 Crvd Iris
Scissors



QTY
1 Adson
Forceps



QTY
1 Gauze
Pad, 3 x 3



QTY
1 Alcohol
Prep Pad



QTY
1 Iodine
Prep Pad

ITEM #

ZZ-0665

SPECIFICATIONS

STERILE

Pkg: L 6 in. x W 3.5 in. x D 0.38 in.
Weight: 2 oz

SUTURE REMOVAL KIT #3



QTY
1 KC 100
Kimguard



QTY
1 Str Iris
Scissors



QTY
1 Adson
Forceps



QTY
1 Gauze
Pad, 3 x 3



QTY
1 Alcohol
Prep Pad

ITEM #

ZZ-0666

SPECIFICATIONS

STERILE

Pkg: L 6 in. x W 3.5 in. x D 0.38 in.
Weight: 2 oz

SUTURE REMOVAL KIT #4



QTY
1 KC 100
Kimguard



QTY
1 Littauer
Scissors



QTY
1 Adson
Forceps



QTY
1 Alcohol
Prep Pad



QTY
1 Iodine
Prep Pad

ITEM #

ZZ-0667

SPECIFICATIONS

STERILE

Pkg: L 6 in. x W 3.5 in. x D 0.38 in.
Weight: 2 oz

SURGICAL INSTRUMENTS

Precision standard surgical instruments available individually

STERILE

ADSON FORCEPS

FP, 1x2 • 4.75 in.



ITEM #

ZZ-0705

HALSTEAD MOSQUITO

Curved • 5 in.



ITEM #

ZZ-0695

HALSTEAD MOSQUITO

Straight • 5 in.



ITEM #

ZZ-0694

IRIS SCISSORS

Curved • 4.5 in.




ITEM #

ZZ-0697

IRIS SCISSORS

Straight • 4.5 in.



ITEM #

ZZ-0696

KELLY HEMOSTAT

Curved • 5.5 in.




ITEM #

ZZ-0701

KELLY HEMOSTAT

Straight • 5.5 in.



ITEM #

ZZ-0700

LISTER SCISSORS

5.5 in.



ITEM #

ZZ-0698

O.R. SCISSORS

SB Straight 5.5 in.



ITEM #

ZZ-0699

ROCHESTER PEAN

Hemostat, Curved 6.25 in.




ITEM #

ZZ-0702

WEBSTER NEEDLE

Holder, 5 in.




ITEM #

ZZ-0703



FACIAL LACERATION TRAY

- Cost effective, sterile, single use kit for facial lacerations with a 3 year shelf life
- Conveniently packaged in a peel-to-open pouch for easy access



ITEM #	NSN#
ZZ-0683	6515-01-661-1802

SPECIFICATIONS

Pkg: L 9 in. x W 5 in. x D 2.75 in.
Weight: 12 oz

STERILE



MINOR LACERATION TRAY

- Cost effective, sterile, single use kit for facial lacerations with a 3 year shelf life
- Conveniently packaged in a peel-to-open pouch for easy access



ITEM #	NSN#
ZZ-0684	6515-01-661-1786

SPECIFICATIONS

Pkg: L 9 in. x W 5 in. x D 2.75 in.
Weight: 12 oz

STERILE



- Cost effective, sterile, single use kit for facial lacerations with a 3 year shelf life
- Conveniently packaged in a peel-to-open pouch for easy access to components


SPECIFICATIONS

Pkg: L 9 in. x W 5 in. x D 2.75 in.
Weight: 12 oz

STERILE

ITEM

ZZ-0685



Custom Procedure Trays/Kits available • Contact North American Rescue for information





RESCUE

EQUIPMENT & ACCESSORIES

see pages 74-85



PEGASUS ROPE®

WARHORSE™ SERIES

100% Innegra™ fiber braided core


- High quality, multipurpose rope
- Static rope provides reliability, high performance and great handling
- Excellent knot-ability with high strength-to-weight ratio
- Low stretch/sheath slippage – wet, dry or in any other climatic condition
- Ultra-lightweight, strong Innegra™ core for reliable strength, durability & flotation
- Properties remain unchanged even when wet, dry or contaminated with chemicals
- Low dielectric and UV stable

SPECIFICATIONS

Diameter	9.1 mm	9.5 mm
Weight:	4.88 lb/100 ft	3.34 lb/100 ft
Tensile Strength:	5,670 lbf (22.7 kN)	5,100 lbf (22.7 kN)
Fibers (core/sheath):	Innegra™/Polyester	Innegra™/Polyester
Sheath Ratio:	63% Polyester	28% Polyester
Ratings:	NFPA 1983/2012 Throwline (Pending)	NFPA 1983/2012 Throwline (Pending)

ORDERING INFORMATION

DIAMETER	COLOR	RUNNER	ITEM #
9.1 mm	Black	Gray	90-0207
9.5 mm	Black	Gray Reflective	90-0206
9.5 mm	Yellow	Blue Reflective	90-0205
9.5 mm	White	Blue	90-0208





OLYMPUS™ SERIES

Engineered for rough terrain

- Ideal for use in challenging and abusive environments
- Durable, lightweight and great handling offers reliability in all weather conditions
- Technora sheath adds strength and abrasion resistance
- Buoyant even with a 62% Technora sheath (9.5mm)

ORDERING INFORMATION

DIAMETER	COLOR	RUNNER	ITEM #
9.5 mm	Gold	Black/Red	90-0401
11 mm	Gold	Black/Blue	90-0400

SPECIFICATIONS

Diameter	9.5 mm	11 mm
Weight:	3.82 lb/100 ft	4.99 lb/100 ft
Tensile Strength:	5,495 lbf (24.4 kN)	6,834 lbf (30.4 kN)
Fibers (core/sheath):	Innegra™/Technora	Innegra™/Technora
Sheath Ratio:	62% Technora	60% Technora



STALLION™ SERIES

Low visibility - high performance

- 100 ft of this 9.5 mm rope offers 5,100 lb of strength under 3 lb – wet or dry
- Chemically resistant to acids, bases & common industrial chemicals
- Low sheath slippage, excellent strength-to-weight ratio and hydrophobic

ORDERING INFORMATION

DIAMETER	COLOR	RUNNER	ITEM #
9.5 mm	Black	Blue/Reflective	90-0103
9.5 mm	Black	Blue	90-0104

SPECIFICATIONS

Diameter	9.5 mm
Weight:	2.9 lb/100 ft
Tensile Strength:	5,100 lbf (22.7 kN)
Fibers (core/sheath):	Innegra™/Innegra™
Sheath Ratio:	24% Innegra™

TRITON™ SERIES

Specially designed for water rescue

- 10% lighter & 250% stronger than polypropylene alone
- 100% Innegra™ fiber braided core, the lightest structure fiber available, provides additional buoyancy, low elongation and consistent strength

SPECIFICATIONS

Diameter	7.5 mm
Weight:	2 lb/100 ft
Tensile Strength:	3,050 lbf (13.6 kN)
Fibers (core/sheath):	Innegra™/Polypropylene
Sheath Ratio:	46% Polypropylene
Ratings:	NFPA 1983/2012
	Throwline (Pending)

ORDERING INFORMATION

DIAMETER	COLOR	RUNNER	ITEM #
7.5 mm	Yellow	Blue Reflective	90-0303
7.5 mm	Red	Yellow	90-0305



RESCUE | Rope Kits

R.A.C.E. SYSTEM Rapid Access Casualty Extraction

- Lightweight, multiple use casualty evacuation system enabling first responders/medical providers to ascend, descend, hoist & lower both rescuer & casualty to safety
- Combines 150 ft of lightweight durable 9.5 mm Olympus Pegasus Rope® with the latest in lightweight multipurpose hardware
- Capable of rigging 2:1, 3:1, 4:1, 5:1, 6:1, and 9:1 systems in different configurations for a wide variety of rescuer options with fully illustrated instruction booklet
- Front pocket allows for optional storage
- Carriage options include 4 Point, D-ring configurable modified side and/or backpack carry-style adjustable strap



R.A.C.E. SYSTEM:



ORDERING INFORMATION

COLOR	ITEM #
Multi-Cam	90-0063
Coyote	90-0077

SPECIFICATIONS

H 15 in. x W 10 in. x D 7.5 in.
Weight: 10 lb 14 oz

Water Rescue BELT THROW BAG

- Multipurpose water rescue solution designed for easy deployment and retrieval for rescuers
- Quick release reverse tapered throw bag for rapid packing of Triton rope and a slim profile for users
- 2 in. adjustable belt with quick release strap, 1.5 in. plastic cam buckle, and red toggle ball for easy ID



ORDERING INFO

ITEM #
90-0065

SPECIFICATIONS

Rope Bag: L 11 in. x W 4.5 in. x D 4 in.
Gear Pouch: L 5.5 in. x W 5 in. x D 2 in.
Total Weight: 2 lb 2 oz

50 Foot THROW BAG

- Must-have necessity for water rescue
- Durable, easy-to-use rescue device with highly visible orange nylon bag and mesh panels that drain quickly for use in multiple rescue operations
- Packed with 50 ft of 7.5 mm Triton Pegasus Rope® rated at 3,050 lbf
- Triton rope is strong enough for multiple uses: as a throw line, tension diagonal, tow line and for pinned kayak recovery



ORDERING INFO

ITEM #
90-0064

SPECIFICATIONS

Rope Bag: L 11 in. x W 4.5 in. x D 4 in.
Gear Pouch: L 5.5 in. x W 5 in. x D 2 in.
Total Weight: 2 lb 2 oz

HELO TAGLINE KIT

Help control litter movement during hoist

- Help control litter movement during hoist operations – a requirement on all Army MEDEVAC helicopters
- Contains 250 ft of 7.5 mm Triton Pegasus Rope® featuring a highly visible polypropylene sheath with a super strong Innegra™ core
- Includes one locking carabiner, ten weak links (designed to break away if the litter is snagged during the hoist) and one stainless steel quick link

SPECIFICATIONS

L 11 in. x W 7 in. x D 7 in.
Weight: 5 lb 13 oz

ORDERING INFO

ITEM

90-0069

HELO TAGLINE KIT:



QTY
1

250 ft
Triton Rope



QTY
1

rockD
Carabiner



QTY
1

Quick
Link, SS



QTY
10

Weak
Links



HIGH THREAT EXTRACTION KIT

Essential tools for rescue ops in high threat environments

- Facilitates immediate casualty extraction so that Care Under Fire casualty treatment may be performed as soon as possible
- Self-contained kit with all of the essential tools for tactical operators and rescue personnel
- Includes a C-A-T® tourniquet to immediately address extremity hemorrhage at the point of wounding
- Rugged nylon bag attaches to the duty belt as a drop down leg rig or to the assault vest/body armor with PALS/MOLLE-style connectors

SPECIFICATIONS

L 8.5 in. x W 8.75 in. x D 4.5 in.
Weight: 3 lb 6 oz

HIGH THREAT EXTRACTION KIT:




QTY
1

C-A-T®
Tourniquet



QTY
1

RAT™
Strap



QTY
1

Dragon
Handle System



QTY
1

Hasty
Harness



QTY
1

Shield
Carabiner



QTY
1

Aero
Carabiner



QTY
1

Revolver
Carabiner



QTY
2

Red
ChemLight™

ORDERING INFORMATION

COLOR	ITEM #
Coyote	80-0141
OD Green	80-0140
Black	80-0139
Digital Camo	80-0142





RESCUE | High Threat Extraction

RAT STRAP

Rescue Assault Tether

- For head/feet first linear drags during shield-assisted, open-field/indoor-corridor rescues or single-rescuer extractions
- Innovative Daisy Chain configuration offers multiple reinforced attachment points

SPECIFICATIONS

Pkg: H 1 in. x W 4 in. x D 5.5 in.
Extended: L 62 in.
Minimum Break Strength: 2,702 lb
Weight: 12.7 oz

ORDERING INFO

ITEM #	NSN #
90-0002	6530-01-632-7872



EX-TRACT STRAP

Point of Wounding Extraction


- Point of wounding extraction
- Back pack carry
- Hasty harness for traditional and confined space lowering and hauling

SPECIFICATIONS

Pkg: H 18 in. x W 3.5 in. x D 1 in.
Extended: L 58 in.
Weight: 9.3 oz

ORDERING INFO

ITEM #
90-0082



DRAGON HARNESS

Under the vest harness

- Built-in handle system for immediate high threat extraction
- Operator harness to aid in the rescue of downed personnel without the need to holster or sling weapons
- Prevents the "ride up" of the casualty's ballistic vest, reducing exposed, unprotected body regions

SPECIFICATIONS

Height: 20 in. from waist to non-deployed handle
Length: 14 in. from harness attachment point to top of handle
Weight: 1 lb 1.6 oz.
Waist Size: 30 in. to 45 in.

ORDERING INFO

ITEM #
90-0003



DRAGON HANDLE

Tool offering head/feet first drag options


- Maintain an anatomically correct posture for rapidly extracting greater weight
- Decreases manpower needs in "zeroed in" locations
- Black Omega, non-locking "D" carabiner rated for 31kN (6,968.8 lb)

SPECIFICATIONS

H 6.3 in. x W 4.1 in. x D 2.25 in.
Weight: 11.3 oz | Device Length: L 42 in.
US Patent(s): Pending

ORDERING INFO

ITEM #	NSN #
90-0004	6530-01-598-9595



TACTICAL BAIL-OUT KIT

Rapid extraction of individuals or teams

- Integral for multiple hasty rescues including rapid casualty lowering & hoisting, rescuer throw bag and rappel line
- Facilitates a variety of anchor locations & rescue options

ITEM


90-0020

SPECIFICATIONS

L 10 in. x W 6 in. x D 4 in.
Weight: 2 lb 14 oz



TACTICAL BAIL-OUT KIT:



HASTY HARNESS

Multi-purpose harness

- Resourceful, multi-purpose harness for unlimited applications & maximum versatility
- Indispensable gear for every webbing equipment toolbox
- Constructed of 22 ft loop of 1 in. Tech Tape (MBS: 4,400 lb) with opposing color thread to aid in visual safety inspection

SPECIFICATIONS

Loop:
L 10 ft 11 in. x W 1 in.
Circumference:
21 ft 10 in.
Min. Break Strength:
Tech Tape (4,400 lb)
Sewn Sling (6,000 lb)

ORDERING INFO NATO NSN*

COLOR	ITEM #	NSN #
Red	90-0044	----
Black	90-0001	6530-01-674-9774*





ROCKD CARABINER

by Rock Exotica®

Narrow asymmetric "D" shape machined frame with snag resistant keylock nose and large, ergonomic angled gate opening

SPECS

Gate Opening: .24 mm
Weight: 2.6 oz (73 g)

ITEM

ZZ-0543



BASIC ASCENDER

by Petzl®

Compact ascender with excellent grip for rope ascent on a fixed line or as progress capture in a hauling system

SPECS

Rope Size: 8-11 mm
Weight: 3 oz (85 g)

ITEM

ZZ-0554



GRIGRI 2 BELAY DEVICE

by Petzl®

Compact, ultra-light belay device allowing excellent control during the descent with assisted braking capability

SPECS

Rope Size: 8-11 mm
Weight: 6 oz (170 g)

ITEM

ZZ-0545



MINI MACHINED PULLEY

by Rock Exotica®

Mini pulley with recessed axle screw, machined sideplates and integrated axle for a stronger, compact, snag-free profile

SPECS

Rope Size: up to 11 mm
Weight: 2.9 oz (83 g)

ITEM

ZZ-0548



TRIRIG RIGGING PLATE

by Rock Exotica®

Lightweight, high performance 3 hole rigplate with numbered holes to assist with detailed rigging

SPECS

Break Strg: 36 kN
Weight: 1.8 oz (51 g)

ITEM

ZZ-0553



REVERSO 4

Belay/Rappel device by Petzl®

Regulates the amount of braking friction on the rope according to diameter and condition

SPECS

Weight: 2.08 oz (59 g)

ITEM

ZZ-0546



TIBLOC

Ascender/Rope Grab by Petzl®

Ultra-light Emergency Ascender/Rope Grab that can be used to build hauling systems or as a Prusik-replacement in self-rescue situations

SPECS

Rope Size: 8-11 mm
Weight: 3 oz (85 g)

ITEM

ZZ-0547



FOOTAPE

Adjustable Foot Loop by Petzl®


Adjustable nylon foot loop for attachment to the BASIC rope clamp for rope ascents

SPECS

Weight: 2.3 oz (65 g)

ITEM

ZZ-0551



ANNEAU SLINGS

by Petzl®

Polyester sling to set up a belay or an anchor

SPECS

W 19 mm
L 60 or 120 cm
Wgt: 60 cm - 2 oz
120 cm - 3.5 oz

LENGTH ITEM

60 cm ZZ-0549

120 cm ZZ-0550



SHIELD CARABINER

Hot forged "I-Beam" construction with patented anti-snag twisted wirelock gate

SPECS

Length: 100 mm
Weight: 1.1 oz

ITEM

90-0005



AERO SCREWGATE

Offset locking D carabiner with keylock nose that is the new standard for entry level lockers

SPECS

L 4 in x 2.46 in
Weight: 2.1 oz

ITEM

90-0006



REVOLVER CARABINER

Unique roller insert performs as an emergency pulley

SPECS

Gate Clear: 24 mm
Weight: 1.6 oz

ITEM

90-0007





RAPTOR™ MULTI-TOOL by Leatherman®

High-performance multi-tool equipped with a variety of mission critical tools essential to emergency medical practitioners in the field

1. 420HC Stainless Steel Folding Medical Shears
2. Strap Cutter
3. Oxygen Tank Wrench
4. Ring Cutter
5. Carbide Glass Breaker
6. 5 cm Ruler
7. Holster with Clip

SPECIFICATIONS

Shear Length: 1.9 in.
Closed Length: 5 in.
Weight: 5.8 oz

ORDERING INFO

ITEM #	NSN #
70-0014	----

See our Treatment Support Section
on page 87 for the compatible Scissor Leash

SUPERTOOL 300 by Leatherman®

Heavy-duty, all purpose multi-tool featuring 19 devices in one easy to carry configuration

SPECIFICATIONS

Closed: L 4.5 in.
Primary Blade: L 3.2 in..
Weight: 9.6 oz

ORDERING INFO

ITEM #
ZZ-0542

1. Needlenose Pliers
2. Regular Pliers
3. Replaceable Wire Cutters
4. Hard-Wire Replaceable Cutters
5. Stranded-Wire Cutters
6. Electrical Crimper
7. Wire Stripper
8. 420HC Knife
9. 420HC Serrated Knife



10. Saw
11. Awl with Thread loop
12. 9 in. Ruler
13. Can Opener
14. Bottle Opener
15. Wood/Metal File
16. Phillips Screwdriver
17. Large Screwdriver
18. Medium Screwdriver
19. Small Screwdriver



MUT® by Leatherman®

Tactical and practical multi-tool for military, law enforcement and civilian shooters - features 16 devices

1. Needlenose Pliers
2. Regular Pliers
3. Replaceable Wire Cutters
4. Hard-Wire Replac. Cutters
5. Stranded-Wire Cutters
6. Electrical Crimper
7. 420HC Knife
8. Saw
9. Replaceable Cutting Hook
10. Hammer
11. Bolt Override Tool
12. Bronze Carbon Scraper
13. #8-32 Cleaning Rod/Brush Adapter
14. Firearm Disassembly Punch
15. Carabiner/Bottle Opener Punch
16. Large Bit Drive

SPECIFICATIONS

Closed: L 5 in.
Primary Blade: L 3 in.
Weight: 11.2 oz

ORDERING INFO

ITEM #
ZZ-0539



BENCHMADE® KNIVES & CUTTERS



1 HOOK SAFETY CUTTERS by Benchmade®

- Hook cutters designed to quickly and reliably cut through a wide array of materials
- Extensively tested and widely used by the US Military for medical and egress purposes

ITEM #	STYLE	NSN #
ZZ-0652	7 Hook	----
ZZ-0250	8 Hook	----

SPECS
7 Hook: L 4.3 in.
8 Hook: 6.5 in.
Weight: 1.5 oz

2 SOCP RESCUE TOOL by Benchmade®

- Fixed 440C high-chromium stainless steel blade hook for rescue performance and less-than-lethal self-defense
- Includes carbide glass breaker, oxygen wrench, single-piece handle with finger loop & molded sheath
- Designed by Special Operations Combatives Program (SOCP)

SPECS
Blade: 2.5 in. Overall: 6.75 in.
Weight: 1.88 oz

ITEM #	NSN #
ZZ-0760	----

3 H₂O 112 SBK by Benchmade®

- N680 ultra stainless steel blade with opposing fixed bevel blade style, blunt tip, & integrated rescue hook
- Boltaron® sheath designed to fit current Navy SAR kits
- 90-0078 (Orange handle) includes carbide glass breaker & is approved for Navy SAR - "So Others May Live"

ITEM #	HANDLE	NSN #
90-0078	Orange	1095-01-677-1013
♦ 90-0097	Black	----

SPECS
Blade: 3.5 in.
Overall: 8.02 in.
Weight: 3.2 oz

4 TRIAGE 916 SBK by Benchmade®

- Dynamic triple utility tool featuring three essential rescue tools in a single compact design
- N680 high corrosion resistant, partially serrated, opposing bevel blade, 440C folding hook & carbide glass break

SPECS
Blade: 3.5 in
Overall: 8.2 in.
Weight: 5.1 oz

ITEM #	HANDLE	NSN #
90-0096	Orange	----
90-0060	Black	----

5 AUTOTRIAGE 9160 SBK by Benchmade®

SPECS
Blade: 3.58 in. Overall: 8.31 in.
Weight: 5.75 oz

ITEM #	NSN #
90-0061	----

- Packaged in high-strength aluminum with AXIS® auto pull release and integrated safety
- Auto opening N680 highly corrosion resistant blade with modified drop-point

6 TRIAGE 917 BK -or- SBK by Benchmade®

- Black drop point style CPM-S30V steel blade with option of a smooth (BK) or partially serrated edge (SBK)
- 440C folding hook & carbide glass breaker

SPECS
Blade: 3.5 in Overall: 8.3 in.
Weight: 5.24 oz

ITEM #	STYLE	NSN #
90-0092	BK	----
90-0093	SBK	----



ORDERING INFO

ITEM

ZZ-0552

SPECIFICATIONS



Light Output: Min. 5 Lumens
Max: 120 Lumens
(160 Lumens on Boost mode)
Weight: (with batteries) 3 oz

TACTIKKA + RGB

Headlamp by Petzl®

- Ambient Mode: Wide beam with reduced brightness, avoids blinding others in a group situation and provides a long burn time
- Close Range Vision Mode - Mixed beam with wide, uniform flood or 3 colored LED lighting modes (red, green & blue) to remain discreet
- Boost Mode for temporary access to maximum power, 160 lumens (60 m)
- Strobe mode for signaling one's presence; accessible at any time
- Reserve Mode- Burn time management automatically switches on by signaling with blinking light when batteries are almost drained. Red lighting activates to preserve a minimal amount of usage

NAR HEADLAMP

- 100 Lumen LED Light
- Dual color (white and red)
- High and low lighting modes (White)
- Sturdy and simple to rotate single arm bracket
- Large push button for operation
- Splash and water resistant (IPX4)
- Comes with 2 AAA batteries
- 58 hour burn time

ORDERING INFO

ITEM #	NSN #
ZZ-0106	6515-01-542-7937

SPECIFICATIONS

Size: W 2.25 in. x D 1.5 in.
Weight: 2.25 oz



Streamlight HEADLAMP

- Ring Sealed construction – Waterproof to 1M for 1 Hr
- Up to 6 hrs, runtime on high up to 24 hrs. on low
- Impact resistant ABS housing
- Includes (2) 1.5V "AAA" batteries
- Solid State Power Regulation provides maximum light output throughout battery life
- Shatter resistant polycarbonate lens



ORDERING INFO

ITEM

ZZ-1090

SPECIFICATIONS

Faceshield Diameter: 1.3 in.
Maximum headlight: W 2.25 in. x D 1.6 in.
Weight: (with batteries & headband) 2.75 oz

GECKO GRIP®

Multi-purpose Tape

Laminated polyethylene cloth tape featuring an aggressive natural rubber based adhesive formulation and perfect for heavy-duty applications

SPECIFICATIONS

Pack of 6:
L 12.25 in. x W 0.9 in. x D 0.7 in.
Individual Roll: W 2 in. x L 5 yd
Weight: Roll: 2.7 oz
Pack: 1 lb 8.7 oz

ORDERING INFO

UNIT	ITEM #
6 Pack	80-0131





TOMAHAWK NIGHT VISION by First Light®

- Right angle tactical light with rotating head and 3 button control panel for one-handed thumb actuation
- Provides White, IR, and Red+Green LED mix for highlighting blood during casualty treatment
- Featuring Integrated high intensity IR signaling beacon
- Includes momentary and three intensity levels of constant-on (all output colors)
- Can be used handheld or hands-free with MOLLE TRS Blade Mount



ORDERING INFORMATION

CASE COLOR	TASK LIGHT	ITEM #
Coyote	UV Task	70-0015
OD Green	UV Task	70-0112
Coyote	Blue Task	70-0111
OD Green	Blue Task	70-0110

SPECIFICATIONS

Output Wavelengths:
White, IR, & Red+Green (Blood Highlighting),
IR Signaling Beacon
Maximum White Light Output: 425 lumens
Max Output Runtime: 2 hours (ANSI Tested)
Power Source: CR123 Batteries (2) | Material: Aerospace-Grade Aluminum
With mount: H 3.4 in. x L 2.4 in. x W 1.5 in. | Weight: (with mount & batteries) 8 oz

TOMAHAWK MULTI-COLOR by First Light®

- Right angle tactical light with rotating head and 3 button control panel
- Provides White, Red, Green, and Red+Green LED mix for highlighting blood during casualty treatment
- Includes momentary and three intensity levels of constant-on

SPECIFICATIONS

Output Wavelengths: White, Red, Green,
& Red+Green (Blood Highlighting)
Maximum White Light: 425 lumens
Max Output Runtime: 2 hours (ANSI Tested)
Power Source: CR123 Batteries (2)
Material: Aerospace-Grade Aluminum
With mount: H 3.4 in. x L 2.4 in. x W 1.5 in.
Weight (with mount & batteries): 7.2 oz

ORDERING INFORMATION

CASE COLOR	TASK LIGHT	ITEM #
Coyote	None	70-0016
OD Green	None	70-0117
Coyote	UV Task	70-0116
OD Green	UV Task	70-0115
Coyote	Blue Task	70-0114
OD Green	Blue Task	70-0113

TORQ by First Light®

- Right angle tactical light, rotating head and 3 button control
- White, Red, Green, and Red+Green LED mix for highlighting blood
- Includes momentary and three intensity levels of constant-on (all output colors)
- Features integrated high intensity white tactical strobe and red/white/green flashing safety beacon
- 10 hour Runtime on high output

ORDERING INFO

ITEM #	NSN #
70-0118	6515-01-646-2636

SPECIFICATIONS

Wavelengths: White, Red, Green, &
Red+Green (Blood Highlighting)
Max Runtime: 10 hours (ANSI Tested)
With mount: H 3.9 in. x L 2.6 in. x W 1.3 in.
Weight (with mount & batteries): 7.1 oz







TREATMENT

see pages 88-96

SUPPORT • PPE • TRIAGE • ACCESSORIES



TREATMENT SUPPORT | Hyperthermia

POLAR SKIN™

- The latest cooling technology for effective Human Thermo-regulation
- Polar Skin Mini-Mesh,™ a unique smart textile with a hydrofunctional polymer that allows for liquid absorptive capacity 350% greater than base weight
- ArcticAqua,™ a proprietary solution that prevents rigidity at temperatures well below 0°F while providing the requisite moisture to support evaporative, convective and conductive cooling
- Adheres to current U.S. Army and National Athletic Trainers' Association protocols for heat illness treatment
- Protracted direct cooling of high blood flow regions

COMPREHENSIVE COOLING SYSTEM

- Highly-effective portable cooling system for the treatment of heat-related illnesses
- Scientifically proven to reduce core temperature by .09°F/minute in 100°F heat when used according to protocol



ORDERING INFO

ITEM

ZZ-0674

SPECIFICATIONS

L 11 in. x W 10 in. x H 4 in.

Weight: 8 lb 5.3 oz

COMPLETE CARE PACK



- Portable, highly insulated option
- Contains one Polar Skin Comprehensive Cooling System and one heavy-duty, leak & sweat proof cooler that keeps ice from melting for 24 hours in +100°F heat

ORDERING INFO

ITEM

80-0451

SPECIFICATIONS

W 14 in. x L 12 in. x D 7 in.

Weight: 10 lb 3.2 oz

COMPLETE CARE BACKPACK

Available in either Single Patient (one comprehensive cooling system) or Double Patient (two cooling systems) configurations

ORDERING INFO

SIZE	ITEM #
Single	80-0543
Double	80-0544

SPECIFICATIONS

W 19.5 in. x L 15 in. x D 7.5 in.

Single Weight: 12 lb 13 oz

Double Weight: 21 lb 13 oz



Available in Standard, Large and Backpack COOLER

Standard size accommodates one Polar Skin™ Comprehensive Cooling System and Large or Backpack Size accommodates two cooling systems

ORDERING INFORMATION

COLOR	SIZE	ITEM #
Red	Standard	ZZ-0682
Black	Standard	ZZ-0759
Black	Large	ZZ-0756
Black	Backpack	ZZ-0794



SPECIFICATIONS

Standard: L 14 in. x W 12 in. x H 7 in.
Weight: 1 lb 3.2 oz

Large: L 18.5 in. x W 13 in. x H 9.5 in.
Weight: 2 lb 13 oz

Backpack: L 19.5 in. x W 15 in. x D 7.5 in.
Weight: 3 lb 13 oz



ICE SHEETS *Standard*

- Principle cooling component of Polar Skin™ cooling systems
- Adheres to current National Athletic Trainers' Assoc and U.S. Army protocols for heat illness treatment
- Sheet Size: W 36 in. x L 30 in.



ORDERING INFORMATION

UNIT	ITEM #	NSN #
Single Pack	ZZ-0668	6532-01-653-6054
Double Pack	ZZ-0669	----
Triple Pack	ZZ-0755	----

SPECIFICATIONS

Single: L 4.5 in. x W 7.75 in. x H 0.5 in. Wgt: 9.6 oz
 Double: L 9.5 in. x W 8 in. x H 0.5 in. Wgt: 1 lb 6.4 oz
 Triple: L 12.5 in. x W 7.5 in. x H 0.5 in. Wgt: 2 lb 3.2 oz

ICE SHEETS *Large*

Same characteristics as the standard size Ice Sheets, but in a larger 36 in. x 60 in. size

ORDERING INFO

ITEM #	SPECIFICATIONS
ZZ-0758	L 8.5 in. x W 5.5 in. x H 1.5 in. Weight: 1 lb 8 oz

C2E ICE PACK BUNDLE

- Four Comfort Gel Packs
- Facilitates Conductive, Convective, and Evaporative (C2E) cooling of high blood flow (large artery) regions of the body



ORDERING INFO SPECIFICATIONS

ITEM #	Indiv. Gel Pk: W 6 in. x L 10 in.
ZZ-0675	Pkg: L 8 in. x W 9.5 in. x H 1.25 in. Weight: 4 lb 4 oz

COOLING COLLAR



- Two gel packs encased in integrated pouches
- Designed to align with the carotid artery

ORDERING INFO SPECIFICATIONS

ITEM #	Indiv. Gel Pk: W 6 in. x L 10 in.
ZZ-0670	Pkg: L 8 in. x W 9.5 in. x H 1.3 in. Weight: 4 lb 4 oz



Splinting

TC SPLINT *Telescoping Collapsible Splint*

- Multi-capable device for splinting a wide variety of major and minor fractures
- Lightweight and compact for easy field carry with 6 easily adjustable interconnected telescoping poles and 4 wide elastic bands with hook & loop securing
- Includes clearly marked directions on the device as well as an easy-to-use, step-by-step instructional insert

SPECIFICATIONS

Pkgd: L 8 in. x W 5 in. x D 3 in.
 Open: W 9.5 in. x L 34 in.
 Weight: 1 lb 4 oz

ORDERING INFO

ITEM #
50-0039



online video available
 < scan to view



COMBAT READY™ *Tactical Body Tape*

Average of **60% pain reduction** observed in over 1400+ patient encounters

- Easy-to-use solution for reduced pain and swelling
- Precut strips designed for specific body area, can be worn for multiple days
- Allows the user to stay in the fight and complete the mission



MultiPack: One each of all 5 Combat Ready® Tape

ORDERING INFORMATION

STYLE	ITEM #	NSN #
Alpha	ZZ-0630	6510-01-659-7999
Bravo	ZZ-0631	6510-01-659-7998
Charlie	ZZ-0632	6510-01-659-8001
Delta	ZZ-0633	6510-01-659-8353
Echo	ZZ-0634	6510-01-659-8349
MultiPack	ZZ-0635	6510-01-659-3942

SPECIFICATIONS

Alpha, Bravo, Charlie, Delta, Echo Packaged:
 H 8 in. x W 5 in. x D 0.2 in.
 Weight: 1.5 oz

MultiPack Packaged:
 H 8 in. x W 5 in. x D 1 in.
 Weight: 7.5 oz



TREATMENT SUPPORT | Splinting

TACTICAL TRACTION SPLINT (TTS) Compact Mechanical Traction Device

SPECIFICATIONS

Packaged:
L 11 in. x W 5.5 in. x D 3 in.
Weight: 1 lb 5 oz

CE

ORDERING INFO

COLOR	ITEM #	NSN #
Black	50-0015	6515-01-346-9186
Coyote	50-0033	----
Digital Camo	50-0034	----
OD Green	50-0032	----

- Relieve pain, increase patient comfort and reduce the risk of vascular injury or nerve damage in femur fractures

- Lightweight and portable kit with clearly marked elastic straps for the proper fit to the thigh, knee & shin

TACTICAL TRACTION SPLINT (TTS):



QTY 1 Traction Pole



QTY 1 Ankle Hitch



QTY 1 Shin Strap



QTY 1 Knee Strap



QTY 1 Thigh Strap



QTY 1 Thigh System

SAM® SPLINT II

- Versatile, malleable splint easily molds to fit most fractures
- Lightweight, comfortable & strong

- Waterproof, Reusable

SPECIFICATIONS

H 9.25 in. x W 4.25 in. x D 1 in.
Unfolded: L 36 in. x W 4.25 in.
Weight: 4.5 oz

ORDERING INFO

ITEM #	NSN #
50-1005	6515-01-494-1951

SAM® PELVIC SLING for Pelvic Ring Fractures

- Force-controlled, ideal compression every time without the guesswork
- Rapid, simple, three-step application, usually in less than one minute



ORDERING INFO

ITEM #	NSN #
50-0004	6515-01-509-6866

SPECIFICATIONS

H 7.25 in. x W 7.75 in. x D 4 in.
Hip Circumference: 32 in. to 50 in.
Weight: 9.8 oz

SAM® SOFT SHELL SPLINT

SPECIFICATIONS

L 15 in. x W 4.25 in. x D 0.5 in.
Weight: 3.2 oz

ORDERING INFO

UNIT	ITEM #	NSN #
Each	50-0008	6515-01-523-5440
Box of 20	50-0009	6515-01-507-8074

SCISSOR LEASH™

SPECIFICATIONS

H 6 in. x W 16 in. x D 0.6 in.
Cable Length - 36 in.
Weight: 16 oz

ORDERING INFO

ITEM #	NSN #
70-0005	6515-01-540-7226

TRAUMA SHEARS

SPECIFICATIONS

SMALL: L 6.25 in. x W 2.75 in. x D 0.3 in.
Weight: 1.5 oz
LARGE: L 7.25 in. x W 3.7 in. x D 0.4 in.
Weight: 2 oz

ORDERING INFO

* NATO NSN

SIZE	ITEM #	NSN #
S 6.25 in.	ZZ-0064	6510-01-622-1624*
L 7.25 in.	ZZ-0063	6515-01-538-9276



6 oz MOLLE Retractable GEAR LEASH

- MOLLE Mount compatible / 360° rotating
- High break strength. salt water resistant self-flushing
- 36 in. Extension / Spectra Nylon Line

ORDERING INFO

ITEM #
♦ 70-0155



SPECIFICATIONS

Retracted:
H 3.5 in. x W 1.6 in. x D 1 in.
Cable Length - 36 in.
Weight: 2.0 oz



NITRILE GLOVES

Nitrile gloves available in Black or Tan colors and come either as pre-packaged pairs (kits) or in bulk (box)

BEAR CLAW®



- Form-fitting, ambidextrous design with textured fingertips
- Thick palm ensures durability
- Textured for improved wet and dry gripping
- Longer and thicker than standard nitrile gloves

BLACK TALON®



SPECIFICATIONS

Kit: L 4.75 in. x W 6 in. x D 4 in.
Weight: 1 lb 1 oz
Box: L 11 in. x W 5 in. x D 3.75 in.
Weight: 2 lb 4 oz

Bear Claw ORDERING INFORMATION

SIZE	PKG QTY	ITEM #	NSN #
S	Kit 25 pair	ZZ-0246	6515-01-619-5350
S	Box 50 pair	70-0009	6515-01-600-6059*
M	Kit 25 pair	ZZ-0247	6515-01-619-5689
M	Box 50 pair	70-0010	6515-01-600-6057
L	Kit 25 pair	ZZ-0248	6515-01-600-0262
L	Box 50 pair	70-0011	6515-01-600-6039
XL	Kit 25 pair	ZZ-0249	6515-01-619-5367
XL	Box 45 pair	70-0012	6515-01-600-6052*

Black Talon ORDERING INFORMATION

SIZE	PKG QTY	ITEM #	NSN #
S	Kit 25 pair	ZZ-0218	----
S	Box 50 pair	70-0001	6515-01-515-0195
M	Kit 25 pair	ZZ-0217	6515-01-672-7977*
M	Box 50 pair	70-0002	6515-01-515-0200
L	Kit 25 pair	ZZ-0216	6515-01-632-0059
L	Box 50 pair	70-0003	6515-01-515-0197
XL	Kit 25 pair	ZZ-0215	6515-01-672-7985*
XL	Box 45 pair	70-0004	6515-01-515-0190

NAR RESPIRATOR

- NIOSH N95 certified particulate barrier
- Adjustable nose bridge for a secure, comfortable fit

SPECIFICATIONS

Packaged:
L 5 in. x W 4.25 in. x H 0.75 in.
Mask size: Adult
Weight: 0.9 oz

ORDERING INFO

ITEM #	NSN #
70-0006	6515-01-537-1531



NAR PPE KIT

Certified NIOSH N95 disposable particulate respirator & 2 pair Nitrile Gloves

SPECIFICATIONS

L 3.5 in. x W 4 in. x H 1.5 in.
Weight: 2 oz

ORDERING INFO

ITEM #	NSN #
70-0007	6515-01-537-1534



BIOLOGICAL PERSONAL PROTECTION KIT

Designed to maximize personal protection when handling and managing detainees, refugees, etc.

SPECIFICATIONS

H 6 in. x W 8 in. x D 1 in. | Weight: 4.6 oz

BIOLOGICAL PERSONAL PROTECTION KIT:



ORDERING INFO

ITEM #	NSN #
80-0488	6515-01-388-1351



TREATMENT SUPPORT | Eye Care & Protection

PES Polycarbonate Eye Shield

- Meets CoTCCC treatment protocols for penetrating eye trauma
- Convex eye shield accommodating either eye
- Translucent orange shell features perforations to allow air circulation



SPECIFICATIONS

3 in. x 2.38 in.
Weight: 0.2 oz



ORDERING INFO

UNIT	ITEM #	NSN #
Single	30-0142	6515-01-647-5456
Box of 100	30-0143	6515-01-630-5551

FOX Eye Shield

- Convex rigid eye shield made of lightweight aluminum with perforations allowing air circulation
- Meets CoTCCC treatment protocols for penetrating eye trauma
- Available with or without garter



SPECIFICATIONS

3 in. x 2.38 in.
Weight: 0.2 oz

ORDERING INFO

STYLE	UNIT	ITEM #	NSN #
No Garter	Box of 12	30-0144	6515-01-449-1016
Garter	Single	30-0055	----
Garter	Box of 50	30-5550	6515-01-598-1877

NAR IPRO® Tactical Goggle System

High performance eye protection with Hexo-vents for maximum ventilation, fine particle filtration and fog-free viewing

SPECIFICATIONS

Packaged:
L 2.6 in. x W 7 in. x D 3 in.
Weight: 8 oz

ORDERING INFO

ITEM #
ZZ-0109



AEGIS ECHO II Ballistic Glasses & Case

- Superior ballistic impact protection per US MIL-PRF-31013 standard clause 3.5.11 and 100% protection from UVA/UVB rays
- Includes clear & gray interchangeable lenses

SPECIFICATIONS

Glasses: W 5.5 in. x H 2 in. x D 4.75 in.
Case: L 6.5 in. x W 4.3 in. x H 3 in.
Weight: 9 oz

ORDERING INFO

ITEM #
90-0067

OTW GOGGLES Outside the Wire - with Case

- 100% UVA/B/C protection with maximum ballistic protection
- Tapered, neutral tint gray lenses with anti-fog and scratch resistant coating, reduced glare and correct distortion

ORDERING INFO

ITEM #
90-0066

SPECIFICATIONS

Goggles: W 7.5 in. x H 2.75 in.
Case: L 8.75 in. x W 3.5 in. x H 3 in.
Weight: 10.1 oz



GRAY MAN EyePro

Low visibility EyePro option with classic, subtle styling

ORDERING INFO

ITEM #
90-0068



SPECIFICATIONS

Glasses: W 5 in. x H 2 in. x D 5 in.
Case: L 6.625 in. x W 3 in. x H 2.5 in.
Weight: 9 oz



BURNTEC® Hydrogel Dressings

- Sterile non-stick hydrogel dressings for treating 1st, 2nd & 3rd-degree burns (thermal, chemical, electric)
- Also treats other skin injuries, such as sunburns, abrasions, skin irritations, bruises and swelling (insect bites)
- Applied dressings provide a cool, soothing effect for up to 24 hours
- Easy to handle without spilling or damaging the hydrogel
- Permeable to water vapor and oxygen but impermeable to bacteria
- Biocompatible and hypo-allergenic; Reduces risk of hypertrophic scars

ORDERING INFORMATION

CE

SIZE	ITEM #	NSN #
2 in. x 2 in.	30-0114	6510-01-631-1913
2.5 in. x 5 in.	30-0115	6510-01-631-2080
4 in. x 4 in.	30-0110	6510-01-631-1883
5 in. x 5 in.	30-0116	6510-01-631-2138
5 in. x 10 in.	30-0117	6510-01-631-1907
8 in. x 8 in.	30-0118	6510-01-631-1897
8.5 in. x 11 in.	30-0119	6510-01-631-1890
4 in. x 16 in.	30-0111	6510-01-631-2124
8 in. x 16 in.	30-0112	6510-01-626-2990
16 in. x 24 in.	30-0120	6510-01-631-1886
10 in. x 10 in. Face	30-0121	6510-01-627-1834
12 in. x 16 in. Face	30-0113	6510-01-627-1855
16 in. x 24 in. Face	30-0122	----

DRY STERILE BURN DRESSING

- Dry sterile, non adherent burn dressing & bandage
- Includes two No. 3 closed safety pins and rubber band individually packaged



SPECIFICATIONS

Bandage: 45 in. x 45 in. x 63 in.
Pkg: W 4.5 in. x H 4 in. x D 0.5 in.
Weight: 1 oz

ORDERING INFORMATION

ITEM #	NSN #
30-0162	6510-01-587-6579



BurnTec® BURN DRESSING KIT

- Convenient array of BurnTec® dressings in a handy carrying case
- Kit contains (1) each of the following dressings:

2 in. x 2 in.	5 in. x 10 in.	8 in. x 16 in.
2.5 in. x 5 in.	8 in. x 8 in.	10 in. x 10 in. (face)
4 in. x 4 in.	4 in. x 16 in.	12 in. x 16 in. (face)
5 in. x 5 in.	8.5 in. x 11 in.	

ORDERING INFORMATION

ITEM #	NSN #
80-0361	6515-01-630-9177

SPECIFICATIONS

CE

L 12 in. x W 15.5 in. x D 3.5 in.
Weight: 6 lb 4 oz



online video
available

Λ scan to view

BurnTec® MINOR BURN DRESSING KIT

- Variety of BurnTec® dressings in a handy carrying case
- Kit contains the following:
 - Two 4 in. x 4 in. dressings
 - Two 4 in. x 16 in. dressings
 - One 8.5 in. x 11 in. dressing

ORDERING INFORMATION

ITEM #	NSN #
80-0415	6510-01-644-5427

SPECIFICATIONS

CE

L 14 in. x W 8.5 in. x D 2.5 in.
Weight: 2 lb 13 oz





TREATMENT SUPPORT | Triage System

T2® TAG NAR


- Triage & treatment tag utilizing the S.A.L.T. Mass Casualty Triage Standard
- To aid in documenting critical patient information throughout the entire mass casualty response process
- Liquid-resistant in all settings; made of rugged materials contained in a plastic sheath with latex-free band

ORDERING INFORMATION


UNIT	ITEM #	NSN #
20 Pack	20-0031	6515-01-624-7148
1000 Pack	20-0032	6545-01-592-3261

SPECIFICATIONS

L 5.69 x W 3.88 in.
Weight: 0.7 oz



T2 TAG:



TACTICAL TRIAGE KIT T2®

- Compact incident command system
- Streamlines medical triage
- Contained in a lightweight, rugged vacuum-sealed pouch


SPECIFICATIONS

L 5.25 x W 7 in. x D 1.25 in.
Weight: 4.9 oz

ORDERING INFORMATION

ITEM #	NSN #
20-0040	6545-01-591-2134

TACTICAL TRIAGE KIT:



COMBAT CASUALTY CARD

- Streamlined triage card that meets all TCCC standards for casualty management
- State-of-the-art ruggedized paper is tear resistant, waterproof & designed to fit into the IFAK
- Includes grommet reinforced punch with latex-free band for attachment
- Currently in use by Special Operations Forces

SPECIFICATIONS

H 6 in. x W 4.5 in. each
Weight: 0.22 oz ea | 5.5 oz 25-Pack

ORDERING INFORMATION

UNIT	ITEM #	NSN #
25 Pack	ZZ-0262	6515-01-587-1700



COMBAT CASUALTY REFERENCE CARD

- Rugged reference card for emergency management featuring the current TCCC Guidelines
- Clear, concise instructions for medics and warfighters
- Information divided into four sections:
 - Care Under Fire
 - Tactical Field Care
 - TACEVAC Protocols
 - 9-line MEDEVAC request/notes
- Pocket sized for convenience
- Twenty-two pages

SPECIFICATIONS

H 6 in. x W 5 in. x D 0.1 in.
Weight: 2.3 oz

ORDERING INFORMATION

ITEM#	NSN #
ZZ-0057	6515-01-537-4161






TRIAGE OFFICER KIT T2®



T2 TRIAGE OFFICER KIT:



- Mid-range solution providing initial triage of up to 20 casualties
- Compact bag with mesh storage pockets equipped with T2/SALT triage cards
- Includes a laminated S.A.L.T. Triage Flowchart & Casualty Count Card

ORDERING INFORMATION

COLOR	ITEM #	NSN #
OD Green	20-0023	6545-01-537-4189
Digital Camo	20-0075	----
Red	20-0076	----

SPECIFICATIONS

Pkg: H 10 in. x W 8 in. x D 5 in.
Weight: 1 lb 14 oz

COMMAND KIT T2®

- Organized, rapid response system for Incident Commanders during a mass casualty/tactical triage incident
- 3 laminated, detachable Command Boards
- Pad of 50 Patient-Tracking Boards
- 3-ring, clear storage sheets for bar codes & mesh pockets to store extra accessories

SPECIFICATIONS

Case Closed: L 21 in. x W 17 in. x D 6 in.
Weight: 4 lb 14 oz

T2 COMMAND KIT:



ORDERING INFORMATION

COLOR	ITEM #	NSN #
OD Green	20-0025	6545-01-537-4197
Digital Camo	20-0074	----
Red	20-0026	----



Mass Casualty Incident Management MCIM KIT T2®

- Comprehensive triage system ideal in disaster response, field hospitals & military troop concentrations.
- Lightweight, modular, man-portable triage kit.
- includes user-friendly tools for every aspect of the uniformed S.A.L.T. model for disaster triage.

SPECIFICATIONS

H 21 in. x W 17 in. x D 9 in. Weight: 15 lb 2 oz

ORDERING INFORMATION

COLOR	ITEM #	NSN #
OD Green	20-0029	6545-01-537-3138
Digital Camo	20-0073	----
Red	20-0077	----

T2 MCIM:





TREATMENT SUPPORT | Apparel & Accessories

OPERATORS CAP NAR

- 100% cotton
- Adjustable strap
- 3.5 in. by 2.5 in. cut loop on the front designed for attaching patches (sold separately)
- 1 in. square of cut loop on the top for attaching an IR square

SPECIFICATIONS

One size fits most
Weight: 3.2 oz

ORDERING INFO

ITEM#

ZZ-0092



EMBROIDERED LOGO CAP NAR

- All American Cap
- Made in the USA
- Comfortable, attractive Fit-4-You™ design fits most wearers
- American flag on the front, NAR logo on back

SPECIFICATIONS

Fit-4-You™ S/M
Weight: 3.2 oz

ORDERING INFO

ITEM#

ZZ-0176



MEDIC CROSS PATCHES PVC

- Soft flexible PVC rubber patch
- Backed with hook material for easy attachment and removal
- Ideal for use on medical kits such as the M-FAK and T.O.R.K.
- 1 in. x 1 in. available at NARescue.com

ORDERING INFORMATION

COLORS	ITEM #
Coyote & Black	ZZ-0797
Black & Gray	ZZ-0798
OD Green & Black	ZZ-0800
Black & Orange	ZZ-0799



SPECIFICATIONS

L 2 in. x W 2 in. x D 0.125 in.
Weight: 0.6 oz

IR FLAGS NAR

- Clearly visible under IR light
- Hook-and-loop attachment
- Sold in Packs of 2

SPECIFICATIONS

L 2 in. x W 2 in. x D 0.125 in.
Weight: 0.6 oz

ORDERING INFORMATION

COLOR	STARS	ITEM #
Full Color	L	ZZ-0093
Coyote	L	ZZ-0094
Coyote	R	ZZ-0235
OD Green	L	ZZ-0107



MED ID PATCH LUMINOUS

Glow in the Dark

- Soft flexible custom rubber PVC
- Ultimate glow-in-the-dark identification technology
- 1 in. block letters
- 2 medical patches per package

SPECIFICATIONS

Each: L 3 in. x W 2 in. x D 0.125 in.
Weight: 0.4 oz

ORDERING INFO

ITEM#

ZZ-0121



IR SQUARES NAR

- Infrared patch printed on honeycomb reflective material
- Black hook fasteners on the back for easy attachment
- Clearly visible under IR light illumination
- Six squares per pack

SPECIFICATIONS

Each: 1 in. x 1 in.
Weight: 0.2 oz.

ORDERING INFO

ITEM#

ZZ-0108





Discover a broad range of educational opportunities for first responders and professional emergency services providers. Earn and maintain your critical skills certifications for EMT, take exam prep courses, or prepare for mass casualty situations in the civilian world by completing bleeding control courses.

The NAR Training instructors offer solid content that is based on proven medical science. We are proud to offer course series that cover approved material from The American Heart Association®, National Registry of EMTs, and National Association of EMTs (NAEMT). See the next few pages and our website for course offerings and other details.



NORTH AMERICAN RESCUE™

Education and Training

Explore the options at
[NARescue.com/](http://NARescue.com/training-education)
training-education

Wherever the
MISSION
takes you...





TRAINING MENU:

Public Access Series:

• NAR Equipment Training • Bleeding Control

Distance Learning Series:

• Online CME • Distance Education

NREMT Series:

• Paramedic • AEMT • EMT • Test Prep Workshops

American Heart Association Series:

• BLS • ACLS • PALS

NAEMT Series:

• TCCC • TECC • PHTLS

Wherever the
MISSION
takes you...

PUBLIC SAFETY



Discover a broad range of educational opportunities for first responders and professional emergency services providers. Earn and maintain your critical skills certifications for EMT, take exam prep courses, or prepare for mass casualty situations in the civilian world by completing bleeding control courses.

The NAR Training instructors offer solid content that is based on proven medical science. We are proud to offer course series that cover approved material from The American Heart Association®, National Registry of Paramedics, and NAEMT. Be sure to visit our website for the latest course schedules so that you can reserve your spot.



OPERATIONAL TACTICAL

Tapping into the in-depth knowledge developed through a commitment to knowledge and real-world know-how, the instructors from NAR Training provide educational experiences that other training organizations can't compare to. Understanding the tactical environment brings a whole new dimension to medical training.

NAR Training delivers a broad range of curricula that will frequently include hands-on skills stations and scenarios designed to develop competence and confidence for handling injuries in the field. From Tactical Combat Casualty Care (TCCC) and TECC to Fed Med and SOF Med, course offerings cover the broad spectrum of considerations for urban, rural and other austere environments while delivering lifesaving measures.

MISSION SPECIFIC

Knowledge and information will only prepare you so much for the rigors of real-world stress experienced in combat. In order to be better prepared for the problem-solving challenges of your operational environment, NAR Training delivers customized training that captures the requirements for ensuring optimal outcomes for you and your entire team.

With a proud history of working with Tier One teams to help sharpen their skills and prepare them for the physical, mental, and emotional challenges that they will experience in the field, these courses cover the range of instructor-led classroom lecture all of the way to fully immersive scenario-driven experiential training events that empower students to develop superior skills and tactics. Contact our team today to discuss your specific needs and prepare for optimal outcomes.

COURSES:

The RAPTOR

Advanced Pre-Deployment Combat Medic Workshop / 11 Days

SOF MED

Special Operations Forces Medical Refresher / 8 Days

FED MED:

US Federal Agent EMS Refresher & TECC Provider / 5 Days

EGT & TCCC

Emergency Ground Tactics & TCCC Provider / 5 Days

The PREDATOR

Point of Injury & Prolonged Field Care Workshop / 5 Days

SUPPORT SOLUTIONS:

CASSIM & SFX

Professional Actors/ Amputees with Special Effects Teams to Enhance Training Events

NAR Simulation

TOMManikin, CASEVAC Simulators, Battlefield Effects brought to your location



TRAINING | TOMManikin™ Training Simulators

TOMMANIKIN™ Series

State-of-the-art training simulators, the TOMManikin - Tactical Operation Medical Manikin provide the ability to enhance training in Tactical Combat Casualty Care (TCCC) from Point of Injury to transfer of higher care

- All-in-one, hyper-realistic trauma simulators enabling instructors to build multiple realistic traumatic casualty scenarios with reactive patients - ideal for use in classroom settings, trauma lanes and field training scenarios
- Multiple, interchangeable extremities present a variety of wounding patterns such as gunshot, blast and burns to provide full mission profiles in combat scenarios, disaster re-creations and mass casualty events (See website for additional extremity options not shown here)
- Constructed of silicone and welded steel
- Weighs more than 160 pounds with extremities that can withstand up to 6,000 pounds of pressure
- Rugged remote controlled simulators the ideal choice for confined space, building collapse, care under fire, and entrapment training scenarios
- Intervention capabilities on all manikins include tourniquet application, airway management (basic to advanced), needle decompression, chest tube, IV/IO, and wound packing

Every TOMManikin™ includes:

- Equipped with Advanced Remote Tablet and Software enabling wireless control of up to 12 TOMManikin units with capabilities of bleeding, breathing, multiple voice tracks, pulse rates, pulse pressure and realistic responses to responder interventions
- Presets, scenario builder software and exportable student evaluation checklists included
- Onboard power via a 12v / 5 amp & 10 amp standard rechargeable nickel metal hydride (NiMH) batteries with smart chargers included
- 1.2 gallon onboard blood reservoir facilitating massive bleeding from multiple locations
- 1 x 1400 Black Pelican Accessories Case
- 1 x T handle tool (for extremity changing) and various air & fluid replacement connectors
- 1 x TOMManikin Operator's Manual
- 1 x 1-year warranty document

BASIC TOMManikin™

Intervention capabilities:

- Tourniquets (proximal, junctional, and abdominal tourniquets)
- Active bleeding controllable by direct pressure on pressure points
- Wound packing with realistic bleeding control
- Nasopharyngeal and oral airways, supraglottic airway insertion and endotracheal intubation
- Surgical airway with realistic bleeding and tissue dissection | 35 operative sites included
- Remote controlled breathing rates: Slow, Normal, and Fast
- Bilateral Needle Decompression at either 2nd/3rd intercostal (mid-clavicular), or 5th/6th intercostal (mid-axillary) locations
- Realistic chest seal application & adhesion
- Chest Tube capable (bilateral)
- Radial pulses in uninjured upper extremities & bilateral carotid pulses
- Variable pulse rate & pressure (Normal, bounding, thready)
- IO (sternal & humoral) and antecubital IV access on uninjured upper extremities

BASIC TOMManikin™:

- 1 Partial amputation to arm (tourniquet-responsive)
- 2 Full amputation to leg (tourniquet-responsive)

SPECIFICATIONS

Pelican Accessories Case:
L 13.74 in. x W 11.63 in. x D 5.9 in.
Weight: Manikin - 180 lb
Total package: 240 lb

1

2

ORDERING INFORMATION

ITEM

◆ 93-0035

Intervention capabilities:

- Tourniquets (proximal, junctional, and abdominal tourniquets)
- Active bleeding controllable by direct pressure on pressure points
- Wound packing with realistic bleeding control
- Nasopharyngeal and oral airways, supraglottic airway insertion and endotracheal intubation
- Surgical airway with realistic bleeding and tissue dissection - 35 operative sites included
- Remote controlled breathing rates: Slow, Normal, and Fast
- Bilateral Needle Decompression at either 2nd/3rd intercostal (mid-clavicular), or 5th/6th intercostal (mid-axillary) locations
- Realistic chest seal application & adhesion
- Chest Tube capable (bilateral)
- Radial pulses in uninjured upper extremities and bilateral carotid pulses
- Variable pulse rate & pressure (Normal, bounding, thready)
- IO (sternal & humoral) and antecubital IV access on uninjured upper extremities
- Packable abdominal evisceration

BLAST TOMManikin™:

- 1 *Blast or burned head skin options*
- 2 *Packable abdominal evisceration*
- 3 *Shrapnel wound right arm*
- 4 *Full amputation left leg (tourniquet-responsive)*



SPECIFICATIONS

Pelican Accessories Case:
L 13.74 in. x W 11.63 in. x D 5.9 in.
Weight: Manikin - 180 lb
Total package: 240 lb

ORDERING INFORMATION

ITEM

◆ 93-0036



All TOMManikin Simulators include a remote tablet, rechargeable batteries, instruction manual and additional tools all in a convenient accessory case

◆ New Product



GSW (Gun Shot Wound) TOMManikin™

Intervention capabilities:

- Tourniquets (proximal, junctional, and abdominal tourniquets)
- Active bleeding controllable by direct pressure on pressure points
- Wound packing with realistic bleeding control
- Nasopharyngeal and oral airways, supraglottic airway insertion and endotracheal intubation
- Surgical airway with realistic bleeding and tissue dissection. 35 operative sites included
- Remote controlled breathing rates: Slow, Normal, and Fast
- Bilateral needle decompression; realistic sucking chest wound (with mid-axillary exit wound)
- Realistic chest seal application & adhesion
- Chest Tube capable (bilateral)
- Radial pulses in uninjured upper extremities and bilateral carotid pulses
- Variable pulse rate & pressure (Normal, bounding, thready)
- IO (sternal & humoral) and antecubital IV access on uninjured upper extremities
- Through & through gunshot wound to right leg
- Pelvic & leg injuries responsive to direct pressure and tourniquets
- Packable gunshot wound to groin (with posterior exit wound)

GSW TOMManikin™:

- 1 Sucking chest wound, right (with anterior entry/posterior exit)**
- 2 Junctional tourniquet/packable GSW to left groin (anterior entry/posterior exit)**
- 3 Tourniquet-responsive/packable GSW right leg (through & through)**

SPECIFICATIONS

Pelican Accessories Case:
L 13.74 in. x W 11.63 in. x D 5.9 in.
Weight: Manikin - 180 lb
Total package: 240 lb

ORDERING INFORMATION

ITEM

◆ 93-0040





K9 CANINE MANIKIN *Training Simulator*

Intervention capabilities:

- The look and feel of a real dog with realistic limb movement
- Realistic weight of approximately 60 lb
- Life-like airway allows for manual manipulation and accepts 10-12mm ET tubes
- Bag ventilation will cause chest to rise and fall
- Stomach will rise and fall if ventilated into the esophagus
- IV training site with flash feedback to show success
- Shaveable replacement fur for IV sites
- Bleeding wounds facilitate various bleeding control methods
- Tracheostomy simulator featuring shaveable fur with multi-layered tissue insert
- Bilateral Needle decompression/thoracocentesis site allow release of air
- Adjustable femoral pulses – features variable rate, frequency and pulse pressure
- Lacerated paw for bandage training – foreleg, hind leg
- Five interchangeable limbs with different injuries – complete amputation, open fracture, soft tissue injuries, road rash/burns
- Customized wireless remote control
- Durable fur makes allows easy care and maintenance
- Rugged exterior for use in all training locations and situations
- Operated by a single onboard rechargeable NiMH battery
- Includes pulse ox
- 1-year warranty
- Customizable wounds available for same price as original build
- Two fur style options: Malinois-type or long coat black

SPECIFICATIONS

Pelican Accessories Case:
L 13.74 in. x W 11.63 in. x D 5.9 in.
Weight: Manikin - 60 lb



K9 Canine Manikin 5 Interchangeable Limbs:

- 1 Complete amputation
- 2 Open fracture
- 3 Soft tissue injuries
- 4 Road rash/burns
- 5 Gun shot wound

ORDERING INFORMATION

ITEM #
◆ 93-0027



TCCC SKILLS SYSTEM

All Combatant

- All the necessary tactical medical and casualty extraction/evacuation equipment needed when teaching the TCCC – All Combatant Skills required tasks for completing the Tactical Combat Casualty Care (TCCC) course
- A turnkey solution for training centers and educational institutions designed to equip up to twenty-five (25) students
- Complete system consists of three self-contained Hard Storage Case modules: Skill Stations | IFAKS | Evacuation

SPECIFICATIONS

Hard Storage Case:

External – L 32.5 in. x W 21 in. x D 13 in.

Internal – L 29 in. x W 18 in. x D 10 in.

Weights: TCCC Skills Set: 55 lb

TCCC IFAK Skills Set: 70 lb

TCCC Evacuation Skills Set: 50 lb

ORDERING INFO

UNIT	ITEM #
3 Modules	80-0970

TCCC SKILLS SET

All Combatant

- Train up to twenty-five (25) students at individual skill stations in the treatment of: *extremity hemorrhage, wound packing, chest seal application, needle chest decompression, basic airway adjuncts, splinting, hypothermia prevention & management, and eye injuries*
- Self-contained, wheeled, rigid case with internal dividers to keep the equipment organized and easy to inventory

SPECIFICATIONS

Hard Storage Case:

External – L 32.5 in. x W 21 in. x D 13 in.

Internal – L 29 in. x W 18 in. x D 10 in.

Weight: 55 lb

ORDERING INFORMATION

MODULE	ITEM #
TCCC Skills Set	80-0968



TCCC SKILLS SET MODULE:





All Combatant **TCCC IFAK SKILLS SET**



- Twenty-five (25) pre-packed Training IFAKs for use during trauma lanes or training scenarios
- IFAKs come packaged in a self contained wheeled hard, case that allows for easy transport and storage

SPECIFICATIONS

Hard Storage Case:
External – L 32.5 in. x W 21 in. x D 13 in.
Internal – L 29 in. x W 18 in. x D 10 in.
Weight: 70 lb

ORDERING INFORMATION

MODULE	ITEM #
TCCC IFAK Skills Set	80-0969

25 TRAINING IFAKs
each
containing:



QTY **1** C-A-T®
Tourniquet, Trainer



QTY **1** NARETD™
6 in. Flat



QTY **1** Gauze
S-Rolled



QTY **1** Combat
Gauze® Trainer



QTY **1** 28F Naso
Airway & Lube



QTY **1** HyFin®
Trainer, Twin Pk



QTY **1** ARS®
Needle, 14g



QTY **1** Surgical
Tape, 2 in.



PAIR **1** Nitrile
Gloves, Lg



QTY **1** Permanent
Marker, Large



QTY **1** Combat
Casualty Card

TCCC EVACUATION SKILLS SET



- Various high threat extraction devices and casualty evacuation systems to train students in different options to move casualties during the all three TCCC Phases of Care
- Packaged in a self-contained, wheeled, rigid case that allows for easy transport and storage

SPECIFICATIONS

Hard Storage Case: External –
L 32.5 in. x W 21 in. x D 13 in.
Internal – L 29 in. x W 18 in. x D 10 in.
Weight: 50 lb

ORDERING INFORMATION

MODULE	ITEM #
TCCC Evacuation Skills Set	80-0971

TCCC EVACUATION SKILL SET MODULE:



QTY **1** HPMK®



QTY **1** Talon®
90C Litter



QTY **1** R.E.E.L.
Litter Black



QTY **1** Hasty
Harness



QTY **1** Rat™
Strap



QTY **1** Dragon
Handle System



QTY **1** Shield
Carabiner



TRAINING | TCCC All Combatant Training Components

TCCC IFAK SKILLS KIT *All Combatant*

- The necessary training equipment for initial TCCC care using the Module TCCC IFAK Skills Set - All Combatant
- Easy to attach modular carrying bag is compact and designed to be worn on the assault vest or body armor via MOLLE-style attachment points

TCCC IFAK SKILLS KIT:




SPECIFICATIONS

H 7 in. x W 4.5 in. x D 2.75 in.
Weight: 2 lb 8 oz

ORDERING INFO

OPTION	ITEM #
TCCC IFAK Skills Kit	85-1958



TRAINING | Reference Handbooks

A portion of the proceeds from these handbooks is donated to wounded warrior foundations supporting Special Operations soldiers & their families



FLIGHT PARAMEDIC *Certification & Review Manual*

154 page textbook, thoroughly covering ALL required areas of knowledge necessary to be successful on the Flight Paramedic Certification (FP-C) Exam

SPECIFICATIONS

H 11 in. x W 8.5 in. x D 0.5 in.
Weight: 1 lb 2 oz
154 pages

ORDERING INFO

ITEM #
ZZ-0467

TACTICAL PARAMEDIC *Certification & Practical Application*

182 page training guide covering all materials required for achieving Tactical Paramedic Certification (TP-C) by the Board For Critical Care Transport Paramedic Certification (BCCTPC)

SPECIFICATIONS

H 10.9 in. x W 8.4 in. x D 0.4 in.
Weight: 1 lb 5 oz
182 pages

ORDERING INFO


ITEM #
ZZ-0617



C-A-T® TRAINER

Combat Application Tourniquet

- Fully functioning C-A-T® tourniquet colored blue to identify it as a training device
- Enables instructors and trainees hands-on training & tourniquet application practice
- Providing a distinct visual means to distinguish training devices from personal/operational life-saving devices



ORDERING INFO

ITEM #	NSN #
30-0033	6910-01-560-2972

SPECIFICATIONS

H 6.5 in. x W 2.4 in. x D 1.5 in.
Weight: 2.7 oz

See page 6 for Tactical and Rescue versions of the C-A-T®

HYFIN® TRAINER

Twin Pack

- Unique adhesive gel formulation does not damage mannequins or leave behind excessive residue
- Packaged with two training chest seals in one unique perforated package allowing students to open only one seal at a time



SPECIFICATIONS

Packaged:
L 7.5 in. x W 4.5 in. x D 0.25 in.
Chest Seal:
H 6 in. x W 6 in. x D 0.1 in.
Weight: 2.5 oz

ORDERING INFO

ITEM #
10-0032

See page 12 for other HyFin® Chest Seal products

Responder

TRAINING ETDs

- Economic, low cost solution for training in the application of immediate direct pressure
- Non-sterile, elastic pressure bandages

SPECIFICATIONS

4 inch ETD:

Bandage: L 60 in. x W 4 in.
Pad: H 6 in. x W 4 in.
Packaged:
L 4.5 in. x W 2.5 in. x D 1.75 in.
Weight: 2.0 oz

6 inch ETD:

Bandage: L 60 in. x W 6 in.
Pad: H 6 in. x W 6 in.
Packaged:
L 6 in. x W 2.5 in. x D 2.1 in.
Wt: 2.9 oz

ORDERING INFO

SIZE	ITEM #
4 in.	30-0145
6 in.	30-0146



COMBAT GAUZE®

Mouflage Trainer

- Thoroughly instruct students on the proper application of Z-Fold Combat Gauze®
- Completely inert (no hemostatic properties) - training purposes only
- Re-sealable packaging for multiple uses

SPECIFICATIONS

Packaged:
H 6 in. x W 4.5 in. x D 1.25 in.
Gauze:
L 4 yd x W 3 in.
Weight: 0.75 oz

ORDERING INFO

ITEM #
30-0063





TRAINING | Part Task Trainers

HEMORRHAGE CONTROL TRAINING KIT *with Bio Feedback*

- Training tool designed to teach effective wound packing skills with Immediate results for both trainer and student
- Verifies correct or incorrect technique

SPECIFICATIONS

Case: 23 in. x 8 in. x 15 in.
Weight: 17 lb 3.2 oz

ORDERING INFO

ITEM

93-0019



HEMORRHAGE CONTROL TRAINING KIT *Combat Gauze®*

- Comes with 10 Combat Gauze Trainers and 5 Combat Gauze LE Demo Samples
- Includes a "training leg" that features an entry/exit bullet wound and a deep laceration

SPECIFICATIONS

Case: H 15 in. x W 23 in. x D 8 in.
Weight: 15 lb 8 oz

ORDERING INFO

ITEM

93-0017



G3 POWER DRIVER TRAINER *by EZ-IO®*

- Red color indicates training only device
- Compact design



SPECIFICATIONS

H 4.5 in. x W 6.5 in. x D 2.5 in.
Weight: 15 oz

ORDERING INFO

ITEM

30-0076

EZ-IO® G3 Power Driver TRAINING NEEDLES

Designed for use with any EZ-IO® Power Driver



ORDERING INFO

SIZE	ITEM #
15 mm (Pediatric)	30-0082
25 mm (Adult)	30-0084
45 mm (Adult, large)	30-0086

T.A.L.O.N. MANUAL TRAINING NEEDLE

Train in manual I/O insertion at all seven optional sites, including the sternum



ORDERING INFO

ITEM

30-0133

MANUAL TRAINING NEEDLE *by EZ-IO®*

For training in six different FDA cleared sites for easily accessible insertion



ORDERING INFO

ITEM

30-0079

G3 EDUCATOR KIT *by EZ-IO®*

- G3 Training Power Driver
- EZ-IO® Training Needle Sets: 15mm Pediatric, 25mm Adult & 45mm Large Adult
- 4 Wristbands
- 4 EZ-Connects
- G3 Vascular Access Pack
- 4 EZ-Stabilizers
- Training Bones: Infant, Tibia, Adult Proximal, Humerus & Adult Tibia
- Softside Nylon Case
- Training CD



ORDERING INFO

ITEM

30-0074



TRUECLOT® Hemorrhage Training Kit™

- Introduce students to realistic hemorrhage control for TCCC, TECC, TEMS, PHTLS, Active Shooter Training and general wound packing techniques for advanced hemostatic gauze training
- Highly realistic Wound Packing Trainer features a large caliber gunshot wound with hidden bleed points and a simulated fracture/defect hard plastic bone insert to aid in teaching deep wound packing techniques
- Includes TrueClot® Clotting Blood Simulant and Simulated Hemostatic Gauze which demonstrates the reaction of Hemostatic Gauze in a wound
- Durable design features a sturdy silicone construction with no moving parts, reinforced tubing mount, and a more realistic wound pattern and outer anatomy
- Ideal for the repeated use of a training environment



SPECIFICATIONS

Task Trainer: L 7.25 in. x H 2.5 in.
Storage Case: L 13.5 in. x W 10 in. x D 5.5 in.
Liquid Capacity: 1 liter
Weight: 5 lb

ORDERING INFORMATION

ITEM #

93-0020

TRUECLOT TRAINER:



QTY
1 Lg. GSW/
Task Trainer

QTY
4 TrueClot®
Blood Simulant

QTY
9 TrueClot®
Hemostatic, 12 ft

QTY
1 Bottle,
Delivery
with Tube &
1-way Valve

QTY
1 Bottle,
Mixing & Storage

QTY
1 Hard Case
with Instructions

PRODUCT by Name index

A

Aegis Echo II Ballistic Glasses	92
Aero Screwgate Carabiner	81

AIRWAY MANAGEMENT

AMBU RDIC Resuscitator	12
AMBU SPUR II BVM	12
BAC-Pack Bougie-Aided Cricothyroidotomy Pack	10
Crickit®, Tactical	11
Cyclone® Improved Pocket BVM	12
Laryngoscope Blades,	12
King LTS-D™	11
Laryngoscope Set, Special Ops	12
Nasopharyngeal Airway	11
Tactical Suction Device, NAR	11
Tracheal Hook, NAR	11
Tracheostomy Kit with Bougie Introducer	10
AMBU RDIC Resuscitator	12
AMBU Spur II BVM	12
Anneau Slings, Petzl	81

APPAREL

Embroidered Logo Cap, NAR	96
IR Flags	96
IR Squares	96
Medic Cross Patches, PVC	96
Med ID Patch Luminous	96
Operators Cap	96
Armadillo Medication Case	67
ARS® Needle Decompression Kit	15
Autotriage 9160 SBK, Benchmark	83

B

BAC-Pack Bougie-Aided Cricothyroidotomy Pack	10
Basic Ascender, Petzl	81
Bear Claw® Gloves	91
Belt Throw Bag, Water Rescue	78

BENCHMARK® KNIVES

AUTOTRIAGE 9160 SBK	83
H ₂ O 112 SBK	83
Safety Hook Cutters	83
SOCPR Rescue Tool	83
TRIAGE 916 SBK	83
TRIAGE 917 BK	83
TRIAGE 917 SBK	83
Biological Personal Protection Kit	91
Black Talon® Gloves	91
BOA® IV Constricting Band	16

BOAT RESPONSE AID KITS

Naval Boat Response Aid Kit	39
USCG Boat Response Aid Kit	39

BOOKS/MANUALS

Flight Paramedic Manual	106
Tactical Paramedic Manual	106
Bougie Introducer, Tracheostomy Kit	10
BP/Stethoscope Combo Kit	66
Burn Dressing, Dry Sterile	93
BURNTEC® Hydrogel Dressings	93

C

CANINE (K-9)

K9 Manikin Training Simulator	103
K-9 Handler IFAK	43
K-9 Litter, NAR	58
K-9 Tactical Field Kit	43
CardioPak ECG, Resolute Medical	65
Casualty Restraint Straps, NAR	59
C-A-T® Combat Application Tourniquet	6
C-A-T® TQ Case, Rigid	7

Celox™ Gauze	8
Cervical Collar, ACE	59
ChitogGauze® PRO	8

CIRCULATION

BOA® IV Constricting Band	16
EZ-IO G3 Power Driver	17
EZ-IO G3 Power Driver Needle Sets	17
EZ IO Manual Extremity Needle Set	17
Field Blood Transfusion Kits, FBTK	20
Needleless Saline Lock Kit, NAR	16
Quantum® Blood & Fluid Warming System, Lifewarmer	18-19
RAPTOR IV Securing Device	17
Saline Lock Kit, NAR	16
Sharps Shuttle	17
T.A.L.O.N.™ Manual Needle Set	17

Clarius, Wireless Ultrasound	62
Combat Casualty Card	94
Combat Casualty Reference Card	94

COMBAT CASUALTY

RESPONSE KITS (CCRK®)

Eagle IFAK	26
Expeditionary Casualty Response Kit, ECRK	41
Individual Kit	31
K-9 Handler IFAK	43
K-9 Tactical Field Kit	43
Mariner Kit	33
Maritime Assault Kit	30
Medic Leg Rig	31
M-FAK Mini First Aid Kit	25
Multi-Mission Trauma Pack	44
NAR-4 Aid Kit	35
NAR-4 Chest Pouch	34
NAR-5 SAR Aid Kit	39
Naval Boat Response Aid Box	39
Naval Boat Response Aid Kit	39
OCHO IFAK	28
Range Trauma Aid Kit Hard Case	42
Range Trauma Kit	42
RIG Series Eagle IFAK	27
Squad Kit	32
Squad Responder Kit	40
Tiered Trauma Medical Bag (TTMB)	45
T.O.R.K.	29
USCG Boat Response Aid Kit	39
USCG IFAK	30
USMC CLS Kit	40
WALK Warrior Aid & Litter Kit	37

Combat Gauze®	8
Combat Gauze® Moulage Trainer	107
Combat Ready™ Body Tape	89
Command Kit, T2®	95
Compressed Gauze, NAR	9
Cravat (Triangular Bandage)	9
Crickit® Tactical Cricothyroidotomy Kit	11
Cutter, 7 Hook Safety	83
Cyclone® BVM	12

D

DIAGNOSTIC EQUIPMENT

BP/Stethoscope Combo Kit, NAR	66
CardioPAK, ECG	65
Clarius, Wireless Handheld Ultrasound	62
Field Corpsman Kits	66
MedicECG	64
MightySat Rx Fingertip Pulse Ox	63
Ophthalmoscope, Heine Mini 3000	66
Stethoscope, ADC	66
Vscan Ultrasound Tool from GE	63

DOD BLEEDING CONTROL

RESPONSE KITS - D-BCRK

D-BCRK 8-Pack	49
D-BCRK Individual Kits	48
D-BCRK Wall Stations	49
Dragon® Handle System	80
Dragon® Harness	80

E

Eagle IFAK	26
------------------	----

EDUCATION & TRAINING

SERVICES, NAR

Emergency Evacuation Litter (EEL)	58
Emergency Survival Wrap	21
Emergency Trauma Dressing (ETD™)	9
Expeditionary Casualty Response Kit, ECRK	41
Ex-tract Strap	80
Eye Shield, Fox	92
Eye Shield, Polycarbonate, PES	92
EZ-IO G3® Educator Kit	108
EZ-IO G3® Power Driver Trainer	108
EZ-IO® Training Needles	108
EZ-IO G3 Power Driver	17
EZ IO Manual Extremity Needle Set	17

F

FBTK Field Blood Transfusion Kits

Blood Surveillance	20
FBTK Basic	20
FBTK Complete	20
Field Corpsman Kits	66

FIRST LIGHT® TACTICAL LIGHTS

Tomahawk Multi Color	85
Tomahawk Night Vision	85
TORQ	85
Footape, Petzl	81
Fox Eye Shield	92

G

G3 Educator Kit, EZ-IO	108
G3 Power Driver	17
G3 Power Driver Trainer, EZ IO	108

GAUZE

Compressed Gauze, NAR	9
S-Rolled Gauze, NAR	9

Gear Leash, 6 oz MOLLE Retractable	90
Gecko Grip® Tape, Pack of 6	84
Gloves, Bear Claw®	91
Gloves, BlackTalon®	91
Goggles, Outside the Wire	92
Gray Man EyePro	92
GriGri 2 Belay Device, Petzl	81

H

H ₂ O 112 SBK	83
Hasty Harness	80
Headlamp, NAR	84
Headlamp, Streamlight	84
Headlamp, Tactikka + RGB	84
Heine Mini 3000 Ophthalmoscope	66
HELO Tagline Kit	79
Hemorrhage Control Training, Biofeedback	108
Hemorrhage Control Training, Combat Gauze®	108

HEMOSTATICS

Celox Gauze	8
ChitoGauze Pro	8
Combat Gauze	8
NuStat	8
High Threat Extraction Kit	79
Hook Safety Cutters, Benchmark	83
HPMK®	21
HRS® Heat Reflective Shell	21
HyFin® Chest Seals	13
HyFin® Trainer Twin Pack	107

HYPOTHERMIA PREVENTION & MANAGEMENT

Emergency Survival Wrap	21
HPMK	21

HRS® Heat Reflective Shell	21
Ready Heat™ Blanket	21

I

Individual Kit	31
Injury & Ailment Case, I&A	67
IPOK	24
IPRO® by NAR	92
IV Pole, Adjustable	55

J

JETT®	6
-------------	---

K

K-9 Handler IFAK	43
K-9 Litter, NAR	58
K-9 Tactical Field Kit	43
King LTS-D™	11

L

LACERATION TRAYS

Facial	73
Minor	73
Toenail Removal	73
Laryngoscope Blades, Fiber Optic	12
Laryngoscope Set, Special Ops	12

LEATHERMAN MULTI-TOOLS

MUT®	82
Raptor Multi-Tool	82
SUPERTOOL 300,	82
LED Field Medical Light, PAXLight	68

LITTER ACCESSORIES

Backrest, Talon® or Raven®	55
Casualty Equipment Bag	56
IV Pole, Adjustable	55
Litter Buckle, Replacement	56
Litter Stands, Standard	55
OSL™ Backrest	57
OSL™ Litter Stands	57
Retention Strap, Litter Support	56
RIG Series Assault Litter Carrier	56
Talon® Litter Carrier	56
Wheeled Litter Carrier	55

LITTERS

Emergency Evacuation Litter (EEL)	58
K-9 Litter, NAR	58
MedEVAC Litter	57
OverSized Litter (OSL™)	57
RAVEN® 90C	54
R.E.E.L. Rapid Expeditionary Evacuation Litter	58
Stingray®	58
T3™ Litter	54
TALON® II 90C	54
TED Tactical Extrication Device	58

M

Manual Extremity Needle Set, EZ-IO®	17
Manual Training Needle, EZ-IO®	108
M.A.R.C.H.	4-21
Mariner Kit	33
Maritime Assault Kit	30
Mark IV Field Hospital Bed	68

MASSIVE HEMORRHAGE

C-A-T® Tourniquet	6
C-A-T® Holders	7
Celox Gauze	8
ChitoGauze Pro	8
Combat Gauze	8
Compressed Gauze, NAR	9
Cravat (Triangular Bandage)	9
Emergency Trauma Dressing (ETD™)	9
JETT®	6
NuStat	8
S-Rolled Gauze, NAR	9
Tactical Compression Wrap	9
MCIM Kit	95
MedEVAC Litter	57

MEDICAL/SURGICAL

Armadillo Medication Case	67
BP/Stethoscope Combo Kit, NAR	66
CardioPAK, ECG	65
Clarius, Wireless Handheld	
Ultrasound	62
Field Corpsman Kits	66
Injury & Ailment Case, I&A	67
Laceration Trays	73
MARK IV Field Hospital Bed	68
MedicECG	64
microMend Skin Closure Devices	69
MightySat Rx Fingertip	
Pulse Oximeter	63
Minor Surgery Sets	71
Multi-Use Luggage & Equipment Carrier, M.U.L.E.	67
Ophthalmoscope, Heine Mini 3000	66
PAXLight LED Field Medical Light	68
Penblade Safety Scalpels	71
R-SICK Rapid Surgical Interventio	
Convenience Kits	70
Stethoscope, ADC	66
Surgical Instruments	72
Surgical Suture Kits	71
Suture Removal Kits	72
Toenail Removal Kit	73
Toenail Removal Trays	71
Vscan with Dual Probe Pocket-sized	
Ultrasound Tool from GE	63
Medic ECG, Resolute Medical	64
Medic Leg Rig	31
Medic Trauma Pack	36
Medic Trauma Sheet	46
M-FAK Mini First Aid Kit	25
microMend Skin Closure Devices	69
MightySat Rx Fingertip	
Pulse Oximeter.....	63
Mini Machined Pulley,	
Rock Exotica.....	81
Minor Surgery Sets	71
Multi-Mission Trauma Pack	44
Multi-Use Luggage & Equipment Carrier, M.U.L.E.	67
MUT® by Leatherman®	82

N

NAR-4 Chest Pouch.....	34
NAR-4™ Aid Kit	35
NAR-5 SAR Aid Kit	39
Nasopharyngeal Airway with lubricant	11
Naval Boat Response Aid Box	39
Naval Boat Response Aid Kit	39
Naval First Aid Box Response Kit	38
Naval Gun Response Aid Kit	38
Navy SAR SOML 112 SBK	83
Needleless Saline Lock Kit, NAR	16

NITRILE GLOVES

Bear Claw® Gloves	91
Black Talon® Gloves	91
NuStat , Hemostatic.....	8

O

OCHO BAGS & ACCESSORIES

OCHO Backpack	46
OCHO Briefcase	46
OCHO IFAK Accessory Bag	46
OCHO Quick Detach	
Gear Organizers	46
OCHO Shoulder Bag	46
Ocho IFAK	28
Ophthalmoscope, Heine Mini	66
OSL™ Backrest	57
OSL™ Litter Stands	57
OTW Goggles, Outside the Wire	92
OverSized Litter (OSL™).....	57

P

PART TASK TRAINERS

G3 Educator Kit, EZ-IO	108
Hemorrhage Control Training,	
Biofeedback	108

Hemorrhage Control Training,	
Combat Gauze®	108
TrueClot® Hemorrhage Training Kit	109
PAXLight LED Field Medical Light	68
PAXLight EZ Light Handlest	68
PAXLight Wheeled Base	68

PEGASUS ROPE®

50 Foot Throw Bag	78
Helo Tagline Kit	79
Olympus™ Series	77
Rapid Access Casualty Extraction	
R.A.C.E. Kit	78
Stallion™ Series	77
Triton™ Series	77
Warhorse™ Series	6
Water Rescue Belt Throw Bag	78
Penblade Safety Scalpels	71

POLAR SKIN™

C2E Ice Pack Bundle	89
Complete Care Backpack	88
Complete Care Pack	88
Comprehensive Cooling System	88
Cooler, Standard or Large	88
Cooling Collar	89
Ice Sheets, Standard or Large	89
Polycarbonate Eye Shield, PES	92
Power Driver, G3	17
Power Driver Needle Set	17
PPE Kit, NAR	91

Q

QUANTUM® Blood & Fluid Warming System, Lifewarmer.....	18-19
Quick Detach Panels, NAR	47

R

R.A.C.E. System, Rapid Access Casualty Extraction	78
Range Trauma Aid Kit Hard Case	42
Range Trauma Kit	42
Raptor® IV Securing Device.....	17
Raptor Multi-Tool, Leatherman	82
Ready Heat™ Blanket.....	21
Rescue Assault Tether (RAT) Strap	80

RESCUE EQUIPMENT

50 Foot Throw Bag	78
Aero Screwgate Carabiner	81
Anneau Slings, Petzl	81
Basic Ascender, Petzl	81
Belt Throw Bag, Water Rescue	78
Dragon Handle System	80
Dragon Harness	80
EX-TRACT Strap	80
Footape, Petzl	81
GriGri 2 Belay Device, Petzl	81
Hasty Harness	80
HELO Tagline Kit	79
High Threat Extraction Kit	79
Mini Machined Pulley, Rock Exotica	81
PEGASUS Ropes	76
R.A.C.E. System, Rapid Access Casualty Extraction	78
Rescue Assault Tether, RAT Strap	80
Reverso 4 Belay/Rappel, Petzl	81
Revolver Carabiner	81
ROCKD Carabiner	81
Shield Carabiner	81
Tactical Bail-Out Kit	80
TiBloc Ascender/Rope Grab, Petzl	81
TriRig Rigging Plate, Rock Exotica	81
Water Rescue Belt Throw Bag	78

Rescue Sleeve II, NAR	9
Respirator, NAR	91

RESPIRATORY SUPPORT

ARS® Needle Decompression Kit	15
HyFin Occlusive Individual	13
HyFin Vent Individual	13
SPEAR™ Simplified Pneumothorax Emergency Air Release	14
USMC Combo Pack Chest Seals	13
Retention Strap, Litter Support Bar	56
Reverso 4 Belay/Rappel, Petzl	81

Revolver Carabiner	81
Rigid C-A-T® TQ Case with Shirt Shield	7
RIG Series Assault Litter Carrier	56
RIG Series Eagle IFAK	27
RockD Carabiner, Rock Exotica	81
R-SICK Rapid Surgical Intervention Convenience Kits	70

S

Saline Lock Kit, NAR	16
SAM® Pelvic Sling	90
SAM® Splints	90
Scissor Leash™	90
Sharps Shuttle	17
Shield Carabiner	81
SIRK Gen 1	24
SIRK Gen 2	24
SPEAR™ Simplified Pneumothorax Emergency Air Release	14
Specialty Kitting Capabilities	51
Spine Boards, NAR	59

SPLINTING

SAM® Pelvic Sling	90
SAM® Soft Shell Splints	90
SAM® Splint II	90
Tactical Traction Splint	90
TC SPLINT, Telescoping Collapsible	89
Squad Kit	32
Squad Responder Kit	40
S-Rolled Gauze, NAR	9
Stethoscope, ADC	66
Stethoscope/BP Combo Kit	66
Stingray® Poleless Litter	58
Streamlight Headlamp	84
Supertool 300 by Leatherman®	82
Surgical Instruments	72
Surgical Suture Kits	71
Survival Wrap, Emergency	21
Suture Removal Kits	72

T

T2® Tag	94
---------------	----

TACEVAC

Backrest, Litter	55
Casualty Equipment Bag	56
Casualty Immobilization System, NAR	59
Casualty Restraint Straps, NAR	59
Cervical Collar, ACE	59
Head Immobilizer, NAR	59
IV Pole, Adjustable	55
K-9 Litter	58
Litter Buckle Replacement Set	56
Litter Stands	55
Litter Support Bar Retention Strap	56
MEDEVAC Litter	57
OSL Backrest	57
OSL Litter, Oversize	57
OSL Litter Stands	57
Raven 90C	54
R.E.E.L. Rapid Expeditionary Evacuation Litter	58
Rescue Sleeve II	59
RIG Series Litter Carrier	56
Spine Boards	59
STINGRAY Poleless Litter	58
T3 Litter	54
TALON® II 90C	54
TALON Litter Carrier	56
TED, Tactical Extrication Device	58
Wheeled Litter Carrier	55
Tactical Bail Out Kit	80
Tactical Compression Wrap	9
Tactical Cricket®	11
Tactical Suction Device, NAR	11
Tactical Traction Splint	90
Tactical Triage Kit, T2®	94
Talon® Litter Carrier	56
T.A.L.O.N. Manual Training Needle	108
T.A.L.O.N.™ Manual Needle Set	17
Tape, Gecko Grip®	84

TCCC ALL COMBATANT TRAINING SYSTEM

TCCC Evacuation Skills Set	105
----------------------------------	-----

TCCC IFAK Skills Kit -	
All Combatant	106
TCCC IFAK Skills Set -	
All Combatant	105
TCCC Skills Set - All Combatant	104
TCCC Skills System -	
All Combatant	104
TC Splint, Telescoping Collapsible	89
Throw Bag, 50 Foot	78
TiBloc Ascender/Rope Grab, Petzl	81
Tiered Trauma Medical Bag (TTMB)	45
Toenail Removal Kit	73
Toenail Removal Trays	71

TOMMANIKIN™

TRAINING SIMULATORS

BASIC TOMManikin™	100
BLAST TOMManikin™	101
GSW TOMManikin™	102
K9 Canine Manikin Training Simulator	103
T.O.R.K.	29
Tracheal Hook, NAR	12
Tracheostomy Kit with Bougie Introducer	10
Training ETD, Responder	107
Training Nneedle Sets,	
EZ-IO G3 Power Driver	108

TRANSPORT & STORAGE

ARMADILLO Medication Case	67
I&A CASE, Injury & Ailment	67
M.U.L.E.	67
Trauma Shears, NAR	90
Triage 916 SBK	83
Triage Officer Kit, T2®	95
TriRig Rigging Plate, Rock Exotica	81
TrueClot® Hemorrhage Training Kit	109

U

ULTRASOUND HANDHELD SCANNERS

Clarius, Wireless Ultrasound	62
Vscan with Dual Probe Pocket-sized Ultrasound Tool from GE	63
USCG Boat Response Aid Kit	39
USCG IFAK	30
USMC CLS KIT	40
USMC Combo Pack -	
HyFin Chest Seals	13

V

Vehicle Trauma Kit, Patrol	48
----------------------------------	----

VIDEOS, ONLINE

ARS Needle Decompression Kit	15
BOA® IV Constricting Band	16
Combat Application Tourniquet	6
ETD™ Emergency Trauma Dressing	9
HyFin® Chest Seals	13
JETT	6
SPEAR™ Simplified Pneumothorax Emergency Air Release	14
Vscan with Dual Probe Pocket-sized Ultrasound Tool from GE	63

W

WALK® Warrior Aid & Litter Kit	37
WALK ReSupply Kit	37
Water Rescue Belt Throw Bag	78
Wheeled Litter Carrier	55

PRODUCT index

by NSN Number

1095-01-677-1013 Navy SAR SOML 112 SBK.....83	6510-01-659-3942 Combat Ready™ Tactical Body Tape, Multi-Pack.....89	6515-01-538-9276 Trauma Shears, Large 7.25".....90
6510-00-201-1755 Cravat, Triangular Bandage.....9	6510-01-659-7998 Combat Ready™ Tactical Body Tape, Bravo.....89	6515-01-540-7206 Tactical Suction Device, NAR.....11
6510-01-503-2117 Compressed Gauze, NAR.....9	6510-01-659-7999 Combat Ready™ Tactical Body Tape, Alpha.....89	6515-01-540-7226 Scissor Leash, NAR.....90
6510-01-532-8019 HyFin® Chest Seal, Occlusive Ind.13	6510-01-659-8001 Combat Ready™ Tactical Body Tape, Charlie.....89	6515-01-540-7568 Tactical Crickit.....11
6510-01-541-8121 ETD™ Abdominal/Stump.....9	6510-01-659-8349 Combat Ready™ Tactical Body Tape, Echo.....89	6515-01-540-7579 Medic Trauma Sheet, Digital Camo.....46
6510-01-558-4108 ETD™ 4 inch, Standard.....9	6510-01-659-8353 Combat Ready™ Tactical Body Tape, Delta.....89	6515-01-540-7585 Medic Trauma Sheet, OD Green.....46
6510-01-558-4114 ETD™ 6 inch, Standard.....9	6515-01-311-0360 Surgical Suture Kits, Pack of 20.....71	6515-01-540-7603 Medic Trauma Sheet, Coyote.....46
6510-01-562-3325 Combat Gauze, Standard.....8	6515-01-338-6602 Ambu RDIC, Resuscitator.....12	6515-01-540-7612 Medic Trauma Sheet, Black.....46
6510-01-587-6579 Dry Sterile Burn Dressing.....93	6515-01-346-9186 Tactical Traction Splint, Black.....90	6515-01-541-0635 ARS Needle Decompression Kit 14 ga.....15
6510-01-591-7740 ChitoGauze PRO.....8	6515-01-388-1351 Biological Personal Protection Kit.....91	6515-01-541-8147 Cervical Collar, OD Green.....59
6510-01-600-4705 Combat Gauze, XL.....8	6515-01-449-1016 Fox Eye Shield, Box of 12.....92	6515-01-542-7723 BP/Stethoscope Kit.....66
6510-01-605-2234 S-Rolled Gauze, NAR.....9	6515-01-450-9790 Laryngoscope Set, Special Ops.....12	6515-01-542-7937 Headlamp, NAR.....84
6510-01-614-0938 Tactical Compression Wrap.....9	6515-01-458-6057 Basic Field Corpsman Kit.....66	6515-01-562-9734 Tracheal Hook, NAR.....11
6510-01-619-8636 Celox Gauze, 10 foot.....8	6515-01-458-6178 Deluxe Field Corpsman Kit.....66	6515-01-568-0193 Cyclone Pocket BVM.....12
6510-01-622-1624* Trauma Shears, Small 6.25".....90	6515-01-494-1951 Sam Splint II.....90	6515-01-587-1700 Combat Casualty Card, Pack of 25.....94
6510-01-624-0840 HyFin® Chest Seal, Vent Individual.....13	6515-01-507-8074 Soft Shell SAM Splint 15 in., Box of 20.....90	6515-01-591-0397 King LTS-D, Purple, Size 5.....11
6510-01-626-2990 BurnTec 8 in. x 16 in.....93	6515-01-509-8866 SAM Pelvic Sling.....90	6515-01-591-1578 King LTS-D, Red, Size 4.....11
6510-01-627-1834 BurnTec, Face 10 in. x 10 in.....93	6515-01-515-0190 Black Talon Gloves, Box, 45 pr, XL.....91	6515-01-598-1877 Fox Eye Shield, Box of 50.....92
6510-01-627-1855 BurnTec Face 12 in. x 16 in.....93	6515-01-515-0195 Black Talon Gloves, Box, 50 pr, Small.....91	6515-01-600-0262 Bear Claw Gloves, Kit, 25 pr, Large.....91
6510-01-630-9177 BurnTec Burn Dressing Kit.....93	6515-01-515-0197 Black Talon Gloves, Box, 50 pr, Large.....91	6515-01-600-6039* Bear Claw Gloves, Box, 50 pr, Large.....91
6510-01-631-1883 BurnTec 4 in. x 4 in.....93	6515-01-515-0200 Black Talon Gloves, Box, 50 pr, Medium.....91	6515-01-600-6052* Bear Claw Gloves, Box, 45 pr, XL.....91
6510-01-631-1886 BurnTec 16 in. x 24 in.....93	6515-01-518-8551 Sharps Shuttle, Single Pack.....17	6515-01-600-6057* Bear Claw Gloves, Box, 50 pr, Medium.....91
6510-01-631-1890 BurnTec 8.5 in. x 11 in.....93	6515-01-519-3612 Ambu Spur II BVM.....12	6515-01-600-6059* Bear Claw Gloves, Box, 50 pr, Small.....91
6510-01-631-1897 BurnTec 8 in. x 8 in.....93	6515-01-521-7976 CAT Tourniquet, Black.....6	6515-01-607-1943 Needleless Saline Lock Kit, NAR.....16
6510-01-631-1907 BurnTec 5 in. x 10 in.....93	6515-01-523-3317 Raptor® IV Securing Device.....17	6515-01-612-7129 CAT Tourniquet, Orange.....6
6510-01-631-1913 BurnTec 2 in. x 2 in.....93	6515-01-523-5440 Soft Shell SAM Splint, 15 inch.....90	6515-01-616-5841 JETT.....6
6510-01-631-2080 BurnTec 2.5 in x 5 in.....93	6515-01-529-1187 Nasopharyngeal Airway & Lube, 28F.....11	6515-01-619-5350 Bear Claw Gloves, Kit, 25 pr, Small.....91
6510-01-631-2124 BurnTec 4 in. x 16 in.....93	6515-01-532-8056 HPMK, Hypothermia Prev & Mgt Kit.....21	6515-01-619-5367 Bear Claw Gloves, Kit, 25 pr, XL.....91
6510-01-631-2138 BurnTec 5 in. x 5 in.....93	6515-01-537-1531 NAR Respirator.....91	6515-01-619-5671 Sharps Shuttle, 6 pack.....17
6510-01-642-6212 HyFin® Vent Twin Pack.....13	6515-01-537-1534 PPE Kit, NAR.....91	6515-01-619-5689 Bear Claw Gloves, Kit, 25 pr, Medium.....91
6510-01-644-5427 BurnTec Minor Burn Dressing Kit.....93	6515-01-537-2611 BOA® IV Constricting Band, Standard.....16	6515-01-624-7148 T2® Tag, NAR, Pack of 20.....94
6510-01-654-6787 USMC Combo Pack - Hyfin Chest Seals..13	6515-01-537-4094 Saline Lock Kit, NAR.....16	6515-01-630-5551 Polycarbonate Eye Shield, Box of 100.....92
6510-01-654-6801 ETD™ 8x10 inch.....9	6515-01-537-4161 Combat Casualty Reference Card.....94	6515-01-630-8050 Toenail Removal Trays, Pack of 10.....71
		6515-01-632-0059 Black Talon Gloves, Kit, 25 pr, Large.....91
		6515-01-646-2636 First Light® TORQ.....85
		6515-01-646-2636 TORQ, First Light.....85
		6515-01-647-5456 Polycarbonate Eye Shield, Single.....92

6515-01-650-0036 Cervical Collar, White.....	59	6532-01-653-6054 Polar Skin Ice Sheets, Single Pack Std.....	89	8465-01-666-7169 Naval Gun Response Aid Bag.....	38
6515-01-661-1786 Minor Laceration Tray.....	73	6540-01-491-6651 Heine Mini 3000 Ophthalmoscope.....	66	8465-01-666-7328 Naval First Aid Box Pouch.....	38
6515-01-661-1802 Facial Laceration Tray.....	73	6545-01-529-4185 Maritime Assault Kit, Black.....	30	8465-01-669-5796* NAR-5 SAR Aid Bag, COY.....	39
6515-01-672-7977* Black Talon Gloves, Kit, 25 pr, Medium.....	91	6545-01-529-4187 Maritime Assault Kit, Coyote.....	30		
6515-01-672-7985* Black Talon Gloves, Kit, 25 pr, XLarge.....	91	6545-01-531-6626 Maritime Assault Kit, OD Green.....	30		
6515-01-673-1701 ARS Needle Decompression Kit - 10 ga.....	15	6545-01-532-4962 Warrior Aid & Litter Kit - DUC.....	37		
6530-01-432-5114 Raven® 90C Bifold Litter.....	54	6545-01-537-3138 MCIM Kit, OD Green.....	95		
6530-01-435-7337 Litter Buckle Replacement.....	56	6545-01-537-4189 Triage Officer Kit, T2® OD Green.....	95		
6530-01-477-8525 Rescue Sleeve II, Orange.....	59	6545-01-537-4197 Command Kit, OD Green.....	95		
6530-01-490-2487 Spine Board, NAR, OD Green.....	59	6545-01-543-2956 Squad Kit, Coyote.....	32		
6530-01-492-4079 Stringray Poleless Litter.....	58	6545-01-543-3159 Squad Kit, OD Green.....	32		
6530-01-492-4114 Tactical Extrication Device (TED), ODG.....	58	6545-01-543-3170 Medic Leg Rig, Coyote.....	31		
6530-01-504-9051 Talon® II 90 C Litter.....	54	6545-01-543-3174 Medic Leg Rig, OD Green.....	31		
6530-01-504-9054 Litter IV Pole, Adjustable.....	55	6545-01-543-3185 Medic Leg Rig, Black.....	31		
6530-01-504-9055 Talon® Assault Litter Carrier, OD Green.....	56	6545-01-543-5721 Squad Kit, Black.....	32		
6530-01-523-3314 Litter Stands, Standard.....	55	6545-01-543-5733 Squad Kit, Digital Camo.....	32		
6530-01-529-1134 Individual Kit, Black.....	31	6545-01-543-5770 Mariner Kit.....	33		
6530-01-531-6649 Rescue Sleeve II, Black.....	59	6545-01-571-4470 USMC CLS Kit, COY.....	40		
6530-01-531-6771 Spine Board, NAR, Orange.....	59	6545-01-572-7248 USCG IFAK, Black.....	30		
6530-01-531-6808 Tactical Extrication Device (TED), Black.....	58	6545-01-572-7251 Medic Trauma Pack, Black.....	36		
6530-01-531-6854 Talon® Assault Litter Carrier, Coyote.....	56	6545-01-587-1199 Warrior Aid & Litter Kit - COY.....	37		
6530-01-531-6907 Talon® Assault Litter Carrier, Black.....	56	6545-01-587-8513 USCG IFAK, Coyote.....	30		
6530-01-549-1089 PAXLight LED Field Medical Light.....	68	6545-01-588-1862 Medic Trauma Pack, Coyote.....	36		
6530-01-552-7708 Oversized (OSL™) Litter.....	57	6545-01-591-2134 Tactical Triage Kit.....	94		
6530-01-586-2411 Casualty Immobilization System, NAR.....	59	6545-01-592-3261 T2® Tag, NAR, 1000 Pack.....	94		
6530-01-594-7444 Backrest, Talon® or Raven®.....	55	6545-01-594-5617* Medic Trauma Pack, OD Green.....	36		
6530-01-598-9595 Dragon® Handle System.....	80	6545-01-594-9399* Armadillo Medication Case, Clear Lid.....	67		
6530-01-605-9432 Nylon C-A-T Holder ODGreen.....	7	6545-01-610-5411* Range Trauma Kit.....	42		
6530-01-619-1777 Head Immobilizer, NAR.....	59	6545-01-646-2623 USCG Boat Response Aid Kit.....	39		
6530-01-624-7626 Mark IV Field Hospital Bed.....	68	6545-01-648-2690 Naval Gun Response Aid Kit, OD Green.....	38		
6530-01-632-7872 Rescue Assault Tether (RAT) Strap.....	80	6545-01-648-2705 Naval First Aid Box Response Kit.....	38		
6530-01-643-9134* Talon® Assault Litter Carrier, DUC.....	56	6545-01-648-3228 Naval Boat Response Aid Kit.....	39		
6530-01-643-9138* RIG Series Assault Litter Carrier, COY.....	56	6910-01-560-2972 CAT Tourniquet, Trainer Blue.....	6, 107		
6530-01-647-8571* RIG Series Litter Carrier, MTC.....	56	8115-01-656-8316 Naval Boat Response Aid Box.....	39		
6530-01-663-2308* Wheeled Litter Carrier.....	55	8145-01-573-2533 Armadillo Case, Black.....	67		
6530-01-674-9774* Hasty Harness, Black.....	80	8145-01-617-9075* Half in. Hook Discs for Armadillo Case.....	67		
6532-01-627-7888 HRS Heat Reflective Shield.....	21	8465-01-575-6644 Casualty Equipment Bag.....	56		

* denotes a NATO NSN#

PRODUCT by Item Number index

10-0001.....12	30-0079108	60-000658	80-009332	80-049725	85-180629	ZZ-022267
10-0015.....13	30-008117	60-000855	80-009446	80-049825	85-189924	ZZ-023596
10-0017.....11	30-0082108	60-000956	80-009636	80-049925	85-1958106	ZZ-024691
10-0018.....11	30-008317	60-001056	80-009737	80-050025	85-205748	ZZ-024791
10-0019.....12	30-0084108	60-001156	80-010733	80-050125	85-205848	ZZ-024891
10-0020.....11	30-008517	60-001455	80-013184	80-050225	85-205948	ZZ-024991
10-0025.....11	30-0086108	60-002254	80-013256	80-050325	85-206648	ZZ-025083
10-0026.....11	30-00886	60-002357	80-013979	80-050425	85-206749	ZZ-025867
10-0027.....11	30-00899	60-002457	80-014079	80-050525	85-206849	ZZ-025967
10-0029.....13	30-00939	60-002556	80-014179	80-050927	85-206948	ZZ-026294
10-0032.....107	30-00949	60-002656	80-014279	80-051027	85-207049	ZZ-029815
10-0037.....13	30-011093	60-002755	80-014930	80-051127	85-207149	ZZ-0467106
10-0043.....12	30-011193	60-002957	80-015830	80-051227	85-207249	ZZ-053982
10-0047.....12	30-011293	60-004155	80-015940	80-051327	85-207349	ZZ-054282
10-0048.....10	30-011393	60-004268	80-016724	80-051427	85-207449	ZZ-054381
10-0051.....14	30-011493	60-004568	80-016824	80-051527	85-207549	ZZ-054581
10-0052.....10	30-011593	60-004668	80-017334	80-051627	85-207649	ZZ-054681
20-0001.....66	30-011693	60-005157	80-017534	80-051727	85-207749	ZZ-054781
20-0002.....66	30-011793	60-005257	80-017634	80-051827	85-230029	ZZ-054881
20-0007.....66	30-011893	60-005656	80-018135	80-051927	85-230129	ZZ-054981
20-0023.....95	30-011993	60-005756	80-018335	80-052027	85-230229	ZZ-055081
20-0025.....95	30-012093	60-006355	80-018435	80-054388	90-000180	ZZ-055181
20-0026.....95	30-012193	60-006455	80-021342	80-054488	90-000280	ZZ-055284
20-0029.....95	30-012293	60-0074W58	80-022324	80-087170	90-000380	ZZ-055381
20-0031.....94	30-013217	60-007558	80-022524	80-088841	90-000480	ZZ-055481
20-0032.....94	30-0133108	60-007658	80-024145	80-088941	90-000581	ZZ-0617106
20-0040.....94	30-014292	66-001968	80-024435	80-089347	90-000681	ZZ-063089
20-0057.....66	30-014392	70-000191	80-024530	80-089447	90-000781	ZZ-063189
20-0060.....66	30-014492	70-000291	80-024734	80-089547	90-002080	ZZ-063289
20-0073.....95	30-0145107	70-000391	80-024845	80-089647	90-006083	ZZ-063389
20-0074.....95	30-0146107	70-000491	80-024931	80-089747	90-006183	ZZ-063489
20-0075.....95	30-016293	70-000590	80-025132	80-089847	90-006378	ZZ-063589
20-0076.....95	30-01659	70-000691	80-025646	80-091528	90-006478	ZZ-065283
20-0077.....95	30-01708	70-000791	80-025746	80-091628	90-006578	ZZ-066472
20-0083.....64	30-01847	70-000991	80-025846	80-091728	90-006692	ZZ-066572
20-0084.....64	30-01857	70-001091	80-025946	80-092128	90-006792	ZZ-066672
20-0085.....64	30-01867	70-001191	80-026046	80-092228	90-006892	ZZ-066772
20-0086.....65	30-020369	70-001291	80-026146	80-093240	90-006979	ZZ-066889
20-0087.....65	30-020469	70-001482	80-027145	80-09447	90-007778	ZZ-066989
20-0089.....62	30-02087	70-001585	80-027245	80-09457	90-007883	ZZ-067089
20-0089.....62	30-02127	70-001685	80-027345	80-09467	90-008280	ZZ-067488
20-0091.....62	30-02137	70-011085	80-027445	80-094944	90-009283	ZZ-067589
20-0092.....63	30-023069	70-011185	80-028934	80-095044	90-009383	ZZ-068288
20-0093.....63	30-023169	70-011285	80-029124	80-0968104	90-009683	ZZ-068373
20-0094.....63	30-023269	70-011385	80-029942	80-0969105	90-009783	ZZ-068473
20-0095.....63	30-02368	70-011485	80-030043	80-0970104	90-010377	ZZ-068573
30-0001.....6	30-555092	70-011585	80-030143	80-0971105	90-010477	ZZ-069472
30-0003.....9	35-000118	70-011685	80-030343	80-099127	90-020576	ZZ-069572
30-0006.....17	35-000218	70-011785	80-030835	80-099227	90-020676	ZZ-069672
30-0009.....16	35-000318	70-011885	80-033846	80-099327	90-020776	ZZ-069772
30-0012.....9	35-000418	70-015590	80-033946	80-099725	90-020876	ZZ-069872
30-0013.....16	35-000518	80-000331	80-034046	83-000020	90-030377	ZZ-069972
30-0014.....17	35-000618	80-000431	80-035339	83-001670	90-030577	ZZ-070072
30-0016.....9	50-000490	80-001330	80-035439	83-001770	90-040077	ZZ-070172
30-0023.....6	50-000890	80-001430	80-03556	83-002220	90-040177	ZZ-070272
30-0031.....9	50-000990	80-001530	80-03566	83-002520	93-0017108	ZZ-070372
30-0032.....9	50-001059	80-002721	80-035738	85-000829	93-0019108	ZZ-070572
30-0033.....6, 107	50-001159	80-003732	80-036193	85-000929	93-0020109	ZZ-075589
30-0039.....8	50-001259	80-003832	80-037046	85-001229	93-0027103	ZZ-075688
30-0044.....7	50-001359	80-003932	80-037146	85-002037	93-0035100	ZZ-075889
30-0046.....16	50-001459	80-004631	80-037228	85-004829	93-0036101	ZZ-075988
30-0049.....8	50-001590	80-004731	80-037328	85-004929	93-0040102	ZZ-076083
30-0052.....9	50-002059	80-004831	80-037428	85-005029	ZZ-003256	ZZ-079488
30-0053.....7	50-002159	80-005830	80-037528	85-005229	ZZ-003411	ZZ-079796
30-0055.....92	50-002759	80-005937	80-039939	85-005329	ZZ-005121	ZZ-079896
30-0057.....7	50-003159	80-006037	80-041593	85-005429	ZZ-005221	ZZ-079996
30-0059.....7	50-003290	80-006746	80-041639	85-010430	ZZ-005615	ZZ-080096
30-0063.....107	50-003390	80-006846	80-041738	85-017726	ZZ-005794	ZZ-082249
30-0065.....8	50-003490	80-006946	80-041838	85-017826	ZZ-006390	ZZ-083421
30-0070.....8	50-003989	80-007036	80-041938	85-017926	ZZ-006490	ZZ-087771
30-0071.....16	50-100590	80-007136	80-042039	85-018926	ZZ-008167	ZZ-087871
30-0074.....108	50-100759	80-007236	80-042440	85-018926	ZZ-009296	ZZ-087971
30-0075.....17	60-000154	80-008731	80-043943	85-019026	ZZ-009396	ZZ-095671
30-0076.....108	60-000254	80-008937	80-044043	85-019126	ZZ-009496	ZZ-095771
30-0077.....17	60-000458	80-009131	80-044143	85-019226	ZZ-010684	ZZ-095871
30-0078.....17	60-000558		80-044243	85-019640	ZZ-010796	ZZ-095971
			80-044343	85-019740	ZZ-010896	ZZ-096071
			80-044443	85-088942	ZZ-010992	ZZ-096171
			80-044543	85-104513	ZZ-012196	ZZ-100458
			80-044643	85-173824	ZZ-017696	ZZ-109084
			80-045188	85-174948	ZZ-021267	
			80-048891	85-175049	ZZ-021591	
			80-049425	85-175149	ZZ-021691	
			80-049525	85-176440	ZZ-021791	
			80-049625	85-180529	ZZ-021891	

STAY informed

and KEEP IN TOUCH



Keep up with the latest North American Rescue® news and happenings!

Join us on Facebook:
facebook.com/NARescueLLC



Providing medical solutions and celebrating the life-savers who use them!

Join us on Instagram:
instagram.com/northamericanrescue



Keep up with North American Rescue® with 140 characters or less!

Join us on Twitter:
twitter.com/NARescue



View North American Rescue® instructional and informational videos!

Join us on You Tube:
youtube.com/NARescue01



Visit North American Rescue® online to sign up for our newsletter!

Sign up at:
www.NARescue.com



Office Hours: M-F 8:00am - 7:00pm ET
Call toll-free 888.689.6277 • www.NARescue.com

CUSTOMER SUPPORT

- **COMMITMENT:** We pride ourselves in providing an exceptional customer experience beginning with your initial inquiry up to order placement and beyond. Our customer service staff members are trained by certified instructors on all NAR products and are standing by to help with any questions you may have
- **SUPPORT:** Our in-house technical support professionals consist of an elite group of Nationally Registered Paramedics and Special Operations Medics from military, law enforcement, tactical law enforcement, fire-rescue and EMS backgrounds

ORDERING

- **ONLINE:** www.NARescue.com
SPECIAL PRICING: automatically available ONLINE for all military customers with a .mil e-mail address! Register your account in just minutes by creating an online account using your .mil e-mail address
ESTABLISHED CUSTOMERS: Want to order online with special pricing? Contact Customer Service at 888.689.6277 to have your online account set up
NEW CUSTOMERS: Simply create an online account and you will have access to customer specific pricing – anytime, anywhere you are logged in to NARescue.com
QUOTE REQUESTS: Online accounts have no waiting to receive your custom pricing to place an order or, we will gladly accept your submittal for a pricing quote request
- **CONTRACT:**
Tax ID: 27-1024029
Cage Code: 06ST7
DUNS: 83-24-26782
SAM Registered: Yes
ECAT: dmmonline.dscp.dla.mil
Prime Vendor: dmmonline.dscp.dla.mil/primevendor/primehome
- **EMAIL:** info@narescue.com
- **PHONE:** 888.689.6277
- **ADDRESS:**
35 Tedwall Ct | Greer, SC 29650
- **FAX:** 864.675.9880

SHIPPING

- In-stock merchandise ships within two business days
- Need it ASAP? Ask Customer Service about overnight or second day delivery



PRODUCTS[®] WITH A MISSION[®]



North American Rescue[®] • 35 Tedwall Court • Greer, SC 29650 • 888.689.6277 • NARescue.com

MILITARY CATALOG | Volume 18 | Number 1