



OUTDOOR

EXTREME LIFESAVING GEAR

OUTDOOR | ADVENTURE | SPORTS



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277



COMBAT PROVEN ADVENTURE READY

Be prepared wherever you are...

Whether caused by explosives, bullets, vehicle accidents, boat propellers, surfboard fins, spear guns, shark bites, coral scrapes or bike spills, ***trauma is trauma and must be treated with speed using the proper equipment to stop the bleeding and save a life.***

Nobody understands this better than leading Military medical equipment provider, North American Rescue® (NAR). ***Since 1996, NAR has been dedicated to decreasing preventable death by providing the most effective and highest quality mission critical medical products*** to our military, federal agencies, civilian law enforcement, EMS and pre-hospital life savers. We design, develop and test these products based on the military's Tactical Combat Casualty Care (TCCC) guidelines, Rescue Human Factors® engineering, evidenced based medicine and an in-depth understanding of the requirements of our valued customers.

For the past few years we have been ***committed to the process of better understanding the potential medical risks, needs and criteria for those engaged in active Outdoor Adventure Sports.***

During this investigation, we "talked trauma" with professional and enthusiast adventurers, expeditionary guides and resort activity coordinators that participate in activities/events ranging from diving, big wave surfing, fishing, kayaking and boating to mountain biking, hunting, rock climbing, skiing and snowmobiling. Hearing about, witnessing and even treating the kinds of injuries likely to occur in these scenarios, ***we have focused our experience and knowledge on developing kits that meet the diverse first aid scenarios that are encountered in a wide variety of adventure sports activities and environments.***

At the end of the day, we just want our friends to make it home safe to their loved ones so that they can tell their story. Whether it's major trauma or a minor cut, these kits will cover your needs – ensuring peace of mind on your "adventure trip of a lifetime" or simply a trip to the store. These kits provide the necessary tools to save your own life or a friend's life. Accidents happen. When you're in the middle of the chaos it's too late to plan. Ask yourself today, "What is your life worth?" ***Remember, life is full of the unexpected...***



Understanding your GEAR...

Get to know the innovative medical gear that sets our Trauma Kits apart from other First Aid products on today's market:

C-A-T® *Combat Application Tourniquet*

Quickly controls life-threatening extremity bleeding and is the official Tourniquet of the US Army



QUIKCLOT® *Hemostatics*

The simplicity of gauze impregnated with kaolin, the same advanced hemostatic used by the US Military to control severe bleeding



HYFIN® VENT *Chest Seals*

Advanced adhesive technology designed in a convenient Twin Pack to provide a superior seal when treating entry/exit wounds or multiple penetrating injuries



Scientifically proven as the **Fastest, Safest,** Most **Effective** Prehospital Field Tourniquet

THE 7 ADVANTAGE DIFFERENCE

1

Single Routing Buckle

- Simplified training with single protocol directions for all applications
- Effective slack removal
- Fewer windlass turns
- Faster application
- Decreased blood loss

Free-moving Internal Band

7

- Patented band within a band
- True evenly distributed circumferential pressure

Red Tip Technology®

6

- Red elliptical tip
- Added visual cue for application

RED TIP
TECHNOLOGY®

Windlass Strap

2

- Sonic welded to clip for constant contact
- Color changed to Gray for tactical considerations (30-0001 only)

Windlass Clip

3

- Bilateral entry
- Rapid windlass lock
- Bilateral buttress
- Added strength

Stabilization Plate

4

- Reinforced, beveled contact bar

Windlass Rod

5

- Increased diameter
- Enhanced strength
- Aggressive ribbing
- Improved grip



online video
available

< scan to view

Gen 7 Combat Application Tourniquet **C-A-T®**

C-A-T® COMBAT APPLICATION TOURNIQUET

COLOR	STYLE	ITEM #	NSN #
Black	Tactical	30-0001	6515-01-521-7976
Orange	Rescue	30-0023	6515-01-612-7129
Blue	Trainer	30-0033	6910-01-560-2972

SPECIFICATIONS

Packaged: H 6.5 in. x W 2.4 in. x D 1.5 in.
Device: L 37.5 in. | Weight: 2.7 oz

Tourniquets & Accessories

RIGID C-A-T® TQ CASE with Shirt Shield

- Allows you to carry a Gen 7 C-A-T® Tourniquet on your duty, patrol, or pistol belt so you have life-saving equipment on your person at all times
- Includes a Shirt Shield complete with mounting hardware to protect your uniform shirt from undue wear caused by the friction from the windlass clip
- Blade-Tech TekLok Belt Attachment fits up to a 2.25 in. belt
- Can be mounted either vertically or horizontally or attached to a Safariland Drop Leg Shroud
- Injection Molded Nylon Case stands up to rugged LE/MIL usage
- Cross medical identifier patch included

RIGID TQ CASE with SHIRT SHIELD

COLOR	ITEM #
Black	80-0944
Black, Basket Weave	80-0962
Black, Polished	80-0963
Coyote	80-0945
Ranger Green	80-0946

SPECIFICATIONS

H 5.88 in. x W 2.5 in. x D 2 in.
Weight: 5.36 oz

Explore other TQ carriers also available at **NARescue.com**



Easily attaches to case with just two screws!



RIGID C-A-T® TQ CASE SHIRT SHIELD (Shield Only)

Upgrade your Rigid TQ Case to protect your shirt from undue wear!

SPECIFICATIONS

H 1.75 in. x W 2.5 in. x D 1.5 in.
Weight: 1.44 oz

RIGID SHIRT SHIELD

COLOR	ITEM #
Black	80-0941
Coyote	80-0942
Ranger Green	80-0943

Nylon C-A-T® HOLDER

- Tailored specifically for your C-A-T®, to protect it from exposure to the elements
- PALS/MOLLE-style connectors for versatile attachment to the duty vest, gear or belt (mounted horizontally or vertically) for rapid & easy access
- 1000D IR Signature Reduced Nylon
- Easy-open pull elastic tab with tourniquet ID patch – "TQ"

NOTE: Accommodates both Gen 6 and Gen 7 C-A-T® Tourniquets (sold separately)

NYLON CAT HOLDER

COLOR	ITEM #	NSN #
Black	30-0057	----
Coyote	30-0044	----
Multi-Cam	30-0059	----
Red	30-0128	----
Ranger Green	30-0208	----
RealTree®	30-0207	----
OD Green	30-0053	6530-01-605-9432*

NATO NSN*



SPECIFICATIONS

L 6.25 in. x W 1.75 in. x D 1.5 in.
Weight: 1.26 oz

ANKLE TRAUMA KIT



SPECIFICATIONS

L 17.65 in. x H 6.5 in. x D 1.5 in.
Weight: 8.25 oz

ANKLE TRAUMA KIT

KIT/HOLSTER OPTION	ITEM #
Kit with Combat Gauze LE	80-0870
Kit with S-Rolled Gauze	80-0869
Ankle Holster Only	30-0177T

ANKLE TRAUMA KIT contains:



M-FAK® Mini First Aid Kit

- Contains the essential point-of-wounding equipment needed for treating penetrating, blast and other traumatic injuries
- Super compact, rugged nylon platform that allows attachment both vertically (MOLLE backing) and horizontally (3 in. belt loop)
- Clam shell configuration utilizes (2) two main sleeves that open on both ends for easy access



NEW! ▲
Ranger Green

M-FAK MINI FIRST AID KIT

COLOR	LEVEL	ITEM #
Black	Basic	80-0494
Coyote	Basic	80-0495
Gray	Basic	80-0872
Multi-Cam	Basic	80-0496
OD Green	Basic	80-0497
Ranger Green	Basic	80-0996
Black	Combat Gauze LE	80-0498
Coyote	Combat Gauze LE	80-0499
Gray	Combat Gauze LE	80-0873
Multi-Cam	Combat Gauze LE	80-0500
OD Green	Combat Gauze LE	80-0501
Ranger Green	Combat Gauze LE	80-0998

Every M-FAK contains:



PLUS these Individual Components based upon configuration of choice:



SPECIFICATIONS

H 6 in. x W 3 in. x D 3 in.
Weight: 13 oz



SPECIFICATIONS

H 6 in. x W 3 in. x D 3 in.
Weight: 13 oz

RealTree®

Sportsman also
includes a PVC
MED Cross Patch



Every Sportsman **M-FAK** contains:



SPORTSMAN M-FAK

COLOR	LEVEL	ITEM #
RealTree	Basic	80-0913
RealTree	with Bleeding Control Dressing	80-0914
Red	Basic	80-0890
Red	with Bleeding Control Dressing	80-0891

PLUS these Individual Components
based upon configuration of choice:



Individual First Aid Kits IFAK

T.O.R.K. Tactical Operator Response Kits

- Streamlined, compact IFAK designed to MOLLE to the assault vest or body armor
- Includes essential life-saving equipment to be used at the point of wounding
- For self-aid/buddy-aid injuries caused by penetrating or blast trauma



RealTree® Sportsman ▲
Includes a PVC MED Cross Patch

Every Sportsman T.O.R.K. contains:



PLUS these Individual Components based upon configuration of choice:



Sportsman T.O.R.K.

COLOR	LEVEL	ITEM #
RealTree	Basic	80-1006
RealTree	with Bleeding Control Dressing	80-1007

SPECIFICATIONS

H 7 in. x W 4.5 in. x D 2.75 in.
Weight: 1 lb 5 oz


T.O.R.K. BASIC | BASIC with Combat Gauze

COLOR	LEVEL	ITEM #
Black	Basic	85-0008
Coyote	Basic	85-0009
Digital Camo	Basic	85-2300
Multi-Cam	Basic	85-0012
OD Green	Basic	85-0054
Red	Basic	85-1805
Ranger Green	Basic	85-2301
Black	with Combat Gauze	80-1097
Coyote	with Combat Gauze	80-1098
Digital Camo	with Combat Gauze	80-1099
Multi-Cam	with Combat Gauze	80-1100
OD Green	with Combat Gauze	80-1101
Ranger Green	with Combat Gauze	80-1102
Red	with Combat Gauze	80-1108



T.O.R.K. BASIC / BASIC with Combat Gauze contains:





OCHO IFAK

- Compact, low visibility, point-of-wounding IFAK for your designed duty/pistol belt
- Purposely designed to look like a non-tactical, accessory bag
- Equipped with a 3 in. belt loop attachment to accommodate a duty/pistol belt
- Constructed with rip stop nylon, elastic loop organization al storage and easy open zipper

SPECIFICATIONS

L 7 in. x H 2.25 in. x D 3.25 in. | Weight: 13 oz

OCHO IFAK

COLOR	LEVEL	ITEM #
Black/Gray	Basic	80-0374
Blue/Gray	Basic	80-0375
Red/Gray	Basic	80-0915
Black/Gray	Combat Gauze	80-0921
Blue/Gray	Combat Gauze	80-0922
Red/Gray	Combat Gauze	80-0917

Every **OCHO IFAK** contains: //////////////



PLUS these Individual Components based upon configuration of choice:



MICRO IFAK® Individual First Aid Kit

- Ruggedized, vacuum packed, re-sealable IFAK designed to treat injuries ranging from major bleeding trauma and to minor cuts and abrasions
- First aid for minor to severe injuries for single patient
- Waterproof, Ultra-compact & lightweight
- Provides emergency treatment solutions until definitive care is available

SPECIFICATIONS

L 5.5 in. x W 2.25 in. x H 7 in.
Weight: 0.71 lb

MICRO IFAK

ITEM

82-0002

MICRO IFAK:



Adventure Trauma Kits

Individual First Aid Kit **SOLO IFAK**

- Ruggedized, nylon kit designed to treat injuries ranging from major bleeding trauma and to minor cuts and abrasions
- First aid for minor to severe injuries for single patient
- Ultra-compact & lightweight
- Optimized organization of contents
- Provides emergency treatment solutions until definitive care is available

SPECIFICATIONS

L 5 in. x W 3.5 in. x H 7.5 in.
Weight: 135 lb

SOLO IFAK

ITEM #

82-0000

SOLO IFAK:

 QTY 1 C-A-T® Tourniquet	 QTY 1 ETD™ 6 in. Responder	 QTY 1 Bleeding Control Dressing	 QTY 4 Gauze, Compressed
 QTY 1 Survival Blanket	 QTY 1 BurnTec 4"x4" Dressing	 QTY 1 Bandage Silver 4" x 4.5"	 QTY 1 Bandage Silver 2" x 2"
 QTY 2 Gauze Pads 2"x2" (2pk)	 QTY 2 Gauze Pads 4"x4" (2pk)	 1 Surgical Tape 2" Flatfold	 QTY 10 Flexible Fabric Bandages
 QTY 2 Gauze, Petrolatum	 QTY 1 Tagaderm Film 4" x 4.75"	 QTY 1 Tagaderm Film 2.375" x 2.75"	 QTY 1 (10pk) Wound Closure Strips
 QTY 2 Triangular Bandage	 QTY 1 MicroShield CPR Face Mask	 QTY 1 Syringe, Irrigation	
 QTY 2 PES Eye Shield	 QTY 1 Mini Duct Tape, 2 in.	 QTY 1 Permanent Marker, Small	 QTY 1 MED ID Patch
 PAIR 1 Nitrile Gloves, Lg	 QTY 1 Pen Light	 QTY 1 Forceps, Splinter, 4.5 in.	 QTY 1 Trauma Shears, 6.25 in.

X-PAK Expeditionary Personnel Aid Kit

- Compact nylon kit designed to address the first aid needs of a small group with injuries ranging from major bleeding trauma and fractures to minor cuts and abrasions
- Treat minor to severe injuries – 2 to 3 patients
- Provide emergency treatment solutions until definitive care is available
- Compact and capable – ideal for adventure activities or extended travel

SPECIFICATIONS

L 9 in. x W 5 in. x H 6.5 in.
Weight: 2.72 lb


X-PAK

ITEM

82-0001



X-PAK IFAK:



Adventure Trauma Kits

ATAK

Amphibious Trauma Aid Kit

- Expedition quality first aid kit designed with input from water safety, wilderness medicine and tactical medical experts
- Provides first aid for minor to severe injuries for multiple patients
- Rugged, compact & waterproof case ideal for the adventure sports and outdoor enthusiasts
- Internally compartmentalized into treatment categories for rapid and successful management of an injury
- Provide emergency treatment solutions until definitive care is available

SPECIFICATIONS

L 9.5 in. x W 13.5 in. x H 6 in.
Weight: 9.8 lb


ATAK

ITEM

85-0639



ATAK:





OUTDOOR

Range Trauma Kits

RANGE TRAUMA AID KIT Hard Case

- Water-tight case maintains the integrity of your medical products when you need them most
- Kit contents arranged in a modular system featuring rugged lid storage organizer and removable nylon pouches
- Labeled internal pouches for ease of use

RANGE TRAUMA AID

COLOR	ITEM #
Orange	85-0889

SPECIFICATIONS

H 9.5 in. x W 13.5 in. x D 6 in.
Weight: 8 lb 12 oz



RANGE TRAUMA AID KIT HARD CASE:



RANGE TRAUMA KIT

- Essential equipment for medical and non-medical personnel to address injuries common to firearm ranges and training facilities
- Treatment capabilities include penetrating injuries, chest wounds, burns, eye injuries, minor lacerations, sprains and fractures

RANGE TRAUMA KIT

*NATO NSN#

HEMOSTATIC	COLOR	ITEM #	NSN #
Combat Gauze	Orange	80-0299	----
ChitoGauze	Orange	80-0213	6545-01-610-5411*

SPECIFICATIONS

H 7.75 in. x W 13 in. x D 9 in. | Wgt: 5 lb 10 oz



RANGE TRAUMA KIT:



K-9 HANDLER IFAK

K-9 HANDLER IFAK:



SPECIFICATIONS

H 7 in. x W 4 in. x D 2.75 in
Weight: 1 lb 3 oz

Bag only option also available at NARescue.com

- Designed for immediate lifesaving treatment to a wounded handler, operator or K-9
- Treatment capabilities range from major trauma typical to operating in tactical or patrol environments to K-9 emergencies

K-9 HANDLER IFAK

COLOR	HEMOSTATIC	ITEM #
Black	----	80-0439
Coyote	----	80-0440
Multi-Cam	----	80-0441
OD Green	----	80-0442
Black	Combat Gauze	80-0443
Coyote	Combat Gauze	80-0444
Multi-Cam	Combat Gauze	80-0445
OD Green	Combat Gauze	80-0446

K-9 TACTICAL FIELD KIT




Bag only option also available at NARescue.com

- Designed for the handler or medic to deliver lifesaving treatment to a wounded K-9 at the point of injury
- K-9 treatment capabilities from major trauma typical to operating in tactical or patrol environments to more common K-9 emergencies
- Provides the essential equipment to render care to the injured K9 when veterinary care is not immediately available

SPECIFICATIONS

H 7.75 in. x W 13 in. x D 9 in.
Weight: 5 lb 2 oz

K-9 TACTICAL FIELD KIT:



K-9 TACTICAL FIELD KIT

COLOR	HEMOSTATIC	ITEM #
Black	Combat Gauze	80-0300
Coyote	Combat Gauze	80-0301
OD Green	Combat Gauze	80-0303
3 options online	ChitoGauze	See NARescue.com

Also available: **K-9 Litter** (#60-0076), See p. 25

Patrol **VEHICLE TRAUMA KIT**

- Provides personnel with an easily accessible, quick to deploy trauma kit from their vehicle
- Modular storage system comes organized with an array of medical products ranging from basic first aid to major trauma care
- Equipment is packed in the kit to offer basic first aid and minor injury equipment on one side and major trauma equipment on the other
- Exterior includes MOLLE-style system for extra equipment as needed

SPECIFICATIONS

H 7 in. x W 12 in. x D 5 in.
Weight: 6 lb




Bag only option also available at NARescue.com

PATROL VEHICLE TRAUMA KIT:

VEHICLE TRAUMA KIT

LEVEL	ITEM #
Basic	80-0542
Combat Gauze LE	80-0859



Workplace Injury Kits

Class A **TRAUMA & FIRST AID KIT**

- Designed as a portable First Aid Kit for the workplace **exceeding ANSI/ISEA Z308.1-2015 Standard for Class A kits**
- Able to treat the most common workplace injuries with the added benefit of being able to control major bleeding with the C-A-T Tourniquet - the safest and most effective Tourniquet on the market!
- Packaged in a convenient zippered nylon bag or rugged hard case that is easily transported to the patient's side



TRAUMA & FIRST AID KIT CLASS A

CASE	DRESSING	ITEM #
Soft Case	Compressed Gauze	80-0947
Soft Case	Bleeding Control Dressing	80-1045
Hard Case	Compressed Gauze	80-1033
Hard Case	Bleeding Control Dressing	80-1046

SPECIFICATIONS

Soft Case:

L 14.75 in. x H 8 in. x D 4 in.
Weight: 2.94 lb

Hard Case:

L 13 in. x H 9.5 in. x D 5 in.
Weight: 6 lb

TRAUMA & FIRST AID KIT CLASS A:

<p>QTY 1 C-A-T® Tourniquet</p>	<p>QTY 2 Responder Compressed Gauze</p>	<p>-OR-</p>	<p>QTY 1 Responder Bleeding Control Dressing</p>	<p>QTY 1 Bleeding Control Dressing</p>
<p>QTY 1 ETD™ 4 in. Responder, Flat</p>	<p>QTY 1 MicroShield CPR Face Mask</p>	<p>QTY 1 Elastic Wrap 2 in. Bandage</p>	<p>QTY 1 Elastic Wrap 4 in. Bandage</p>	<p>QTY 1 Gauze Pads 4"x4" (2pk)</p>
<p>QTY 2 Gauze Pad, 5" x 9"</p>	<p>QTY 16 Flexible Fabric Bandages</p>	<p>QTY 10 Burn Gel</p>	<p>QTY 1 BurnTec 4"x4" Dressing</p>	<p>QTY 1 Instant Cold Pack</p>
<p>QTY 1 Surgical Tape, 1 in.</p>	<p>QTY 1 Triangular Bandage</p>	<p>QTY 1 SAM® Splint II</p>	<p>QTY 1 Eyewash Solution 1 oz</p>	<p>QTY 2 PES Eye Shield</p>
<p>QTY 10 Antiseptic Towelette</p>	<p>QTY 10 Antibiotic Ointment</p>	<p>QTY 6 Hand Sanitizer</p>	<p>PAIR 2 Responder Nitrile Gloves, Lg</p>	<p>QTY 1 Responder Trauma Shears, Lg</p>
<p>QTY 1 First Aid Pocket Guide</p>	<p>QTY 1 Just in Time Instructions</p>	<p>QTY 1 Kit Contents Inventory Card</p>	<p>QTY 1 *Card, Inspection</p>	

* Inspection Card is included in the Soft Case version only

Workplace Injury Kits

Class B **TRAUMA & FIRST AID KIT**

- Designed as a portable First Aid Kit for the workplace **exceeding ANSI/ISEA Z308.1-2015 Standard for Class B kits**

- Same supplies as the Class A Kit but with larger quantities to accommodate larger workplaces

- Packaged in a convenient zippered nylon bag or rugged hard case that is easily transported to the patient's side

TRAUMA & FIRST AID KIT CLASS B

CASE	DRESSING	ITEM #
Soft Case	Compressed Gauze	80-0948
Soft Case	Bleeding Control Dressing	80-1047
Hard Case	Compressed Gauze	80-1034
Hard Case	Bleeding Control Dressing	80-1048

SPECIFICATIONS


Soft Case:

L 14.75 in. x H 8 in. x D 7.5 in.
Weight: 4.87 lb

Hard Case:

L 13.5 in. x H 11 in. x D 6 in.
Weight: 9 lb

TRAUMA & FIRST AID KIT CLASS B:



ANSI Kits designed to treat the most common workplace injuries with the added benefit of being able to control major bleeding...



COMMUNITY
PREPAREDNESS

PUBLIC ACCESS

BLEEDING CONTROL



SPECIFICATIONS

H 7 in. x W 5 in. x D 2.75 in.
Weight: 7 oz

STB KIT

- Designed to provide the essential bleeding control equipment needed to help save a life
- Small, compact and easy to use
- Kits contain the same Tourniquet and Hemostatic Dressing technology recommended by the U.S. Military & the American College of Surgeons Committee on Trauma to control bleeding

STB KIT

ITEM #

85-1574

STB KIT contains:



Public Access Bleeding Control Kits

Public Access **INDIVIDUAL BLEEDING CONTROL KIT**

Essential equipment needed for persons to take action as immediate responders in helping to stop life threatening bleeding

Nylon Bag Option:

- Easy to restock in cases of use
- Convenient to carry and easy to open in stressful situations
- Opens flat for ease of use and immediate identification of enclosed devices
- Easily recognizable red nylon bag with attached PVC Bleeding Control Patch

Vacuum Sealed Option:

- Cost-effective alternative to Nylon
- Compact and rugged Vacuum-Sealed pouch
- Large Identifying Card Vacuum-Sealed inside the pouch for ease of identification

SPECIFICATIONS

Nylon Bag: H 7 in. x W 5 in. x D 2.75 in.

Vacuum Sealed: H 7 in. x W 5 in. x D 2.75 in.

PUBLIC ACCESS INDIVIDUAL BLEEDING CONTROL KIT


LEVEL	BAG	WEIGHT	ITEM #
Basic	Red Nylon	15 oz	80-0452
Intermediate	Red Nylon	16.5 oz	80-0453
Advanced BCD	Red Nylon	17.4 oz	80-0523
Advanced	Red Nylon	17.4 oz	80-0454
Basic	Vacuum Sealed	11.7 oz	80-0465
Intermediate	Vacuum Sealed	13.25 oz	80-0466
Advanced BCD	Vacuum Sealed	14 oz	80-0522
Advanced	Vacuum Sealed	14 oz	80-0467



Every **INDIVIDUAL BLEEDING CONTROL KIT** contains:



PLUS these Individual Components based upon configuration of choice:



Public Access Bleeding Control Kits

Public Access Bleeding Control **TWIN PACK**

- Designed to provide essential equipment that empowers the general public to take action as immediate responders in stopping life threatening bleeding
- Provides intuitive and easy-to-use tools that are proven to help to save lives
- Convenient kit that includes the contents of 2 Public Access Individual Bleeding Control Kits
- Packaged in a compact, ruggedized nylon bag designed to fit next to an AED in the cabinet

2-in-1 Kit sized to fit next to an AED in the cabinet!



SPECIFICATIONS

Nylon Bag: L 10 in. x H 9 in. x D 3 in..

BLEEDING CONTROL TWIN PACK

LEVEL	WEIGHT	ITEM #
Basic	166 lb	80-0951
Intermediate	187 lb	80-0952
Advanced BCD	180 lb	80-0953

Every PUBLIC ACCESS BLEEDING CONTROL TWIN PACK contains:



QTY **2** C-A-T® Tourniquet

QTY **2** ETD™ 4 in. Responder, Flat

PAIR **2** Responder Nitrile Gloves, Lg



QTY **2** Just in Time Instructions

QTY **2** Permanent Marker, Small

QTY **2** Survival Blanket

QTY **1** Responder Trauma Shears, Lg

PLUS these Individual Components based upon configuration of choice:



BASIC

QTY **2** Responder Compressed Gauze

INTERMEDIATE

QTY **2** Responder Compressed Gauze

QTY **2** HyFin® Vent Compact, Twin

ADVANCED BCD

QTY **2** Bleeding Control Dressing

QTY **2** HyFin® Vent Compact, Twin



See our complete line of **PUBLIC ACCESS BLEEDING CONTROL SOLUTIONS** including **Wall-Mounted Stations** at www.NARescue.com

Public Access Bleeding Control Trainers

INDIVIDUAL TRAINER Public Access Bleeding Control

- The appropriate equipment needed to teach individuals the necessary life-saving skills of bleeding control
- Corresponds with NAR's single Public Access Bleeding Control Kits in Basic, Intermediate, or Advanced
- Packaged in a compact, blue nylon carry bag with PVC Bleeding Control Patch

INDIVIDUAL BLEEDING CONTROL TRAINER

LEVEL	WEIGHT	ITEM #
Basic	15 oz	80-0901
Intermediate	16.5 oz	80-0902
Advanced	16.5 oz	80-0903




Every **INDIVIDUAL PUBLIC ACCESS BLEEDING CONTROL TRAINER** contains:



SPECIFICATIONS

Nylon Bag:
H 7 in. x W 5 in. x D 2.75 in.

PLUS these Individual Components based upon configuration of choice:



SPECIFICATIONS

H 7 in. x W 5 in. x D 2.5 in.
Weight: 7 oz

STB TRAINER

ITEM #
85-1880

Public Access Bleeding Control **STB TRAINER**

- Designed to provide the essential training equipment to teach how to control bleeding
- Contains training items that correspond with the necessary tools recommended by the American College of Surgeons Committee on Trauma to equip bystanders to help save lives in the instance of life threatening bleeding
- Includes illustrated "Just in Time" pictorial instructions

PUBLIC ACCESS BLEEDING CONTROL STB TRAINER contains:



READY HEAT™ BLANKET



As an individual unit, the self-heating, 4-cell liner provides 10 hours of continuous dry heat after reaching peak temperature approximately 15-20 minutes after opening

READY HEAT

ITEM #
ZZ-0052

SPECIFICATIONS

Open: L 48 in. x W 36 in.
Weight: 1 lb 11 oz

Emergency SURVIVAL WRAP

- Reflects heat back to the body
- Larger size than standard survival blankets for full body coverage
- Reversible high visibility orange and silver color



SURVIVAL WRAP

ITEM #
ZZ-0834

SPECIFICATIONS

Unfolded: L 90 in. X W 60 in.
Weight: 3.7 oz

DRAGON HANDLE

- Maintain an anatomically correct posture for rapidly extracting greater weight
- Decreases manpower needs in "zeroed in" locations
- Black Omega, non-locking "D" carabiner rated for 31kN (6,968.8 lb).

SPECIFICATIONS

H 6.3 in. x W 4.1 in. x D 2.25 in.
Weight: 11.3 oz
Device Length: L 42 in.
US Patent(s): Pending

DRAGON HANDLE

ITEM #	NSN #
90-0004	6530-01-598-9595



NAR QUIKLITTER

LITE, Standard, Responder or XL

- Multi-use, low cost disposable evacuation platforms for safe, effective emergency casualty transport & transfer
- Water-resistant lightweight heavy duty construction with ten large hand-holds
- Super compact LITE is packaged to fit in a cargo pocket
- The XL is 12 in. wider and 6 in. longer featuring 12 hand-holds and two built-in patient retention straps

SPECIFICATIONS

LITE: Payload 750 lb.
Litter Bed: L 74 in. x W 27 in.
Weight: 12 oz

STANDARD/RESPONDER:

Payload 500 lb.
Litter Bed: L 74 in. x W 28 in.
Weight: 13.3 oz

XL: Payload 1,000 lb
Litter Bed: L 80 in. x W 40 in.
Weight: 2.2 lb

QUIKLITTERS

SIZE	ITEM #
LITE	60-0071
Standard	60-0061
Responder	60-0080
XL	60-0065



HASTY HARNESS *Multi-purpose Harness*

- Indispensable gear for every webbing equipment toolbox
- Resourceful, multi-purpose harness for unlimited applications & maximum versatility
- Constructed of 22 ft loop of 1 in. Tech Tape (MBS: 4,400 lb) with opposing color thread to aid in visual safety inspection

SPECIFICATIONS

Loop: L 10 ft 11 in. x W 1 in.
Circumference - 21 ft 10 in.
Min. Break Str: Tech Tape (4,400 lb)
Sewn Sling (6,000 lb)

HASTY HARNESS NATO NSN*

COLOR	ITEM #	NSN #
Red	90-0044	----
Black	90-0001	6530-01-674-9774*





NAR
K-9 LITTER
& Carrying Strap Accessories

- Lightweight, compact litter for the movement of an injured canine
- Litter can be carried in cargo pocket or pouch
- Leg cut outs to safely carry the canine in the litter
- Six large and four small carry handles that allow for one to four person carries
- Optional K-9 Carry Strap & 2 Carabiners for over the shoulder carry of injured K-9 also available

SPECIFICATIONS

Litter: W 35 in. x H 48 in.
Adjustable Carry Strap: L 9 to 54 in. x W 2 in.
Carabiner open gate clearance: 27 mm
Packaged: W 5.5 in. x H 10 in. x D 2.5 in. | Weight: 138 lb

K-9 LITTER

OPTIONS	ITEM#
K-9 Litter Only	60-0074W
Carry Strap & Carabiners	ZZ-1004
K-9 Litter with Carry Strap & Carabiners	60-0076

BENCHMADE® KNIVES & CUTTERS



1 HOOK SAFETY CUTTERS by Benchmade®

- Hook cutters designed to quickly and reliably cut through a wide array of materials
- Extensively tested and widely used by the US Military for medical and egress purposes

ITEM #	STYLE	NSN #
ZZ-0652	7 Hook	----
ZZ-0250	8 Hook	----

SPECS

7 Hook: L 4.3 in.
8 Hook: 6.5 in.
Weight: 1.5 oz

2 SOCP RESCUE TOOL by Benchmade®

- Fixed 440C high-chromium stainless steel blade hook for rescue performance and less-than-lethal self-defense
- Includes carbide glass breaker, oxygen wrench, single-piece handle with finger loop & molded sheath
- Designed by Special Operations Combatives Program (SOCP)

SPECS

Blade: 2.5 in. Overall: 6.75 in.
Weight: 188 oz

ITEM #	NSN #
ZZ-0760	----

3 H₂O 112 SBK by Benchmade®

- N680 ultra stainless steel blade with opposing fixed bevel blade style, blunt tip, & integrated rescue hook
- Boltaron® sheath designed to fit current Navy SAR kits
- 90-0078 (Orange handle) includes carbide glass breaker & is approved for Navy SAR - "So Others May Live"

ITEM #	HANDLE	NSN #
90-0078	Orange	1095-01-677-1013
90-0097	Black	----

SPECS

Blade: 3.5 in.
Overall: 8.02 in.
Weight: 3.2 oz

4 TRIAGE 916 SBK by Benchmade®

- Dynamic triple utility tool featuring three essential rescue tools in a single compact design
- N680 high corrosion resistant, partially serrated, opposing bevel blade, 440C folding hook & carbide glass break

SPECS

Blade: 3.5 in
Overall: 8.2 in.
Weight: 51 oz

ITEM #	HANDLE	NSN #
90-0096	Orange	----
90-0060	Black	----

5 AUTOTRIAGE 9160 SBK by Benchmade®

SPECS

Blade: 3.58 in. Overall: 8.31 in.
Weight: 5.75 oz

ITEM #	NSN #
90-0061	----

- Packaged in high-strength aluminum with AXIS® auto pull release and integrated safety
- Auto opening N680 highly corrosion resistant blade with modified drop-point

6 TRIAGE 917 BK -or- SBK by Benchmade®

- Black drop point style CPM-S30V steel blade with option of a smooth (BK) or partially serrated edge (SBK)
- 440C folding hook & carbide glass breaker

SPECS

Blade: 3.5 in Overall: 8.3 in.
Weight: 5.24 oz

ITEM #	STYLE	NSN #
90-0092	BK	----
90-0093	SBK	----

Leatherman®

RAPTOR™ MULTI-TOOL

High-performance multi-tool equipped with a variety of mission critical tools essential to emergency medical practitioners in the field

1. 420HC Stainless Steel Folding Medical Shears

2. Strap Cutter

3. Oxygen Tank Wrench

4. Ring Cutter

5. Carbide Glass Breaker

6. 5 cm Ruler

7. Holster with Clip

RAPTOR MULTI-TOOL

HANDLE	ITEM #	NSN #
Black	70-0014	5110-01-626-8766
Blue	70-0264	----
Coyote	70-0265	----
Red	70-0266	----

**SPECS**

Shear Length: 19 in.

Closed Length: 5 in.

Weight: 5.8 oz

HANDLE COLORS:

**Leatherman® SUPERTOOL 300**

Heavy-duty, all purpose multi-tool featuring 19 devices in one configuration

SPECIFICATIONS

Closed: L 4.5 in.

Primary Blade: L 3.2 in.

Weight: 9.6 oz

SUPERTOOL 300

ITEM #

ZZ-0542

1. Needlenose Pliers

2. Regular Pliers

3. Replaceable Wire Cutters

4. Hard-Wire Replaceable Cutters

5. Stranded-Wire Cutters

6. Electrical Crimper

7. Wire Stripper

8. 420HC Knife

9. 420HC Serrated Knife

10. Saw

11. Awl with Thread Loop

12. 9 in. Ruler

13. Can Opener

14. Bottle Opener

15. Wood/Metal File

16. Phillips Screwdriver

17. Large Screwdriver

18. Medium Screwdriver

19. Small Screwdriver

**Leatherman® MUT®**

Tactical and practical multi-tool for military, law enforcement and civilian shooters - features 16 devices

1. Needlenose Pliers

2. Regular Pliers

3. Replaceable Wire Cutters

4. Hard-Wire Replac. Cutters

5. Stranded-Wire Cutters

6. Electrical Crimper

7. 420HC Knife

8. Saw

9. Replaceable Cutting Hook

10. Hammer

11. Bolt Override Tool

12. Bronze Carbon Scraper

13. #8-32 Cleaning Rod/Brush Adapter

14. Firearm Disassembly Punch

15. Carabiner/Bottle Opener Punch

16. Large Bit Drive

**SPECIFICATIONS**

Closed: L 5 in.

Primary Blade: L 3 in.

Weight: 11.2 oz

MUT

ITEM #

ZZ-0539



PREPARE *for the* **ADVENTURE**



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

35 TEDWALL COURT | GREER SC 29650

REV071020