



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS WITH A **MISSION®**

PUBLIC SAFETY | COMMUNITY PREPAREDNESS | FEDERAL AGENCIES

Volume 8 | Number 1 | REV062520



PRODUCTS WITH A MISSION®

***Over twenty years of helping to reduce
preventable causes of death
on the battlefield and in our home towns...***



Table of Contents

Volume 8 | Issue 1

4

BALLISTIC PROTECTION

Body Armor Overview | Ballistic Armor | Ballistic Plates | Ballistic Helmets
Ballistic Eye Pro | Accessories

10

M.A.R.C.H. PRODUCTS *for Point of Wounding*

Massive Hemorrhage: Tourniquets/Holders | Bandages & Gauze | Hemostatics
Airway Management: Naso | Cric Kits | BVMs | Supraglottic | Laryngoscopy
Respiratory Support: Chest Seals | Tension Pneumothorax
Circulation/Vascular Access: Saline Lock | Intraosseous | Fluid Warming
Hypothermia Prevention: HPMK | Survival Blankets

27

MEDICAL TRAUMA KITS

Individual First Aid Kits (IFAKs) | Multiple Injury Trauma | Range Trauma
Canine | School Crisis | School Resource Officer | Nylon Bags & Accessories
Workplace Injury: ANSI Class A & B

60

PUBLIC ACCESS BLEEDING CONTROL KITS

Individual (Vacuum-Sealed or Nylon) | Multi Packs | Stations | Accessories

68

TRAINING

Public Access Trainers | PHTLS TCCC/TECC Modular System
Part Task Trainers | Individual Training Components

76

EVACUATION & PATIENT MOVEMENT

High Threat Extraction Gear | Rope Rescue Equipment | Framed Litters
Soft Litters | Litter Accessories

84

TREATMENT SUPPORT ACCESSORIES

Triage Cards | Medical/Surgical Equipment | Field Hospital Equipment
Diagnostic Equipment | Hyperthermia | Splinting | Burn Dressings | Eye Care
Trauma Shears | Rescue Tools | Storage & Transport | Nitrile Gloves | Accessories

100

ORDERING INFORMATION

Quick reference guide to Ordering



TECHNOLOGY®

Items marked with this icon have been specifically branded with NAR's signature Red Tip Technology® packaging. Use these items with confidence, knowing that they have been engineered for less-than-optimal, often hostile conditions that exist when our equipment is needed the most... **Accept no imitations!**

BALLISTIC

M.A.R.C.H.

TRAUMA KITS

PUBLIC ACCESS

TRAINING

EVACUATION

TREATMENT SUPPORT

ORDER

OVERVIEW | Body Armor

The term body armor is generally used to describe vests designed to provide ballistic protection to the vital organs in the torso, armor has typically been associated with tactical medics and rescue task force members whose primary objective is to provide medical care in high threat environments. With the increase of active shooter and mass casualty incidents in today's field of operation, Fire & EMS personnel are realizing the critical need for ballistic protection on a daily basis.


Though Fire & EMS personnel have been limited to adapting body armor traditionally designed for the military & law enforcement, North American Rescue is now confronting this challenge by designing ballistic solutions specifically for the unique needs of the Fire & EMS community.

Soft & Hard | 2 Main Types of Body Armor

SOFT ARMOR offers critical protection against handgun threats along with the capability for extended daily wear. NAR provides an optimal over-the-uniform soft armor vest that features flexible layers of ballistic resistant materials constructed from a synthetic Aramid fiber with a very high strength to weight ratio. These super strong fibers are woven into multiple layers that are stacked to create an armor panel. The number and orientation of each layer influences its overall performance. When a handgun bullet strikes the armor panel, it gets caught in a web of high performance fibers that stretch and disperse the impact energy causing the bullet to deform or "mushroom". This dispersion of energy not only prevents bullet perforation, but also helps minimize the possible injury to underlying tissue.

HARD ARMOR is designed to provide a higher degree of protection from rifle threats. Constructed from a variety of hard ballistic-resistant materials such as ceramic, compressed UHMWPE (ultra-high-molecular-weight-polyethylene), metallic plates, or a composite of multiple materials, each of these rigid materials offer their own unique properties. Ceramic materials, for example, are required for so called "armor piercing" bullets or special threats. Hard armor plates are classified as either Stand-Alone (no soft armor required) or ICW - *in conjunction with* (soft armor backing panel required) and rated as single hit or multi-hit. ▶

How Soft Armor disperses energy from a handgun projectile...



NO armor is 100% bullet proof!

One truth stands out when it comes to protective gear – ballistic protection is NOT created equal. Differing armors are designed to protect against specific types of ammunition, and no body armor will stop every threat. So when choosing a ballistic solution, consider the threats most likely to be faced in your area of operation. At a bare minimum, armor should protect against handguns.



NIJ 0101.06 Ballistic Standards

Special Threats

In addition, the FBI and DEA have developed their own stringent protocols to test body armor that exceeds the NIJ standards.

5.56mm ammunition (M855)

was designed for military use and has a steel core. It is the most common and widely used AR ammo on the market today. To provide requisite protection, a ceramic component must be included in the hard armor plate.



Choose the right ARMOR for YOU...

There are several key considerations in choosing the right body armor: level of protection, weight, comfort, mobility, and cost. It is important to choose an armor system that will provide protection but also enhance your ability to perform your job. Particularly for Fire & EMS personnel, it is critical that the system facilitates your ability to provide the life saving skills and equipment.

North American Rescue has the optimal ballistic protective solution designed as a system specifically for Fire & EMS personnel and the duties they must perform.

The **National Institute of Justice (NIJ)** is the government agency responsible for establishing performance standards and testing certification for body armor. These body armor classifications: levels IIA, II, IIIA, III and IV correspond with the weaponry they protect against. NIJ certified armor ensures ballistic performance and the completion of proper testing standards.

The NIJ conducts two types of ballistic testing: Perforation/BackFace Signature (P-BFS) and Ballistic Limit (V50). No complete or partial perforation of the armor panel is allowed during testing. In addition, no more than 44mm (1.73 inches) of indentation is allowed into the ballistic clay backing utilized during shoot testing. The back face signature translates to the upper tolerable limits of blunt force injury to the user. The ballistic limit (V50) is the velocity at which the projectile is expected to perforate the armor 50% of the time. It is important to ensure that any armor package has been certified by an NIJ approved laboratory.

Level IIA:

Tested to stop 9mm and .40 S&W ammunition fired from short barrel handguns

No rifle ammunition protection

Level II:

Tested to stop 9mm and .357 Magnum ammunition fired from short barrel handguns

No rifle ammunition protection.

Level IIIA:

Tested to stop .357 SIG and .44 Magnum ammunition fired from longer barrel handguns.

No rifle ammunition protection

Level III:

Tested to stop 7.62mm FMJ (full metal jacket) lead core rifle ammunition

Level IV:

Tested to stop .30 cal steel core armor piercing rifle ammunition

NIJ Threat Projectiles & Related Armor Classifications



MKIII

**TACTICAL RESPONDER
BALLISTIC VEST** *with Cummerbund*

- Developed for providing Fire and EMS personnel in active shooter and mass casualty incidents with Level IIIA Soft Armor Ballistic Protection
- Allows for full range of motion while providing superb coverage area
- Includes front and rear NIJ 0101.06 Level IIIA ballistic panels and two integrated plate pockets designed to accommodate additional 10 in. x 12 in. or 11 in. x 14 in. hard armor plates (each sold separately)
- Featuring extra side protection with Level IIIA armor cummerbund equipped with easy hook and loop attachment to the vest for rapid user adjustment
- MOLLE-style PALS webbing on the front, back and around the waist allows for attaching extra gear and for unit customization
- Rapidly adjustable armor system customizes to any team member's frame within seconds
- Can be stored in and deployed from all standard emergency vehicles

**MKIII TACTICAL RESPONDER
VEST with CUMMERBUND**

COLOR	ITEM #
Black	70-0251
Blue	70-0253
Red	70-0254
Coyote	70-0252

DIMENSIONS

Vest: H 19 in. x W 12.5 in. x D 4 in.

Vest with Cummerbund:

H 19 in. x W 15.5 in. x D 4 in.

Maximum Waist Size: 62 in.

Total Weight: 9 lb 6 oz

BALLISTIC SPECIFICATIONS

Protection Type/Threat Level: NIJ 0101.06 IIIA
Special Threat Tested

Construction: 3000D Twaron® Microlam |
Twaron® 2-Ply Unidirectional Fabric
W.R.A.P.™ - Water Repellant Armor Perimeter
(patent pending)

Protective Cover: Brookwood 70D
with Laminated Polyurethane

Thickness: 0.275 in. | AD= 1.17 psf

Model#: Fusion IIIA AFC-T3A-525
(NIJ Compliance Product List)

Testing Certification: NIJ 0101.06
and Special Threat Tested

Special Threats:

.44 Mag V50: 1615 fps

.357 Sig V50: 1735 fps

FN SS197 40 grain VO

Win SXT 9mm+P+ VO

7.62 x 25 Tokarev 85 grain VO

.357 Sig WIN SXT 125 grain VO

Speer GD .357 Mag VO

2 grain RCC V50: 3013 fps

4 grain RCC V50: 2545 fps

16 grain RCC V50: 2225 fps

17 grain RCC V50: 2077 fps

MKIII

RTF TACTICAL RESPONDER BALLISTIC VEST SYSTEM

The practical solution for Fire and EMS personnel in active shooter and mass casualty incidents, with an emphasis on integrating Tactical Emergency Casualty Care...

- NIJ 0101.06 Level IIIA Soft Armor Ballistic Protection equipped with medical supplies that address the leading causes of preventable death in the tactical environment
- Front and back hard armor plate pockets able to accommodate additional 10 in. x 12 in. or 11 in. x 14 in. hard armor plates (each sold separately)
- Includes Cummerbund-style Side Armor Protection with easy hook and loop attachment to the vest for rapid user adjustment - expandable up to 62 in.
- MKIII RTF Vest Tactical Responder Vest System includes:
 - 1 x MKIII Tactical Responder Ballistic Vest with Cummerbund
 - 1 x RTF Chest Pouch Medical Kit with Luminous Med ID Patch
 - 2 x Nylon Tourniquet Holder with C-A-T®
- ID patches are also available for separate purchase and are customizable upon special request.



MKIII RTF TACTICAL RESPONDER BALLISTIC VEST SYSTEM

VEST	ACCESSORIES	ITEM #
Black	Black	80-1092
Blue	Black	80-1093
Red	Red	80-1095
Coyote	Coyote	80-1094

DIMENSIONS

Vest, Cummerbund & Medical Kit:
H 19 in. x W 15.5 in. x D 7 in.

Maximum Waist 62 in.

Total Weight: 12.4 lb

Medical Kit (separately):
H 7.75 in. x W 9.75 in. x D 4.5 in.
Weight: 2 lb 10.3 oz

RTF CHEST POUCH MEDICAL KIT (also available separately)

COLOR	ITEM #
Black	80-0407
Coyote	80-1109
Red	80-0406

on MKIII CUMMERBUND:



RTF CHEST POUCH MEDICAL KIT Contains:



Ballistic Rifle Plates

Comfortable & lightweight...

*The clear choice for an economical, multi-strike
NIJ certified advanced threat armor plate under 5 lb...*

NIJ Level III Plus Enhanced **M7X**

- Protection Type/Threat Level:
NIJ III Plus
- Configuration: Stand-Alone
- Construction: Ceramic
/Polyethylene /Composite
- Protective Cover:
Polyurethane
- Hit Rating: Multi-Strike
- Model#: DM855 (NIJ
Compliance Product List)
- Positive Buoyant
- 7 Year Warranty

- Testing Certification: NIJ 0101.06 &
Special Threat Tested | Tested to DEA
standards against M855* and 7.62 x 39 MSC

Special Threats:

- 5.56mm x 45mm 55gr M193 FMJ (NATO)
- 5.56mm x 45mm 62gr Ball (M855)*
- 7.62mm x 39mm 123gr PS Ball
(Mild Steel Core)
- 7.62mm x 51mm 147gr M80 FMJ



**Tested to stop M855 5.56 "Green Tip" -
the cheapest, most widely available AR15 ammunition*



M7X NIJ LEVEL III+ ENHANCED

STYLE	WEIGHT	ITEM #
Full Size	5.15 lb	70-0255
Shooter	4.65 lb	70-0256

SPECIFICATIONS

W 10 in x H 12 in x D 125 in

Ballistic Helmets

Batlskin® Full Cut **VIPER A3 HELMET** with Modular Suspension System



SPECIFICATIONS

Weight: 3 lb

VIPER A3 with MSS

COLOR	SIZE	ITEM #	COLOR	SIZE	ITEM #
Black	S	70-0081	Coyote	S	70-0247
Black	M	70-0082	Coyote	M	70-0248
Black	L	70-0083	Coyote	L	70-0249
Black	XL	70-0084	Coyote	XL	70-0250

- The cost effective ballistic and head impact protection solution
- Exceeds U.S. Army ACH impact standards and tested with multiple hits and across all temperature ranges
- Stops NIJ IIIA level threats (9mm and 44 mag) (NIJ 0101.06)
- Offers unsurpassed impact protection and ballistic rated hardware tested for ballistic performance from -40° F to 158° F
- Compatible with NVGs, Chem/Bio Masks & communication devices
- **Modular Suspension System (MSS)** features 4-point integrated harness with unique buckle system (no loose straps and Dial Ratchet Mechanism to adjust the fit band with one hand
- Adjustable nape pad ensures helmet stability
- Adjustable crown net allows sweat wicking, comfort and low pressure wear, providing a free flow of air to aid in cooling
- Low profile impact material does not contact the head in normal operation



M3

NIJ Level III Plus

The lightest armor in its category!

M3 NIJ LEVEL III +

STYLE	SIZE	WGT	ITEM #
Full Size	10 x 12	3.15 lb	70-0257
Shooter	10 x 12	2.9 lb	70-0259

STYLE	SIZE	WGT	ITEM #
Full Size	11 x 14	3.97 lb	70-0258
Shooter	11 x 14	3.5 lb	70-0260

- Protection Type/Threat Level: NIJ III Plus
- Configuration: Stand-Alone
- Construction: Polyethylene
- Protective Cover: Polyurethane
- Hit Rating: Multi-Strike
- Model#: DKX MAX III (NIJ Compliant Product List)
- Positive Buoyant
- 7 Year Warranty
- Testing Certification: NIJ 0101.06 & Special Threat Tested

Special Threats:

- .223 Rem 60gr Partition
- 7.62mm x 39mm 123gr Mild Steel Core (MSC)
- 7.62mm x 51mm 147gr M80 FMJ, V50 = 3288 fps
- 5.56mm x 45mm 55gr M193 FMJ (NATO)

Ballistic Eye Protection

NAR IPRO®

Tactical Goggle System

- Fully enclosed, high-impact eye protection compatible with helmet and night-vision goggles
- 2 clear & 2 tinted lenses (neutral gray) with 100% UV shield
- Hexo-vents for fog-free viewing
- Certified to ANSI Z87.1-2003 & EN 166:2001
- Includes cleaning kit & padded zipper case

SPECIFICATIONS

Case: L 2.6 in. x W 7 in. x D 3 in.
Weight: 8 oz

NAR IPRO

ITEM

ZZ-0109



Nylon Accessories

ID PATCH

Rescue Task Force

Nylon Patches with Hook Backing for agency ID on Ballistic Vests

RESCUE TASK FORCE

RESCUE TASK FORCE

Custom patches also available.
Call Customer Service for more details.

RTF ID PATCHES

COLOR	SIZE	ITEM #
Red	3 in. x 11 in.	70-0200
Coyote	3 in. x 11 in.	70-0199
Black	3 in. x 11 in.	70-0198
Red	2 in. x 11 in.	ZZ-0537
Coyote	2 in. x 11 in.	ZZ-0510
Black	2 in. x 11 in.	ZZ-0509

MKIII Side Armor

PLATE POUCHES

PLATE POUCHES

COLOR	ITEM #
Black	70-0175
Coyote	70-0177
Blue	70-0176
Red	70-0178

- Hard armor plate holder facilitates additional side armor protection
- Nylon pouch accommodates up to 6 in. x 8 in. hard armor plates (not included)

SPECIFICATIONS

L 8.5 in. x H 6.5 in. x D 1 in. | Wgt: 2.72 oz



Nylon GEAR BAG

Single compartment bag for maximum capacity - large enough to store essential ballistic gear including vest & helmet



SPECS

L 30 in. x W 14 in.
x D 14 in.
Weight: 1 lb 4 oz

NYLON GEAR BAG

COLOR	ITEM #
Black	80-0547



NAR SOFT LITTER HOLDER

SPECIFICATIONS

H 14 in. x W 14 in. x D 0.25 in.
Weight: 9 oz

LITTER HOLDER

COLOR	ITEM #
Black	80-0556
Blue	80-0557
Red	80-0558

- Carries (2) NAR QuikLitters or (1) NAR QuikLitter XL (each sold separately)
- Allows deployment of litter from the bottom or top of the carrier
- Featuring a MOLLE-style attachment that connects to the rear of protective vest leaving front of vest open for other mission critical equipment



Massive Hemorrhage | Tourniquets

C-A-T® Combat Application Tourniquet Gen 7

- **Single route buckle** decreases application time to help decrease blood loss, allows for slack removal for proper TQ application
- **Applied same way** to both upper and lower extremities
- True **one-handed tourniquet**
- **Scientifically proven** 100% effective in occluding blood flow in both upper and lower extremities by the USAISR and recommended by CoTCCC



online video available

< scan to view

SPECIFICATIONS

Packaged:
H 6.5 in. x W 2.4 in. x D 1.5 in.
Device: L 37.5 in.
Weight: 2.7 oz



**The FASTEST, SAFEST,
Most EFFECTIVE
Prehospital Field Tourniquet!**

COMBAT APPLICATION TOURNIQUET

COLOR	STYLE	ITEM #	NSN #
Black	Tactical	30-0001	6515-01-521-7976
Orange	Rescue	30-0023	6515-01-612-7129
Blue	Trainer	30-0033	6910-01-560-2972

See page 68 for more about the Blue Trainer version of the C-A-T®

TACO KYDEX TQ HOLDERS

- Designed to retain windlass tourniquets on a belt, rig or bag
- Body constructed of aircraft-grade Kydex® a smooth finish for easy insertion and draw (.080 thickness)
- Resistant to moisture and solvents & won't add weight from rain, humidity or sweat
- Will not damage the tourniquet or the attached platform
- Unique shock cord design compatible with most windlass-style tourniquets
- TACO® hole pattern allows for stacking two pouches in the space of one
- Can be mounted to belts up to 2 in. or MOLLE with included two universal clips
- Includes embossed medical cross ID

SPECIFICATIONS

H 4.5 in. x W 2.75 in.
x D 1.5 in.
Weight: 0.21 lb

Textured

Slick

TACO KYDEX TOURNIQUET HOLDERS

COLOR	FINISH	ITEM #
Black	Slick	30-0272
Coyote Brown	Slick	30-0274
LE Blue	Slick	30-0273
Olive Drab	Slick	30-0276
Black Basket Weave	Slick	30-0277
Wolf Gray	Slick	30-0275

COLOR	FINISH	ITEM #
Black	Textured	30-0267
Coyote Brown	Textured	30-0269
LE Blue	Textured	30-0268
Olive Drab	Textured	30-0271
Black Basket Weave	Textured	30-0278
Wolf Gray	Textured	30-0270

TACO KYDEX SHIRT SHIELD

COLOR	ITEM #
Black	30-0279
Coyote Brown	30-0281
LE Blue	30-0280
Olive Drab	30-0283
Wolf Gray	30-0282

TACO KYDEX TQ HOLDER SHIRT SHIELD Accessory

Upgrade your Kydex TACO Kydex TQ Holders to protect your shirt from undue wear caused by the friction from the windlass clip!

SPECIFICATIONS

H 2 in. x W 3.25 in. x D 1 in. | Weight: 0.06 lb

Available in 5 Colors

RIGID C-A-T® TQ CASE *with Shirt Shield*

- Allows you to carry a Gen 7 C-A-T® Tourniquet on your duty, patrol, or pistol belt so you have life-saving equipment on your person at all times
- Includes a Shirt Shield complete with mounting hardware to protect your uniform shirt from undue wear caused by the friction from the windlass clip
- Blade-Tech TekLok Belt Attachment fits up to a 2.25 in. belt
- Can be mounted either vertically or horizontally or attached to a Safariland Drop Leg Shroud
- Injection Molded Nylon Case stands up to rugged LE/MIL usage
- Cross medical identifier patch included

RIGID C-A-T TQ CASE *with Shirt Shield*

COLOR	ITEM #
Black	80-0944
Black, Basket Weave	80-0962
Black, Polished	80-0963
Coyote	80-0945
Ranger Green	80-0946

SPECIFICATIONS

H 5.88 in. x W 2.5 in. x D 2 in.
Weight: 5.36 oz

Cases & Shirt Shields (see below) also available for separate purchase



NOTE: All Tourniquets sold separately

RIGID C-A-T® TQ CASE *(Case Only)*

This rugged Case is also available as a stand alone item, available in multiple color options

SPECIFICATIONS

H 4.75 in. x W 2.25 in. x D 2 in.
Weight: 3.8 oz

RIGID C-A-T TQ CASE

COLOR	ITEM #
Black	80-0479
Black, Basket Weave	80-0560
Black, Polished	80-0562
Blue	80-0480
Coyote	80-0481
Gray	80-0484
Multi-Cam	80-0561
Orange	80-0483
Ranger Green	80-0482



Easily attaches to case with just two screws!



RIGID C-A-T® TQ CASE SHIRT SHIELD *(Shield Only)*

Upgrade your Rigid TQ Case to protect your shirt from undue wear!

SPECIFICATIONS

H 1.75 in. x W 2.5 in. x D 1.5 in.
Weight: 1.44 oz

RIGID SHIRT SHIELD

COLOR	ITEM #
Black	80-0941
Coyote	80-0942
Ranger Green	80-0943

TQ Accessories | Massive Hemorrhage

LEATHER C-A-T® HOLDERS

- Designed to meet the specific duty gear needs of law enforcement officers - for use on daily duty belt to fit within organizational uniform standards
- Securely fastens to your duty belt with easy one-handed access to allow you to take control of extreme medical situations in seconds
- Constructed of black full grain leather in different finishes, attachment options, and snaps colors
- High quality snaps ensure secure placement on your duty belt and easy accessibility to your tourniquet
- Attachment options: Clip, Loop or Hook & Loop Strap attachment
- Ideal for the C-A-T® Gen 7 but also accommodates other popular windlass TQs

SPECIFICATIONS

L 5.75 in. x W 2 in. x D 2 in.
Weight: 3.2 oz

UNIFORM DUTY BELT LEATHER TQ HOLDER

FINISH	ATTACHMENT	SNAP	ITEM #
Basket Weave	Clip	Silver	30-0192
High Gloss	Clip	Silver	30-0194
Full Grain	Clip	Silver	30-0193
Basket Weave	Loop	Silver	30-0195
High Gloss	Loop	Silver	30-0198
Full Grain	Loop	Silver	30-0197
Basket Weave	Strap Attach	Silver	30-0228
High Gloss	Strap Attach	Silver	30-0229
Full Grain	Strap Attach	Silver	30-0196

See Gold & Black snap options online

with Clip ▼

with Strap Attachment ►

with hook and loop allows the holder to be added to a radio strap without removing the radio pouch from the strap

See website for more Uniform Duty Belt Tourniquet Holder options

▲ with Loop

TQ 911 C-A-T® HOLDER

Solid one-piece injection molded design

SPECIFICATIONS

H 4 in. x W 2.25 in.
x D 2.25 in.
Weight: 1.6 oz

TQ 911 HOLDER

ITEM #
30-0160

C-A-T® ELASTIC TQ HOLDER

- PALS/MOLLE-style strap with tuck finish eliminating the need for a button snap
- Easy-open quick release hood
- Extended elastic body for TQ release outwards under pressure
- Accommodates both Gen 6 & 7 C-A-T®s (sold separately)

SPECIFICATIONS

L 7 in. x W 1.75 in. x D 1.5 in.
Weight: 2 oz

ELASTIC TQ HOLDER

COLOR	ITEM #
Black	30-0184
Coyote	30-0185
Multi-Cam	30-0186

Nylon C-A-T® HOLDER

SPECIFICATIONS

L 6.25 in. x W 1.75 in. x D 1.5 in.
Weight: 1.26 oz

NYLON TQ HOLDER

COLOR	ITEM #
Black	30-0057
Coyote	30-0044
Multi-Cam	30-0059
Red	30-0128
Ranger Green	30-0208
RealTree®	30-0207
OD Green	30-0053

Ankle Trauma Holster also available (See page 25)

Inside the Waistband C-A-T® HOLDER

- Inside the waist band comfortable solution for carrying your C-A-T®
- Clips inside your waistband or outside on your pants belt loop or board shorts
- Can MOLLE to your gear

SPECIFICATIONS

L 4 in. x W 2 in. x D 2 in.
Weight: 2.08 oz

INSIDE THE WAISTBAND HOLDER

COLOR	ITEM #
Black	30-0205
Red	30-0206



Cummerbund C-A-T® HOLDER

- Attaches to tactical armor such as assault and plate carriers
- Utilizes 4 in. elastic coupled with hook & loop attachment to external cummerbund flap
- Rapid, easy access with ambidextrous one handed TQ removal
- Compatible with both Gen 6 & Gen 7 C-A-T®s (sold separately)



CUMMERBUND TQ HOLDER

COLOR	ITEM #
Black	30-0212
Coyote	30-0213

SPECIFICATIONS

L 5 in. x W 4 in. x D 0.25 in.
Weight: 0.6 oz

Junctional Tourniquets | Massive Hemorrhage

AAJT-S

Abdominal Aortic
and Junctional
Tourniquet -
STABILIZED



- Easy to use and quick to apply – Lightest, most compact junctional tourniquet on the market
- Stabilizes the pelvis in addition to stopping bleeding
- Large area of effective pressure for greater reliability
- Provides stable and complete occlusion of blood flow to the lower extremities, upper extremities, or pelvis
- Can be applied at junctional sites (inguinal, axillary) or mid-abdomen
- Demonstrated ability to remain in place and effective during patient movement
- Similar hemostatic efficacy to Zone 3 REBOA¹
- Lower training burden compared to REBOA procedure

AAJT-S

ITEM #
30-0248

SPECIFICATIONS

Packaged: L 9 in. x W 8 in. x D 2 in.
Stored: L 7.5 in. x W 6.5 in. x D 2 in.
Weight: 17 oz

JETT® Junctional Emergency Treatment Tool

- Immediately treat high extremity wounds in the thigh/groin area that is not amenable to TQ application
- Unilateral or bilateral treatment solution
- Simultaneous occlusion of blood flow to both lower limbs that will not impede respiration

JETT

ITEM #	NSN #
30-0088	6515-01-616-5841



online video
available
< scan to view

SPECIFICATIONS

Packaged:
H 5.5 in. x W 7 in. x D 3.75 in.
Overall Device: L 56.75 in.
Buttocks Pad:
L 19 in. x W 6.5 in.
Compression Pad:
L 4.85 in. x W 2.5 in. x D 1.1 in.
T-Handle:
H 4.62 in. x W 4.04 in. x D 0.63 in.
Weight: 1 lb 9.6 oz



See NARescue.com for more info about the
JETT Pack: Item# 80-0355 & 80-0356
a nylon pouch designed specifically for carrying the JETT

Massive Hemorrhage | Bandages & Gauze

NAR COMPRESSED GAUZE

- Low-cube, ultra-compact sterile gauze
- Premium rolled 100% cotton gauze for maximum absorption

COMPRESSED GAUZE

STYLE	ITEM #
Responder	30-0182
Standard	30-0052

SPECS

Gauze: W 4.5 in. x L 4.1 yd
Packaged: L 3 in. x W 2.5 in. x D 1 in.
Weight: 12 oz

S-ROLLED GAUZE

- Unique dispenser, allows controlled linear feed application
- Premium 100% cotton

S-ROLLED GAUZE

ITEM #
30-0003

SPECS

Gauze: W 4.5 in. x L 4.1 yd
Packaged: L 3 in. x W 2.5 in. x D 1 in.
Weight: 12 oz

RESPONDER ETD™ Emergency Trauma Dressing

- Sterile pressure bandages that deliver extremely effective pressure & excellent absorption
- Suitable for applying immediate direct pressure
- Equipped with a non-adherent pad
- Easy-to-use hook & loop securing device with no confusing pressure bars or hooks
- Simple application with the Quick grip roll control system which prevents the dressing from unrolling if dropped
- Rugged vacuum-sealed, low-cube packaging with easy-open tear notches

RESPONDER ETDs

STYLE	BANDAGE	PAD	WGT	PACKAGE	ITEM #
4 in. Flat	W 4 in. x L 60 in.	4x6	2 oz	L 4.25 in. x W 2.75 in. x D 0.88 in.	30-0171
4 in. Flat (Mini)	W 4 in. x L 42 in.	4x6	2.4 oz	L 5 in. x W 3.5 in. x D 0.75 in.	30-0179
6 in. Rolled	W 6 in. x L 60 in.	6x6	2.9 oz	L 6 in. x W 2.5 in. x D 2.1 in.	30-0168

NAR WOUND PACKING GAUZE

- Ultra compressed, Z-fold gauze indicated for wound packing
- Soft, white 1-ply gauze with an embedded x-ray detectable line

WOUND PACKING GAUZE

ITEM #
30-0054

SPECS

Gauze: W 3 in. x L 5 yd
Packaged: H 3 in. x W 2 in. x D 1 in.
Weight: 12 oz

ETD™ Emergency Trauma Dressing

- Premium cotton construction for max quality & absorption
- Durable C-clasp securing device

ETD EMERGENCY TRAUMA DRESSINGS

STYLE	BANDAGE	PAD	WGT.	PACKAGE	ITEM #
4 in. Std.	W 4 in. x L 70 in.	4x6	2.6 oz	L 4.5 in. x W 2.5 in. x D 1.75 in.	30-0031
4 in. Flat	W 4 in. x L 70 in.	4x6	2.6 oz	L 4.25 in. x W 2.75 in. x D 0.88 in.	30-0093
6 in. Std.	W 6 in. x L 70 in.	6x6	4 oz	L 6 in. x W 2.5 in. x D 2.1 in.	30-0032
6 in. Flat	W 6 in. x L 70 in.	6x6	2.6 oz	L 6 in. x W 3.5 in. x D 0.88 in.	30-0094
ABD/ Stump	W 6 in. x L 70 in.	12x12	5 oz	L 5.25 in. x W 6 in. x D 1.5 in.	30-0012
8 x 10 in.	W 4 in. x L 70 in.	8x10	4 oz	L 4.125 in. x W 7 in. x D 1 in.	30-0165

CRAVAT (Triangular Bandage)

Bandage: 37 in. x 37 in. x 52 in.
Packaged: 3.5 in. x 2 in. x 0.63 in.
Weight: 16 oz

CRAVAT

ITEM #
30-0089

Tactical COMPRESSION WRAP

Wrap: L 38 in. x W 5 in.
Packaged: 5.4 in. x 2.7 in. x 2 in.
Weight: 6 oz

TACTICAL COMPRESSION WRAP

ITEM #
30-0016

COMBAT GAUZE® | COMBAT GAUZE® LE

- Sterile, non-woven Z-fold gauze impregnated with kaolin for temporary external control of bleeding
- Effective, simple and safe to use on wounds
- Features an embedded x-ray detectable strip
- Does not generate heat

SPECIFICATIONS

Standard & LE Gauze:

L 4 yd x W 3 in.
Pkg: H 5 in. x W 5.5 in. x D 0.625 in.
Weight: 0.8 oz

XL Gauze:

L 4 yd x W 4 in. (2 ply)
Pkg: H 6.5 in. x W 6.5 in. x D 0.75 in.
Weight: 1.75 oz



See page 74
for more info
on the Trainer Version



COMBAT GAUZE Z-FOLD

VERSION	ITEM #
Standard	30-0039
XL	30-0070
LE	30-0109
Trainer	30-0063

QuikClot® BLEEDING CONTROL DRESSING

- Advanced hemostatic agent combined with the simplicity of gauze available without medical device authorization (MDA)
- Utilizes the same inorganic component as QuikClot Combat Gauze® used by the US Military to control severe bleeding



SPECIFICATIONS

Z-fold Gauze: L 4 yd x W 3 in.
Pkg: H 5 in. x W 5.5 in. x D 0.625 in.
Weight: 0.8 oz

Roll Gauze: L 4 ft x W 3 in.
Pkg: H 5 in. x W 3.5 in. x D 0.5 in.
Weight: 0.5 oz

QUIKCLOT BCD

STYLE	ITEM #
Z-Fold	30-0161
Roll	30-0164



SPECS

Standard Gauze: L 48 in. x W 3 in.
Pkg: H 5.75 in. x W 5 in. x D 0.185 in.
Weight: 0.36 oz

XL Gauze: L 4 yd x W 3 in.
Pkg: H 5.75 in. x W 5 in. x D 0.45 in.
Weight: 0.68 oz

CHITOGAUZE® OTC

- Hemostatic chitosan-impregnated gauze for rapid control of minor bleeding
- Z-folded in sterile package for ease of application
- Safe and effective
- Flexible & conformable to tissue surfaces

CHITOGAUZE OTC

SIZE	ITEM #
Standard	30-0250
XL	30-0251

CHITOGAUZE® XR Pro | XR2 Pro

- Non-woven, polyester/rayon blend Z-fold medical gauze coated with chitosan to control severe bleeding
- Includes embedded x-ray detectable lines

SPECIFICATIONS

Gauze: L 4 yd x W 3 in.
Packaged: H 5.75 in. x W 5 in. x D 0.65 in.
Weight: 0.7 oz

CHITOGAUZE XR Pro

STYLE	ITEM #
XR	30-0049
XR2	30-0265



XR2:
Two X-ray lines



XR: One
X-ray line



CELOX™ RAPID

- Stops severe arterial bleeding with only one minute of compression applied
- Compact Z-folded format making it easy to use & carry
- Generates no heat

CELOX RAPID

ITEM #
30-0069

SPECIFICATIONS

Gauze: W 3 in. x 5 ft
Packaged: H 5.25 in. x W 4.25 in. x D 0.5 in.
Weight: 11 oz



BAC-PACK™ Bougie-Aided Cricothyroidotomy Pack

- All-in-one cricothyroidotomy kit for use in any environment
- Pre-inserted bougie in tube allows for a more rapid, easy procedure
- Tube meets TCCC recommendations for bougie-aided open surgical technique
- PenBlade #10 Safety Scalpel included for safety and ease of use
- Small, durable pouch keeps the items organized and ready to use
- Lightweight and ruggedized
- Components configured by order of use
- PenBlade safety scalpel, tube, syringe all provided sterile
- Patent pending

INCLUDED IN THE BAC-PACK:



SPECIFICATIONS

Packaged: L 5.5 in. x H 3.625 in. x D 0.5 in.
 Inner tube diameter: 6.0mm
 Outer tube diameter: 8.0mm
 Tube length: 67.0mm
 Weight: 0.84 oz
 US Patents: Pending

BAC-PACK

ITEM

10-0052

NAR Tracheostomy Kit with BOUGIE INTRODUCER

- Convenient device that gives the user a pre-inserted bougie for introduction of the tube into the trachea and indication of correct placement prior to insertion
- Adjustable for maximum user preference
- Latex free tube features 6.0mm inner diameter and inflatable cuff
- Contains the same preferred flexible tube as NAR's Tactical CricKit®

TRACHEOSTOMY KIT

ITEM

10-0048

TRACHEAL TUBE with FLANGE

SPECIFICATIONS

Tube diameter:
 Inner: 6.0mm | Outer: 8.0mm
 Tube length: 67.0mm
 Packaged:
 L 5.5 in. x H 3.625 in. x D 0.5 in.
 Weight: 0.84 oz

INFLATABLE
DISTAL CUFF

PILOT
BALLOON

PRE-INSERTED
BOUGIE INTRODUCER
with DEPTH GAUGE

a Insert
bougie

b Advance tube
& inflate cuff

c Remove
bougie

TACTICAL CRICKIT® Tactical Surgical Cricothyroidotomy Kit

- Contains all features required by CoTCCC for a surgical airway kit in the tactical setting
- Ideal length & size tracheal tube with securing device
- Configured in the sequence of use, ideal for high-stress, low-light, austere environments
- Lightweight, ruggedized and minimal in cube size
- Included rolled pouch provides a clean working surface

INCLUDED IN THE TACTICAL CRICKIT:



TACTICAL CRICKIT

ITEM

10-0017

SPECIFICATIONS

Packaged: H 6.5 in. x W 2.5 in. x D 125 in.

Unrolled: H 6.5 in. x W 7.5 in. x D 125 in.

Weight: 3.7 oz

NASOPHARYNGEAL AIRWAY with lubricant

Packaged with surgical lubricant attached

L 6 in. x W 16 in. x D 0.5 in.
Weight: 0.5 oz

NASOPHARYNGEAL AIRWAY

SIZE	ITEM #
28F	ZZ-0034

NAR TACTICAL SUCTION DEVICE

Suction Capacity: 1000 mL
Vacuum Force: 100 mmHg
Weight: 7.3 oz
L 9.75 in. x W 3.75 in. x D 3 in.

TACTICAL SUCTION DEVICE

ITEM #
10-0018

CYCLONE® Improved POCKET BVM



AMBU® SPUR II BVM



L 11.63 in. x W 5 in. x H 5 in.
Weight: 11.2 oz

AMBU® SPUR II BVM

ITEM #
10-0043

NAR TRACHEAL HOOK



L 4.2 in. x W 0.5 in. x D 0.8 in.
Weight: 0.5 oz

NAR TRACHEAL HOOK

ITEM #
10-0020

- Unique collapsible design packaged 25% smaller than conventional BVMs
- Requires minimum pressure for compression of the bag to ventilate

Diameter: 5.3 in. x H 2.9 in.
Weight: 17.4 oz

CYCLONE® POCKET BVM

ITEM #
10-0019

Airway Management

KING LTS-D™ Supraglottic Airway Kit



- Unique dual channel design, anatomical shaped distal tip/cuff & open second lumen
- Tapered ramp design provides additional ventilation outlet

SPECIFICATIONS

Folded: H 8 in. x W 11.43 in. x D 1.92 in.
Unfolded: H 8 in. x W 13.36 in. x D 1.22 in.
Weight: 5.7 oz

INCLUDED WITH KING LTS-D™



KING LTS-D

CLR	SIZE	ITEM #	NSN #
YEL	3	10-0025	----
RED	4	10-0026	6515-01-591-1578
PUR	5	10-0027	6515-01-591-0397

AMBU® RDIC Mark III Resuscitator



AMBU RDIC RESUSCITATOR

ITEM #	NSN #
10-0047	6515-01-338-6602

SPECIFICATIONS

L 14 in. x W 7.5 in. x H 7 in.
Weight: 5.4 lb

Special Operations

LARYNGOSCOPE SET




- Perform direct visual laryngoscopy even in the most severe environments
- Lightweight, compact kit to equip medical care providers in the field with tools to secure an airway
- Contains 2 single-use fiber optic blades, textured handle & removable light source all rolled up in a rugged nylon pouch

SPECIFICATIONS

Open: L 8 in. x W 5.5 in. x D 2 in.
Rolled: L 8 in. x W 2 in. x D 2 in.
Weight: 10 oz

INCLUDED IN THE LARYNGOSCOPE SET:



LARYNGOSCOPE SET

ITEM #	NSN #
10-0001	6515-01-450-9790

NAR Fiber Optic LARYNGOSCOPE BLADES

Disposable, safe and convenient

LARYNGOSCOPE BLADES

AVAILABLE SIZES	ITEM #	NSN# available on all sizes offered
Miller #1, #2, #3 & #4	See NARescue.com	See NARescue.com for details
Mac #2, #3 & #4	See NARescue.com	See NARescue.com for details



HYFIN® CHEST SEALS



HYFIN® VENT Twin Pack

- Two individually packaged, vented/valve-type chest seals in one unique package designed for treatment of both entry/exit wounds or multiple penetrating injuries
- Patented 3 channel pressure relief vents
- Ruggedized dual perforated pouch designed to allow the provider to open only one dressing at a time, as needed
- Advanced adhesive technology designed to provide a superior seal



HYFIN VENT TWIN PACK

ITEM #	NSN #
10-0037	6510-01-642-6212

SPECIFICATIONS

Device (excludes tab): 6 in. x 6 in. Wgt: 2.5 oz
 Packaged: H 7.5 in. x W 4.5 in. x D 0.25 in.
 Unfolded: H 7.5 in. x W 9 in. x D 0.13 in.



online video
available

< scan to view

HYFIN® VENT Individual

HYFIN VENT INDIVIDUAL

ITEM #	NSN #
10-0029	6510-01-624-0840

SPECIFICATIONS

Device (excludes tab): 6 in. x 6 in. Wgt: 1.31 oz
 Packaged: H 4.63 in. x W 7.5 in. x D 0.1 in.



HYFIN® Occlusive Individual

HYFIN OCCLUSIVE INDIVIDUAL

ITEM #	NSN #
10-0015	6510-01-532-8019

SPECIFICATIONS

Device (excludes tab) 6 in. x 6 in Wgt: 1.31 oz
 Packaged: H 4.63 in. x W 7.5 in. x D 0.1 in.



HYFIN® VENT COMPACT, Twin Pack

- Two Chest Seals packaged approximately 25% smaller than our standard HyFin®
- 3-channel pressure relief vents
- Two Chest Seals designed for the treatment of entry/exit or multiple penetrating wounds
- Large red tab designed for single step, peel and apply application and for burping the wound if necessary

SPECIFICATIONS


Device (excludes tab):
 4.75 in. x 4.75 in.
 Wgt: 1.55 oz
 Packaged:
 H 3.88 in x W 6.6 in x D 0.13 in.

HYFIN VENT COMPACT

ITEM #
10-0042

Also available at
www.NARescue.com

HyFin® Twin Pack Trainer
 Item# 10-0032
 See page 75 for more information



10 ga x 3.75 inch

SPEAR™ Simplified Pneumothorax Emergency Air Release

- Latest innovation in Needle Decompression technology
- Designed based on known gaps in the management of tension pneumothorax
- Enhanced needle tip for improved insertion control
- User friendly features like spin-lock connection, depth markings, and detachable low-pressure one-way valve
- Fenestrations along catheter for additional air flow
- Familiar, rugged, pen-style case for ease of storage and deployment
- Allows for either lateral decompression or anterior needle thoracostomy
- Catheter made of flexible kink-resistant material with a soft tip
- Single use only. Supplied sterile



online video available

< scan to view

SPECIFICATIONS

Packaged:
L 7.625 in. x Diameter 0.75 in.
Catheter: 10Ga x 3.75 in.
Weight: 1 oz
US Patents: Pending

SPEAR

GAUGE	ITEM #	NSN #
10 ga	10-0051	---



**DETACHABLE
ONE-WAY VALVE**

SPIN LOCK


**10 gauge 3.75 in.
FLEXIBLE CATHETER with
DEPTH GAUGE**

DISTAL FENESTRATIONS


SOFT SAFETY TIP

**NEEDLE with
MODIFIED TIP**

*Versatility that allows for
Lateral Decompression as well as
Anterior Needle Thoracostomy...*



Lateral Approach



Anterior Approach

3.25 inch **ARS® NEEDLE DECOMPRESSION KIT**

- Penetrate the pleural space in 99% of casualties presenting the signs and symptoms of a tension pneumothorax
- Quality, reliable devices with the appropriate length needle/catheter
- Capless flash chamber provides immediate confirmation of correct placement
- Patented design packaged in an easy-to-open, crush-proof protective tube featuring easy-to-identify textured twist top with clip



ARS NEEDLE DECOMPRESSION KIT



GAUGE	ITEM #	NSN #
14 ga	ZZ-0056	6515-01-541-0635
10 ga	ZZ-0298	6515-01-673-1701

U.S. Design Patents: D584,409 S; D595,847S; D584,410S; U.S. Patent 7,874,426B2; © North American Rescue, LLC., 127820, 127822, 127823, 001002372-0001, 001013940-0001, 001013940-0002

SPECIFICATIONS

Pkg: H 5.7 in. x W 0.7 in. x D 0.8 in.
Weight: 0.6 oz



online video available

< scan to view

NAR SIMPLE THORACOSTOMY KIT

KIT CONTENTS:

- 1 x Rochester Pean Curved Hemostat, 6.25 in.
- 1 x PenBlade Scalpel #10
- 2 x Gauze Pads 4 in. x 4 in. (2 Pk)
- 1 x HyFin® Vent Compact Twin Pack
- 1 x Chloraprep® Swabsticks (3 Pk)
- 1 x Permanent Marker, Small - Black



SIMPLE THORACOSTOMY KIT



ITEM #	UNIT	NSN#
83-0027	Single	-----

SPECIFICATIONS

Packaged: L 4.5 in. x W 7.5 in. x D 0.75 in.
Weight: 5.28 oz

CHEST TUBE INSERTION KIT

KIT CONTENTS:

- 1x HyFin® Chest Seal Dressing
- 1x Petrolatum Gauze
- 1x Scalpel, #10 - Protected
- 1x Hemostat, Kelly Curved - 5.5"
- 1x Drape, Sterile 18 X 25 Fenestrated
- 1pr Glove, 8.5 Sterile Surgical Latex
- 1x Chest Drainage Valve
- 1x Chest Tube, 36 Fr Soft
- 1x Chloraprep Swabsticks, 3 Per Pkg
- 1x 4 X 4 Split Sponge
- 1x O Silk Fsl, 30 Inch Suture
- 1x 8" Cvd Roch-Peon, Sterile



CHEST TUBE INSERTION KIT



ITEM #	UNIT	NSN#
83-0007	Single	-----

SPECIFICATIONS

Packaged: L 10 in. x W 6 in. x D 2.25 in.
Weight: 11.2 oz

Circulation & Vascular Access

NAR Needleless SALINE LOCK KIT

IV safety catheter & needle-free valve port to reduce or eliminate accidental needle-stick exposures

NEEDLELESS SALINE LOCK KIT

ITEM #	NSN #
30-0046	6515-01-607-1943



SPECIFICATIONS


Pkg: H 5.5 in. x W 3.25 in. x D 1 in.
Weight: 1.3 oz



NEEDLELESS SALINE LOCK KIT CONTENTS:



NAR SALINE LOCK KIT *Meets all requirements for US Army Task: 081-831-1011*



Essential components for establishing a saline lock for safer, quicker and more effective IV cannulation in the most demanding environments

SALINE LOCK KIT

ITEM #	NSN #
30-0013	6515-01-537-4094

SPECIFICATIONS

Pkg: H 5.5 in. x W 3.25 in. x D 1 in.
Weight: 1.3 oz



SALINE LOCK KIT CONTENTS:



BOA® IV Constricting Band

- Direct peripheral IV access utilizing a revolutionary rolling action for transfusing venous blood to the distal veins while still allowing arterial flow
- Quick Click connector & two-finger Quick Release improves "time on task"
- XL version stretches approximately 50% longer than the standard BOA®



SPECIFICATIONS

Standard: L 7.5 in. to 18 in.
XL: L 10 in. to 27 in.
Packaged: L 4 in. x W 3 in. x D 0.9 in.
Weight: 15 oz
US Patents: 5,304,202 & 5,607,448

BOA IV CONSTRICTING BAND

SIZE	ITEM #	NSN #
Standard	30-0009	6515-01-537-2611
XL	30-0071	----

SHARPS SHUTTLE

SPECIFICATIONS

L 6.38 in. x Diameter 1.11 in.
Weight: 15 oz

SHARPS SHUTTLE

UNIT	ITEM #	NSN #
Single	30-0014	6515-01-518-8551
6 pack	30-00146	6515-01-619-5671



RAPTOR® IV Securing Device

SPECIFICATIONS

Extended: L 19 in. x W 1 in.
Weight: 19 oz

RAPTOR IV BAND

UNIT	ITEM #	NSN #
6 pack	30-0006	6515-01-523-3317



SAM[®] IO Intraosseous Access System

The **SAM IO system** is designed to be intuitive for current IO users, featuring color coding and depth indicators which match industry norms. The system components for effective IO placement are sold separately

SAM IO DRIVER

A multi-use manually operated intraosseous access driver

Passive ratcheting mechanism rotates clockwise when resistance is encountered.
Provides controlled insertion of the IO needle with a tactile feel of cortical penetration.

SAM IO NEEDLE SETS

A stylet and catheter assembly secured by a snap-fit connection. 15 mm, 25 mm, and 45 mm lengths are intended to address a broad range of patient demographics and multiple application sites

15 mm

25 mm

45 mm

EXTENSION TUBE*

Connects to catheter for flexible infusion

*Extension Tube & NeedleVISE® included in each needle set

NEEDLEVISE®*


Single-use sharps containment for stylet

INSERTION SITES:

Proximal Humerus (left or right)
Proximal Tibia (left or right)
Distal Tibia (left or right)
Distal Femur (pediatric patients)

SAM STABILIZER

Cylinder design intended to protect the catheter from dislodgement



SAM IO NEEDLE ADAPTOR

The SAM IO Adaptor is designed to accommodate alternative IO needle designs with the SAM IO Driver

REV022420

R_X ONLY

SAM
MEDICAL

SAM IO INTRAOSSEOUS ACCESS SYSTEM

PRODUCT	UNIT	ITEM #
Driver	Single	30-0255
Needle 15mm	Single	30-0257
Needle 15mm	5-Pack	30-0258
Needle 25mm	Single	30-0259
Needle 25mm	5-Pack	30-0260
Needle 45mm	Single	30-0261
Needle 45mm	5-Pack	30-0262
Stabilizer	Single	30-0263
Stabilizer	5-Pack	30-0264
Adaptor	Single	30-0253
Training Kit	Kit	30-0256
Storage Case	Single	30-0254


SAM IO FIELD STORAGE CASE

Custom field storage case designed for the SAM IO with elastic loops for storing needles and the SAM IO Driver



LIFEWARMER™ Quantum™ Blood & Fluid Warming System

- **Lightweight** The Quantum™ has the lightest operational weight of any portable blood & fluid warming system (Total Operational Weight: 22 oz)
- **Powerful** Warms two units of 39°F refrigerated whole blood at 100mL/min+ on a single charge
- **Compact** Fewest components to carry (three)
- **Fast** Warm-up time of only 24 seconds
- **Intuitive** The thermal transfusion/infusion set looks, handles, and sets up/primes like any conventional IV set
- **Safe** Tested safe for blood, blood products & IV solutions
- **Effective** The system is plug-and-play for ease of setup under stressful conditions such as tactical field care (only two connection points, with no plug orientation or set-up order requirements – system status constantly visible at a glance)
- **Complete** Does not require additional standard transfusion or infusion set – all other warmers do
- **Tactical** Only system FDA cleared allowing easy disable of audible alerts to protect operators in combat



OPERATIONAL CAPABILITIES

QUANTUM COMPLETE SYSTEM

ITEM #
35-0006

BLOOD or BLOOD PRODUCTS

From a **single battery charge**, capable of warming **2 units** of **39°F Whole Blood** to an output temperature of **100.4°F (38°/- 2°C)**. Warm up time: **24 seconds**. Delivery Flow rate: **100mL/min**.


COMPLETE QUANTUM SYSTEM CONTENTS:

1 x Quantum Charger

1 x Quantum Battery

1 x Quantum Controller

2 x TTS-B Thermal Transfusion Set - Blood
(Blood, blood products & IV solutions)



SPECIFICATIONS

Total Operational Weight: 22 oz
(Quantum™ Controller, Battery & TTS-B)


Battery: 16 oz

Controller: 3 oz

Battery Charger: 22 oz

Controller Cover: 0.23 oz

Thermal Transfusion Set (TTS-B): 2.7 oz



IV FLUIDS

From a **single battery charge**, capable of warming **1700 mL** of **68°F IV fluids** to an output temperature of **100.4°F (38°/- 2°C)**. Warm up time: **24 seconds**. Delivery Flow rate: **200mL/min**.


REV021220

TTS-B Thermal Transfusion Set

- Compatible with whole blood, component therapy (including FDP), colloids, and isotonics
- Similar size and weight to traditional IV tubing
- Intelligent sensing and heating elements integrated into the walls of the tubing set
- No disruption to the fluid path or laminar flow of the fluid
- Easy bi-directional connection to the Quantum Controller

SPECIFICATIONS

TTS-B package (each): 6.5 in. x 11 in.
Tubing Length: 80 in.
TTS-B Weight (each): 2.6 oz
Packaged Weight (2 TTS-B & box): 11 oz



TTS-B (2 Pack)

ITEM

35-0005

TIS Thermal Infusion Set

- Similar size and weight to traditional IV tubing
- Intelligent sensing and heating elements integrated into the walls of the tubing set
- No disruption to the fluid path or laminar flow of the fluid
- Easy bi-directional connection to the Quantum Controller
- Compatible with standard IV fluids such as colloids and isotonics

SPECIFICATIONS

TIS packaged (each): 6.5 in. x 11 in.
Tubing Length: 80 in.
Weight: TIS (each): 2.1 oz
Packaged (2 TIS & box): 10 oz

TIS (2 Pack)

ITEM

35-0004



Quantum BATTERY

- Lightweight at only 16 oz
- Designed for the tactical environment- water resistant, subdued color, and the option to disable the audible alarm
- 44V Lithium-Polymer Battery
- Intelligent, state-of-the-art battery technology
- Status button to check charge state at any time

QUANTUM BATTERY

ITEM

35-0001

SPECIFICATIONS

H 4.5 in. x W 3.25 in. x D 1.25 in.
Weight: 16 oz



Quantum HARDCASE

SPECS

L 10.63 in.
x W 9.63 in.
x D 6.73 in.
Weight: 3.22 lb



- Rugged hard case designed to protect your Quantum from dust, dirt, and moisture
- Features pressure relief valve for automatic pressure equalization
- Custom laser cut high density foam holds components securely in place
- Flexible storage capacity for multiple batteries, controllers and thermal tubing sets

QUANTUM HARDCASE

ITEM

35-0008

Quantum CONTROLLER

- Moderates the energy flow from the Quantum Battery to the Thermal Tubing
- Designed for ease of connecting and durability
- Colored LED lights constantly show during administration of blood/fluid to show status

Replacement cover is also available



- Extremely lightweight and durable
- IP53 rated for dust and water resistance with jacket
- Minimum service life - 1000 insertions/removals

SPECIFICATIONS

Controller Body: W 1.25 in. x H 3 in. x D 0.75 in.
Controller Total Length: 51.25 in.
Weight (Controller only): 2.5 oz
Packaged Weight: 5 oz

QUANTUM CONTROLLER

OPTION	ITEM #
Controller	35-0003
Cover (Replacement)	35-0010

- Fully charge the Quantum Battery in <90 minutes
- Smart, safe charging
- Includes 6 ft AC power cord
- Colored LED lights to inform the user of charging and battery status
- IP 22 rated for dust and water resistance
- 100V to 264V AC input; 50Hz, 60Hz
- Minimum service life- 2500 insertions/removals

SPECIFICATIONS


Charger: H 6 in. x W 2.25 in. x D 1.75 in.
AC power cord length: 6 feet
Charger Weight: 1 lb 4 oz
Weight with cord: 1 lb 10 oz

Quantum CHARGER

QUANTUM CHARGER

ITEM

35-0002



Hypothermia Prevention & Management

HPMK[®] Hypothermia Prevention & Management Kit

- Includes reinforced Heat Reflective Shell (HRS[™]) of 4-ply, composite fabric with a protected non-conductive reflective layer
- Built-in hood, integrated fluid absorption pad & top-to-bottom taper providing excellent thermal insulation
- Continuous 1.5 in. Velcro[®] closures provide rapid 360° access to the casualty with minimal exposure to the elements
- Self-heating, 4-cell liner (Ready-Heat[™] Blanket) provides 10 hours of continuous dry heat up to 125°F (52°C)

SPECIFICATIONS

Shell Packaged:
H 6.75 in. x W 10.5 in. x D 5.5 in.
Shell Open: L 78 in. x W 43 in.
Weight: 3 lb 8 oz

US Patent: 7,766,950 B2

HPMK[®] CONTENTS:



HYPOTHERMIA PREVENTION & MANAGEMENT KIT - HPMK

ITEM #	NSN #
80-0027	6515-01-532-8056



NOTE: Litter is sold separately

HRS[®] Heat Reflective Shell



As an individual unit, reinforced 4-ply, composite shell with built-in hood, weight tested to support 250+ lb.

HEAT REFLECTIVE SHELL

ITEM #	NSN #
ZZ-0051	6532-01-627-7888

SPECIFICATIONS

Shell Packaged: L 10 in. x W 6.75 in. x D 2.5 in.
Open: L 78 in. x W 43 in. US Patent: 7,766,950 B2
Weight: 1 lb 13 oz

NAR SURVIVAL BLANKET

Aluminized, reflective PET to help retain up to 90% of your body heat

SURVIVAL BLANKET

ITEM #
ZZ-0038

SPECS

Unfolded:
L 84 in. x W 52 in..
Weight: 17 oz



READY HEAT[™] BLANKET



As an individual unit, the self-heating, 4-cell liner provides 10 hours of continuous dry heat after reaching peak temperature approximately 15-20 minutes after opening

READY HEAT BLANKET

ITEM #
ZZ-0052

SPECIFICATIONS

Open: L 48 in. x W 36 in.
Weight: 1 lb 11 oz

Emergency SURVIVAL WRAP

Reflects heat back to the body

SURVIVAL WRAP

ITEM #
ZZ-0834

SPECIFICATIONS

Unfolded: L 90 in. X W 60 in,
Weight: 3.7 oz



ANKLE TRAUMA KIT

- Very lightweight, low profile, versatile carrier that allows for carry of essential trauma equipment
- Designed to accommodate a C-A-T® Tourniquet, Emergency Trauma Dressing and Gauze/Hemostatic
- Advanced breathable materials that provide user comfort and long-term wearability
- Protective, water resistant pouches for each carried component with strong hook and loop closures
- Quick, single-handed, ambidextrous access to components
- One size fits most with convenient hook and loop adjustability

ANKLE TRAUMA KIT

KIT/HOLSTER OPTION	ITEM #
Kit with Combat Gauze LE	80-0870
Kit with S-Rolled Gauze	80-0869
Ankle Holster Only	30-0177T

ANKLE TRAUMA with S-ROLLED:



with COMBAT GAUZE LE:



SPECIFICATIONS

L 17.65 in. x H 6.5 in. x D 1.5 in.
Weight: 8.25 oz

BELT TRAUMA KIT (BTK)

- Rugged black pouch fits a belt up to 2.25 in. and can be worn rotated 90° to make room for other gear
- Mesh pockets allow quick and easy contents identification

SPECIFICATIONS

BASIC: H 4.5 in. x W 6 in. x D 2.5 in

Weight: 6.6 oz


COMPLETE: H 4.5 in. x W 6.75 in. x D 3.75 in.

Weight: 10.5 oz


BELT TRAUMA KIT

LEVEL	ITEM#
Basic	80-0426
Complete	80-0427

BASIC BELT TRAUMA KIT:



COMPLETE BELT TRAUMA KIT:



IFAK Individual First Aid Kits

M-FAK® Mini First Aid Kit - LE

- Contains the essential point-of-wounding equipment needed for treating penetrating, blast and other traumatic injuries
- Super compact, rugged nylon platform that allows attachment both vertically (MOLLE backing) and horizontally (3 in. belt loop)
- Clam shell configuration utilizes (2) two main sleeves that open on both ends for easy access

M-FAK MINI FIRST AID KIT - LE

COLOR	LEVEL	ITEM #
Black	Basic	80-0494
Coyote	Basic	80-0495
Gray	Basic	80-0872
Multi-Cam	Basic	80-0496
OD Green	Basic	80-0497
Ranger Green	Basic	80-0996
Black	Combat Gauze LE	80-0498
Coyote	Combat Gauze LE	80-0499
Gray	Combat Gauze LE	80-0873
Multi-Cam	Combat Gauze LE	80-0500
OD Green	Combat Gauze LE	80-0501
Ranger Green	Combat Gauze LE	80-0998
Black	Advanced	80-0502
Coyote	Advanced	80-0503
Gray	Advanced	80-0874
Multi-Cam	Advanced	80-0504
OD Green	Advanced	80-0505
Ranger Green	Advanced	80-0997



NEW!
Ranger Green

See page 100
for the compatible
2 in. Medic Cross Patch
(sold separately)




SPECIFICATIONS

H 6 in. x W 3 in. x D 3 in.
Weight: 13 oz

LEM-FAK | ROO M-FAK BASIC or BASIC with Combat Gauze:



LEM-FAK | ROO M-FAK ADVANCED:



Bag only option also
available at NARescue.com

Individual First Aid Kits IFAK

ROO M-FAK *Mini First Aid Kit*

Rapid one handed access to your tourniquet without ever having to open the full kit!

Contains the same packout options as the M-FAK:
Basic, Basic with Combat Gauze LE and **Advanced**
(see page 28)

ROO M-FAK MINI FIRST AID KIT

COLOR	LEVEL	ITEM #
Black	Basic	80-1049
Coyote	Basic	80-1050
OD Green	Basic	80-1051
Ranger Green	Basic	80-1052
Black	Combat Gauze LE	80-1053
Coyote	Combat Gauze LE	80-1054
OD Green	Combat Gauze LE	80-1055
Ranger Green	Combat Gauze LE	80-1056
Black	Advanced	80-1057
Coyote	Advanced	80-1058
OD Green	Advanced	80-1059
Ranger Green	Advanced	80-1060

Bag only option also
available at NARescue.com



NEW!
ROO bag
includes a
Mini 1 in.
PVC ID
Patch

Sportsman M-FAK[®] *Mini First Aid Kit*

- The same super compact, M-FAK nylon platform designed with the Sportsman in mind
- Essential point-of-wounding medical gear for treating penetrating, blast and other traumatic injuries in the field



RealTree[®] ▶
Sportsman also includes
2 in. PVC MED Cross Patch

M-FAK MINI FIRST AID KIT - SPORTSMAN

COLOR	LEVEL	ITEM #
RealTree	Basic	80-0913
RealTree	with BCD	80-0914
Red	Basic	80-0890
Red	with BCD	80-0891
Red	Advanced	80-0892

SPECIFICATIONS

H 6 in. x W 3 in. x D 3 in.
Weight: 13 oz

Sportsman M-FAK BASIC or BASIC with BCD:



Sportsman M-FAK ADVANCED:



IFAK Individual First Aid Kits

EAGLE IFAK

- Can connect to the assault vest/body armor or to a duty/pistol belt
- Contains the essential equipment needed for self-aid/buddy-aid to address the leading causes of potentially preventable death in a tactical environment
- Quick Pull Tab easily opens the kit with just one hand
- Unique, detachable shock-cord inner panel that can be adjusted to your components
- Inner panel is removable for easy access to your kit

SPECIFICATIONS

H 8 in. x W 4 in. x D 3 in.
Weight: 1 lb 2 oz

EAGLE IFAK

COLOR	LEVEL	HEMOSTATIC	ITEM #
Black	Basic	None	85-0177
Coyote	Basic	None	85-0178
Multi-Cam	Basic	None	85-0179
OD Green	Basic	None	85-0180
Black	Advanced	Combat Gauze	85-0189
Coyote	Advanced	Combat Gauze	85-0190
Multi-Cam	Advanced	Combat Gauze	85-0191
OD Green	Advanced	Combat Gauze	85-0192
4 options online	Advanced	ChitoGauze	See NARescue.com

Bag only option also available at NARescue.com

EAGLE IFAK BASIC:



EAGLE IFAK ADVANCED:



RIG SERIES EAGLE IFAK®

Bag only option
also available at
NARescue.com

- Tourniquet attached for quick deployment by banding it to inside of front opening flap
- Protective hood with quick pull latch release helps protect IFAK components and includes sewn in embroidered medical cross identifier
- Detachable (hook & loop backed) shock-cord inner panel adjusts to secure internal components and is easily removable for rapid access

RIG SERIES EAGLE IFAK

COLOR	LEVEL	ITEM #
Black	Basic	80-0509
Coyote	Basic	80-0510
Multi-Cam	Basic	80-0511
OD Green	Basic	80-0512
Ranger Green	Basic	80-0991
Black	Combat Gauze LE	80-0517
Coyote	Combat Gauze LE	80-0518
Multi-Cam	Combat Gauze LE	80-0519
OD Green	Combat Gauze LE	80-0520
Ranger Green	Combat Gauze LE	80-0993
Black	Advanced	80-0513
Coyote	Advanced	80-0514
Multi-Cam	Advanced	80-0515
OD Green	Advanced	80-0516
Ranger Green	Advanced	80-0992

SPECIFICATIONS

H 8 in. x W 4 in. x D 3 in.
Kit Packed Weight: 1.18 lb
Empty Bag Weight: 7 oz


RIG SERIES EAGLE BASIC:



RIG SERIES EAGLE with COMBAT GAUZE LE:



RIG SERIES EAGLE ADVANCED:



OCHO IFAK

- Compact, low visibility, point-of-wounding IFAK designed for your duty/pistol belt
- Purposely designed to look like a non-tactical, accessory bag
- Equipped with a 3 in. belt loop attachment to accommodate a duty/pistol belt
- Constructed with rip stop nylon, elastic loop organizational storage and easy open zipper

SPECIFICATIONS

L 7 in. x H 2.25 in. x D 3.25 in.

Weight: 13 oz

OCHO IFAK

COLOR	LEVEL	ITEM #
Black/Gray	Basic	80-0374
Blue/Gray	Basic	80-0375
Red/Gray	Basic	80-0915
Black/Gray	Combat Gauze	80-0921
Blue/Gray	Combat Gauze	80-0922
Red/Gray	Combat Gauze	80-0917
Black/Gray	Advanced	80-0372
Blue/Gray	Advanced	80-0373
Red/Gray	Advanced	80-0916

See OCHO Accessory Bag only option on page 50



OCHO IFAK BASIC:



OCHO IFAK with COMBAT GAUZE:



OCHO IFAK ADVANCED:



Individual First Aid Kits IFAK

Individual Patrol Officer Kit IPOK



- Individual hemorrhage control kit
- Compact, vacuum-sealed & packaged to fit into a pocket, vest pouch or individual bag
- Rugged, durable pouch features NAR's signature Red Tip Technology® ensuring quality & performance

SPECIFICATIONS


H 4 in. x W 6.5 in. x D 2.75 in.
Weight: 7.4 oz



IPOK - INDIVIDUAL PATROL OFFICER KIT

HEMOSTATIC/GAUZE	ITEM#
S-Rolled Gauze	80-0167
Combat Gauze, LE	80-0168
ChitoGauze	80-0223

IPOK:



IPOK HOLDER

- Protects your medical kit from the elements
- Enclosure can be secured by Fastex snap buckle, hook-and-loop or both

IPOK HOLDER

COLOR	ITEM#
Black	80-0205
Coyote	80-0207
Digital Camo	80-0208
OD Green	80-0254
Red	80-0358



SPECIFICATIONS

H 6.5 in. x W 6 in. x D 3.5 in.
Weight: 4.1 oz

VISOR MED POUCH

- Securely mounts to a vehicle visor
- Orange webbing handle to grab for easy deployment
- Pouch with easy-open buckle



ORDERING INFO

COLOR	ITEM#
Black	80-0269
Coyote	80-0270
Red	80-0360

SPECIFICATIONS

H 6.5 in. x W 5.5 in. x D 3 in.
Weight: 5.4 oz

Compact Officer Response Emergency Kit CORE KIT

- Compact design with the essential life-saving equipment needed to address moderate and severe traumatic bleeding
- Contains the #1 CoTCCC Recommended Prehospital Tourniquet and Hemostatic Gauze



SPECIFICATIONS

H 7 in. x W 3.5 in. x D 2.25 in.
Weight: 6.4 oz

CORE KIT

ITEM #
85-0834

CORE KIT:



IFAK Individual First Aid Kits

T.O.R.K. TACTICAL OPERATOR RESPONSE KIT

COLOR	LEVEL/HEMOSTATIC	ITEM #
Black	Basic	85-0008
Coyote	Basic	85-0009
Digital Camo	Basic	85-2300
Multi-Cam	Basic	85-0012
OD Green	Basic	85-0054
Red	Basic	85-1805
Ranger Green	Basic	85-2301
Black	Basic w/Combat Gauze	80-1097
Coyote	Basic w/Combat Gauze	80-1098
Digital Camo	Basic w/Combat Gauze	80-1099
Multi-Cam	Basic w/Combat Gauze	80-1100
OD Green	Basic w/Combat Gauze	80-1101
Ranger Green	Basic w/Combat Gauze	80-1102
Red	Basic w/Combat Gauze	80-1108
Black	Advanced w/Combat Gauze	85-0048
Coyote	Advanced w/Combat Gauze	85-0049
Digital Camo	Advanced w/Combat Gauze	85-0050
Multi-Cam	Advanced w/Combat Gauze	85-0052
OD Green	Advanced w/Combat Gauze	85-0053
Ranger Green	Advanced w/Combat Gauze	85-2302
Red	Advanced w/Combat Gauze	85-1806
7 options online	Advanced w/ChitoGauze	NARescue.com

T.O.R.K. Tactical Operator Response Kit

- Streamlined, compact IFAK designed to MOLLE to the assault vest or body armor
- Includes essential life-saving equipment to be used at the point of wounding
- For self-aid/buddy-aid injuries caused by penetrating or blast trauma



Bag only options
also available at
NARescue.com

T.O.R.K. BASIC/BASIC with Combat Gauze:



SPECIFICATIONS

H 7 in. x W 4.5 in. x D 2.75 in.
Weight: 1 lb 5 oz

T.O.R.K. ADVANCED:



Individual First Aid Kits IFAK

T.O.R.K. Tactical Operator Response Kit - Sportsman





RealTree® Sportsman Includes a PVC MED Cross Patch

The same compact TORK design with modified contents and RealTree® nylon bag with PVC Medic Cross Patch to suit the needs of sportsmen addressing traumatic injuries at the point of wounding.

REALTREE / SPORTSMAN TORK

COLOR	LEVEL	ITEM #
RealTree	Basic	80-1006
RealTree	with Bleeding Control Dressing	80-1007

Sportsman T.O.R.K.:



SPECIFICATIONS

H 7 in. x W 4.5 in. x D 2.75 in.
Weight: 1 lb 5 oz


NAR-4 CHEST POUCH

- Streamlined pouch that provides quick access to essential equipment from the chest
- Designed to be secured to the front of the body armor using a MOLLE-style attachment system
- Ambidextrous trauma shear accessibility

SPECIFICATIONS

H 7.75 in. x W 9.75 in. x D 4.5 in
Weight: 2 lb 10.3 oz

NAR-4 CHEST POUCH:



Bag only option also available at NARescue.com

NAR-4 CHEST POUCH

COLOR	HEMOSTATIC	ITEM #
Black	Combat Gauze	80-0173
Coyote	Combat Gauze	80-0175
Digital Camo	Combat Gauze	80-0176
Multi-Camo	Combat Gauze	80-0247
OD Green	Combat Gauze	80-0289
5 options online	ChitoGauze	See NARescue.com

IFAK *Individual First Aid Kits*

OPERATOR BLS IFAK *Basic Life Support*

- Essential life-saving equipment that can be used for self-aid/buddy-aid at the point of wounding
- Attaches directly to the duty belt as a drop down leg rig or to the assault vest or body armor with MOLLE/PALS-style connectors

OPERATOR BLS IFAK

COLOR	ITEM #
Black	80-0117
OD Green	80-0118

SPECIFICATIONS

H 8.5 in. x W 8.75 in. x D 4.5 in.
Weight: 1 lb 12 oz

OPERATOR BLS IFAK:



Bag only option
also available at
NARescue.com

S.T.O.R.M. OPERATOR IFAK *Specialized Tactics for Operational Rescue & Medicine*

- Provides the necessary equipment to administer immediate self-aid/buddy-aid to casualties in hostile environments
- Attaches directly to the duty belt as a drop down leg rig or to the assault vest/body armor with MOLLE/PALS-style connectors

S.T.O.R.M. OPERATOR IFAK

COLOR	ITEM #
Black	80-0133
Coyote	80-0135
Digital Camo	80-0136
OD Green	80-0134

SPECIFICATIONS

H 7.4 in. x W 6.9 in. x D 5 in.
Weight: 1 lb 12 oz

S.T.O.R.M. OPERATOR IFAK:



Bag only option
also available at
NARescue.com

Multiple Injury Trauma Kits

High Risk Warrant

CASUALTY KIT

- Designed with the capability of both casualty extraction & medical treatment in a high threat situation
- Provides the capability to address the leading causes of preventable death in the tactical environment.
- Equipped with a Dragon Handle System to perform rapid assault rescue
- Multiple carriage capabilities (fanny pack, over the shoulder or attachment to packs with MOLLE-style webbing)

SPECIFICATIONS

H 7.75 in. x W 13 in. x D 9 in.
Weight: 4 lb 12.8 oz



HIGH RISK WARRANT CASUALTY KIT:



HIGH RISK WARRANT

COLOR	ITEM #
Black	80-0121
OD Green	80-0122

Bag only option also available at NARescue.com

EMS RAPID DEPLOYMENT KIT

- Highly visible, compact kit for pre-hospital healthcare providers
- Provides necessary equipment to treat several casualties

SPECIFICATIONS

H 7.75 in. x W 13 in. x D 9 in.
Weight: 4 lb 6.5 oz

EMS RAPID DEPLOYMENT

COLOR	ITEM #
Orange	80-0151



EMS RAPID DEPLOYMENT KIT:



Bag only option also available at NARescue.com

Multiple Injury Trauma Kits

Patrol **VEHICLE TRAUMA KIT**

- Provides personnel with an easily accessible, quick to deploy trauma kit from their vehicle
- Modular storage system comes organized with an array of medical products ranging from basic first aid to major trauma care
- Equipment is packed in the kit to offer basic first aid and minor injury equipment on one side and major trauma equipment on the other
- Exterior includes MOLLE-style system for extra equipment as needed



SPECIFICATIONS

H 7 in. x W 12 in. x D 5 in.
Weight: 6 lb

PATROL VEHICLE TRAUMA KIT:



PATROL VEHICLE TRAUMA

LEVEL	ITEM #
Basic	80-0542
Combat Gauze LE	80-0859

Tactical **RAPID DEPLOYMENT KIT**

- Medium size go-bag with necessary equipment to treat critical casualties
- Kit and its contents are compact, durable and allow reliable and versatile interventions wherever the need arises

SPECIFICATIONS

H 7.75 in. x W 13 in. x D 9 in.
Weight: 4 lb 12.8 oz

TACTICAL RAPID DEPLOYMENT KIT:



Bag only option also available at NARescue.com

TACTICAL RAPID DEPLOYMENT

COLOR	ITEM #
Black	80-0201
Coyote	80-0203
Digital Camo	80-0204
OD Green	80-0202

RUF KIT

- Slim profile that can be worn with single bandolier style sling strap and with a detachable waist strap
- Two tear away internal panels designed to hold standard IFAK contents
- Two additional tear away clear accessory pouches for additional equipment
- Internal shock cord area to hold compact, soft litter or additional equipment
- External MOLLE on sides of bag for external tourniquet pouches
- Large panel on exterior of kit for ID patches

RUF KIT


LEVEL	ITEM #
Basic	80-0974
Combat Gauze LE	80-0975
Advanced	80-0976

Bag only option
also available at
NARescue.com

SPECIFICATIONS

H 18 in. x W 8 in. x D 2 in.
Weight: 6.3 lb

Every RUF KIT includes these core items:




PLUS+ BASIC

also includes:





COMBAT GAUZE LE

also includes:



ADVANCED

also includes:



Can be worn with single bandolier style sling strap as well as with a detachable waist strap



MINI MEDIC KIT

- Compact, multi-purpose trauma aid bag with vital medical equipment for point of wounding treatment of multiple injuries
- Two external pockets with two external tourniquet holders
- One external pouch with additional drop-down pouch for use with breaching tools or storage for a hydration pack
- Expansion zipper to increase load capacity
- External MOLLE on sides and bottom of bag / Internal MOLLE rows on inside back of bag
- Three clear tear away pouches on the inside main compartment
- Two internal mesh pouches and one larger internal pouch on inner lid
- Internal tape roll strap


Every MINI MEDIC KIT includes these core items:




PLUS+ BASIC
also includes:



COMBAT GAUZE
also includes:



ADVANCED also includes:



Inside Back Compartment for stowing extra equipment

SPECIFICATIONS

H 14 in. x W 11 in. x D 5 in.
(Depth expandable to 7.25 in.)
Weight: 12.5 lb

*Bag only option
also available at
NARescue.com*

MINI MEDIC KIT

COLOR	LEVEL	ITEM #
Black	Basic	80-0979
Coyote	Basic	80-0980
OD Green	Basic	80-0981
Ranger Green	Basic	80-1003
Black	Combat Gauze	80-0985
Coyote	Combat Gauze	80-0986
OD Green	Combat Gauze	80-0987
Ranger Green	Combat Gauze	80-1005
Black	Advanced	80-0982
Coyote	Advanced	80-0983
OD Green	Advanced	80-0984
Ranger Green	Advanced	80-1004

RAPID RESPONSE KIT LE Edition

- Modular, personalized organization of mission essential gear and medical equipment
- Includes 2 C-A-T® tourniquets with holders, double mag pouch and nylon pouch to store additional equipment

RAPID RESPONSE - LE

COLOR	ITEM #
Black	85-0715

SPECIFICATIONS

H 12 in. x W 12 in. x D 6 in.
Shoulder strap: 28 in. to 50 in.
Weight: 3 lb 12.5 oz

RAPID RESPONSE KIT - LE EDITION:



Gun Magazines are not included in mag pouch

RAPID RESPONSE KIT RTF Edition

- Equipped to provide initial triage of a mass casualty incident and provide treatment for 4 or more patients
- Includes 4 C-A-T® tourniquets with holders and nylon dispenser with 4 rolls of color-coded triage tape

RAPID RESPONSE -RTF

COLOR	ITEM #
Black	85-0716
Red	85-0717

SPECIFICATIONS

H 12 in. x W 12 in. x D 6 in.
Shoulder strap: 28 in. to 50 in.
Weight: 6 lb 4.5 oz

RAPID RESPONSE KIT - RTF EDITION:



Red Rapid Response RTF Kit comes standard with orange C-A-T® Tourniquets, red holders and a red tape dispenser

Triage Tape, Set of 4
(Item# ZZ-0591) also available



Multiple Injury Trauma Kits

THROW KIT

INDIVIDUAL KIT WITH COMBAT GAUZE

KIT CONTENTS:

- 1x C-A-T® Tourniquet, Black
- 1x Combat Gauze with Hemostatic Bandage
- 4 in. Flat Emergency Trauma Dressing
- HyFin Vent Chest Seal, Compact, Twin Pack
- 1x Blue Responder Gloves, Large

PN: 80-1105

INDIVIDUAL THROW KITS

- Compact vacuum-sealed pouch with essential point of wounding equipment provided for rapid deployment to on-scene responders during a Mass Casualty Incident
- Designed to treat bleeding from penetrating and other traumatic injuries
- Small cube space fits into a BDU pocket, vest pouch, or individual bag

SPECIFICATIONS

H 7 in. x W 5 in. x D 2.5 in.
Weight: 8 oz

INDIVIDUAL THROW KITS

OPTIONS	ITEM #
Basic	80-1104
with Combat Gauze LE	80-1105

HASTY HARNESS

Multi-purpose Harness

Resourceful, multi-purpose harness for unlimited applications & maximum versatility... a must for every webbing equipment toolbox

SPECIFICATIONS

Loop: L 10 ft 11 in. x W 1 in.
Circumference - 21 ft 10 in.
Min. Break Str: Tech Tape (4,400 lb)
Sewn Sling (6,000 lb)

HASTY HARNESS

COLOR	ITEM #
Red	90-0044
Black	90-0001

INDIVIDUAL THROW KITS:



MCI RAPID RESPONSE THROW KITS

Hands Free, Point of Wounding Solution with the enhanced ability of rapid distribution of medical gear to treat a large number of injured patients in a mass casualty situation...

- Rapid response sling bag designed to take the guesswork out of the management of a critical incident by providing the rapid distribution of necessary life saving equipment over a large scene
- Includes **6 Individual Throw Kits** (see contents above) with essential life saving bleeding control equipment designed for rapid deployment to on-scene responders
- Nylon sling bag features exterior MOLLE-style attachments so the user can customize the bag according to their needs
- Provides ample storage for extra gear such as a Rescue Harness for casualty evacuation (sold separately)

SPECIFICATIONS

H 12 in. x W 12 in. x D 6 in.
Shoulder strap: 28 in. to 50 in.
Weight: 6 lb 4.5 oz

MCI RAPID RESPONSE THROW KITS

OPTIONS	ITEM #
Basic	80-1106
with Combat Gauze LE	80-1107



Multiple Injury Trauma Kits

Mass Casualty Incident **MCI WALK**

- All-in-one trauma treatment & evacuation platform for treatment of multiple trauma patients
- Choose from two different litter configurations: 8 non-rigid, disposable QuikLitters or 1 full-size, quad-fold, lightweight MedEVAC Litter
- Pocket system holds self- or buddy-aid point-of-wounding medical supplies
- Bag features back pack straps for easy transport from the vehicle to the point of need
- Ideal for use as a vehicle trauma kit

(80-0565/80-0458) Red MCI WALKs contain orange C-A-T® Tourniquets

(#80-0458) Red MCI WALK with MedEVAC Litter contains an orange MedEVAC Litter

MCI WALK

COLOR	LITTER OPTION	C-A-T	ITEM #
Black	8 QuikLitters	Black	80-0564
Red	8 QuikLitters	Orange	80-0565
Black	MedEVAC Litter (Black)	Black	80-0457
Red	MedEVAC Litter (Orange)	Orange	80-0458

SPECIFICATIONS

H 22.75 in. x W 16.5 in. x D 11.5 in.

Weight with 8 QuikLitters: 18 lb 12 oz

Weight with MedEVAC Litter: 29 lb 15 oz

MCI WALK with QUIKLITTERS:

 QTY 8 C-A-T® Tourniquet	 QTY 8 ETD™ 6 in. Dressing	 QTY 5 Combat Gauze® LE	 QTY 8 Gauze, Compressed	 QTY 5 28F Naso Airway & Lube
 QTY 8 HyFin® Vent. Twin Pk	 QTY 8 NAR QuikLitter	 PAIR 8 Nitrile Gloves, Lg	 QTY 20 Combat Casualty Card	 QTY 5 Permanent Marker, Large

MCI WALK with MEDEVAC LITTER:

 QTY 8 C-A-T® Tourniquet	 QTY 8 ETD™ 6 in. Dressing	 QTY 5 Combat Gauze® LE	 QTY 8 Gauze, Compressed	 QTY 5 28F Naso Airway & Lube
 QTY 8 HyFin® Vent. Twin Pk	 QTY 1 Rigid Litter MedEVAC	 PAIR 8 Nitrile Gloves, Lg	 QTY 20 Combat Casualty Card	 QTY 5 Permanent Marker, Large



Bag only option also available
at NARescue.com



online video
available

< scan to view

Multiple Injury Trauma Kits

MEDIC TRAUMA PACK

- 9 removable, color-coded pull bags and 2 removable side compartments provide modular organization
- Over 1,000 cubic inches of gear capacity and MOLLE-style attachment system compatible with removable side pouches

SPECIFICATIONS

H 19 in. x W 19 in. x D 10 in

Weight: 18 lb 1 oz

MEDIC TRAUMA PACK

**NATO NSN

COLOR	ITEM #	NSN #
Black	80-0070	6545-01-572-7251
Coyote	80-0072	6545-01-588-1862
Digital Camo	80-0096	----
OD Green	80-0071	6545-01-594-5617**



MEDIC TRAUMA PACK:

 QTY 2 C-A-T® Tourniquet	 QTY 6 ETD™ 6 in. Dressing	 QTY 2 ETD™ Abdominal	 QTY 4 Gauze S-Rolled	 QTY 2 28F Naso Airway & Lube
 QTY 1 King LT-D Size 4	 QTY 1 Cyclone Pocket BVM	 QTY 1 Tactical Suction Device	 QTY 2 Gauze, Petrolatum	 QTY 2 ARS® Needle, 14g
 QTY 1 BOA® IV Band	 QTY 4 Saline Lock Kit	 QTY 1 Sharps Shuttle	 QTY 1 HRS, Heat Reflective Shell	 QTY 1 Tactical Traction Splint
 QTY 4 SAM® Splint II	 QTY 1 Trauma Shears, Lg	 QTY 2 Surgical Tape, 2 in.	 QTY 1 Armadillo Med Case	 QTY 1 Head Lamp
 PAIR 5 Nitrile Gloves, Lg	 QTY 2 PPE Kit	 QTY 6 PES Eye Shield	 QTY 4 T2® Tag Triage Card	 QTY 1 Casualty Reference Card



Bag only option
also available at
NARescue.com

NAR-4 AID KIT

- Low-profile, load carriage for Tactical Healthcare Professionals
- Modular system with 5 segregated removable equipment pouches to carry equipment consistent with TCCC Guidelines
- Designed to attach directly to the body armor with four rapid attachment straps or carried using shoulder straps

SPECIFICATIONS

H 19 in. x W 15 in. x D 7 in.
Weight: 12 lb 8.3 oz

NAR-4 AID KIT

COLOR	HEMOSTATIC	ITEM #
Black	Combat Gauze	80-0181
Coyote	Combat Gauze	80-0183
Digital Camo	Combat Gauze	80-0184
Multi-Cam	Combat Gauze	80-0244
OD Green	Combat Gauze	80-0308
5 options online	ChitoGauze	See NARescue.com

NAR-4 AID KIT:

 QTY 2 C-A-T® Tourniquet	 QTY 2 ETD™ 6 in. Dressing	 QTY 1 ETD™ Abdominal	 QTY 4 Gauze S-Rolled	 QTY 2 Choice of Hemostatic
 QTY 1 Tactical Compression Wrap	 QTY 2 28F Naso Airway & Lube	 QTY 1 King LTS-D	 QTY 1 Cyclone Pocket BVM	 QTY 1 Tactical Suction Device
 QTY 1 Tactical CricKit®	 QTY 1 HyFin® Vent. Twin Pk	 QTY 2 ARS® Needle, 14g	 QTY 1 BOA® IV Band	 QTY 4 Saline Lock Kit
 QTY 1 Sharps Shuttle	 QTY 1 Tactical Traction Splint	 QTY 2 SAM® Splint II	 QTY 1 Trauma Shears, Lg	 QTY 2 Tape, Mini 2" Multi-Purpose
 QTY 1 Armadillo Med Case	 PAIR 8 Nitrile Gloves, Lg	 QTY 6 PES Eye Shield	 QTY 4 Combat Casualty Card	 QTY 1 Casualty Reference Card
			 QTY 1 Permanent Marker, Large	

Bag only option
also available at
NARescue.com

Multiple Injury Trauma Kits

WALK Warrior Aid & Litter Kit

- Official vehicle life-saving kit of the U.S. Army
- For treatment of multiple trauma casualties and the immediate evacuation of one patient
- Includes full-size Talon® 90C quad-fold litter
- Pocket system holds self or buddy-aid medical supplies
- Back pack straps provide easy transport from the vehicle to the point of need
- Includes Hypothermia Prevention & Management Kit™ to help stabilize patients during transport
- Interfaces with any standard or non-standard TACEVAC asset

















WARRIOR AID & LITTER KIT

COLOR	ITEM #	NSN #
Black	80-0097	----
Coyote	80-0089	6545-01-587-1199
Digital Camo	80-0059	6545-01-532-4962
OD Green	80-0060	----

SPECIFICATIONS

H 22.75 in. x W 16.5 in. x D 11.5 in.
Weight: 29 lb 15 oz

WARRIOR AID & LITTER KIT:

 QTY 2 C-A-T® Tourniquet	 QTY 6 ETD™ 6 in. Dressing	 QTY 1 ETD™ Abdominal	 QTY 4 Gauze S-Rolled	 QTY 2 28F Naso Airway & Lube
 QTY 2 Hyfin® Occlusive, Sgl.	 QTY 2 ARS® Needle, 14g	 QTY 1 HPMK® Kit	 QTY 2 SAM® Splint II	 QTY 1 Talon 90C Litter
 QTY 4 Litter Ret. Strap	 QTY 1 Casualty Equip. Bag	 QTY 1 Aviation Panel	 QTY 1 Trauma Shears, Lg	 PAIR 5 Nitrile Gloves, Lg
 QTY 1 Tape, 2 in. Gecko-Grip	 QTY 1 Surgical Tape, 2 in.	 QTY 6 PES Eye Shield	 QTY 2 Combat Casualty Card	 QTY 1 Casualty Reference Card

WALK ReSupply Kit (Item# 85-0020):

Contains WALK components for resupply - excluding nylon bag & litter



Bag only option
also available at
NARescue.com

Range Trauma Kits

RANGE TRAUMA KIT

- Essential equipment for medical and non-medical personnel to address injuries common to firearm ranges and training facilities
- Treatment capabilities include penetrating injuries, chest wounds, burns, eye injuries, minor lacerations, sprains and fractures

RANGE TRAUMA KIT

*NATO NSN#

HEMOSTATIC	COLOR	ITEM #	NSN #
Combat Gauze	Orange	80-0299	----
ChitoGauze	Orange	80-0213	6545-01-610-5411*

SPECIFICATIONS

H 7.75 in. x W 13 in. x D 9 in
Weight: 5 lb 10 oz



RANGE TRAUMA KIT:



Bag only option also available at NARescue.com



RANGE TRAUMA AID KIT Hard Case

- Water-tight case maintains the integrity of your medical products when you need them most
- Kit contents arranged in a modular system featuring rugged lid storage organizer and removable nylon pouches
- Labeled internal pouches for ease of use

RANGE TRAUMA AID

COLOR	ITEM #
Orange	85-0889

SPECIFICATIONS

H 9.5 in. x W 13.5 in. x D 6 in.
Weight: 8 lb 12 oz



RANGE TRAUMA AID KIT HARD CASE:



Canine Trauma Kits

K-9 HANDLER ROO M-FAK

- One of NAR's most compact, versatile Mini First Aid Kits featuring rapid one handed access to your tourniquet from an easy-to-open front pouch
- Contains the First Responder's most requested critical point-of-wounding medical equipment for treating penetrating, blast or other traumatic injuries in the line of duty
- Super compact, rugged nylon platform that allows attachment both vertically (MOLLE backing) and horizontally (3 in. belt loop)
- Multiple elastic loops for secure gear organization
- Includes 1 in. Mini PVC Med Cross ID Patch

K-9 HANDLER ROO M-FAK

COLOR	LEVEL	ITEM #
Black	Basic	80-1067
Coyote	Basic	80-1068
OD Green	Basic	80-1069
Ranger Green	Basic	80-1070
Black	Combat Gauze LE	80-1071
Coyote	Combat Gauze LE	80-1072
OD Green	Combat Gauze LE	80-1073
Ranger Green	Combat Gauze LE	80-1074

Bag only option also available at NARescue.com

SPECIFICATIONS

L 7 in. x W 3.75 in. x D 4 in.
Weight: 13 oz



Rapid one handed access to your tourniquet without ever having to open the full kit!

K-9 HANDLER ROO M-FAK:



NAR K-9 LITTER & Carrying Strap Accessories

SPECIFICATIONS

Litter: W 35 in. x H 48 in.
Adjustable Carry Strap:
L 9 to 54 in. x W 2 in.
Carabiner open gate
clearance: 27 mm
Packaged:
W 5.5 in. x H 10 in. x D 2.5 in.
Weight: 1.38 lb

- Lightweight, compact litter for the movement of an injured canine
- Litter can be carried in cargo pocket or pouch
- Leg cut outs to safely carry the canine in the litter
- Six large and four small carry handles that allow for one to four person carries
- Optional K-9 Carry Strap & 2 Carabiners for over the shoulder carry of injured K-9 also available

K-9 LITTER

OPTIONS	ITEM#
K-9 Litter Only	60-0074W
Carry Strap & Carabiners	ZZ-1004
K-9 Litter with Carry Strap & Carabiners	60-0076



K-9 HANDLER IFAK

K-9 HANDLER IFAK:



SPECIFICATIONS

H 7 in. x W 4 in. x D 2.75 in
Weight: 1 lb 3 oz

Bag only option also available at NARescue.com

- Designed for immediate lifesaving treatment to a wounded handler, operator or K-9
- Treatment capabilities range from major trauma typical to operating in tactical or patrol environments to K-9 emergencies

K-9 HANDLER IFAK

COLOR	HEMOSTATIC	ITEM #
Black	----	80-0439
Coyote	----	80-0440
Multi-Cam	----	80-0441
OD Green	----	80-0442
Black	Combat Gauze	80-0443
Coyote	Combat Gauze	80-0444
Multi-Cam	Combat Gauze	80-0445
OD Green	Combat Gauze	80-0446

K-9 TACTICAL FIELD KIT



- Designed for the handler or medic to deliver lifesaving treatment to a wounded K-9 at the point of injury
- K-9 treatment capabilities from major trauma typical to operating in tactical or patrol environments to more common K-9 emergencies
- Provides the essential equipment to render care to the injured K9 when veterinary care is not immediately available

SPECIFICATIONS


H 7.75 in. x W 13 in. x D 9 in.
Weight: 5 lb 2 oz

Bag only option also available at NARescue.com

K-9 TACTICAL FIELD KIT

COLOR	HEMOSTATIC	ITEM #
Black	Combat Gauze	80-0300
Coyote	Combat Gauze	80-0301
OD Green	Combat Gauze	80-0303
3 options online	ChitoGauze	See NARescue.com

K-9 TACTICAL FIELD KIT:



K-9 MINI MEDIC KIT

- Compact, multi-purpose trauma aid bag with life-saving trauma equipment and K9 specific gear for point of wounding treatment of multiple traumatic injuries
- K9 specific gear enables handlers to render care to an injured canine when veterinary care is not immediately available
- Versatile platform features two external pockets with two external tourniquet holders, one external pouch with additional drop-down pouch for extra storage
- External MOLLE on sides and bottom of bag, internal MOLLE rows on inside back of bag, three clear tear away pouches | two internal mesh pouches on the inside main compartment, and one larger internal pouch on inner lid
- Internal tape roll strap also included

Inside Back
Compartment for
stowing extra
equipment

K-9 MINI MEDIC KIT:

 QTY 4 C-A-T® Tourniquet	 QTY 2 C-A-T® TQ Holder	 QTY 2 Combat Gauze Z-fold	 QTY 2 Gauze S-Rolled
 QTY 4 ETD™ 4 in. Flat	 QTY 2 28F Naso Airway & Lube	 QTY 2 HyFin® Vent. Twin Pk	 QTY 4 Gauze, Petrolatum
 QTY 2 ARS® Needle, 14g	 QTY 2 Survival Blanket, Org	 QTY 1 Tactical Muzzle	 QTY 1 NAR K-9 Litter
 QTY 4 Green ChemLight™	 QTY 2 (6 Pk) IR Square	 QTY 1 Head Lamp	 QTY 10 Flexible Fabric Bandages
 QTY 1 Surgical Tape, 2 in.	 QTY 1 Trauma Shears, Lg	 QTY 1 Veterinary Thermometer	 PAIR 5 Nitrile Gloves, Lg
 QTY 2 Permanent Marker, Large	 QTY 1 Medical Cross Patch		

SPECIFICATIONS

H 14 in. x W 11 in. x D 5 in.
(Depth expandable to 7.25 in.)
Weight: 12.5 lb

K-9 MINI MEDIC KIT

COLOR	ITEM #
Black	80-1075
Coyote	80-1076
OD Green	80-1077
Ranger Green	80-1078

Bag only option also available at NARescue.com

K-9 RUF KIT

- Slim profile trauma aid pack designed for the K9 Handler or Tracking Team member
- Carries a minimum of two complete sets of IFAK contents along with K9 medical supplies
- Carriage options include single bandolier style sling strap as well as detachable waist strap
- Two tear away internal panels hold standard IFAK contents
- Two additional tear away clear accessory pouches to store additional equipment
- Internal shock cord area to hold K9 Litter
- External MOLLE on sides of bag for external tourniquet pouches
- Large exterior loop panel for ID patches


Bag only option also available at NARescue.com

K-9 RUF KIT:



K-9 RUF KIT

OPTION	ITEM #
S-Rolled Gauze	80-1027
with Combat Gauze LE	80-1026



SPECIFICATIONS

H 18 in. x W 8 in. x D 2 in.
Weight: 6.3 lb



Can be worn with single bandolier style sling strap as well as with a detachable waist strap

CRISIS INCIDENT RESPONSE KIT

- Designed to treat penetrating trauma, major bleeding, rescue breathing and basic first aid on multiple injured people
- Rugged carrying bag constructed of durable nylon with multiple carry options

SPECIFICATIONS

H 7.75 in. x W 13 in. x D 9 in.
Weight: 5 lb

CRISIS INCIDENT RESPONSE

COLOR	ITEM #
Black	85-1338
Orange	85-0410

CRISIS INCIDENT RESPONSE KIT:



Black Kit (#85-1338)
contains Black TQs
Orange Kit (#85-0410)
contains Orange TQs

Bag only option also available at NARescue.com

INDIVIDUAL AID KIT

- Designed for non-medical trained personnel to help control bleeding until medically trained responders arrive on scene
- Vacuum packed & easy-to-open

INDIVIDUAL AID KIT:



SPECIFICATIONS

H 5.75 in. x W 3.5 in. x D 2 in.
Weight: 6.6 oz

INDIVIDUAL AID KIT

ITEM

85-0404

MASS CRISIS INCIDENT KIT

- Designed to store & transport Individual Aid Kits to classrooms, cafeterias, auditoriums or other areas where a large mass of injured persons are gathered
- Contains 40 Individual Aid Kits (see above for Individual Aid Kits contents) packaged in a medium-sized red duffle-style bag

MASS CRISIS INCIDENT

ITEM

85-0411

SPECIFICATIONS

H 10 in. x W 20 in. x D 10 in.
Weight: 16.5 lb



SRO INDIVIDUAL RESPONSE KIT

- Designed to treat penetrating trauma and severe bleeding for immediate care to injuries
- Attaches to the duty belt as a drop down leg rig
- Includes an extraction device for casualty evacuation

SPECIFICATIONS


H 8.5 in. x W 8.75 in. x D 4.5 in.
Weight: 2 lb 4.8 oz

SRO INDIVIDUAL RESPONSE KIT

ITEM

85-0163

SRO INDIVIDUAL RESPONSE KIT:



Bag only option also available at NARescue.com

SRO CRISIS RESPONSE KIT

- Equip School Resource Officers with the essential life-saving tools to begin treating injured victims until medical personnel arrive
- Empowers an SRO (School Resource Officer) to administer self-aid in the event he or she is injured or provide immediate care to injured casualties during a school incident
- Multiple carry options

SRO CRISIS RESPONSE KIT

ITEM

85-0417


SPECIFICATIONS

H 7.75 in. x W 13 in. x D 9 in.
Weight: 5 lb 16 oz

SRO CRISIS RESPONSE KIT:



Bag only option
also available at NARescue.com



Bag only option
also available at NARescue.com

Nylon Bags

OCHO Lo-Visibility Nylon Bags

- Low visibility aid bags purposely designed to resemble standard off-the-shelf bags to minimize the threat of theft and maximize low visibility discretion
- Featuring hidden compartments with interchangeable pouches backed with hook and loop fields
- Compatible with CrossBreed Holsters (not included)

OCHO ACCESSORY BAG

SPECIFICATIONS

H 7 in. x W 2.25 in. x D 3.25 in.
Weight: 2.75 oz

COLOR	ITEM #
Blue/Gray	80-0371
Gray/Black	80-0370
Red/Gray	80-0571

OCHO Accessory Bag does NOT include Mesh Aid Pouches



OCHO BACKPACK



COLOR	ITEM #
Blue/Gray	80-0256
Red/Gray	80-0338
Green/Gray	80-0339
Gray/Black	80-0340

SPECIFICATIONS

H 19 in. x W 14 in. x D 9 in.
Weight: 3 lb 14 oz

OCHO MESH AID POUCH

SPECIFICATIONS

H 7 in. x W 4.5 in. x D 3 in.
Weight: 1.87 oz

COLOR	ITEM #
Black/ODG	80-0259



OCHO SMALL LOOP PANEL

SPECIFICATIONS

H 7 in. x W 5.5 in. x D 0.25 in.
Weight: 1.83 oz

COLOR	ITEM #
Black	80-0261



OCHO SHOCK CORD PANEL

SPECIFICATIONS

H 6 in. x W 3 in. x D 0.25 in.
Weight: 1.16 oz

COLOR	ITEM #
Black	80-0260

Gear sold separately



MEDIC TRAUMA SHEET

- Versatile, rugged, lightweight equipment access platform for managing casualties during TACEVAC operations
- Oversized pockets & large grommet attachment points facilitate rapid set-up and equipment retrieval
- Quick attach and release straps allow for secure attachment to the inside wall of an aircraft or vehicle

MEDIC TRAUMA SHEET

COLOR	ITEM #	NSN #
Black	80-0067	6515-01-540-7612
Coyote	80-0069	6515-01-540-7603
Digital Camo	80-0094	6515-01-540-7579
OD Green	80-0068	6515-01-540-7585

SPECIFICATIONS

Open: H 41.5 in. x W 21.5 in. x D 2.5 in.
Folded: H 20.25 in. x W 21.5 in. x D 5 in.
Weight: 4 lb 12 oz

Bag only; contents are sold separately



Horizontal & Vertical Mounting capability available with the 4.5 x 6.5 in. Panel



NAR QUICK DETACH PANELS

- Allows the operator the fastest access to their kit during the critical moments after wounding
- Allows any Molle/PALS kit to be detachable
- Gives user faster access to Molle/PALS equipped kits
- Mil-Spec hook & loop and heavy duty snaps to ensure kit security
- Designed to be compatible with all of NAR IFAK's, including a **NEW horizontal mount option** in the medium version (W 4.5 x H 6.5 in.)

Small QUICK DETACH PANEL (3 in. x 6 in.)

- Utilizes only 3 rows by 2 columns on user vest/bag, vertical mount configuration
- Designed specifically for NAR's most compact M-FAK kits

SPECIFICATIONS

W 3 in. x H 6 in. x D 0.5 in.
Weight: 2.8 oz

SMALL PANEL

SIZE	COLOR	ITEM #
3 x 6 in.	Black	80-0893
3 x 6 in.	Coyote	80-0894



Medium DUAL DIRECTION QUICK DETACH PANEL (4.5 in. x 6.5 in.)


- Utilizes 3 rows by 3 columns vertical mount or 2 rows by 4 columns horizontal mount
- Designed for compatibility with all NAR IFAKs such as horizontal mount for M-FAK and vertical or horizontal mount of the TORK

SPECIFICATIONS

W 4.5 in. x H 6.5 in. x D 0.5 in.
Weight: 4.5 oz

MEDIUM DUAL DIRECTION PANEL

SIZE	COLOR	ITEM #
4.5 x 6.5 in.	Black	80-0895
4.5 x 6.5 in.	Coyote	80-0896



Large QUICK DETACH PANEL (9 in. x 6.5 in.)

- Utilizes 3 rows by 7 columns on user vest/bag, vertical mount configuration
- Designed specifically for use with NAR 4 Chest Pouch & Emergency Evacuation Litter (EEL)

SPECIFICATIONS

W 9 in. x H 6.5 in. x D 0.5 in..
Weight: 5.6 oz

LARGE PANEL

SIZE	COLOR	ITEM #
9 x 6.5 in.	Black	80-0897
9 x 6.5 in.	Coyote	80-0898



Workplace Injury Kits

Class A **TRAUMA & FIRST AID KIT**

- Designed as a portable First Aid Kit for the workplace **exceeding ANSI/ISEA Z308.1-2015 Standard for Class A kits**
- Able to treat the most common workplace injuries with the added benefit of being able to control major bleeding with the C-A-T Tourniquet - the safest and most effective Tourniquet on the market!
- Packaged in a convenient zippered nylon bag or rugged hard case that is easily transported to the patient's side



TRAUMA & FIRST AID KIT - CLASS A

CASE	DRESSING	ITEM #
Soft Case	Compressed Gauze	80-0947
Soft Case	Bleeding Control Dressing	80-1045
Hard Case	Compressed Gauze	80-1033
Hard Case	Bleeding Control Dressing	80-1046

SPECIFICATIONS

Soft Case:

L 14.75 in. x H 8 in. x D 4 in.
Weight: 2.94 lb

Hard Case:

L 13 in. x H 9.5 in. x D 5 in.
Weight: 6 lb

Case only option (#80-0113) also available at NARescue.com

TRAUMA & FIRST AID KIT CLASS A:

<p>QTY 1 C-A-T® Tourniquet</p>	<p>QTY 2 Responder Compressed Gauze</p>	-OR-	<p>QTY 1 Responder Compressed Gauze</p>	<p>QTY 1 Bleeding Control Dressing</p>
<p>QTY 1 ETD™ 4 in. Responder, Flat</p>	<p>QTY 1 MicroShield CPR Face Mask</p>	<p>QTY 1 Elastic Wrap 2 in. Bandage</p>	<p>QTY 1 Elastic Wrap 4 in. Bandage</p>	<p>QTY 1 Gauze Pads 4"x4" (2pk)</p>
<p>QTY 2 Gauze Pad, 5" x 9"</p>	<p>QTY 16 Flexible Fabric Bandages</p>	<p>QTY 10 Burn Gel</p>	<p>QTY 1 BurnTec 4"x4" Dressing</p>	<p>QTY 1 Instant Cold Pack</p>
<p>QTY 1 Surgical Tape, 1 in.</p>	<p>QTY 1 Triangular Bandage</p>	<p>QTY 1 SAM® Splint II</p>	<p>QTY 1 Eyewash Solution 1 oz</p>	<p>QTY 2 PES Eye Shield</p>
<p>QTY 10 Antiseptic Towelette</p>	<p>QTY 10 Antibiotic Ointment</p>	<p>QTY 6 Hand Sanitizer</p>	<p>PAIR 2 Responder Nitrile Gloves, Lg</p>	<p>QTY 1 Responder Trauma Shears, Lg</p>
<p>QTY 1 First Aid Pocket Guide</p>	<p>QTY 1 Instructions for Use Public Access Bleeding Control</p>	<p>QTY 1 Kit Contents Inventory Card</p>	<p>QTY 1 *Card, Inspection</p>	

* Inspection Card is included in the Soft Case version only

Workplace Injury Kits

Class B **TRAUMA & FIRST AID KIT**

- Designed as a portable First Aid Kit for the workplace **exceeding ANSI/ISEA Z308.1-2015 Standard for Class B kits**
- Same supplies as the Class A Kit but with larger quantities to accommodate larger workplaces
- Packaged in a convenient zippered nylon bag or rugged hard case that is easily transported to the patient's side

SPECIFICATIONS

Soft Case:

L 14.75 in. x H 8 in. x D 7.5 in.
Weight: 4.87 lb

Hard Case:

L 13.5 in. x H 11 in. x D 6 in.
Weight: 9 lb

TRAUMA & FIRST AID KIT - CLASS B

CASE	DRESSING	ITEM #
Soft Case	Compressed Gauze	80-0948
Soft Case	Bleeding Control Dressing	80-1047
Hard Case	Compressed Gauze	80-1034
Hard Case	Bleeding Control Dressing	80-1048

Case only option (#80-0114) also available at NARescue.com

TRAUMA & FIRST AID KIT CLASS B:

 QTY 1 C-A-T® Tourniquet	 QTY 2 Responder Compressed Gauze	-OR-	 QTY 1 Responder Bleeding Control Dressing	 QTY 1 Bleeding Control Dressing	 QTY 1 ETD™ 4 in. Responder, Flat
 QTY 1 MicroShield CPR Face Mask	 QTY 2 Elastic Wrap 2 in. Bandage	 QTY 1 Elastic Wrap 4 in. Bandage	 QTY 2 Gauze Pads 4"x4" (2pk)	 QTY 4 Gauze Pad, 5" x 9"	 QTY 50 Flexible Fabric Bandages
 QTY 25 Burn Gel	 QTY 2 BurnTec 4"x4" Dressing	 QTY 2 Instant Cold Pack	 QTY 1 Surgical Tape, 1 in.	 QTY 1 Surgical Tape, 2" Mini	 QTY 2 Triangular Bandage
 QTY 1 SAM® Splint II	 QTY 1 Eyewash Solution 4 oz	 QTY 2 PES Eye Shield	 QTY 50 Antiseptic Towelette	 QTY 25 Antibiotic Ointment	 QTY 10 Hand Sanitizer
 PAIR 4 Responder Nitrile Gloves, Lg	 QTY 1 Responder Trauma Shears, Lg	 QTY 1 First Aid Pocket Guide	 QTY 1 Just in Time Instructions	 QTY 1 Kit Contents Inventory Card	 QTY 1 *Card, Inspection

ANSI Kits designed to treat the most common workplace injuries with the added benefit of being able to control major bleeding...



**COMMUNITY
PREPAREDNESS**



**3 to 5
MINUTES**

The Difference between Life and Death

- ▶ A victim who is bleeding can quickly die BEFORE professional emergency responders can even arrive on the scene
- ▶ Serious bleeding accounts for 35% of prehospital deaths from traumatic injuries
- ▶ Life-threatening bleeding warrants immediate intervention
- ▶ Basic bleeding control provided by bystanders can often be the difference between life and death...

Every Second Counts...



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

Trauma is a national health problem and after a traumatic injury, bleeding is responsible for over 35% of pre-hospital deaths. Examples of trauma include traffic & industrial accidents, natural disasters, falls, accidental injuries and intentional acts of violence such as shootings, stabbings and bombings.

These types of injuries often result in uncontrolled bleeding. Without immediate intervention these injuries can result in death within minutes – even before emergency medical services or other first responders can arrive.

North American Rescue

is leading the way to equip bystanders and initial first responders with the essential medical gear needed to take immediate action in life threatening emergencies. These lifesaving modular aid components are designed for individuals or groups in the workplace, public buildings, and populated public areas to administer immediate bleeding control at the point of wounding.

***Death from
Uncontrolled
BLEEDING
can occur
in as little as
3 to 5 Minutes***

H.

TRAUMA KITS

PUBLIC ACCESS

TRAINING

EVACUATION

TREATMENT SUPPORT

ORDER

Public Access **INDIVIDUAL BLEEDING CONTROL KITS**

Nylon Bag Option:

- Easy to restock in cases of use
- Convenient to carry and easy to open in stressful situations
- Opens flat for ease of use and immediate identification of enclosed devices
- Easily recognizable red nylon bag with attached PVC Bleeding Control ID Patch

Vacuum Sealed Option:

- Cost-effective alternative to Nylon
- Compact and rugged Vacuum-Sealed pouch
- Large Identifying Card Vacuum-Sealed inside the pouch for ease of identification

SPECIFICATIONS

Nylon Bag: H 7 in. x W 5 in. x D 2.5 in.

Vacuum Sealed: H 7 in. x W 5 in. x D 2.75 in.


*Essential equipment needed
for persons to take action
as immediate responders in
helping to stop life threatening
bleeding in Public Areas...*



Every **INDIVIDUAL BLEEDING CONTROL KIT** contains:



PLUS these Individual Components based upon YOUR configuration choice:



PUBLIC ACCESS INDIVIDUAL BLEEDING CONTROL KITS

LEVEL	BAG	WEIGHT	ITEM #
Basic	Red Nylon	15 oz	80-0452
Intermediate	Red Nylon	16.5 oz	80-0453
Advanced BCD	Red Nylon	17.4 oz	80-0523
Advanced	Red Nylon	17.4 oz	80-0454
Basic	Vacuum Sealed	11.7 oz	80-0465
Intermediate	Vacuum Sealed	13.25 oz	80-0466
Advanced BCD	Vacuum Sealed	14 oz	80-0522
Advanced	Vacuum Sealed	14 oz	80-0467

PVC Bleeding Control Patch (#ZZ-0796) is also available for separate purchase (See page 99 for more info)

**BLEEDING CONTROL
MULTI-PACKS****5-PACKS** *Vacuum*

- Five (5) Vacuum-sealed Individual Bleeding Control Kits inside an easy-to-open nylon tote
- Includes One (1) Responder QuikLitter™ for moving patients
- Available in Basic, Intermediate, Advanced BCD & Advanced

SPECIFICATIONS**5-Pack:**

H 7.5 in. x W 13 in. x D 8 in.

Each Individual

Bleeding Control Kit:

H 7 in. x W 5 in. x D 2.5 in

PUBLIC ACCESS BLEEDING CONTROL 5-PACKS

LEVEL	BAG	WEIGHT	ITEM #
Basic	Vacuum Sealed	5.7 lb	80-0573
Intermediate	Vacuum Sealed	5.75 lb	80-0574
Advanced BCD	Vacuum Sealed	6.5 lb	80-0576
Advanced	Vacuum Sealed	6.5 lb	80-0575

**TWIN
PACKS**

- Contents of 2 Individual Kits
- Compact Nylon Bag designed to fit next to an AED in the cabinet

SPECIFICATIONS

Nylon Bag:

H 10 in. x W 9 in. x D 3 in.

TWIN PACKS

LEVEL	BAG	WEIGHT	ITEM #
Basic	Red Nylon	1.66 lb	80-0951
Intermediate	Red Nylon	1.8 lb	80-0952
Advanced BCD	Red Nylon	1.8 lb	80-0953

8-PACKS *Nylon | Vacuum Sealed*

- Eight (8) Nylon or Vacuum-sealed Individual Bleeding Control Kits inside an easy-to-open nylon carrying case
- Includes One (1) Responder QuikLitter™ for moving patients
- Available in Basic, Intermediate, Advanced & Advanced BCD

SPECIFICATIONS**8-Pack:**

H 14 in. x W 11 in. x D 8 in.

Each Individual Bleeding Control Kit:

H 7 in. x W 5 in. x D 2.5 in.

PUBLIC ACCESS BLEEDING CONTROL 8-PACKS

LEVEL	BAG	WEIGHT	ITEM #
Basic	Red Nylon	8.5 lb	80-0460
Intermediate	Red Nylon	9.3 lb	80-0461
Advanced BCD	Red Nylon	9.7 lb	80-0525
Advanced	Red Nylon	9.7 lb	80-0462
Basic	Vacuum Sealed	6.8 lb	80-0468
Intermediate	Vacuum Sealed	7.6 lb	80-0469
Advanced BCD	Vacuum Sealed	8 lb	80-0524
Advanced	Vacuum Sealed	8 lb	80-0470

Every **BLEEDING CONTROL TWIN PACK** contains:

QTY
2 C-A-T®
Tourniquet



QTY
2 ETD 4in.
Responder Flat



PAIR
2 Responder
Nitrile Gloves, Lg



QTY
2 Survival
Blanket



QTY
1 Responder
Trauma Shears, Lg




QTY
2 Just in Time
Instructions



QTY
2 Permanent
Marker, Small

PLUS these
Individual
Components
based
upon YOUR
configuration
choice:



BASIC
QTY
2 Responder
Compressed Gauze



INTERMEDIATE
QTY
2 Responder
Compressed Gauze



ADVANCED BCD
QTY
2 HyFin®
Vent Compact, Twin



ADVANCED BCD
QTY
2 Bleeding
Control Dressing



ADVANCED BCD
QTY
2 HyFin®
Vent Compact, Twin

Clear Polycarbonate WALL STATIONS

- Includes Individual Bleeding Control Kits packaged in red nylon bags or Vacuum-Sealed pouches.
- The 5 & 8 Pack Stations also include an easy to-open nylon carrier and one (1) Responder QuikLitter™ for moving injured casualties.

SPECIFICATIONS

5 & 8 Pack Stations: H 18 in. x W 15.5 in. X D 10 in.
Empty Case Wt: 15.0 lbs (includes rack)
Large Case Item# ZZ-0822

Individual Stations: H 10 in. x W 10 in. X D 5.25 in.
Empty Case Wt: 2.39 lb
Small Case Item# ZZ-1000

8-PACK CLEAR POLYCARBONATE WALL STATIONS

LEVEL	#KITS	BAG	WGT	ITEM #
Basic	8	Vacuum	21.8 lb	80-0471
Intermediate	8	Vacuum	22.6 lb	80-0472
Advanced BCD	8	Vacuum	23 lb	80-0526
Advanced	8	Vacuum	23 lb	80-0473
Basic	8	Nylon	23.5 lb	80-0474
Intermediate	8	Nylon	24.3 lb	80-0475
Advanced BCD	8	Nylon	27.4 lb	80-0527
Advanced	8	Nylon	27.4 lb	80-0476

5-PACK CLEAR POLYCARBONATE WALL STATIONS

LEVEL	#KITS	BAG	WGT	ITEM #
Basic	5	Vacuum	20.7 lb	80-0581
Intermediate	5	Vacuum	20.7 lb	80-0582
Advanced BCD	5	Vacuum	21.5 lb	80-0584
Advanced	5	Vacuum	21.5 lb	80-0583

IND. CLEAR POLYCARBONATE WALL STATIONS

LEVEL	#KITS	BAG	WGT	ITEM #
Basic	1	Vacuum	3.3 lb	80-0938
Intermediate	1	Vacuum	3.42 lb	80-0939
Advanced BCD	1	Vacuum	3.48 lb	80-0940
Basic	1	Nylon	3.3 lb	80-0935
Intermediate	1	Nylon	3.6 lb	80-0936
Advanced BCD	1	Nylon	3.7 lb	80-0937

Alarms (#ZZ-1075) is also available for the Clear Polycarbonate Wall Cabinet (See page 67)



8 Pack Station



Individual Kit Station



5 Pack Station



Small Individual Clear Case: ZZ-1000



Large Multi-Pack Clear Case: ZZ-0822

- Durable, clear polycarbonate cabinet with rounded corners provides both physical and environmental protection for medical equipment
- Large case features wire shelf to support the contents
- External eyelet allows a secondary safety seal to easily identify if the station has been opened or tampered with
- External mounting tabs and hardware enable easy installation on any wall
- Clear design allows visibility at any angle

Polycarbonate
CLEAR CASE
features

Designed to be mounted on a wall at a designated location or to be co-located with Automatic External Defibrillators (AEDs) to provide publically accessible medical supplies to control major bleeding



5 Pack Station



8 Pack Station

Semi-Recessed
METAL CASE
features

- Mounts permanently inside the wall with a low profile 1 in. trim
- White powder painted corrosive resistant galvanized steel cabinet with clear plexiglass window
- Stainless steel continuous hinge & chrome plated handle designed to open with one quick pull

Semi-Recessed Metal WALL STATIONS

- Includes Individual Bleeding Control Kits packaged in either red nylon bags or Vacuum-Sealed pouches stored within an easy to-open nylon carrier and one (1) Responder QuikLitter™ for moving injured patients
- Outside Dimensions: H 17.5 in. x W 17.5 in. x D 1 in.
Inside Dimensions: H 14.13 in. x W 14.13 in. x D 6.88 in.
Empty Case Wt: 10 lb (Empty Case Item# ZZ-0921)

8-PACK SEMI-RECESSED METAL WALL STATIONS

LEVEL	#KITS	BAG	WGT	ITEM #
Basic	8	Vacuum	16.8 lb	80-0861
Intermediate	8	Vacuum	17.6 lb	80-0862
Advanced BCD	8	Vacuum	18 lb	80-0863
Advanced	8	Vacuum	18 lb	80-0864
Basic	8	Nylon	18.5 lb	80-0865
Intermediate	8	Nylon	19.3 lb	80-0866
Advanced BCD	8	Nylon	15.7 lb	80-0867
Advanced	8	Nylon	15.7 lb	80-0868

5-PACK SEMI-RECESSED METAL WALL STATIONS

LEVEL	#KITS	BAG	WGT	ITEM #
Basic	5	Vacuum	15.7 lb	80-0577
Intermediate	5	Vacuum	15.7 lb	80-0578
Advanced BCD	5	Vacuum	16.5 lb	80-0580
Advanced	5	Vacuum	16.5 lb	80-0579

Semi-Recessed Metal Case: ZZ-0921



*All Multi-Pack Wall Stations feature:
Highly visible "Bleeding Control Kit" letters for fast ID.
Instructions for Use are posted on the interior back panel*



8 Pack Nylon Station



Low Profile Metal Cabinet: ZZ-1199

Designed to be mounted on a wall with a slim 4 in. clearance and featuring easy pull-to-open design for rapid public access to major bleeding control supplies

Low Profile Metal WALL STATIONS

- Low profile metal cabinet mounts on the wall with a slim 4 in. clearance
- Designed to safely store 8 Individual Bleeding Control Kits in public areas
- Kits are available in red nylon bags or Vacuum-Sealed pouches and stored in Two (2) easily removable Mesh Panels each equipped with 4 pockets and nylon handle
- Features highly visible "Bleeding Control Kit" lettering for fast identification and convenient Instructions for Use on the interior door panel

SPECIFICATIONS

H 15.5 in. x W 22.5 in. X D 4 in.

Empty Cabinet Weight: 9.365 lb

Low Pro Cabinet with ID Stickers:
Item# ZZ-1199



8 Pack Vacuum Sealed Station

8-PACK LOW PROFILE METAL WALL STATIONS

LEVEL	#KITS	BAG	WGT	ITEM #
Basic	8	Vacuum	17.49 lb	80-1083
Intermediate	8	Vacuum	18.19 lb	80-1084
Advanced BCD	8	Vacuum	18.69 lb	80-1085
Basic	8	Nylon	17.49 lb	80-1080
Intermediate	8	Nylon	18.19 lb	80-1081
Advanced BCD	8	Nylon	18.69 lb	80-1082

Metal Low-Pro CABINET features

- Low profile metal cabinet with 4 in. wall clearance
- White powder painted corrosive resistant galvanized steel cabinet with stainless steel continuous hinge
- Provides convenient slimline storage for medical equipment
- Magnetic closure provides quick easy access with one pull to open - No twist and pull



Alarms (#ZZ-1075) is also available for the Low Profile Metal Cabinet (See page 67)

- Small, compact and easy to use
- Kits contain the same Tourniquet and Hemostatic Dressing technology recommended by the U.S. Military & the American College of Surgeons Committee on Trauma to control bleeding
- Available in Vacuum-sealed Pouch or Red Nylon Bag

STB KITS

Nylon or Vacuum-Sealed



VAC STB KIT

ITEM

85-1574

SPECIFICATIONS

H 7 in. x W 5 in. x D 2.5 in.

Weight: 8 oz

VACUUM-SEALED STB KIT contains:



SPECIFICATIONS

H 7 in. x W 5 in. x D 2.5 in.

Weight: 15 oz

NYLON STB KIT

ITEM

85-2995

NYLON STB KIT contains:



Cabinet Security ALARMS

- Alarm compatible with Clear Polycarbonate and Low Profile Metal Cabinets
- Loud 120 decibel alarm sounds when magnet separates from base more than 0.5 in.
- Easy and convenient to install
- Requires no wiring
- Mounting hardware and instructions included
- Designed for indoor use only



SPECIFICATIONS

W 163 in. x H 2.41 in. x D 0.72 in.

Weight: 6.4 oz each

CABINET SECURITY ALARMS

ITEM#	UNIT
ZZ-1075	Pack of 2

Accessories

3-WAY SIGN

Bleeding Control



- Highly visible for locating lifesaving bleeding control kits immediately
- 3 way "V" shape projects from wall for easy identification from multiple directions

SPECIFICATIONS

H 7 in. x W 16.5 in. x D 6.75 in.

Weight: 0.35 oz

3-WAY SIGN

ITEM

ZZ-1117

Public Access **INDIVIDUAL BLEEDING CONTROL TRAINER KITS**

- The appropriate equipment to teach individuals the necessary life-saving skills of bleeding control
- Corresponds with NAR's single Public Access Bleeding Control Kits in Basic, Intermediate, or Advanced
- Packaged in a compact, blue nylon bag with PVC Bleeding Control Patch


Every **INDIVIDUAL TRAINER KIT** contains:



SPECIFICATIONS
H 7 in. x W 5 in. x D 2.5 in.



PLUS these Individual Components based upon **YOUR** configuration choice:



BLEEDING CONTROL INDIVIDUAL TRAINER KITS

LEVEL	WGT	ITEM #
Basic	15 oz	80-0901
Intermediate	16.5 oz	80-0902
Advanced	16.5 oz	80-0903



STB TRAINER


- Designed to provide the essential training equipment to teach how to control bleeding
- Contains training items corresponding with the necessary tools recommended by the American College of Surgeons Committee on Trauma to equip bystanders to help save lives in the instance of life threatening bleeding

SPECIFICATIONS
H 7 in. x W 5 in. x D 2.5 in.
Weight: 8 oz

STB TRAINER

ITEM #
85-1880


STB TRAINER contains:




BASIC BLEEDING CONTROL TRAINING SKILLS KIT CONTENTS:



INTERMEDIATE BLEEDING CONTROL TRAINING SKILLS KIT CONTENTS:



ADVANCED BLEEDING CONTROL TRAINING SKILLS KIT CONTENTS:



BLEEDING CONTROL TRAINING SKILLS KITS

- Provides the appropriate equipment for up to 20 students to learn the necessary lifesaving skills of bleeding control
- Packaged in a self-contained, wheeled, rigid case with internal dividers for easy organization and inventory
- Skill set training addresses extremity bleeding, wound packing, pressure bandages, chest injuries and patient transport
- Corresponds with NAR's Public Access Bleeding Control Kit Levels: Basic, Intermediate & Advanced
- Available with or without Wound Mgt. Simulator

SPECIFICATIONS

Hard Storage Case:
W 32 in.x H 21 in. x D 12.5 in.

BLEEDING CONTROL TRAINING SKILLS KITS

LEVEL	WGT	ITEM #
Basic	40 lb	80-0881
Intermediate	42 lb	80-0879
Advanced	42 lb	80-0877
Basic Wound Mgt. Simulator	44 lb	80-0880
Intermediate Wound Mgt. Simulator	48 lb	80-0878
Advanced Wound Mgt. Simulator	48 lb	80-0876



TCCC/TECC TRAINING SYSTEM 3 Module Set

- Turnkey solution for teaching the PHTLS TCCC/TECC required tasks for completing the Tactical Combat Casualty Care (TCCC/TECC) course
- All the necessary tactical medical and casualty extraction/evacuation equipment needed for training centers and educational institutions
- Designed to equip up to 20 students
- Complete system consists of three self-contained hard storage case modules: Skill Stations Module, IFAKS Module and Extraction Module

SPECIFICATIONS

Each Case: W 30.75 in. x H 20.5 in. x D 11.6 in.
Weight: Module 1- 53 lb • Module 2 - 59 lb
Module 3 - 45 lb • 3 Module Set:: 157 lb

3 MODULE SET

UNIT	ITEM #
3 Modules	85-0550



Skill Stations Module TCCC/TECC MODULE 1



- Train up to twenty (20) students at individual skill stations in the treatment of: *extremity hemorrhage, wound packing, chest seal application, needle chest decompression, junctional hemorrhage, basic airway adjuncts, surgical airway, splinting, hypothermia prevention/management, IV skills and eye injuries*
- Self-contained, wheeled, rigid case with internal dividers to keep the equipment organized and easy to inventory

MODULE ONE

UNIT	ITEM #
Module 1 in Case	85-0519
ReSupply*	85-1900

SPECIFICATIONS

W 30.75 in. x H 20.5 in. x D 11.6 in.
Weight: 53 lb

SKILL STATIONS TRAINING MODULE 1:

 QTY 20 C-A-T® Blue Trainer	 QTY 20 ETD® 6 in. Dressing	 QTY 20 Gauze S-Rolled	 QTY 20 Combat Gauze® Trainer	 QTY 1 JETT® Junctional Device	 QTY 5 28F Naso Airway & Lube	 QTY 1 Cyclone Pocket BVM
 QTY 2 Tactical CricKit®	 QTY 10 HyFin® Trainer, Twin Pk	 QTY 5 ARS® Needle, 14g	 QTY 20 Saline Lock Kit	 QTY 5 Sharps Shuttle	 QTY 1 HPMK® Kit	 QTY 5 Triangular Bandage
			 QTY 5 SAM® Splint II	 QTY 5 PES Eye Shield		

*ReSupply Kit available with all contents of Module 1, excluding Rigid Storage Case

PHTLS TCCC/TECC Modular Training System

IFAK
Module

TCCC/TECC MODULE 2

- Twenty pre-packed Training IFAKs for use during trauma lanes or training scenarios
- IFAKs come packaged in a self contained wheeled hard, case that allows for easy transport and storage

MODULE TWO

UNIT	ITEM #
Module 2 in Case	85-0520

SPECIFICATIONS

W 30.75 in. x H 20.5 in. x D 11.6 in.
Weight: 59 lb

IFAK TRAINING MODULE 2: 20 TRAINING IFAKs each containing:



Extraction
Module

TCCC/TECC MODULE 3

High-threat extraction devices and casualty evacuation platforms packaged in a self-contained, wheeled, rigid case that allows for easy transport and storage


MODULE THREE

UNIT	ITEM #
Module 3 in Case	85-0521

SPECIFICATIONS

W 30.75 in. x H 20.5 in. x D 11.6 in.
Weight: 45 lb

EXTRACTION TRAINING MODULE 3:



SKILL STATIONS/IFAK Combo

- Train up to ten (10) students at individual skill stations as well as with 10 pre-packed TCCC Training IFAKs for use during trauma lanes or training scenarios
- Includes self-contained, wheeled, rigid case with internal dividers for easy equipment inventory and organization

SPECIFICATIONS

W 30.75 in. x H 20.5 in. x D 11.6 in.
Weight: 60 lb

COMBO MODULE

UNIT	ITEM #
SKILLS/IFAK Combo	85-0740

SKILL STATIONS/IFAK COMBO:



IFAK TRAINER Individual First Aid Training Kit

- Pre-packaged training kit with all the essential equipment for training initial care
- The same point-of-wounding training kit found in the TCCC/TECC Training Module #2 (IFAKS Module) & Skill Stations/IFAKs Combo, sold individually
- Easy-to-attach modular carrying bag can be worn as a drop-down leg rig, attached to a belt with a quick-detach buckle system or MOLLE'd to a vest
- ReSupply Kit (contents only, no nylon bag) also available in a clear plastic pouch

IFAK TRAINER | TCCC-TECC

UNIT	ITEM #
Complete Kit	85-0493
ReSupply	85-0518

SPECIFICATIONS

H 8.25 in. x W 7.1 in. x D 5.5 in.
Weight: 2 lb 10 oz

IFAK TRAINER:



Part Task Trainers

WOUND MANAGEMENT SIMULATOR

- Rugged part-task trainer for learning essential lifesaving skills of tourniquet application, wound packing and applying pressure bandages.
- Multiple wounding patterns to enhance training objectives
- Simulates severe laceration, penetrating trauma (entrance and exit wounds), along with the ability to apply tourniquets
- Proprietary silicon that is rugged and durable and will not break down due to environmental exposure
- Can be used in the classroom or in field training scenarios


WOUND MGT SIMULATOR

ITEM

93-0020

SPECIFICATIONS

L 15 in. x Diameter 6 in.
Weight: 9.3 lb



Portable training aid to help learn wound management techniques to control bleeding

HEMORRHAGE CONTROL TRAINING KIT *with Bio Feedback*

- Training tool designed to teach effective wound packing skills with Immediate results for both trainer and student
- Verifies correct or incorrect technique

SPECIFICATIONS

Case: 23 in. x 8 in. x 15 in.
Weight: 17 lb 3.2 oz

BIO FEEDBACK

ITEM

93-0019



HEMORRHAGE CONTROL TRAINING KIT *Combat Gauze®*

- Comes with 10 Combat Gauze Trainers and 5 Combat Gauze LE Demo Samples
- Includes a "training leg" that features an entry/exit bullet wound and a deep laceration

SPECIFICATIONS

Case: H 15 in. x
W 23 in. x D 8 in.
Weight:
15 lb 8 oz

COMBAT GAUZE

ITEM

93-0017



TRUECLOT® Hemorrhage Training Kit™

- Highly realistic Wound Packing Trainer
- Large caliber gunshot wound with hidden bleed points and a simulated fracture/defect hard plastic bone insert to aid in teaching deep wound packing techniques
- Sturdy silicone construction with no moving parts, reinforced tubing mount, and a more realistic wound pattern and outer anatomy
- Ideal for the repeated use of a training environment



SPECIFICATIONS

Task Trainer:
L 7.25 in. x H 2.5 in.
Storage Case:
L 13.5 in. x W 10 in.
x D 5.5 in.
Liquid Capacity:
1 liter
Weight: 5 lb

TRUECLOT

ITEM

93-0020



SAM IO TRAINING KIT *Single Bone*



- All of the equipment needed to familiarize trainees with the SAM IO Intraosseous Access hardware and usage
- Includes single training humerus bone with tissue block
- Rugged nylon carrying case features elastic loops for storing needles, along with attachment points for the driver once you've removed these items from their original packaging
- Equipment is for training use only

SAM IO TRAINING KIT - SINGLE BONE:



SPECIFICATIONS

H 7 in. x W 11.25 in. x D 5.5 in.
Weight: 2.32 lb

SAM IO SINGLE BONE

ITEM

30-0256

SAM IO FIELD STORAGE CASE

- Custom field storage case designed for the SAM IO Intraosseous Access System
- Elastic loops for storing needles and the SAM IO Driver
- Two compartments separated by a divider with convenient mesh pouch

SPECIFICATIONS

H 7 in. x W 11.25 in. x D 5.5 in.
Weight: 1 lb

SAM IO FIELD STORAGE CASE

ITEM

30-0254



*Contents shown
for display purposes
Case is sold separately*

C-A-T® TRAINER

Combat Application Tourniquet

- Fully functioning C-A-T® tourniquet colored blue to identify it as a training device
- Enables instructors and trainees hands-on training & tourniquet application practice
- Providing a distinct visual means to distinguish training devices from personal/operational life-saving devices



C-A-T TRAINER

ITEM #	NSN #
30-0033	6910-01-560-2972

SPECIFICATIONS

H 6.5 in. x W 2.4 in. x D 1.5 in.
Weight: 2.7 oz

See page 11 for Tactical and Rescue versions of the C-A-T®

HYFIN® TRAINER

Twin Pack

- Unique adhesive gel formulation does not damage mannequins or leave behind excessive residue
- Packaged with two training chest seals in one unique perforated package allowing students to open only one seal at a time



SPECIFICATIONS

Packaged:
L 7.5 in. x W 4.5 in. x D 0.25 in.
Chest Seal:
H 6 in. x W 6 in. x D 0.1 in.
Weight: 2.5 oz

HYFIN TRAINER

ITEM #
10-0032

See page 19 for other HyFin® Chest Seal products

COMBAT GAUZE®

Moulage Trainer

- Thoroughly instruct students on the proper application of Z-Fold Combat Gauze®
- Completely inert (no hemostatic properties) for training purposes only
- Re-sealable packaging for multiple uses



SPECIFICATIONS

Packaged:
H 6 in. x W 4.5 in. x D 1.25 in.
Gauze: L 4 yd x W 3 in.
Weight: 0.75 oz

COMBAT GAUZE MOULAGE TRAINER

ITEM #
30-0063

Responder

TRAINING ETDs

- Economic, low cost solution for training in the application of immediate direct pressure
- Non-sterile, elastic pressure bandages

ETD TRAINER

SIZE	ITEM #
4 in.	30-0145



SPECS

Bandage:
L 60 in. x W 4 in.
Pad: H 6 in. x W 4 in.
Packaged:
L 4.5 in. x W 2.5 in. x D 1.75 in.
Weight: 2.0 oz

High Threat Extraction Gear

HIGH THREAT EXTRACTION KIT

Essential tools for rescue ops in high threat environments



- Facilitates immediate casualty extraction so that Care Under Fire casualty treatment may be performed as soon as possible
- Self-contained kit with all of the essential tools for tactical operators and rescue personnel
- Includes a C-A-T® tourniquet to immediately address extremity hemorrhage at the point of wounding
- Rugged nylon bag attaches to the duty belt as a drop down leg rig or to the assault vest/body armor with PALS/MOLLE-style connectors


SPECIFICATIONS

L 8.5 in. x W 8.75 in. x D 4.5 in.
Weight: 3 lb 6 oz

HIGH THREAT EXTRACTION

COLOR	ITEM #
Black	80-0139
Coyote	80-0141
Digital Camo	80-0142
OD Green	80-0140

HIGH THREAT EXTRACTION KIT CONTENTS:



TACTICAL BAIL-OUT KIT

Rapid extraction of individuals or teams

- Low profile, multi-purpose kit for rapid exfil of an individual or team
- Integral for multiple hasty rescues including rapid casualty lowering & hoisting, rescuer throw bag and rappel line
- Attach to belt via straps or mounted via 2 non-locking carabiners included in the kit
- Featuring Auto-locking Kong Tango carabiner

SPECIFICATIONS


L 10 in. x W 6 in. x D 4 in.
Weight: 2 lb 14 oz

TACTICAL BAILOUT

ITEM #
90-0020



TACTICAL BAILOUT CONTENTS:



High Threat Extraction Gear

RAT STRAP

Rescue Assault Tether

- For head/feet first linear drags during shield-assisted, open-field/indoor-corridor rescues or single-rescuer extractions
- Innovative Daisy Chain configuration offers multiple reinforced attachment points

SPECIFICATIONS

Pkg: H 1 in. x W 4 in. x D 5.5 in.
Extended: L 62 in.
Minimum Break Strength: 2,702 lb
Weight: 12.7 oz

RAT STRAP

ITEM #	NSN #
90-0002	6530-01-632-7872



EX-TRACT STRAP

Point of Wounding Extraction


- Point of wounding extraction
- Back pack carry
- Hasty harness for traditional and confined space lowering and hauling

SPECIFICATIONS

Pkg: H 18 in. x W 3.5 in. x D 1 in.
Extended: L 58 in.
Weight: 9.3 oz

EX-TRACT STRAP

ITEM #
90-0082



DRAGON HARNESS

Under the Vest Harness


- Built-in handle system for immediate high threat extraction
- Operator harness to aid in the rescue of downed personnel without the need to holster or sling weapons
- Prevents the "ride up" of the casualty's ballistic vest, reducing exposed, unprotected body regions

SPECIFICATIONS

Height: 20 in. from waist to non-deployed handle
Length: 14 in. from harness attachment point to top of handle
Weight: 1 lb 1.6 oz.
Waist Size: 30 in. to 45 in.

DRAGON HARNESS

ITEM #
90-0003



DRAGON HANDLE

Tool offering head/feet first drags


- Maintain an anatomically correct posture for rapidly extracting greater weight
- Decreases manpower needs in "zeroed in" locations
- Black Omega, non-locking "D" carabiner rated for 31kN (6,968.8 lb).

SPECIFICATIONS

H 6.3 in. x W 4.1 in. x D 2.25 in.
Weight: 11.3 oz
Device Length: L 42 in.
US Patent(s): Pending

DRAGON HANDLE

ITEM #	NSN #
90-0004	6530-01-598-9595



HASTY HARNESS *Multi-purpose Harness*

- Indispensable gear for every webbing equipment toolbox
- Resourceful, multi-purpose harness for unlimited applications & maximum versatility
- Constructed of 22 ft loop of 1 in. Tech Tape (MBS: 4,400 lb) with opposing color thread to aid in visual safety inspection

SPECIFICATIONS

Loop: L 10 ft 11 in. x W 1 in.
Circumference - 21 ft 10 in.
Min. Break Str: Tech Tape (4,400 lb)
Sewn Sling (6,000 lb)



HASTY HARNESS NATO NSN*


COLOR	ITEM #	NSN #
Red	90-0044	----
Black	90-0001	6530-01-674-9774*

Rescue Rope Equipment


R.A.C.E. SYSTEM

Rapid Access Casualty Extraction

- Lightweight, multiple use casualty evacuation system enabling first responders/medical providers to ascend, descend, hoist & lower both rescuer & casualty to safety
- Combines 150 ft of lightweight durable 9.5 mm Olympus Pegasus Rope® with the latest in lightweight multipurpose hardware
- Capable of rigging 2:1, 3:1, 4:1, 5:1, 6:1, and 9:1 systems in different configurations for a wide variety of rescuer options with fully illustrated instruction booklet
- Front pocket allows for optional storage
- Carriage options include 4 Point, D-ring configurable modified side and/or backpack carry-style adjustable strap



R.A.C.E SYSTEM CONTENTS:



R.A.C.E. SYSTEM

COLOR	ITEM #
Multi-Cam	90-0063
Coyote	90-0077

SPECIFICATIONS

H 15 in. x W 10 in. x D 7.5 in.
Weight: 10 lb 14 oz

Water Rescue BELT THROW BAG

- Multipurpose water rescue solution designed for easy deployment and retrieval for rescuers
- Quick release reverse tapered throw bag for rapid packing of Triton rope and a slim profile for users
- 2 in. adjustable belt with quick release strap, 1.5 in. plastic cam buckle, and red toggle ball for easy ID



BELT THROW

ITEM #
90-0065

SPECIFICATIONS

Rope Bag: L 11 in. x W 4.5 in. x D 4 in.
Gear Pouch: L 5.5 in. x W 5 in. x D 2 in.
Total Weight: 2 lb 2 oz



50 Foot THROW BAG

- Must-have necessity for water rescue
- Durable, easy-to-use rescue device with highly visible orange nylon bag and mesh panels that drain quickly for use in multiple rescue operations
- Packed with 50 ft of 7.5 mm Triton Pegasus Rope® rated at 3,050 lbf
- Triton rope is strong enough for multiple uses: as a throw line, tension diagonal, tow line and for pinned kayak recovery



50 FT THROW

ITEM #
90-0064

SPECIFICATIONS

Rope Bag: L 11 in. x W 4.5 in. x D 4 in.
Gear Pouch: L 5.5 in. x W 5 in. x D 2 in.
Total Weight: 2 lb 2 oz



HELO TAGLINE KIT



SPECIFICATIONS

L 11 in. x W 7 in. x D 7 in.
Weight: 5 lb 13 oz

HELO TAGLINE

ITEM

90-0069

- **Help control litter movement during hoist operations** - a requirement on all Army MEDEVAC helicopters

- Contains 250 ft of 7.5 mm Triton Pegasus Rope® featuring a highly visible polypropylene sheath with an Innegra™ core

HELO TAGLINE KIT CONTENTS:



QTY **1** 250 ft Triton Rope



QTY **1** rockD Carabiner



QTY **1** Quick Link, SS



QTY **10** Weak Links

A-VRR Adaptable Vertical Response Rack Kit


The Adaptable Vertical Response Rack was designed meet and exceed full-spectrum operational vertical requirements for urban through mountain / canyon access, rescue, and evacuation variables. The kit capabilities were driven utilizing NFPA 1006 (Operations & Technician Level) JPR's and the KPP's (Key Performance Parameters) within the USSOCOM SOF CASEVAC Program framework. The A-VRR Kit's organic versatility allows it to be scaled up to meet all relevant job performance requirements (JPR's) for complex team access and/or rescue or be scaled down for a "fast & light" mission specific assault configuration.

P-A-C-E methodology was engineered throughout the kit design, enabling multiple options for progressive capture and rope grab capabilities depending on the environmental pathology.

- Haul / Mechanical Advantage Capability – Multiple simple, compound and complex configurations ranging from a 2:1 through 9:1
- Full spectrum Friction Configurations for Rappels, Bail-outs, And One or Two Person Lowering Systems
- Traditional and MA-assisted Ascending Techniques
- Multiple Basic and Anchor Configurations to Include Contingency and Releasable Systems
- Load Transfer Capabilities to include Knot Passing on Rappel and within a system and casualty Pick-Offs

SPECIFICATIONS

Packaged: 12 in. x 15 in.
Weight: 3.25 lb



A-VRR KIT


ITEM

90-5000

A-VRR KIT CONTENTS:




QTY **1** Totem Decent Control Device



QTY **3** HMS Auto-lock Carabiner



QTY **1** Oval Auto-lock Carabiner



QTY **1** Oval Auto-lock Carabiner w/Wire Eye



QTY **1** MagnaPulley Small Pulley Pair 7/16 in. w&w/o Beckett



QTY **1** Pulley with Rope Grab




QTY **1** H-Frame Auto-lock Carabiner w/Pulley



QTY **1** Ascender / Rope Clamp



QTY **1** 16ft 6mm XTEC Cord



QTY **1** 8mm VT Prusik



QTY **1** 120cm Sling w/Aramid Core



QTY **1** 7mm Screw Link

MEDEVAC LITTER *Quadfold Litter*

- The ultimate, cost efficient patient evacuation stretcher for First Responders operating in Active Shooter, Mass Casualty and Disaster incidents
- Compact, collapsible, quad-fold litter is easy to deploy and carry
- Two patient securing buckles, six IV-attachment points, ergonomically designed handles & patented joint configuration with locking leg crossbar

SPECIFICATIONS

Folded: H 21.5 in. x W 8.5 in. x D 7.5 in.
Open: H 84.5 in. x W 22.5 in. x D 4 in.
Weight: 13 lb 8 oz
Static Payload: 1,200 lb (544.3 kg)
Ground Clearance: 1.5 in.

MEDEVAC LITTER

BED COLOR	FRAME COLOR	ITEM #
Black	OD Green	60-0051
Orange	Black	60-0052



SPECIFICATIONS

Folded, handles collapsed:
H 20.5 in. x W 8.5 in. x D 9.5 in.
Open, handles collapsed:
L 78 in. x W 22.5 in. x D 6 in.
Open, handles extended:
L 90 in. x W 22.5 in. x D 6 in.
Weight: 15 lb 6 oz (6.97 kg)
Static Payload: 1,800 lb (816 kg)
Ground Clearance: 1.5 in.
U.S. Patents: 6,842,923 & 5,598,592,
Patent(s) Pending

TALON 90C LITTER

ITEM#	NSN#
60-0002	6530-01-504-9051

TALON® II 90C

Quadfold Litter with Collapsible Handle



- Standard 90 in. NATO-compatible evacuation platform
- Developed to U.S. Army requirements
- Advanced, innovative combat evacuation system featuring patented portable quad-fold design with retractable handles
- Slip-resistant, chemical-resistant rip-stop fabric tested and approved by Federal & private labs using live warfare agents



WHEELED LITTER CARRIER

Lightweight litter cart with case

- Lightest wheeled litter cart on the market for single person transport to and from MedEVAC vehicles
- Locking clamp system keeps the litter firmly in place and allows for different litter widths
- Compatible with the Talon® II 90C, Raven® 90C, MedEvac & T3, plus all NATO litters
- Quick & easy to deploy or disassemble



SPECIFICATIONS

Deployed: L 56.29 in. x W 27.5 in. x H 36.2 in.

Case: L 48 in. x W 29.5 in. x H 12.5 in.

Weight: Carrier: 42 lb Carrier with Case: 98 lb

Payload: 550 lb

WHEELED LITTER CARRIER

ITEM#

60-0063

OSL™ LITTER *Oversized Litter*

- Safely evacuate larger patients with 30% larger surface area compared to standard-sized evacuation platforms
- Includes 6 stirrups, 6 padded carrying loops & 3 patient securing straps

SPECIFICATIONS

Folded: H 45 in. x W 9 in. x D 8 in.

Open: L 90 in. x W 31 in. x D 6 in.

Weight: 22 lb 6 oz (10.2 kg) Ground Clearance: 1.5 in.

Static Payload: 2,000 lb (909 kg)

OSL LITTER

ITEM#	NSN#
60-0023	6530-01-552-7708



RAVEN® 90C *Bifold Litter*

- The official litter of the U.S. Navy and certified safe-to-fly by the U.S. Air Force
- Slip-resistant, chemical-resistant, rip-stop fabric

SPECIFICATIONS

Folded: H 45.5 in. x W 7.5 in. x D 8 in.

Open: 90 in. x W 22.5 in. x D 6 in.

Weight: 16 lb 8 oz (7.5 kg)


Static Payload: 2,000 lb (909 kg)

Ground Clearance: 1.5 in.

U.S. Patents: 6,842,923 5,598,592

RAVEN 90C

ITEM#	NSN#
60-0001	6530-01-432-5114



T3™ LITTER *Lightweight, Assault*

- Ultra light, ultra strong litter weighing less than 12 lb
- Specifically designed for Special Operations Forces
- Overall length enables floor loading on an UH-60 helicopter

SPECIFICATIONS

Folded: H 21 in. x W 8.75 in. x D 7.25 in.

Open: L 81 in. x W 22.5 in. x D 6 in.

Weight: 11 lb 13 oz (5.34 kg)

Static Payload: 1,200 lb (544.3 kg)

Ground Clearance: 1.5 in.

T3 LITTER

ITEM#

60-0022



STINGRAY® Poleless Litter

- Ideal for decontamination roller systems & moving wet patients
- Non-skid, flame-retardant monofilament polypropylene fabric that is resistant to mildew, acids, alkalis and penetrating liquids
- Tested by Federal & private labs using live warfare agents
- Includes six heavy-duty carrying handles & 2 casualty securing straps
- Easy to store and clean

SPECIFICATIONS

Folded: H 20 in. x W 12 in. x D 2 in.
Open: L 78 in. x W 24 in.
Weight: 13 oz (0.37 kg)
Payload: 1,000 lb (453.6 kg)

STINGRAY

ITEM #	NSN #
60-0004	6530-01-492-4079



TED Tactical Extrication Device

- Designed to allow two people to carry a patient to safety
- Lighter, more compact non-rigid evacuation platform designed specially for use in high-threat environments and limited spaces
- Constructed of heavy duty vinyl with 4 sturdy carrying handles
- Packaged in a sturdy, transport friendly bag with adjustable shoulder strap

SPECIFICATIONS

Folded: H 5 in. x W 10.25 in. x D 3.5 in.
Open: L 43.5 in. x W 19 in.
Weight: 1 lb 8 oz (0.68 kg)
Payload: 600 lb (272 kg)

TACTICAL EXTRICATION DEVICE

COLOR	ITEM #	NSN #
Black	60-0005	6530-01-531-6808
OD Green	60-0006	6530-01-492-4114



NAR QUIKLITTER LITE, Standard, Responder or XL

- Multi-use, low cost disposable evacuation platforms for safe, effective emergency patient transport & transfer
- Water-resistant lightweight heavy duty construction with ten large hand-holds
- Super compact LITE is packaged to fit in a cargo pocket
- The XL is 12 in. wider and 6 in. longer featuring 12 hand-holds and two built-in patient retention straps

SPECIFICATIONS

LITE: Payload 750 lb.
Litter Bed: L 74 in. x W 27 in.
Weight: 12 oz

STANDARD/RESPONDER:
Payload 500 lb.
Litter Bed: L 74 in. x W 28 in.
Weight: 13.3 oz

XL: Payload 1,000 lb
Litter Bed: L 80 in. x W 40 in.
Weight: 2.2 lb

QUIKLITTER

SIZE	ITEM #
LITE	60-0071
Standard	60-0061
Responder	60-0080
XL	60-0065



EEL Emergency Evacuation Litter

- Cost-effective, non-rigid evacuation platform designed for fast, effective patient evacuation in the tactical environment
- Rapid deploy handle allows quick removal from storage pouch
- Rugged, lightweight black nylon with 8 carry handles
- Compact, small-cube storage pouch attaches to the assault vest with PALS/MOLLE

SPECIFICATIONS

Pouch: L 9.75 in. x W 7.5 in. x D 2.5 in.
Litter (Open): L 78 in. x W 30 in.
Weight: 1 lb 14 oz (0.57 kg)
Payload: 700 lb (317.5 kg)

EEL

COLOR	ITEM #
Black	60-0047
Coyote	60-0048
Multi-Cam	60-0050
OD Green	60-0049



LITTER STANDS

33 in. Standard

- Set of 2 lightweight stanchions for standard sized litters
- Constructed of a lightweight aluminum frame that conveniently folds flat for ease in transport & storage
- Quickly locks open to safely secure the litter in place

SPECIFICATIONS

Folded: H 34 in. x W 25.5 in. x D 6.92 in.
Extended: 33 in.
Base Width: 23.5 in.
Weight: 8 lb 8 oz

33 IN. LITTER STANDS

UNIT	ITEM#	NSN#
Set of 2	60-0014	6530-01-523-3314

OSL™ LITTER STANDS

for the OverSized Litter™

- Set of 3 lightweight aluminum stanchions accommodate our oversized OSL™ Litter
- Quickly opens and securely locks the litter in place
- Easily folds for transport and storage

SPECIFICATIONS

Folded: H 34 in. x W 25 in. x D 2.5 in.
Extended Height: 32 in.
Base Width: 23.5 in.
Weight: 28 lb 8 oz (12.93 kg)

OSL LITTER STANDS

UNIT	ITEM#
Set of 3	60-0029

OSL™ BACKREST

for the OverSized Litter™

- Position patient at 45° & 90° angles for maximized treatment options & increased comfort
- Constructed of a lightweight aluminum frame with durable Nomex/Kevlar-blend fabric cover
- Compatible with our OSL™ litter

SPECIFICATIONS

H 31.5 in. x W 30.25 in. x D 1.6 in.
Weight: 2 lb 15 oz (1.33 kg)

OSL BACKREST

ITEM#
60-0024

LITTER BACKREST

Talon® or Raven®

- Allows patients on a Talon® 90C or Raven® 90C litter to be placed at 90° or 45° angles
- Constructed of a lightweight aluminum frame with Nomex/Kevlar-blend fabric covering

LITTER BACKREST

ITEM#	NSN#
60-0027	6530-01-594-7444

SPECIFICATIONS

H 31.5 in. x W 29.25 in. x D 1.6 in.
Weight: 2 lb 13 oz

Talon® Assault

LITTER CARRIER

- Transport NAR's quad-fold litters (Talon, MedEvac & T3)
- MOLLE-style attachment points, thin back plates, drain hole, non-slip, adjustable shoulder straps & top-carrying handle

SPECIFICATIONS

Collapsed: H 20 in. x W 9 in. x D 2 in.
Expanded: H 20 in. x W 9 in. x D 10 in.
Weight: 1 lb 7 oz

LITTER CARRIER

*NATO NSN

COLOR	ITEM #	NSN #
Black	60-0009	6530-01-531-6907
Coyote	60-0011	6530-01-531-6854
Digital Camo	60-0025	6530-01-643-9134*
OD Green	60-0010	6530-01-504-9055

Adjustable IV POLE

Exclusively designed to fit the pre-molded holes of our framed litters

SPECIFICATIONS

H 19-33 in. x W 8.75 in. x D 0.6 in. | Wgt: 4.8 oz

ADJUSTABLE IV POLE

ITEM#	NSN#
60-0008	6530-01-504-9054

Triage Cards

T2® TAG NAR

- Triage & treatment tag utilizing the S.A.L.T. Mass Casualty Triage Standard
- To aid in documenting critical patient information throughout the entire mass casualty response process
- Liquid-resistant in all settings; made of rugged materials contained in a plastic sheath with latex-free band


T2 TRIAGE TAG

UNIT	ITEM #	NSN #
20 Pack	20-0031	6515-01-624-7148
1000 Pack	20-0032	6545-01-592-3261

SPECIFICATIONS

L 5.69 x W 3.88 in.
Weight: 0.7 oz

INCLUDED WITH THE T2 TAG:



TACTICAL TRIAGE KIT T2®


- Compact incident command system
- Streamlines medical triage
- Contained in a lightweight, rugged vacuum-sealed pouch

SPECIFICATIONS

L 5.25 x W 7 in. x D 1.25 in.
Weight: 4.9 oz

TACTICAL TRIAGE KIT

ITEM #	NSN #
20-0040	6545-01-591-2134



COMBAT CASUALTY CARD

- Streamlined triage card that meets all TCCC standards for casualty management
- State-of-the-art ruggedized paper is tear resistant, waterproof & designed to fit into the IFAK
- Includes grommet reinforced punch with latex-free band for attachment
- Currently in use by Special Operations Forces

SPECIFICATIONS

H 6 in. x W 4.5 in. each
Weight: 0.22 oz ea | 5.5 oz 25-Pack

COMBAT CASUALTY CARD

UNIT	ITEM #	NSN #
25 Pack	ZZ-0262	6515-01-587-1700



COMBAT CASUALTY REFERENCE CARD

- Rugged reference card for emergency management featuring the current TCCC Guidelines
- Clear, concise instructions for medics and warfighters
- Information divided into four sections:
 - Care Under Fire
 - Tactical Field Care
 - TACEVAC Protocols
 - 9-line MEDEVAC request/notes
- Pocket sized for convenience
- Twenty-two pages

SPECIFICATIONS

H 6 in. x W 5 in. x D 0.1 in.
Weight: 2.3 oz

REFERENCE CARD

ITEM#	NSN #
ZZ-0057	6515-01-537-4161



R-SICK *Rapid Surgical Intervention Convenience Kits*

- Flexible austere and prolonged field care surgical instrument storage and recovery system tray
- For both major and minor surgical procedures in less than ideal or highly mobile environments
- Sequentially set up for rapid "incision, spreading, damage control, and closure"
- Ideal for complex, emergent procedures such as surgical airways, hemorrhage control, field amputation, and REBOA

R-SICK

INSTRUMENTS	ITEM #
Floor Grade	83-0016
OR Grade	83-0017
Bag Only	80-0871

SPECIFICATIONS

Packaged: H 8 in. x W 7 in. x D 3.5 in.
Fully deployed (flat):
L 24 in. x W 43 in. x D 0.1 in.
Weight: Kits 3 lb • Bag Only 11 oz



Instruments are available in Floor & OR Grade

R-SICK CONTENTS:

2 x Weitlaner Retractor, 3 x 4 Blunt (5.5 in.)
1 x Trousseau Tracheal Dilator (5.25 in.)
2 x Adson Forceps, 1 x 2 (4.75 in.)
1 x Ewald Forceps 1 x 2
1 x Plain Splinter Forceps (4.5 in.) Serr. Handle
1 x Kelly Forcep Heavy
1 x Converse Needle Holder (4.24 in.) Long
1 x Olsen-Hegar Needle Holder (7.25 in.) Serr
1 x Stevens Tenotomy Scissors Trbl/Bl

1 x Stevens Tenotomy Scissors Cvdbl/Bl
1 x Scissors, Mayo (5.5 in.) Str
1 x Mayo Dissecting Scissors, Crv (5.5 in.)
1 x Knife Handle (#3)
2 x Gigli Saw Blade (20 in.)
2 x Gigli Saw Handle (3.5 in.) Loop
1 x Hiupp Tracheal Hook, (6.5 in.)
2 x Forcep, Rochester-Pean (7.25 in.) Curv
6 x Scalpel #10 w/Handles
6 x Scalpel #11

6 x Scalpel #15
6 x Blade, #10 Pluss/S(Tefl)
10 x Blade, #15 Pluss/S(Tefl)
2 x Suture Ethilon 4-0 Blk Mono Fs-218
2 x Suture 3-0 Silk
2 x Black Braided Silk Suture, 2-0/2-60 in.
2 x 2-0 Suture Ct-2
4 x 4-0 Silk Suture Blk (30 in.)
1 x Sharps Shuttle
1 x Farabeuf Retractor

FBTK *Field Blood Transfusion Kits*



- Comprehensive kit that provides all the needed supplies to allow the user to rapidly perform a field blood transfusion
- Equipped to allow for verifying the blood types of the donor and recipient and collecting blood samples from the donor for additional testing

FBTK Basic:

Donor Module:

5 x Alcohol Prep Pad
2 x 2-pack Gauze Pad - 4 in. x 4 in.
1 x IV Constricting Band - 1 in. x 18 in.
1 x Clear Adhesive Dressing - 4 in. x 4.75 in.
1 x Blood Bag Label

Recipient Module:

5 x Alcohol Prep Pad
1 x IV Catheter - 18 ga x 1.25 in.
2 x 2-pack Gauze Pad - 4 in. x 4 in.
1 x IV Constricting Band - 1 in. x 18 in.
1 x Clear Adhesive Dressing - 4 in. x 4.75 in.
1 x IV Catheter - 16 ga x 1.25 in.

Additional Supplies:

1 x roll 1 in. Surgical Tape
2 x pair Black Nitrile Gloves, Large
1 x Black Permanent Marker
1 x Kelly Curved Hemostat - 5.5 in.
1 x Single Blood Bag (450ml)
1 x 550 Cord - 9.5 in.
1 x Beaded Security Tie - 11 in.
1 x IV Set Y-Type with Filter 170-260 UM
1 x Questionnaire Card - DD572 Modified

FBTK Blood Surveillance Module:

Donor Module:

1 x IV Hypodermic Needle - 18 ga x 1.5 in.
1 x IV Syringe, 10cc
1 x Eldon Card
1 x Interlink Injection Site Needle
1 x Blood Transfer Device
2 x Vacutainer, Gold Top (5ml)
1 x Vacutainer, Purple Top (5ml)
2 x Vacutainer, Pearl Top (5ml)

Recipient Module:

1 x 3 Way Stopcock
1 x Eldon Card
1 x Interlink Injection Site Needle

*FBTK Complete:

Contains all components of *FBTK Basic* & *Blood Surveillance Module*

FIELD BLOOD TRANSFUSION

OPTIONS	ITEM #
FBTK Basic	83-0022
Blood Surveillance	83-0025
FBTK Complete*	83-0000

For more information about the FBTK, see NARescue.com

PENBLADE

Safety Scalpels

- Pen-like activation that uses existing muscle memory
- Blade retraction without repositioning of the hand
- Suture trimming capability when blade is retracted
- High-quality materials and construction

SPECIFICATIONS

Packaged: L 7.5 in. x W 1.7 in. x D 0.45 in.
Unpackaged: L 6 in. x W 0.7 in. x D 0.4 in.
Weight: 0.4 oz

PENBLADE SAFETY SCALPEL

BLADE	COLOR	ITEM #
#10	Blue	ZZ-0877
#11	Orange	ZZ-0878
#15	Green	ZZ-0879

MINOR SURGERY SET



MINOR SURGERY CONTENTS:

- 1 x Tray, Shallow
- 2 x Kelly Forceps, Straight 5.5 in.
- 2 x Kelly Forcep, Curved 5.5 in.
- 1 x Mayo Hegar Needle Holder, 6 in.
- 1 x Adson Forceps, Serr 4.75 in.
- 1 x Forceps Tissue, 5 in.

- 1 x Collier Needle Holder, Straight 5 in.
- 1 x Handle Surgical Knife #3
- 1 x Probe W. Eye 5 in.
- 1 x O.R. Scissors, Sb St 5.5 in.
- 1 x Wrap, Kimguard 24 x 24 in.


MINOR SURGERY SET

ITEM #	UNIT	NSN#
ZZ-0960	Single	-----
ZZ-0959	Pack of 20	-----

SPECIFICATIONS

Packaged: L 9 in. x W 5 in. x D 1 in.
Weight: 5 oz

SURGICAL SUTURE KIT



SURGICAL SUTURE CONTENTS:


- 1 x Tray, Shallow
- 1 x Adson Forceps w/Teeth
- 1 x Webster Needle Holder
- 1 x Iris Scissors, Straight
- 1 x Mosquito Hemostat, Curved
- 1 x Poly Bag w/Twist Tie
- 1 x Poly-backed Towel
- 1 x Fenestrated Drape
- 1 x Alcohol Prep Swab
- 1 x PVP Prep Swab
- 1 x Cotton Filled Sponge.

SPECIFICATIONS

Packaged: L 9 in. x W 5 in. x D 1 in.
Weight: 6 oz

SURGICAL SUTURE SET

ITEM #	UNIT	NSN#
ZZ-0958	Single	-----
ZZ-0961	Pack of 20	6515-01-311-0360



SUTURE REMOVAL KIT #1



ITEM #

ZZ-0664

SPECIFICATIONS

Pkg: L 6 in. x W 3.5 in. x D 0.38 in.
Weight: 2 oz

STERILE

SUTURE REMOVAL KIT #2



ITEM #

ZZ-0665

SPECIFICATIONS

Pkg: L 6 in. x W 3.5 in. x D 0.38 in.
Weight: 2 oz

STERILE

SUTURE REMOVAL KIT #3



ITEM #

ZZ-0666

SPECIFICATIONS

Pkg: L 6 in. x W 3.5 in. x D 0.38 in.
Weight: 2 oz

STERILE

SUTURE REMOVAL KIT #4



ITEM #

ZZ-0667

SPECIFICATIONS

Pkg: L 6 in. x W 3.5 in. x D 0.38 in.
Weight: 2 oz

STERILE

SURGICAL INSTRUMENTS

Precision standard surgical instruments available individually

STERILE

ADSON FORCEPS

FP, 1x2 • 4.75 in.




ITEM #

ZZ-0705

HALSTEAD MOSQUITO

Curved • 5 in.




ITEM #

ZZ-0695

HALSTEAD MOSQUITO

Straight • 5 in.



ITEM #

ZZ-0694

IRIS SCISSORS

Curved • 4.5 in.



ITEM #

ZZ-0697

IRIS SCISSORS

Straight • 4.5 in.



ITEM #

ZZ-0696

KELLY HEMOSTAT

Curved • 5.5 in.




ITEM #

ZZ-0701

KELLY HEMOSTAT

Straight • 5.5 in.



ITEM #

ZZ-0700

LISTER SCISSORS

5.5 in.



ITEM #

ZZ-0698

O.R. SCISSORS

SB Straight
5.5 in.




ITEM #

ZZ-0699

ROCHESTER PEAN

Hemostat,
Curved
6.25 in.




ITEM #

ZZ-0702

WEBSTER NEEDLE

Holder, 5 in.



ITEM #

ZZ-0703

PAXLIGHT® LED FIELD MEDICAL LIGHT

- Compact, lightweight precision lighting solution with the rugged engineering required for use in field hospitals.
- LED Technology with intense light fifteen times brighter than traditional halogen field lights.
- Patented 3-pod head eliminates shadow interference common with single point light sources.
- Energy efficient, cool white color with an exceptionally long life.
- Stays cool to the touch, projecting no heat on the patient.
- Locking joints make it fully adaptable and secure with adjustable mounting clamp for attachment to hospital bed, litter or litter stand.
- Includes a rugged storage case with wheels for easy transport.

SPECIFICATIONS

Light in case: L 47.2 in. x W 16.5 in. x D 6.7 in.
Light Head: 12.8 in. Diameter
Light Arm/Post: L 78 in. x 2.8 in. Diameter
Weight: Case & Light • 37 lb 8 oz
Light with Arm/Post • 11 lb 8 oz
Light Output: 50,000 lux (4590 foot candles)

PAXLIGHT

ITEM #	NSN #
60-0042	6530-01-549-1089



WHEELED BASE for PAXLIGHT®

- Independent 5 caster wheel base with low center of gravity
- 13.5 lb counter weight for stabilization
- Includes convenient attachment point for power cord & AC adapter



SPECIFICATIONS


Base: W 38 in.
Assembled: H 38 in.
Weight: Base: 13.5 lb
Assembled: 22 lb

WHEELED BASE

ITEM #
60-0045

EZ LIGHT HANDLES for PAXLIGHT®

- Sterile, disposable one-piece light handles to eliminate the time and labor of reprocessing
- Conveniently applies to the PAXLight handle via screw-on attachment for a solid feel



EZ LIGHT HANDLES

UNIT	ITEM #
Pack of 10	60-0046

SPECIFICATIONS

L 4.5 in. x W 1 in. x D 1 in.
Shield Diameter: 4.25 in.



SPECIFICATIONS

Folded: L 42 in. x W 17 in. x D 32 in.
 Opened: L 79.25 in. x W 33 in. x H 35.25 in.
 Weight: 54 lb (Mattress and Frame)
 Weight: 68 lb
 (Mattress, Frame and Accessories)
 U.S. Patents Pending

*Extremely compact
 for easy transport
 and storage*

MARK IV Field Hospital Bed

- Lightweight solution designed for patient comfort in the combat support or field hospital/treatment tent
- Durable & easy-to-open, supporting up to 850 lb of static weight
- Designed to be rapidly deployed with minimal training in less than 60 seconds by a single person
- Constructed of a bi-fold aluminum frame that is adjustable to 6 positions
- 3 in. polyurethane foam pad with water-resistant decontaminable cover easily cleans with 10% bleach solution or other standard hospital cleaning products
- Includes IV Pole and a removable table platform that rests on its adjustable bed rails

MARK IV BED

ITEM #	NSN #
66-0019	6530-01-624-7626

Diagnostic Equipment

FIELD CORPSMAN KITS

Portable ENT diagnostics

- Lightweight, durable ENT diagnostic sets for portable field care of the eyes, ears & nose
- A wide selection of high intensity lights & magnifiers for enhancing both day and night examinations



Additional items in Deluxe Kit only ▼



SPECIFICATIONS

H 8 in. x W 4.25 in. x D 2.25 in.
 Weight: Basic Kit: 6.2 oz
 Deluxe Kit: 11.6 oz
 Ophthalmoscope only: 3.7 oz

FIELD CORPSMAN KITS

KIT OPTIONS	ITEM #	NSN #
Basic	20-0001	6515-01-458-6057
Deluxe (with Heine)	20-0002	6515-01-458-6178
Ophthalmoscope only	20-0007	6540-01-491-6651

Hyperthermia

Introducing the latest
Cooling Technology for effective
Human Thermo-Regulation

POLAR SKIN™ Cooling Systems

- Polar Skin Mini-Mesh,™ a unique smart textile with a hydrofunctional polymer that allows for liquid absorptive capacity 350% greater than base weight
- ArcticAqua,™ a proprietary solution that prevents rigidity at temperatures well below 0°F while providing the requisite moisture to support evaporative, convective and conductive cooling
- Adheres to current U.S. Army and National Athletic Trainers' Association protocols for heat illness treatment
- Protracted direct cooling of high blood flow regions

ICE SHEETS *Standard*

- Principle cooling component of Polar Skin™ cooling systems
- Adheres to current National Athletic Trainers' Assoc and U.S. Army protocols for heat illness treatment
- Sheet Size: W 36 in. x L 30 in.



STANDARD ICE SHEETS

UNIT	ITEM #	NSN #
Single Pack	ZZ-0668	6532-01-653-6054
Double Pack	ZZ-0669	----
Triple Pack	ZZ-0755	----

SPECIFICATIONS

Single: L 4.5 in. x W 77.5 in. x H 0.5 in.	Wgt: 9.6 oz
Double: L 9.5 in. x W 8 in. x H 0.5 in.	Wgt: 1 lb 6.4 oz
Triple: L 12.5 in. x W 7.5 in. x H 0.5 in.	Wgt: 2 lb 3.2 oz

ICE SHEETS *Large*

Same characteristics as the standard size Ice Sheets, but in a larger 36 in. x 60 in. size

SPECIFICATIONS

L 8.5 in. x W 5.5 in. x H 1.5 in.
Weight: 1 lb 8 oz

LARGE ICE SHEETS

ITEM #
ZZ-0758



C2E ICE PACK BUNDLE

- Four Comfort Gel Packs
- Facilitates Conductive, Convective, and Evaporative (C2E) cooling of high blood flow (large artery) regions of the body

SPECIFICATIONS

Indiv. Gel Pk: W 6 in. x L 10 in.
Pkg: L 8 in. x W 9.5 in. x H 1.25 in.
Weight: 4 lb 4 oz

C2E ICE PACK

ITEM

ZZ-0675

COMPREHENSIVE COOLING SYSTEM

- Highly-effective portable cooling system for the treatment of heat-related illnesses
- Scientifically shown to reduce core temperature by .09°F/minute in 100°F heat when used according to protocol

SPECIFICATIONS

L 11 in. x W 10 in. x H 4 in.
Weight: 8 lb 5.3 oz

COMPREHENSIVE

ITEM

ZZ-0674

COMPLETE CARE PACK

- Portable, highly insulated option
- Contains one Polar Skin Comprehensive Cooling System and one heavy-duty, leak & sweat proof cooler that keeps ice from melting for 24 hours in +100°F heat

SPECIFICATIONS

W 14 in. x L 12 in. x D 7 in.
Weight: 10 lb 3.2 oz

COMPLETE CARE

ITEM

80-0451

COMPLETE CARE BACKPACK

Available in either Single Patient (one comprehensive cooling system) or Double Patient (two cooling systems) configurations

BACKPACK

SIZE	ITEM #
Single	80-0543
Double	80-0544

SPECIFICATIONS

W 19.5 in. x L 15 in. x D 7.5 in.
Single Weight: 12 lb 13 oz
Double Weight: 21 lb 13 oz



COOLING COLLAR

- Two gel packs encased in integrated pouches
- Designed to align with the carotid artery

SPECIFICATIONS

Indiv. Gel Pk: W 6 in. x L 10 in.
Pkg: W 8.5 in. x L 5.25 in. x D 2 in.
Weight: 2.5 lb

COOLING COLLAR

ITEM

ZZ-0670



Standard size cooler accommodates one Polar Skin™ Comprehensive Cooling System.

Large cooler and Backpack accommodates two cooling systems

Available in Standard, Large & Backpack COOLERS

SPECIFICATIONS

Standard:

L 14 in. x W 12 in. x H 7 in. • Weight: 1 lb 3.2 oz

Large:

L 18.5 in. x W 13 in. x H 9.5 in. • Weight: 2 lb 13 oz

Backpack:

L 19.5 in. x W 15 in. x D 7.5 in. • Weight: 3 lb 13 oz

POLAR SKIN COOLERS

COLOR	SIZE	ITEM #
Red	Standard	ZZ-0682
Black	Standard	ZZ-0759
Black	Large	ZZ-0756
Black	Backpack	ZZ-0794

TACTICAL TRACTION SPLINT

(TTS) Compact Mechanical Traction Device

SPECIFICATIONS

Packaged:
L 11 in. x W 5.5 in. x D 3 in.
Weight: 1 lb 5 oz

TACTICAL TRACTION SPLINT

BAG COLOR	ITEM #	NSN #
Black	50-0015	6515-01-346-9186
Coyote	50-0033	----
Digital Camo	50-0034	----
OD Green	50-0032	----

- Relieve pain, increase patient comfort and reduce the risk of vascular injury or nerve damage in femur fractures
- Lightweight and portable kit with clearly marked elastic straps for the proper fit to the thigh, knee & shin

TACTICAL TRACTION SPLINT CONTENTS:



SAM[®] SPLINT II

Malleable Splint

- Versatile, malleable splint easily molds to fit most fractures
- Lightweight, comfortable and strong
- Waterproof, Reusable
- Easy to clean, Compact
- Illustrated step-by-step instructions included

SAM SPLINT

ITEM #	NSN #
50-1005	6515-01-494-1951

SPECIFICATIONS

Packaged: H 9.25 in. x W 4.25 in. x D 1 in.
Unfolded: L 36 in. x W 4.25 in.
Weight: 4.5 oz

NAR Elastic Wrap BANDAGES

- Latex-free Elastic Wrap Bandage with Hook Closure
- Length: 5 yd

ELASTIC WRAP

WIDTH	ITEM #	NSN #
6 in.	30-0007	6510-01-532-6522
4 in.	30-0026	6510-01-171-5130
3 in.	30-0028	6510-01-230-8702
2 in.	30-0027	6510-00-935-5820

TC SPLINT

Telescoping Collapsible Splint

- Multi-capable device for splinting a wide variety of major and minor fractures
- Lightweight and compact for easy field carry with 6 easily adjustable interconnected telescoping poles and 4 wide elastic bands with hook and loop securing
- Includes clearly marked directions on the device as well as an easy-to-use, step-by-step instructional insert



online video
available

< scan to view

TC SPLINT

ITEM #
50-0039

SPECIFICATIONS

Pkgd: L 8 in. x W 5 in. x D 3 in.
Open: W 9.5 in. x L 34 in.
Weight: 1 lb 4 oz

SAM[®] PELVIC SLING

for Pelvic Ring Fractures

- Force-controlled, ideal compression every time without the guesswork
- Rapid, simple, three-step application, usually in less than one minute

SPECIFICATIONS

H 7.25 in. x W 7.75 in. x D 4 in.
Hip Circum: 32 in. to 50 in.
Weight: 9.8 oz



SAM PELVIC SLING

ITEM #	NSN #
50-0004	6515-01-509-6866



BURNTEC® Hydrogel Dressings

- Sterile non-stick hydrogel dressings for treating 1st, 2nd & 3rd-degree burns (thermal, chemical, electric)
- Also treats other skin injuries, such as sunburns, abrasions, skin irritations, bruises and swelling (insect bites)
- Applied dressings provide a cool, soothing effect for up to 24 hours
- Easy to handle without spilling or damaging the hydrogel
- Permeable to water vapor and oxygen but impermeable to bacteria
- Biocompatible and hypo-allergenic; Reduces risk of hypertrophic scars

BURNTEC BURN DRESSINGS

SIZE	ITEM #	NSN #
2 in. x 2 in.	30-0114	6510-01-631-1913
2.5 in. x 5 in.	30-0115	6510-01-631-2080
4 in. x 4 in.	30-0110	6510-01-631-1883
5 in. x 5 in.	30-0116	6510-01-631-2138
5 in. x 10 in.	30-0117	6510-01-631-1907
8 in. x 8 in.	30-0118	6510-01-631-1897
8.5 in. x 11 in.	30-0119	6510-01-631-1890
4 in. x 16 in.	30-0111	6510-01-631-2124
8 in. x 16 in.	30-0112	6510-01-626-2990
16 in. x 24 in.	30-0120	6510-01-631-1886
10 in. x 10 in. Face	30-0121	6510-01-627-1834
12 in. x 16 in. Face	30-0113	6510-01-627-1855
16 in. x 24 in. Face	30-0122	----

Dry Sterile BURN DRESSING

- Dry sterile, non adherent burn dressing & bandage
- Includes two No. 3 closed safety pins and rubber band individually packaged



SPECIFICATIONS

Bandage: 45 in. x 45 in. x 63 in.
Pkg: W 4.5 in. x H 4 in. x D 0.5 in.
Weight: 1 oz

DRY STERILE BURN

ITEM #	NSN #
30-0162	6510-01-587-6579



BurnTec® BURN DRESSING KIT

- Convenient array of BurnTec® dressings in a handy carrying case
- Kit contains (1) each of the following dressings:

2 in. x 2 in.	5 in. x 10 in.	8 in. x 16 in.
2.5 in. x 5 in.	8 in. x 8 in.	10 in. x 10 in. (face)
4 in. x 4 in.	4 in. x 16 in.	12 in. x 16 in. (face)
5 in. x 5 in.	8.5 in. x 11 in.	

BURNTEC BURN KIT

ITEM #	NSN #
80-0361	6515-01-630-9177

SPECIFICATIONS

L 12 in. x W 15.5 in. x D 3.5 in.
Weight: 6 lb 4 oz



online video available
A scan to view

BurnTec® MINOR BURN DRESSING KIT

- Variety of BurnTec® dressings in a handy carrying case
- Kit contains the following:
 - Two 4 in. x 4 in. dressings
 - Two 4 in. x 16 in. dressings
 - One 8.5 in. x 11 in. dressing

BURNTEC MINOR BURN KIT

ITEM #	NSN #
80-0415	6510-01-644-5427

SPECIFICATIONS

L 14 in. x W 8.5 in. x D 2.5 in.
Weight: 2 lb 13 oz



Eye Care

PES Polycarbonate Eye Shield

- Meets CoTCCC treatment protocols for penetrating eye trauma
- Convex eye shield to accommodate either eye
- Translucent orange shell features perforations to allow air circulation

SPECIFICATIONS

3 in. x 2.38 in.
Weight: 0.2 oz

PES EYESHIELD

UNIT	ITEM #	NSN #
Single	30-0142	6515-01-647-5456
Box of 100	30-0143	6515-01-630-5551



FOX Eye Shield

- Convex rigid eye shield made of lightweight aluminum with perforations allowing air circulation
- Meets CoTCCC treatment protocols for penetrating eye trauma

SPECIFICATIONS

3 in. x 2.38 in.
Weight: 0.2 oz

FOX EYESHIELD

OPTION	UNIT	ITEM #	NSN #
No Garter	Box of 12	30-0144	6515-01-449-1016
Garter	Single	30-0055	----
Garter	Box of 50	30-5550	6515-01-598-1877



Trauma Shears & Accessories

Responder 7.25 in. TRAUMA SHEARS

Large, premium-grade Trauma Shears designed for safety and durability while cutting through tough materials

SPECIFICATIONS

L 7.25 in. x W 3.7 in. x D 0.4 in.
Weight: 2 oz

RESPONDER SHEARS

ITEM

ZZ-1019



NAR TRAUMA SHEARS Black

- Premium-grade stainless steel with precision sharpness
- Unique all-over black color

SPECIFICATIONS

SMALL: L 6.25 in. x W 2.75 in. x D 0.3 in.
Weight: 1.5 oz

LARGE: L 7.25 in. x W 3.7 in. x D 0.4 in.
Weight: 2 oz

TRAUMA SHEARS

SIZE	ITEM #	NSN #
S 6.25 in.	ZZ-0064	6510-01-622-1624*
L 7.25 in.	ZZ-0063	6515-01-538-9276



*NATO NSN

6 oz MOLLE Retractable GEAR LEASH

- MOLLE Mount compatible with 360°rotating
- High break strength
- Salt water resistant
- Self-Flushing
- 36 in. Extension / Spectra Nylon Line
- Stainless Steel Spring & Hardware
- Quick Connect/Disconnect

SPECIFICATIONS

Retracted:
H 3.5 in. x W 1.6 in. x D 1.0 in.
Cable Length - 36 in.
Weight: 20 oz

GEAR LEASH

ITEM

70-0155



SCISSOR LEASH™

- Salt water resistant
- Self cleaning
- Helps keep track of your shears

SCISSOR LEASH

ITEM #	NSN #
70-0005	6515-01-540-7226

SPECIFICATIONS

H 6 in. x W 1.6 in. x D 0.6 in.
Cable Length - 36 in.
Weight: 16 oz



Leatherman®

RAPTOR™ MULTI-TOOL

High-performance multi-tool equipped with a variety of mission critical tools essential to emergency medical practitioners in the field

1. 420HC Stainless Steel Folding Medical Shears

2. Strap Cutter

3. Oxygen Tank Wrench

4. Ring Cutter

5. Carbide Glass Breaker

6. 5 cm Ruler

7. Holster with Clip

RAPTOR MULTI-TOOL

HANDLE	ITEM #	NSN #
Black	70-0014	5110-01-626-8766
Blue	70-0264	----
Coyote	70-0265	----
Red	70-0266	----



SPECS

Shear Length: 19 in.
Closed Length: 5 in.
Weight: 5.8 oz

HANDLE COLORS:



Leatherman® SUPERTOOL 300

Heavy-duty, all purpose multi-tool featuring 19 devices in one configuration

SPECIFICATIONS

Closed: L 4.5 in.
Primary Blade: L 3.2 in..
Weight: 9.6 oz

SUPERTOOL 300

ITEM

ZZ-0542

- 1. Needlenose Pliers
- 2. Regular Pliers
- 3. Replaceable Wire Cutters
- 4. Hard-Wire Replaceable Cutters
- 5. Stranded-Wire Cutters
- 6. Electrical Crimper
- 7. Wire Stripper
- 8. 420HC Knife
- 9. 420HC Serrated Knife
- 10. Saw
- 11. Awl with Thread loop
- 12. 9 in. Ruler
- 13. Can Opener
- 14. Bottle Opener

- 15. Wood/Metal File
- 16. Phillips Screwdriver
- 17. Large Screwdriver
- 18. Medium Screwdriver
- 19. Small Screwdriver



Leatherman® MUT®

Tactical and practical multi-tool for military, law enforcement and civilian shooters - features 16 devices



- 1. Needlenose Pliers
- 2. Regular Pliers
- 3. Replaceable Wire Cutters
- 4. Hard-Wire Replac. Cutters
- 5. Stranded-Wire Cutters
- 6. Electrical Crimper
- 7. 420HC Knife
- 8. Saw
- 9. Replaceable Cutting Hook
- 10. Hammer
- 11. Bolt Override Tool
- 12. Bronze Carbon Scraper
- 13. #8-32 Cleaning Rod/Brush Adapter
- 14. Firearm Disassembly Punch
- 15. Carabiner/Bottle Opener Punch
- 16. Large Bit Drive

SPECIFICATIONS

Closed: L 5 in.
Primary Blade: L 3 in.
Weight: 11.2 oz

MUT

ITEM

ZZ-0539



BENCHMADE® KNIVES & CUTTERS



1 HOOK SAFETY CUTTERS by Benchmade®

- Hook cutters designed to quickly and reliably cut through a wide array of materials
- Extensively tested and widely used by the US Military for medical and egress purposes

ITEM #	STYLE	NSN #
ZZ-0652	7 Hook	----
ZZ-0250	8 Hook	----

SPECS

7 Hook: L 4.3 in.
8 Hook: 6.5 in.
Weight: 1.5 oz

2 SOCP RESCUE TOOL by Benchmade®

- Fixed 440C high-chromium stainless steel blade hook for rescue performance and less-than-lethal self-defense
- Includes carbide glass breaker, oxygen wrench, single-piece handle with finger loop & molded sheath
- Designed by Special Operations Combatives Program (SOCP)

SPECS

Blade: 2.5 in. Overall: 6.75 in.
Weight: 188 oz

ITEM #	NSN #
ZZ-0760	----

3 H₂O 112 SBK by Benchmade®

- N680 ultra stainless steel blade with opposing fixed bevel blade style, blunt tip, & integrated rescue hook
- Boltaron® sheath designed to fit current Navy SAR kits
- 90-0078 (Orange handle) includes carbide glass breaker & is approved for Navy SAR - "So Others May Live"

ITEM #	HANDLE	NSN #
90-0078	Orange	1095-01-677-1013
90-0097	Black	----

SPECS

Blade: 3.5 in.
Overall: 8.02 in.
Weight: 3.2 oz

4 TRIAGE 916 SBK by Benchmade®

- Dynamic triple utility tool featuring three essential rescue tools in a single compact design
- N680 high corrosion resistant, partially serrated, opposing bevel blade, 440C folding hook & carbide glass break

SPECS

Blade: 3.5 in
Overall: 8.2 in.
Weight: 5.1 oz

ITEM #	HANDLE	NSN #
90-0096	Orange	----
90-0060	Black	----

5 AUTOTRIAGE 9160 SBK by Benchmade®

SPECS

Blade: 3.58 in. Overall: 8.31 in.
Weight: 5.75 oz

ITEM #	NSN #
90-0061	----

- Packaged in high-strength aluminum with AXIS® auto pull release and integrated safety
- Auto opening N680 highly corrosion resistant blade with modified drop-point

6 TRIAGE 917 BK -or- SBK by Benchmade®

- Black drop point style CPM-S30V steel blade with option of a smooth (BK) or partially serrated edge (SBK)
- 440C folding hook & carbide glass breaker

SPECS

Blade: 3.5 in Overall: 8.3 in.
Weight: 5.24 oz

ITEM #	STYLE	NSN #
90-0092	BK	----
90-0093	SBK	----

ARMADILLO Medication & Storage Case

- Compact solution for keeping critical medications safe, secure and easily accessible
- Features Automatic Pressure Equalization Valve
- Hardware constructed of Marine grade 316 SS
- Case is water-resistant and crush-proof
- Hook & loop attachment system allows storage compartment to be tailored to hundreds of custom configurations

SPECIFICATIONS

Ext: L 7.95 in. x W 5.1 in. x D 2.14 in.
Int: L 6.45 in. x W 3.87 in. x D 1.75 in.
Weight: 10.5 oz

ARMADILLO MEDICATION CASE * NATO NSN

LID STYLE	ITEM #	NSN #
Black	ZZ-0081	8145-01-573-2533
Clear	ZZ-0259	6545-01-594-9399*

also available: Half Inch Hook Discs with pressure sensitive adhesive
Set of 250 • Item# ZZ-0258 • NSN# 8145-01-617-9075*

Clear Lid
ZZ-0259 only
includes Loop Liner
on the bottom

I&A CASE Injury & Ailment Case

- Superbly engineered medical equipment organization system
- Lid's clear pockets and hook & loop fasteners allow secure storage as well as superior inventory management
- Polymer O-ring provides dust and water-proof seal

- Automatic Gore-TEX® pressure equalization valve stops moisture from entering the case – also prevents vacuum-lock
- Double-throw latches engineered to tolerate up to 400 lb of pressure
- All-weather, sure-grip handle and reinforced padlock protectors

SPECIFICATIONS

H 8.5 in. x W 20.5 in. x D 17 in.
Wgt: 13 lb 12 oz

I&A CASE

ITEM

ZZ-0222

M.U.L.E. Multi-Use Luggage & Equipment Carrier

- Heavy-duty, water-resistant nylon platform with large compartment and mesh pockets on the inside lid for maximum storage of critical items
- Outside end pockets feature open channels for expanded storage capability
- Suitcase-style carrying handles on top and side; Telescoping handle and wheels for easy transport
- Hidden rucksack shoulder pads

SPECIFICATIONS

H 30.5 in. x W 18.5 in. x D 11 in.
Wgt: 11 lb 10 oz | Payload: 70 lb

M.U.L.E.

ITEM

ZZ-0212

M.U.L.E. II Multi-Use Luggage & Equipment Carrier II

- Extra large, heavy-duty, water-resistant 1000 denier nylon platform featuring 9,000 cubic in. of storage with large reinforced open compartment and multiple zippered mesh pockets for maximum storage
- Top & side web carry handles and 3 durable 3.75 in. diameter wheels for easy transport

Available in Black & Coyote

M.U.L.E. II

COLOR	ITEM #
Black	ZZ-1157T
Coyote	ZZ-1158T

SPECS

Exterior: H 40 in.
x W 21 in. x D 15 in.
Interior: H 40 in. x W 15 in. x D 15 in.
Weight: 15.4 lb
Capacity: 9,000 cubic in.

Nitrile Gloves

Bulk Box | Pre-packaged Pairs NITRILE GLOVES

- Form-fitting, ambidextrous design with textured fingertips
- Thick palm ensures durability
- Textured for improved wet and dry gripping
- Longer and thicker than standard nitrile gloves
- Nitrile gloves available in Black, Tan and Blue colors and come either as pre-packaged pairs (kits) or in bulk (box)

SPECIFICATIONS

Kit (25 pair): L 4.75 in. x W 6 in. x D 4 in.
Weight: 1 lb 1 oz

Box (Bulk): L 11 in. x W 5 in. x D 3.75 in.
Weight: 2 lb 4 oz

BEAR CLAW®



COLOR:
Tan/Sand

Bear Claw BULK

*NATO NSN

SIZE	PKG QTY	ITEM #	NSN #
S	Box 50 pair	70-0009	6515-01-600-6059*
M	Box 50 pair	70-0010	6515-01-600-6057
L	Box 50 pair	70-0011	6515-01-600-6039
XL	Box 45 pair	70-0012	6515-01-600-6052*

Bear Claw PAIRS

*NATO NSN

SIZE	PKG QTY	ITEM #	NSN #
S	Kit 25 pair	ZZ-0246	6515-01-619-5350
M	Kit 25 pair	ZZ-0247	6515-01-619-5689
L	Kit 25 pair	ZZ-0248	6515-01-600-0262
XL	Kit 25 pair	ZZ-0249	6515-01-619-5367

COLOR:
Tactical
Black



BLACK TALON®

Black Talon BULK

*NATO NSN

SIZE	PKG QTY	ITEM #	NSN #
S	Box 50 pair	70-0001	6515-01-515-0195
M	Box 50 pair	70-0002	6515-01-515-0200
L	Box 50 pair	70-0003	6515-01-515-0197
XL	Box 45 pair	70-0004	6515-01-515-0190

Black Talon PAIRS

*NATO NSN

SIZE	PKG QTY	ITEM #	NSN #
S	Kit 25 pair	ZZ-0218	----
M	Kit 25 pair	ZZ-0217	6515-01-672-7977*
L	Kit 25 pair	ZZ-0216	6515-01-632-0059
XL	Kit 25 pair	ZZ-0215	6515-01-672-7985*

Blue RESPONDER

COLOR:
Blue



Blue Responder PAIRS

SIZE	PKG QTY	ITEM #
S	Kit 25 pair	ZZ-1059
M	Kit 25 pair	ZZ-1060
L	Kit 25 pair	ZZ-1043
XL	Kit 25 pair	ZZ-1061

Third-party tested compliant with the EN ISO 374-5:2016 Standard for **Virus Protection** against **bacteria, fungi and viruses** including, but not limited to:

- ▶ **COVID-19** (Coronavirus)
 - ▶ **HIV** virus
 - ▶ **Ebola** virus
- ▶ **Avian influenza** virus (Bird Flu)
 - ▶ **Hantavirus** [causes Hantavirus Pulmonary Syndrome (HPS)]
 - ▶ **Marburg** virus
 - ▶ **Rabies** virus
 - ▶ **Smallpox** virus
- ▶ **Influenza** viruses (A, B and C)
 - ▶ **Rotavirus**
 - ▶ **HPV** virus

Pre-packaged Pairs

BLACK RESPONDER NITRILE GLOVES

- 25 individually wrapped pairs for convenient rapid deployment
- High performance personal protection compliant with the EN ISO 374-5:2016 Standard for Virus Protection
- Third-party tested for effective personal protection from Fentanyl and Viruses including COVID-19 Coronavirus
- Latex- and Powder- Free with decreased risk of allergies
- Biodegradable and disposable with second-skin fit & feel
- Available in 6 sizes: XS, S, M, L, XL, & XXL

SPECIFICATIONS

Gloves: Length: 9.5 in.
Thickness (mils): 4.0
Cuff: Rolled
Packaged (25 pairs):
L 4.75 in. x W 6 in. x D 4 in.
Individual pair:
L 2.5 in. x W 1 in. x D .8 in.
Weight: 16.1 oz

Black Responder PAIRS

SIZE	PKG QTY	ITEM #
XS	Kit 25 pair	ZZ-1285
S	Kit 25 pair	ZZ-1286
M	Kit 25 pair	ZZ-1287
L	Kit 25 pair	ZZ-1288
XL	Kit 25 pair	ZZ-1289
XXL	Kit 25 pair	ZZ-1290

COLOR:
Tactical
Black

Personal
Protection from
Fentanyl
and
COVID-19



Accessories

MEDIC CROSS ID PATCHES PVC

- Soft flexible PVC patch
- Backed with hook material for easy attachment and removal
- 1 in. x 1 in. available online

MEDIC CROSS PATCHES

COLORS	ITEM #
Coyote & Black	ZZ-0797
Black & Gray	ZZ-0798
OD Green & Black	ZZ-0800
Black & Orange	ZZ-0799

SPECIFICATIONS

L 2 in. x W 2 in. x D 0.125 in.
Weight: 0.6 oz

IR FLAGS 2 Pack

- Clearly visible under IR light
- Black hook fasteners on the back for easy hook-and-loop attachment and removal
- 2 flags per pack

IR FLAGS

COLOR	Stars	ITEM #
Full Color	L	ZZ-0093
Coyote	L	ZZ-0094
Coyote	R	ZZ-0235
OD Green	L	ZZ-0107

SPECIFICATIONS

L 2 in. x W 2 in. x D 0.125 in.
Weight: 0.6 oz

6-Pack of 1 in. IR Squares: #ZZ-0108



BLEEDING CONTROL PATCH

ITEM#

ZZ-0796

SPECIFICATIONS

L 4.375 in. x H 1.75 in. x D 0.125 in. | Weight: 0.275 oz

BLEEDING CONTROL PATCH PVC

- Soft flexible PVC rubber patch for labeling your medical gear
- Ideal for the Public Access Bleeding Control Kit or any surface with an area for hook-&-loop attachment

Luminous MED ID PATCH

- Soft flexible custom PVC with hoop and loop backing
- Ultimate glow-in-the-dark ID
- 1 in. block letters
- 2 pack *Glow in the Dark*

MED ID PATCHES

ITEM#

ZZ-0121

SPECIFICATIONS

Each:
L 3 in. x W 2 in. x D 0.125 in.
Weight: 0.4 oz



STAY informed

and KEEP IN TOUCH



Keep up with the latest North American Rescue® news and happenings!

Join us on Facebook:
facebook.com/NARescueLLC



Providing medical solutions and celebrating the life-savers who use them!

Join us on Instagram:
instagram.com/northamericanrescue



Keep up with North American Rescue® with 140 characters or less!

Join us on Twitter:
twitter.com/NARescue



View North American Rescue® instructional and informational videos!

Join us on You Tube:
youtube.com/NARescue01



Visit North American Rescue® online to sign up to Join the MISSION!

Sign up at:
www.NARescue.com



Office Hours: M-F 8:00am - 7:00pm ET
Call toll-free 888.689.6277 • www.NARescue.com

CUSTOMER SUPPORT

- **COMMITMENT:** We pride ourselves in providing an exceptional customer experience beginning with your initial inquiry up to order placement and beyond. Our customer service staff members are trained by certified instructors on all NAR products and are standing by to help with any questions you may have.
- **SUPPORT:** Our in-house technical support professionals consist of an elite group of Nationally Registered Paramedics and Special Operations Medics from military, law enforcement, tactical law enforcement, fire-rescue and EMS backgrounds.

ORDERING

- **ONLINE:** www.NARescue.com

SPECIAL PRICING: Military customers will automatically be eligible for discounted pricing when you create an online account with your .mil email address. Register and respond to the verification email and you'll immediately be ready to go!

Distributors or federal, state, and local agencies may also qualify for discounts off of retail prices. Submit the NAR Pricing Application online or call Customer Service at 888.689.6277 to shop online with your qualified pricing discount.

ESTABLISHED CUSTOMERS: Already have an account, but want to shop online with your special pricing? Contact Customer Service at 888.689.6277 to get your account enabled for online convenience. View your invoices, order history, order status and more!

NEW CUSTOMERS: Simply create an online account TODAY and contact us to see if you qualify for a pricing discount. If so, you will have access to customer specific pricing – anytime, anywhere you are logged in to NARescue.com

QUOTE REQUESTS: Online accounts will not have to wait for custom pricing to place an order. Shop like normal and then convert your cart to a Request For Quote with the click of a button! View our FAQ for details.

- **CONTRACT:**

Tax ID: 27-1024029	ECAT: dmmonline.
Cage Code: 06ST7	dscp.dla.mil
DUNS: 83-24-26782	Prime Vendor:
CCR Registered: Yes	dmmonline.dscp.dla.
ORCA Reg: Yes	mil/primevendor/
	primehome
- **EMAIL:** info@narescue.com
- **PHONE:** 888.689.6277
- **ADDRESS:** 35 Tedwall Ct | Greer, SC
- **FAX:** 864.675.9880

SHIPPING

Need it ASAP? Ask Customer Service about overnight or second day delivery.

Customizable Solutions built to meet YOUR requirements

If you haven't found exactly what you need, our team of subject matter experts at North American Rescue are standing by to guide you to the ideal lifesaving solution tailored to meet all your specifications.

Whether it's the essential gear needed for individuals or customized kits for multiple casualties, an Unconventional Medic® Specialty Kit is your answer to the precise aid kit you and your organization need to increase survivability should injury occur.

Call or Email us TODAY to explore your options for the tailored lifesaving solution YOU can live with!



High
Quality
Tactical
Nylon

-OR-

Vacuum
Bag



LIFE SAVING ► MISSION CRITICAL ► CUSTOM BUILT



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
— WITH A —
MISSION®

35 TEDWALL COURT | GREER SC 29650 | www.NARESCUE.COM