

OSL™ OVERSIZED LITTER *Specialty Litter*

REV012219



Specialty litter designed for safely evacuating larger patients...

Designed to safely evacuate larger patients with both comfort and ease. This design has a 30% increase in patient surface area compared to standard-sized evacuation platforms, extra padded carrying handles, flame-retardant cover, and additional stirrups offering firmer patient support as well as better weight distribution on the floor.

Like the Raven 90C, the OSL has a bi-fold configuration. The cover material is a very durable Nomex/Kevlar blend.

Also, a supplemental Back Rest (not shown) to support the patient at both 90° and 45° angles may be purchased separately.

Specifications:

- Folded: H 45 in. x W 9 in. x D 8 in.
- Open: L 90 in. x W 31 in. x D 6 in.
- Weight: 22 lb 6 oz (10.2 kg)
- Static Payload: 2000 lb (909 kg)
- Ground Clearance: 1.5 in.

Features:

- Certified safe-to-fly by the U.S. Air Force
- Flame-retardant cover made of a durable Nomex/Kevlar-blend fabric
- 6 stirrups for firmer patient support and better weight distribution
- Ergonomically designed handles
- Additional carrying loops
- 6 IV attachment points
- 3 patient securing straps
- Attached storage pouch, Bifold configuration for easy storage
- U.S. Patent 5,598,592, Patent(s) Pending



*Designed with
6 stirrups for firmer
patient support and
better weight distribution*

OSL OverSized Litter



Item#	NSN#
60-0023	6530-01-552-7708



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®