

NAR TRAUMA SHEARS



6.25 in. Small

7.25 in. Large
featuring O₂ wrench

*Unique all-over black
Trauma Shears ideal for the
tactical environment*

Premium-grade stainless steel with precision sharpness...

Built for safety and efficiency, the NAR Trauma Shears from North American Rescue let tactical healthcare professionals get the job done, and get it done right.

These paramedic shears are made from premium-grade stainless steel. Their serrated lower edge cuts through the toughest materials. Their soft-grip, non-slip contoured handles provide maximum control, even when wet. Their clothing lift lip and rounded tip provide the ultimate in patient safety.

NAR Trauma Shears are made to meet the challenges of casualty care with precision and dependability, their unique all-over black color makes them the shears of choice in tactical environments.

Our medical scissors come in two convenient lengths, 6.25 or 7.25 inches.

Features:

- Premium-grade stainless steel with precision sharpness
- Enhanced serrated lower edge cuts through the toughest materials
- Clothing-lift lip and rounded tip for ultimate patient safety
- Soft-grip, non-slip contoured handles for max control even when wet
- Unique all-over black color

Dimensions:

- **6.25 in Shears:**
L 6.25 in.
W 2.75 in.
D 0.3 in.
Weight: 1.5 oz
- **7.25 in Shears:**
L 7.25 in.
W 3.7 in.
D 0.4 in.
Weight: 2 oz

NAR Trauma Shears

Size	Item#	NSN#
7.25	ZZ-0063	6515-01-538-9276
6.25	ZZ-0064	6510-01-622-1624*

* NATO NSN#



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®