



MOUNTAIN LITTER



In tactical situations, evacuating a wounded casualty from mountainous or unlevel terrains can be both challenging and awkward. Packing limitations, when traversing difficult terrain, require lighter and more compact equipment. Taking these constraints into consideration, the NAR Mountain Litter (patents pending) was developed to specifically aid the front line soldier with a solution for rapid packaging and transport of casualties in these extreme tactical environments.

- Bungee is coated to prevent fraying or dry rotting
- Has no more than 1" of deflection while weight is exerted on the litter
- Bed made of lightweight Ripstop 400D
- Slip resistant, ergonomic handles

TALON® ASSAULT LITTER CARRIER



- Transport the Talon® Litter into battle with maximum comfort for the rescuer
- For use by itself or as part of a load carriage system
- Equipment attachment points (MOLLE-style) sewn on exterior to maximize ability to carry additional kit
- Thin back plates & drain hole
- Non-slip, adjustable shoulder straps & top-carrying handle
- 500D CORDURA® Solution-Dyed Nylon with built in Near Infrared Signature Reduction; GhillieTEX™ IR signature reduction hardware



NOTE: Litters sold separately