


OCHO IFAK

REV050622



Available Colors:



Purposely designed to look like a Non-Tactical, Accessory-Style Bag...

The Ocho IFAK from North American Rescue is a very compact, lo-vis individual first aid kit equipped with a 3 in. belt loop for convenient attachment to your duty/pistol belt.

Though the exterior of this small rip-stop nylon case is purposely designed to look like a non-tactical, accessory-style bag, the Ocho IFAK comes fully loaded with critical medical equipment for point-of-wounding treatment of penetrating, blast or other traumatic injuries.

Life support equipment options include Basic, Basic with Combat Gauze, and Advanced configurations. The clam shell design of the platform features rip stop nylon, interior elastic loop organization and easy to open zipper for rapid access to lifesaving equipment. Ocho IFAK bags are also available for individual purchase.

Features:

- Very compact, lo-vis individual first aid kit equipped with a 3 in. belt loop attachment to accommodate a duty/pistol belt
- Contains critical medical equipment for point-of-wounding treatment of penetrating, blast or other traumatic injuries
- Purposely designed to look like a non-tactical, accessory-style bag
- Constructed with rip stop nylon, elastic loop organizational storage and easy open zipper

Dimensions:

- L 7 in. x H 2.25 in. x D 3.25 in.
- Weight: 13 oz

OCHO IFAK

LEVEL	COLOR	ITEM #
Basic	BLK/GRY	80-0374
Basic	BLU/GRY	80-0375
w/Combat Gauze	BLK/GRY	80-0921
w/Combat Gauze	BLU/GRY	80-0922
Advanced	BLK/GRY	80-0372
Advanced	BLU/GRY	80-0373



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®

OCHO IFAK

REV050622

Ocho IFAK Basic



QTY
1 C-A-T®
Tourniquet

QTY
1 Gauze
S-Rolled

QTY
1 ETD™ 4 in.
Mini Responder

QTY
1 HyFin®
Vent Compact, Twin

PAIR
1 Nitrile
Gloves, Lg

Ocho IFAK Basic with Combat Gauze



QTY
1 C-A-T®
Tourniquet

QTY
1 Combat
Gauze Z-fold

QTY
1 ETD™ 4 in.
Mini Responder

QTY
1 HyFin®
Vent Compact, Twin

PAIR
1 Nitrile
Gloves, Lg

Ocho IFAK Advanced



QTY
1 C-A-T®
Tourniquet

QTY
1 Combat
Gauze Z-fold

QTY
1 ETD™ 4 in.
Mini Responder

QTY
1 28F Naso
Airway & Lube

QTY
1 HyFin®
Vent Compact, Twin

QTY
1 ARS®
Needle, 14g

PAIR
1 Nitrile
Gloves, Lg



*Constructed with rip stop nylon,
elastic loop organizational
storage and easy open zipper*



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®