

MARK IV FIELD HOSPITAL BED

MED/SURGE



Designed for Rapid Deployment by a Single Person...

With minimal training the Mark IV Field Hospital Bed can be set up into the operational configuration in less the 60 seconds.

Constructed of a bi-fold 6061 aluminum frame (OD Green), the bed includes a polyurethane foam pad with a water resistant, decontaminable cover.

The bed is manufactured to support a resting force of 850 lb without material fatigue, and is extremely compact for easy transportation and storage.

Cleaning of the Field Hospital Bed is easily done with either a 10% bleach solution or other standard hospital cleaning products. The bed frame and mattress were designed to be extremely functional while taking patient comfort into consideration. The bed frame was designed to evenly distribute the patient's weight across the entire surface area and reduce pressure points.

The 3 in. thick polyurethane mattress provides patient comfort and prevents patient contact with the bed frame to further reduce the possibility of pressure sores over the use of the bed. Both the head and the foot of the bed frame are adjustable to 6 positions.

*Both the head and the foot
of the bed frame are
adjustable to 6 positions.*



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

Mark IV Field Hospital Bed

Item#	NSN#
66-0019	6530-01-549-1089

*With minimal training the Mark IV Field
Hospital Bed can be set up into operational
configuration in less the 60 seconds.*

PRODUCTS WITH A
MISSION®

Features:

- Durable, easy-to-open, Field Hospital Bed that supports up to 850 lb of static weight
- Designed to be rapidly deployed by a single person
- Constructed of a bi-fold 6061 aluminum frame & 3 in. polyurethane foam pad with water-resistant decontaminable cover
- Includes IV Pole and a removable table platform that rests on its adjustable bed rails
- Easily cleaned with 10% bleach solution or other standard hospital cleaning products
- Can be set up to Operational Status in approximately 60 seconds
- Supports a force of 850 lb without material fatigue
- Extremely compact for ease of transport and storage
- Both the head and foot of the bed frame are adjustable to 6 positions
- 3 in. thick polyurethane mattress provides patient comfort

Specifications:

- Folded: L 42 in. x W 17 in. x D 32 in.
- Opened: L 79.25 in. x W 33 in. x H 35.25 in.
- Weight: 54 lb (Mattress and Frame)

The field deployable, lightweight solution for patient comfort in the combat support hospital.

Extremely compact for easy transport and storage...



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS WITH A
MISSION®