

TALON® BACKREST *with Nomex/Kevlar Blend*

TACEVAC



Efficient and Reliable Solution...

For appropriate treatment and improved comfort, supine is not always the best position for patients. When casualties can be better served by being placed at a 45° or 90° angle, the Talon® Series Back Rest from North American Rescue offers an efficient and reliable solution.

This litter backrest attaches to the 2 in. poles of the Talon® system, the premiere litter system on the market. The backrest is also engineered to fit the Raven® Model 90C.

The backrest cover is made of a Nomex/Kevlar blend: Nomex for flame-resistance and heat-energy absorption; Kevlar for abrasion-resistance and long life. The fabric can be repeatedly washed or dry-cleaned without losing these critical properties.

When the patient's condition requires a 45° or 90° angle rely on the Talon® Series Backrest.

Specifications:

- H 31.5 in. x W 22.5 in. x D 1.6 in.
- Weight: 2 lb 13 oz

Features:

- Allows patients to be placed at 90° or 45° angles
- Made of Nomex/Kevlar-blend fabric
- Attaches to the 2 in. poles of the Talon® litter system or the Raven® Model 90C
- Firm, comfortable, and secure

Compatible with the Talon, Raven and Medevac litters.

Talon® Backrest

Item#	NSN#
60-0027	6530-01-594-7444

Seamlessly attaches to the 2 inch poles of the Talon® litter system or the Raven® Model 90C.

PRODUCTS WITH A
MISSION®

NOTE: All accessories are sold separately. Litters are not included.



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277