

# SCISSOR LEASH



# TREATMENT

## Retractable Lanyard...

The Scissor Leash is the ultimate solution for maintaining positive control of your trauma shears in a tactical environment.

The Scissor Leash is lightweight and highly reliable. The retractable lanyard snaps back immediately after use allowing the operator to know where his scissors are at all times.

### Features:

- High break strength
- Salt water resistant
- Self cleaning
- Extends from a retracted length of 6 in. to 42 in.

### Dimensions:

- Retracted: H 6 in. x W 1.6 in. x D 0.6 in.
- Cable Length - 36 in.
- Weight: 1.6 oz

## Scissor Leash

Item#	NSN#
70-0005	6515-01-540-7226

*Lightweight, reliable leash with a retractable snap back lanyard.*

PRODUCTS WITH A  
**MISSION**<sup>®</sup>



**NORTH AMERICAN RESCUE<sup>®</sup>**

[www.NARescue.com](http://www.NARescue.com) • 888.689.6277