


# OSL™ LITTER STANDS *Set of 3*

REV121919



*OverSized Litter (OSL) Stands are easy to open, quick to lock into place and are easily folded for transport and storage.*

## ***Added stability for the patient, better access for the medic...***

Lightweight and compact aluminum folding litter stands, designed to lock into place and provide a stable platform for the patient. The OSL litter stands come in a set of 3 and accommodate the OverSized OSL litter.

Equipped with litter stirrups to restrict the litter from sliding front to back. Once the stands are opened, the cross members easily lock into place.

### **Features:**

- Set of 3 lightweight aluminum stanchions to accommodate our oversized OSL™ Litter
- Added stability for the patient
- Better access for the medic
- Easy to open; quickly locks to secure litter
- Easily folded for transport & storage

### **Specifications:**

- Folded: H 34 in. x W 25 in. x D 2.5 in.
- Extended Height: 32 in.
- Base Width: 23.5 in.
- Weight: 28.5 lb
- Payload: 1600 lbs (set of 3)
- Unit of Issue: Set of 3

### ***OSL Litter Stands***

CE

Unit	Item#	NSN#
Set of 3	60-0029	

*NOTE: All accessories are sold separately. Litters are not included.*



**NORTH AMERICAN RESCUE®**

[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

**PRODUCTS**  
WITH A  
**MISSION®**