

TACTICAL OPERATOR RESPONSE KIT™ TORK™



The Tactical Operator Response Kit™, from North American Rescue, is our streamlined, point-of-wounding IFAK designed for the assault vest or body armor. The TORK™ meets the need for a compact individual first aid kit to carry all the necessary equipment to address the leading causes of preventable death in a tactical environment. This kit contains essential life-saving equipment that can be used at the point of wounding for self-aid/buddy-aid to address injuries caused by penetrating or blast trauma.



- Compact individual first aid kit providing Advanced Life Support for immediate Self-Aid/Buddy-Aid of penetrating or blast trauma
- Critical life-saving tools that address the top leading causes of potentially preventable death in the tactical environment
- Berry Amendment Compliant Bag with easy MOLLE/PALS style attachment to the assault vest or body armor
- Also attaches to NAR's Quick Detach Panel (sold separately)



TORK™ Tactical Operator Response Kit™	
BLK:	Item# 80-0292 (Celox Gauze)
	Item# 85-0013 (ChitoGauze)
	Item# 85-0048 (Combat Gauze)
FOL:	Item# 80-0295 (Celox Gauze)
	Item# 85-0016 (ChitoGauze)
	Item# 85-0051 (Combat Gauze)
COY:	Item# 85-0014 (ChitoGauze)
	Item# 85-0049 (Combat Gauze)
	Item# 80-0293 (Celox Gauze)
DUC:	Item# 80-0294 (Celox Gauze)
	Item# 85-0015 (ChitoGauze)
	Item# 85-0050 (Combat Gauze)
MTC:	Item# 80-0296 (Celox Gauze)
	Item# 85-0017 (ChitoGauze)
	Item# 85-0052 (Combat Gauze)
ODG:	Item# 80-0297 (Celox Gauze)
	Item# 85-0055 (ChitoGauze)
	Item# 85-0053 (Combat Gauze)
Kit Dimensions: Case Closed - H 7 in. x W 4 in. x D 2.75 in.	
Weight: 1 lb 5 oz	

TORK™ (Bag Only)
Item# 80-0227 (Black)
Item# 80-0230 (Foliage)
Item# 80-0228 (Coyote)
Item# 80-0229 (DUC)
Item# 80-0231 (MTC)
Item# 80-0237 (ODGreen)



CHOICE OF HEMOSTATIC



(1)
Celox™ Gauze
Z-fold Hemostatic (5 ft)
Item# 30-0064



(1)
Z-Fold Combat Gauze™ Hemostatic
Item# 30-0039

-OR-



(1)
ChitoGauze® PRO
Z-fold Hemostatic
Item# 30-0049
NSN# 6510-01-591-7740

-OR-



(1)
Nasopharyngeal Airway 28F
with Lubricant
Item# ZZ-0034
NSN# 6515-01-529-1187

◀ 3 hemostatic dressing options available – help control bleeding in hemorrhage not amenable to tourniquet application

*NOTE: The Item & NSN numbers depicted are for the purpose of restocking individual components & may differ from those individual items packed in the kit. Kit contents subject to change dependent upon DoD requirements.



(1)
C-A-T® Combat Application Tourniquet®
Item# 30-0001
NSN# 6515-01-521-7976



(1)
ETD 6 in.
Emergency Trauma Dressing
Item# 30-0032
NSN# 6510-01-558-4114



(1)
S-Rolled Gauze
4.5 in. x 4.1 yd.
Item# 30-0003
NSN# 6510-01-605-2234

▲ MASSIVE HEMORRHAGE



(1)
ARS®
14 gauge x 3.25 in.
Needle Decompression Kit
Item# ZZ-0056
NSN# 6515-01-541-0635



(1)
HyFin® Twin Pack
Chest Seals
Item# 10-0031

▲ RESPIRATION




(2 pair)
Bear Claw®*
Nitrile Trauma Gloves, lg
Item# ZZ-0248
NSN# 6515-01-600-0262



(1 pair)
Trauma Shears
(7.25 in.)
Item# ZZ-0063
NSN# 6515-01-538-9276

▲ TREATMENT SUPPORT



*NOTE: The Item & NSN numbers depicted are for the purpose of restocking individual components & may differ from those individual items packed in the kit. Kit contents subject to change dependent upon DoD requirements.

BLS TORK™ CONTENTS



(1) **C-A-T® Combat Application Tourniquet®**
Item# 30-0001
NSN# 6515-01-521-7976



(1) **ETD 6 in.**
Emergency Trauma Dressing
Item# 30-0032
NSN# 6510-01-558-4114



(2) **S-Rolled Gauze**
4.5 in. x 4.1 yd.
Item# 30-0003
NSN# 6510-01-605-2234

▲ MASSIVE HEMORRHAGE

- Compact individual first aid kit providing basic essential life-saving tools to address injuries caused by penetrating or blast trauma
- Sleek, compact size allows for easy attachment to the assault vest/body armor using MOLLE/PALS style connectors
- Berry Amendment Compliant Bag with easy MOLLE/PALS style attachment to the assault vest or body armor
- Also attaches to NAR's Quick Detach Panel (sold separately)



BASIC LIFE SUPPORT TACTICAL OPERATOR RESPONSE KIT™

BLS TORK™

Basic life support IFAK to provide immediate Self-Aid/Buddy-Aid at the point of injury



BLS TORK™ CONTENTS



(1) **Nasopharyngeal Airway 28F**
with Lubricant
Item# ZZ-0034
NSN# 6515-01-529-1187

▲ AIRWAY



(2) **HyFin® Vent** Chest Seal
Item# 10-0029

▲ RESPIRATION



(2 pair) **Bear Claw®**
Nitrile Trauma Gloves, lg
Item# ZZ-0248
NSN# 6515-01-600-0262

▲ TREATMENT SUPPORT



(1 pair) **Trauma Shears**
(7.25 in.)
Item# ZZ-0063
NSN# 6515-01-538-9276

NAR QUICK DETACH PANEL

Allows the operator the fastest access to their kit during the critical moments after wounding



NAR Quick Detach Panel

- Item# 80-0232 (Black)
- Item# 80-0235 (Foliage)
- Item# 80-0233 (Coyote)
- Item# 80-0234 (DUC)
- Item# 80-0236 (MultiCam)
- Item# 80-0238 (ODGreen)

Kit Dimensions: Case Closed - L 7.25 in. x W 4.5 in. x D 0.5 in.

Weight: 4 oz

- Constructed of Mil-spec Velcro® & heavy-duty metal snaps to ensure that the detachable front panel remains securely fastened
- Designed to compliment the TORK™ IFAK
- Compatible with any Mil-spec standard MOLLE/PALS system
- Can be used to support any NAR IFAK, leg rig or chest pouch up to four pounds