

MARITIME ASSAULT KIT™

REV120518

Complete With All the Tools to Help Mitigate the Three Leading Causes of Preventable Combat Death

North American Rescue's Maritime Assault Kit™ meets the unique and specific requirements of the Navy Expeditionary Combat Command for a small, efficient IFAK.

Complete with all the tools to mitigate the three leading causes of preventable combat death (60% extremity bleeding / 33% tension pneumothorax / 6% airway obstruction).



QTY
1 C-A-T®
Tourniquet



QTY
1 ETD™ 6 in.
Responder



QTY
1 Gauze
S-Rolled



QTY
1 Combat
Gauze



QTY
1 28F Naso
Airway & Lube



QTY
1 HyFin®
Vent. Twin Pk



QTY
1 ARS®
Needle, 14g



PAIR
2 Nitrile
Gloves, Lg



QTY
1 PES
Eye Shield



QTY
1 T2® Tag
Triage Card



QTY
1 Casualty
Reference Card



The Maritime Assault Kit Modular system can be attached to the operator in three ways:

1. Using the split-tale connector with height adjustment and Quick-Grip™ leg strap
2. PALS/MOLLE style webbing for attachment to the vest/body armor
3. Fasten to the web/gun belt using the quick weave method. Belt and leg attachment straps can be removed and reattached as needed to allow maximum flexibility during load carriage configuration

Specifications:

- Case Closed: H 7.4 in. x W 6.9 in. x D 5 in.
- Weight: 1 lb 13 oz

Color	Item#	NSN#
BLK	80-0013	6545-01-529-4185
ODG	80-0014	6545-01-531-6626
COY	80-0015	6545-01-529-4187
DUC	80-0058	
MTC	80-0245	

Available Colors:



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®

MARITIME ASSAULT KIT™

REV120518

Features:

- 500D CORDURA® Solution-Dyed Nylon with Near Infrared Signature Reduction
- Ghillie-Tex IR signature reduction hardware
- Sand & water resistant, self-repairing zipper with silent 550 cord pulls
- MOLLE-style Attachment System Compatible
- Removable Belt and Leg Straps
- Leg Strap with Quick-Grip™
- Available in NAR's 11 Point Berry Amendment Nylon Standard

Kit Contents:

- 1 x Bag (CCRK®—Maritime)
- 1 x C-A-T® (Combat Application Tourniquet®)
- 1 x ETD™ 6 in. Emergency Trauma Dressing
- 1 x Z-Fold Vacuum Packed QuikClot Combat Gauze
- 1 x S-Rolled Gauze™ (4.5 in. x 4.1 yd)
- 1 x Nasopharyngeal Airway 28F with Lubricant
- 1 x HyFin® Vent Chest Seal, Twin Pack
- 1 x ARS® Needle Decompression Kit (14 G x 3.25 in.)
- 1 x Polycarbonate Eye Shield
- 2 x pair Black Talon® Nitrile Trauma Gloves
- 1 x Combat Casualty Reference Card
- 1 x Triage Card (NAR T2 Tag®)



PRODUCTS
— WITH A —
MISSION®



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277