

# INDIVIDUAL KIT™

Combat Casualty Response Kits® **CCRK®**

## Designed with all the Tools for the Warfighter/Operator to Mitigate Extremity Bleeding, Tension Pneumothorax and Airway Obstruction.

The Individual Kit from North American Rescue has been designed for the Warfighter/Operator and packed to provide all of the tools to mitigate the three leading causes of preventable combat death (60% extremity bleeding / 33% tension pneumothorax / 6% airway obstruction).

Our modular system can be attached to the operator using the split-tale connector, PALS/MOLLE style webbing or by fastening to the web/ gun belt using the quick weave method. Belt and leg attachment straps can be removed and reattached as needed to allow maximum flexibility during load carriage configuration.



*Attachment with split-tail  
connector and Quick-Grip  
leg strap or PALS/MOLLE  
style webbing.*

### Available Colors:



Color	Item#	NSN#
BLK	80-0003	6530-01-529-1134
ODG	80-0004	6545-01-531-6626
COY	80-0005	6545-01-529-4187
DUC	80-0091	



**NORTH AMERICAN RESCUE®**  
www.NARescue.com • 888.689.6277

PRODUCTS WITH A  
**MISSION®**

## Features:

- 500D CORDURA® Solution-Dyed Nylon with Near Infrared Signature Reduction
- Ghillie-Tex IR signature reduction hardware
- Sand & water resistant, self-repairing zipper with silent 550 cord pulls
- MOLLE Style Attachment System Compatible
- Removable Belt and Leg Straps
- Leg Strap with Quick-Grip™
- Available in NAR's 11 Point Berry Amendment Nylon Standard

## Specifications:

- Case Closed: H 8.5 in. x W 8.75 in. x D 4.5 in.
- Weight: 1 lb 12 oz

## Kit Contents:

- 1 x Bag (CCRK®—Individual)
- 4 x Black Talon® Nitrile Trauma Gloves (2 pr.)
- 1 x Nasopharyngeal Airway 28F with Lubricant
- 1 x HyFin® Chest Seal
- 1 x ARS® Needle Decompression Kit (14 G x 3.25 in.)
- 1 x C-A-T® (Combat Application Tourniquet®)
- 1 x ETD™ 6 in. Emergency Trauma Dressing
- 2 x S-Rolled Gauze™ (4.5 in. x 4.1 yd)
- 1 x Polycarbonate Eye Shield
- 1 x Combat Casualty Reference Card
- 1 x Triage Card (NAR T2 Tag®)



**NORTH AMERICAN RESCUE®**

[www.NARescue.com](http://www.NARescue.com) • 888.689.6277