


# NAR® TRACHEAL HOOK

## Meets the unique needs of Tactical Healthcare Professionals...

The NAR® Tracheal Hook is designed to meet the specific and unique needs of the Tactical Healthcare Professional when performing surgical cricothyroidotomy. Unlike other tracheal hooks the NAR Tracheal Hook has a unique finger grip that allows you to stabilize and anteriorly displace the cricoid cartilage while providing maximum access and visualization.

### FEATURES/SPECIFICATIONS

- For maximum access & visualization during an emergency cricothyroidotomy
- Stabilizes the cricoid cartilage for anterior displacement
- Unique finger grip design



### NAR TRACHEAL HOOK

ITEM #	NSN
10-0020	6515-01-562-9734

### DIMENSIONS

- Device: L 5.25 in. x W 0.5 in. x D 0.8 in.
- Weight: 0.5 oz



Learn more about this and other *Products with a Mission®* at [www.NARescue.com](http://www.NARescue.com)



**NORTH AMERICAN RESCUE®**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277