



COMBAT CASUALTY RESPONSE KIT® NAR-4™ AID KIT

NOTE: IR Flag sold separately

NAR-4™ Aid (Bag Only)

Item# 80-0177 (Black)
Item# 80-0178 (Foliage)
Item# 80-0179 (Coyote)
Item# 80-0180 (DUC)
Item# 80-0241 (MultiCam)
Item# 80-0322 (ODGreen)



The constant challenge of striking a balance between the increasing amount of equipment carried into battle and the need to optimize performance is the impetus behind this low-profile, load carriage platform for Tactical Healthcare Professionals. The NAR-4™ Aid Kit is an ergo-dynamic, durable, functional solution that supports Tactical Combat Casualty Care equipment standards. This kit is designed to attach directly to the body armor using four rapid attachment straps or carried using shoulder straps which are stowed away until needed.

Contact Customer Service to
order replacement pouches ▶



◀ Features a modular system with 5 removable
nylon bags that can be configured to carry airway
tools, IV equipment & hemorrhage control gear



▶ Can be pre-rigged
to attach directly to
body armor



(2) **Nasopharyngeal Airway 28F**
with Lubricant
Item# ZZ-0034
NSN# 6515-01-529-1187



(1) **Tactical Suction Device**
Item# 10-0018
NSN# 6515-01-540-7206



(1) **King LTS-D™** Size 4
Supraglottic Airway Kit
Item# 10-0026
NSN# 6515-01-591-1578



(1) **Tactical CricKit®**
Cricothyroidotomy Kit
Item# 10-0017
NSN# 6515-01-540-7568



(1) **Cyclone®** New & Improved
Pocket BVM®
Item# 10-0019
NSN# 6515-01-568-0193



(2) **HyFin®** Chest Seal
Item# 10-0015
NSN# 6515-01-532-8019



(2) **ARS®**
14 gauge x 3.25 in.
Needle Decompression Kit
Item# ZZ-0056
NSN# 6515-01-541-0635



(1) **BOA®**
IV Constricting Band
Item# 30-0009
NSN# 6515-01-537-2111



(4) **NAR Saline Lock Kit**
Item# 30-0013
NSN# 6515-01-537-4094

▲ AIRWAY

▲ RESPIRATION

▲ CIRCULATION



▲ MASSIVE HEMORRHAGE

..... CHOICE OF HEMOSTATIC


- Configurable according to individual skill set & mission requirements
- Rapid access & organization is enhanced by segregated removable equipment pouches
- 500D CORDURA® Berry Amendment Nylon Bag with GillieTEX™ hardware for a reduced IR signature
- Sand & water resistant, self-repairing zippers with silent 550 cord pulls
- PALS/MOLLE-style attachment system on front & side of the bag
- Packing list is consistent with TCCC Guidelines

Slimline, bag can also be carried as a low-profile backpack ▶



NAR-4™ Aid Kit (CCRK®)	
BLK:	Item# 80-0314 (Celox Gauze)
	Item# 80-0309 (ChitoGauze)
	Item# 80-0181 (Combat Gauze)
FOL:	Item# 80-0317 (Celox Gauze)
	Item# 80-0312 (ChitoGauze)
	Item# 80-0182 (Combat Gauze)
COY:	Item# 80-0315 (Celox Gauze)
	Item# 80-0310 (ChitoGauze)
	Item# 80-0183 (Combat Gauze)
DUC:	Item# 80-0316 (Celox Gauze)
	Item# 80-0311 (ChitoGauze)
	Item# 80-0184 (Combat Gauze)
MTC:	Item# 80-0318 (Celox Gauze)
	Item# 80-0313 (ChitoGauze)
	Item# 80-0244 (Combat Gauze)
ODG:	Item# 80-0319 (Celox Gauze)
	Item# 80-0320 (ChitoGauze)
	Item# 80-0308 (Combat Gauze)
Kit Dimensions: Case Closed - H 19 in. x W 15 in. x D 7 in.	
Weight: 8 lb 10 oz	

*NOTE: The Item & NSN numbers depicted are for the purpose of restocking individual components & may differ from those individual items packed in the kit. Kit contents subject to change dependent upon DoD requirements.



▲ SPLINTING

▲ TREATMENT SUPPORT