

# WALK® Warrior Aid & Litter Kit

Combat Casualty Response Kits®

CCCRK®

REV032218

## **Official Vehicle Life-Saving Kit for the Department of Defense**

***Allows first-aid treatment and immediate evacuation of one patient in a non-traditional vehicle, plus it allows multiple other trauma patients to be treated until additional transport arrives***

From numerous observations and after-action reviews compiled from Operations Enduring and Iraqi Freedom, we learned that a new focus was needed for vehicle first-aid kits. From Warfighters and Tactical Healthcare Professionals, we learned the focus ought to be decreasing preventable combat death at the point of wounding.

Our engineers took the state-of-the-art Talon 90C (a quad-fold rigid litter with retractable handle system) and designed a package that holds not only the litter, but also a streamlined pocket-system containing self- and buddy-aid medical supplies. The WALK gives true definition to the term aid and litter kit. The WALK also contains our Hypothermia Prevention & Management Kit (HPMK) to help stabilize patients during transport. Light and compact, the WALK® reduces cubic weight and, thereby, evacuation time.

The WALK® interfaces with any standard or non-standard casualty evacuation asset and is ideal as part of the pre-positioning of medical supplies or cache for disaster response. The WALK® has adjustable shoulder straps for occasions when the kit is used to provide treatment and evacuation on target.



**PRODUCTS**  
— WITH A —  
**MISSION®**



**NORTH AMERICAN RESCUE®**

**www.NARescue.com • 888.689.6277**

# WALK® Warrior Aid & Litter Kit


Combat Casualty Response Kits®

CCCRK®

REV032218

*For the treatment of multiple casualties and the immediate evacuation of one patient*

## WARRIOR AID & LITTER KIT CONTENTS:



## Kit Contents:

- 1 x Bag (CCCRK®-WALK®)
- 5 x pair Black Talon® Nitrile Trauma Gloves
- 2 x Nasopharyngeal Airway 28F with Lubricant
- 1 x Casualty Equipment Bag
- 2 x HyFin® Chest Seal
- 2 x ARS® Needle Decompression Kit (14 G x 3.25 in.)
- 2 x C-A-T® (Combat Application Tourniquet®)
- 6 x Trauma Dressing, 6 in.
- 4 x S-Rolled Gauze™ (4.5 in. x 4.1 yd)
- 1 x ETD™ Abdominal Emergency Trauma Dressing
- 2 x SAM® Splint II
- 1 x Trauma Shears (7.25 in.)
- 1 x Surgical Tape (2 in.)
- 1 x GeckoGrip™ Multipurpose Tape
- 6 x Polycarbonate Eye Shield
- 1 x Combat Casualty Reference Card
- 2 x Combat Casualty Card (Triage)
- 1 x Aviation Panel (Recognition, Orange)
- 1 x Talon 90C® Collapsible Handle Litter
- 1 x Hypothermia Prevention & Management Kit™ (HPMK™)
- 4 x Tie Down Strap (Universal Litter)

**WALK ReSupply Kit** (Item# 85-0020):

Contains WALK components for resupply - *excluding nylon bag & litter*



**NORTH AMERICAN RESCUE®**

www.NARescue.com • 888.689.6277

**PRODUCTS**  
WITH A  
**MISSION®**



# WALK® Warrior Aid & Litter Kit

Combat Casualty Response Kits®

CCCRK®

REV032218



*Includes Talon®II 90C quad-fold litter*

## Features:

- Official US Army vehicle life-saving kit
- For treatment of multiple trauma patients and immediate evacuation of one patient
- Includes Talon®II 90C quad-fold litter
- Hypothermia Prevention & Management Kit™ helps stabilize patients during transport
- Interfaces with any standard or non-standard TACEVAC asset
- Accepts both: Talon®II 81C (Retired) & Talon®II 90C Collapsible Handle

## Specifications:

- Case Closed: H 22.75 in. x W 16.5 in. x D 11.5 in.
- Weight: 29 lb 15 oz

## Available Colors:

BLK

ODG

COY

DUC

Color	Item#	NSN#
BLK	80-0097	
ODG	80-0060	
COY	80-0089	6545-01-587-1199
DUC	80-0059	6545-01-532-4962
ReSupply*	85-0020	

\*Contains all WALK components excluding nylon bag & litter



**NORTH AMERICAN RESCUE®**

[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

**PRODUCTS**  
— WITH A —  
**MISSION®**