




# FIBER OPTIC LARYNGOSCOPE BLADES



◀ Miller #0



◀ Miller #1



◀ Miller #2



◀ Miller #3

## Fiber Optic Laryngoscope Blades, Miller:

Single-use fiber optic laryngoscope blades with enhanced polymer construction offer exceptional strength and rigidity while reducing the risk of damage to the patient's teeth or soft tissues.

Designed for safety and convenience, as well as the superior illumination of fiber optic technology. Compatible with any standard (green spec) fiber optic handle.

### Features:

- Single-use laryngoscope blades compatible with any standard fiber optic handle
- Enhanced polymer construction for less risk of damage to patient's teeth or soft tissues
- Exceptional strength and rigidity plus superior illumination
- Designed for safety and convenience

### Miller #0:

- L 3.2 in. x W 1.1 in. x D 0.9 in.
- Weight: 0.4 oz

### Miller #1:

- L 4.1 in. x W 1.1 in. x D 0.9 in.
- Weight: 0.6 oz

### Miller #2:

- L 6.2 in. x W 1.2 in. x D 0.9 in.
- Weight: 1 ozz

### Miller #3:

- L 7.7 in. x W 1.2 in. x D 0.9 in.
- Weight: 1.2 oz



Size	Item#	SEL#	NSN#
Miller #0	10-0008	09ME-02-AWMG	6515-01-460-4718
Miller #1	10-0009	09ME-02-AWMG	6515-01-460-4724
Miller #2	10-0010	09ME-02-AWMG	6515-01-460-4725
Miller #3	10-0011	09ME-02-AWMG	6515-01-460-4681



**NORTH AMERICAN RESCUE®**

www.NARescue.com • 888.689.6277

PRODUCTS WITH A  
**MISSION®**