

The **Fastest, Safest,**  
Most **Effective** Prehospital Field Tourniquet

THE 7 ADVANTAGE DIFFERENCE



**1** Single Routing  
Buckle

- ▶ **Faster application.**
- ▶ **Effective slack removal.**
- ▶ **Decreased blood loss.**
- ▶ **Fewer windlass turns.**
- ▶ **Simplified training with single protocol directions for all applications.**

**2** Windlass  
Rod

- ▶ **Increased diameter.**
- ▶ **Enhanced strength.**
- ▶ **Aggressive ribbing.**
- ▶ **Improved grip.**

Introducing the **NEW Gen 7**

# CAT<sup>®</sup> COMBAT APPLICATION TOURNIQUET<sup>®</sup>



## 3 Windlass Clip

- ▶ **Bilateral entry.**
- ▶ **Rapid windlass lock.**
- ▶ **Bilateral buttress.**
- ▶ **Added strength.**

## 4 Windlass Strap

- ▶ **Sonic welded to clip for constant contact.**
- ▶ **Color changed to Gray for tactical considerations**  
(Black C-A-T<sup>®</sup> only).

## 5 Stabilization Plate

- ▶ **Reinforced, beveled contact bar.**

## 6 Red Tip Technology<sup>®</sup>

- ▶ **Red elliptical tip.**
- ▶ **Added visual cue for application.**

## 7 Free-moving Internal Band

- ▶ **Patented band within a band.**
- ▶ **True evenly distributed circumferential pressure.**