

# NAR SHARPS SHUTTLE



*Convenient Point  
of Use Containment*

## Makes the Disposal of Sharps Safe:

The NAR Sharps Shuttle is clear for immediate fill-level monitoring. The unit's tapered body (conveniently sized at only 1.25 in. diameter and 6 in. length) holds up to six needles.

The vertical "point first" design is easy to store and fits perfectly into our Medic Leg-Rig or Trauma Pack for convenient point-of-use containment of used needles and small syringes. The product's quick-tab locking top secures contaminants until disposal.

The product's unique single-use construction makes the disposal of sharps safe.

## Features:

- Clear container for fill-level monitoring
- Vertical "point first" design for easy storage
- Holds up to 6 needles
- Tapered body (1.25 in. diameter by 6 in. length) fits snugly into our Medic Leg-Rig™ or Trauma Pack
- Quick-tab locking top secures contaminants until disposal
- Available as a single unit or in a 6-pack

## Specifications:

- Individual: L 6.5 in. x Diameter 1.25 in.
- Weight: 1.5 oz



Unit	Item#	NSN#
Single	30-0014	6515-01-518-8551
Pack of 6	30-00146	6515-01-619-5671



**NORTH AMERICAN RESCUE®**  
www.NARescue.com • 888.689.6277

PRODUCTS WITH A  
**MISSION®**