

TALON® II 90C *Collapsible Handle Litter*

REV070822



Ultimate Evacuation Platform

The upgraded Talon II 90C is the most versatile and compact evacuation and treatment litter available today!

Developed to meet the U.S. Army's requirement to provide urgent casualty evacuation from the battlefield. The new design allows patient transport in restricted compartments such as the up-armored M1114 Humvee and other non-conventional commercial vehicles. By simply extending the ergonomically designed handles it transitions to a standard NATO compatible evacuation platform.

Features:

- Advanced, innovative combat evacuation system featuring patented portable quad-fold design with retractable handles
- Developed to U.S. Army requirements
- Retractable, ergonomically designed handles
- Automatic locking hinges
- Slip-resistant, chemical-resistant ripstop fabric tested & approved by Federal & private labs using live warfare agents
- 6 IV-attachment points, 2 patient-securing straps
- U.S. Patents 6,842,923 & 5,598,592, Patent(s) Pending

Specifications:

- Folded with handles collapsed: H 20.5 in. x W 8.5 in. x D 9.5 in.
- Open with handles collapsed: L 78 in. x W 22.5 in. x D 6 in.
- Open with handles extended: L 90 in. x W 22.5 in. x D 6 in.
- Weight: 16 lb 6 oz (7.4 kg)
- Static Payload: 1800 lb (816 kg)
- Ground Clearance: 1.5 in.

Talon® II 90C



Item#	NSN#
60-0002	6530-01-504-9051
CE 60-0002-EU	



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®