

STINGRAY® POLELESS LITTER

NON-RIGID LITTERS

Tested and Approved Poleless Litter...

Weighing less than 2 pounds with a payload of 1000 pounds, the Stingray® Litter is strong and lightweight. Constructed with a non-skid, flame-retardant monofilament polypropylene material. Two web straps secure the litter around the body; six heavy-duty handles make it easy to carry.

This design allows rescuers to insert a 2 in. pole, such as pike poles, through the litter cover or to use debris through the outer handles. The Stingray® is ideal for decontamination roller systems and moving wet patients.


Stingray® is tested and approved by both Federal and private laboratories using live warfare agents.

FEATURES/SPECIFICATIONS

- Constructed of non-skid, flame-retardant monofilament polypropylene fabric
- Resistant to mildew, acids, alkalis and penetration liquids
- Tested by Federal & private labs using live warfare agents
- 6 heavy-duty carrying handles
- 2 casualty securing straps
- Made for decontamination rollers
- Great for long term storage

DIMENSIONS

- Folded: H 20 in. x W 12 in. x D 2 in.
- Open: L 78 in. x W 24 in.
- Weight: 1 lb 13 oz (0.7 kg)
- Payload: 1000 lb (453.6 kg)



**Resistant to mildew,
acids, alkalis and
penetration liquids.**

STINGRAY® POLELESS LITTER



ITEM #	NSN
60-0004	6530-01-492-4079
60-0004-EU	



Learn more about this and other **Products with a Mission®** at www.NARescue.com



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277