

Learn more about this and other **Products with a Mission®** at www.NARescue.com

TALON® BACKREST

WITH NOMEX/KEVLAR BLEND

Compatible with the Talon, Raven, and Medevac litters

For appropriate treatment and improved comfort, supine is not always the best position for patients. When casualties can be better served by being placed at a 45° or 90° angle, the Talon® Series Back Rest from North American Rescue offers an efficient and reliable solution.

This litter backrest attaches to the 1.5 in. poles of the Talon® system, the premiere litter system on the market. The backrest is also engineered to fit the Raven® Model 90C.

The backrest cover is made of a Nomex/Kevlar blend: Nomex for flame-resistance and heat-energy absorption; Kevlar for abrasion-resistance and long life. The fabric can be repeatedly washed or dry-cleaned without losing these critical properties

When the patient's condition requires a 45° or 90° angle rely on the Talon® Series Backrest.



NOTE: All accessories are sold separately. Litters are not included.

FEATURES/SPECIFICATIONS

- Allows patients to be placed at 90° or 45° angles
- Made of Nomex/Kevlar-blend fabric

- Compatible for use with the Talon®, Raven® and Medevac litters
- Firm, comfortable, and secure

TALON® BACKREST

ITEM #	NSN
60-0014	6530-01-523-3314



DIMENSIONS

- H 31.5 in. x W 22.5 in. x D 1.6 in.
- Weight: 2 lb 13 oz



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277