

CELOX™ HEMOSTATIC GRANULES



CELOX™ Granular Hemostatic Agent:

CELOX,™ a clinically-proven family of chitosan-based products certified safe-to-use and successful in stopping life threatening hemorrhage on the battlefield, in serious emergencies and in the workplace. From the original granular form to the various CELOX™ gauzes, NAR offers the full spectrum of this trusted hemostatic to meet every medic's need in the pursuit of saving lives from severe trauma.

The original version of CELOX is composed of hemostatic granules that can be easily poured directly into a wound. The granular presentation is the most versatile of the CELOX™ products providing a simple, quick way to control bleeding and is particularly useful in treating difficult and deep wounds.

Features:

- Temporary traumatic wound treatment for life threatening emergency bleeding
- Promotes rapid coagulation
- Chitosan-based, works independently of the body's normal clotting mechanisms
- Published success in the field
- Generates no heat in use

Specifications:

- Packaged: H 7.5 in. x W 5 in. x D 0.25 in.
- 35 Grams of Granules
- Weight: 1.6 oz

Item#	NSN#
30-0068	



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS WITH A
MISSION®