

NAR TACTICAL SUCTION DEVICE

REV121420



Compact and Powerful:

Compact and powerful, the NAR Tactical Suction Device is designed to effectively evacuate the oropharynx in emergency situations. The device is designed to provide a vacuum force of approx 100 mmHg generated by the rebound of the bulb after each squeeze. Evacuated contents are expelled into a clear collection bag attached to the outlet port. A unique double-valved design allows for continuous operation with one hand. The other hand can then direct the contoured suction tip to prevent tissue damage.

The Tactical Suction Device from North American Rescue is a reliable, simple, rugged, and inexpensive solution for clearing airway obstruction due to bloody fluids or other secretions.

Features:

- Double-valved design for one-hand operation
- Vacuum force of approximately 100 mmHg
- Clear collection bag & bulb to effectively evacuate the oropharynx
- Contoured suction tip prevents tissue damage
- Reliable, simple, rugged & inexpensive

Specifications:

- Packaged: L 9.75 in. x W 3.75 in. x D 3 in.
- Vacuum Force: 100mmHG
- Weight: 7.3 oz

TACTICAL SUCTION DEVICE

ITEM#	NSN #
10-0018	6515-01-540-7206



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®