

TOMMANIKIN® UNINJURED CHEST SKIN

REV072522

The uninjured TOMManikin® chest skin supports chest seal application and pleural needle decompression training scenarios. The uninjured chest skin is secured by a large section of hook and loop. This consumable item can be ordered to be bilateral chest tube compatible.

Specifications:

- Chest Seal application and needle decompression training skin
- No visible exterior wounds or injuries
- Compatible with all male TOMManikin® simulators

Dimensions:

- L 20 in. x W 12 in. x D 6 in.
- Weight: 9 lbs



TOMMANIKIN® Uninjured Chest Skin

ITEM #	NSN #
93-0148	----



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®