

# UNIVERSAL SEAT BACK HOLDER *Shoulder Straps*

The North American Rescue Universal Seat Back Holder (USB)-Shoulder Straps allow for hasty deployment of life saving medical gear further into the field from a vehicle.

The USB shoulder straps are a set of super low-profile laser cut laminate nylon straps. They easily connect and disconnect with metal G-hook buckles directly to the back of the USB. The straps are designed to hide seamlessly against the back of the Universal Seat Back Holder out of the way until user needs to quickly respond to an incident with their entire USB load out. The USB can be removed, and shoulder straps donned in under 10 seconds. Strap can be also connected to each other creating an adjustable bandolier strap as needed.

Be sure to check out the optional Universal Seat Back Holder Malleable Bar as well.

## **Specifications:**

- Durable Laser Cut Laminate (LCL) nylon construction
- Secured with easily added, removed, moved metal G-hooks
- Allows quick deployment of Universal Seat Back Holder
- Adjustable G-hook design for shoulder straps or connected for bandolier style carry

## **Dimensions:**

- H 13 in. x W 2 in. x D 0.125 in.
- Adjustable Strap Size: 14.5 in. to 29 in.
- Adjustable Strap Size (Bandolier Style): 30 in. to 58 in.
- Weight: 0.1 lbs



### **Universal Seat Back Holder Shoulder Straps**

ITEM #

ZZ-1585

*\* Shoulder Straps only - Universal Seat Back Holder sold separately.*

*\*Universal Seat Back Holder sold separately. Check out [NARescue.com](http://NARescue.com) for quick detach panels, kits, gear, and accessories.*



**NORTH AMERICAN RESCUE®**

[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

**PRODUCTS**  
WITH  **MISSION®**