

LITTER STANDS *33 inch Standard*

REV062722



Our 33" stands work with Talon, Raven, T3, and Medevac Litters

Provides a Stable Casualty Platform for Ease of Patient Access...

Lightweight and compact aluminum folding litter stands, designed to lock into place and provide a stable platform for the patient.

Accommodates the Talon II and the Raven litter stirrups to restrict the litter from sliding front to back. Once the stand is opened, the cross member easily locks into place.

Constructed of lightweight aluminum, these 33 inch Standard Litter Stands conveniently fold flat for ease of transport and storage.

Features:

- Set of 2 lightweight, foldable stands for standard size litters
- Aluminum structure for lightweight strength
- Easily opened to stand 33 in. tall
- Quickly secures the litter to increase stability and access to the patient
- Easily folded for transport and storage
- Quickly locks open to safely secure the litter in place.

ODG

Standard Litter Stands

Unit	Item#	NSN#
Set of 2	60-0014	6530-01-523-3314
Set of 2	60-0014-EU	

CE

Lightweight aluminum structure for stability and ease of transport, these Standard Litter Stands weigh less than 9 pounds.

NOTE: All accessories are sold separately. Litters are not included.

Specs:

- Folded: H 34 in. x W 25.75 in. x D 5 in.
- Extended Height: 33 in.
- Base Width: 23.5 in.
- Weight: 8 lb 8 oz
- Payload: 1200lbs. (set of 2)
- Unit of Issue: Set of 2



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS

WITH A

MISSION®