

NASOPHARYNGEAL AIRWAY

REV051718



Nasopharyngeal Airway with Lubricant

This product shows North American Rescue's attention to details. This NPA was designed to be dependable, efficient, comfortable, and secure.

Dependable. Remains steady during transport. It is not easily dislodged.

Efficient. The product includes a conveniently attached sealed packet of surgical lubricant to facilitate easy insertion.

Comfortable. More compatible with human anatomy than the standard oropharyngeal airway, this product is recommended for use as a nasal airway adjunct in semi-conscious or unconscious patients with an intact gag reflex.

Features:

- Not easily dislodged during transport
- Better tolerated than an oropharyngeal airway by semi-conscious or unconscious patients with an intact gag reflex
- Packaged with surgical lubricant attached
- Multiple sizes available



Specifications:

- Packaged:
L 6 in. x W 1.6 in. x D 0.5 in.
- Device Sizes available:
22 French - 36 French
- Weight: 0.5 oz

Size	Item#	NSN#
22F	ZZ-0424	
24F	ZZ-0422	
26F	ZZ-0387	
28F	ZZ-0034	6515-01-529-1187
30F	ZZ-0423	
32F	ZZ-0388	
34F	ZZ-0776	
36F	ZZ-0777	



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®