

CASUALTY RESTRAINT STRAPS

The Casualty Restraint Strap with 2-piece plastic, or metal, buckle with or without a pair of metal speed clips is constructed of heavy-duty, 2 in. nylon webbing for fast and effective immobilization. The webbing is available in 2 colors: black and orange. Designed for use with standard backboards and basket stretchers.

CASUALTY RESTRAINT STRAPS

OPTION	COLOR	ITEM #
Plastic Buckle	Black	50-0023
Plastic Buckle	Orange	50-0024
Plastic Buckle w/ Speed Clip	Black	50-0026
Plastic Buckle w/ Speed Clip	Orange	50-0025
Metal Buckle	Black	50-0027
Metal Buckle	Orange	50-0028
Metal Buckle w/ Speed Clip	Black	50-0029
Metal Buckle w/ Speed Clip	Orange	50-0030

DIMENSIONS

With Plastic Buckle:

- L 5 ft x W 2 in.
- Weight: 4 oz.

With Plastic Buckle/Speed Clip:

- L 5 ft x W 2 in.
- Weight: 78 oz.

With Metal Buckle:

- L 5 ft x W 2 in.
- Weight: 8 oz.

With Metal Buckle/Speed Clip:

- L 5 ft x W 2 in.
- Weight: 11.7 oz.



Learn more about this and other **Products with a Mission®** at www.NARescue.com



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277