

PUBLIC ACCESS INDIVIDUAL BLEEDING CONTROL KIT *Nylon*

REV061722

Control the bleed with these easy-to-use tools!

Life threatening bleeding can occur from everyday emergencies such as vehicle crashes, incidents involving farming equipment or public transportation, accidental injuries from natural disasters like tornadoes, and inflicted injuries from man-made disasters such as shootings and bombings. Public Access Individual Bleeding Control Kits by North American Rescue are designed to provide essential equipment that empowers the general public to take action as immediate responders in stopping life threatening bleeding.

Contained in a nylon carry bag, this life saving kit provides intuitive and easy-to-use tools that are proven to help to save lives such as the fastest, safest and most effective prehospital tourniquet in the world, the C-A-T®, quality gauze, shears, and pressure dressings/bandages. Included are step-by-step illustrated “just in time” pictorial instructions, so even untrained providers can help save a life.

NAR Individual Bleeding Control Kits are available in Basic, Intermediate, Advanced, and Advanced BCD configurations.

BASIC:

The NAR Individual Bleeding Control Kit provides life saving bleeding control equipment such as tourniquets, pressure dressings and gauze bandages and puts them into the hands of the public to help control bleeding and save lives.

INTERMEDIATE:

Includes all equipment found in the Basic Individual Bleeding Control Kit with the added capability to treat penetrating chest wounds with the HyFin® Vent Compact Chest Seal Twin Pack (for entrance and exit wounds).

ADVANCED / ADVANCED BCD:

Includes all equipment found in the Intermediate Individual Bleeding Control Kit with the added capability to treat major bleeding with a hemostatic dressing that enhances blood clotting abilities—for faster bleeding control.



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277



Features:

- Designed to provide essential equipment that empowers the general public to take action as immediate responders in stopping life threatening bleeding
- Provides intuitive and easy-to-use tools that are proven to help to save lives
- Kits are available in Basic, Intermediate, Advanced, and Advanced BCD configurations
- Packaged in a compact, red nylon carry bag

Specifications:

- H 7 in. x W 5 in. x D 2.5 in.
- Weights:
 - Basic: 15 oz
 - Intermediate: 16.5 oz
 - Advanced: 17.4 oz
 - Advanced with BCD: 17.4 oz
 - Bag Only: 3.41 oz

Individual Bleeding Control Kit

Level	Item#
Basic	80-0452
Intermediate	80-0453
Advanced	80-0454
Advanced with BCD	80-0523
Bag Only	80-0459W

PUBLIC ACCESS INDIVIDUAL BLEEDING CONTROL KIT *Nylon*

REV061722

Every PUBLIC ACCESS INDIVIDUAL BLEEDING CONTROL KIT

includes these core items:



QTY
1 C-A-T®
Tourniquet



QTY
1 Responder
6 in. ETD



QTY
2 Wound
Packing Gauze




PAIR
2 Responder
Nitrile Gloves, Lg



QTY
1 Permanent
Marker, Small



QTY
1 Responder
Trauma Shears, Lg



QTY
1 Just in Time
Instructions



QTY
1 Survival
Blanket

Individual Bleeding Control Kit

Level	Item#
Basic	80-0452
Intermediate	80-0453
Advanced	80-0454
Advanced with BCD	80-0523
Bag Only	80-0459W

PLUS these Individual Components based upon your configuration of choice:



All these kits include:

- 1 x C-A-T® Tourniquet
- 1 x 6 in. NAR Responder ETD (Emergency Trauma Dressing)
- 2 x NAR Wound Packing Gauze
- 2 x Pair, Responder Nitrile Gloves, Lg
- 1 x Responder Trauma Shears, 7.25 in.
- 1 x NAR Survival Blanket
- 1 x Permanent Marker, Small
- 1 x Just in Time Instruction Card
- 1 x PVC Bleeding Control Patch
- 1 x Red Nylon Bag

Intermediate Kit also includes:

- 1 x HyFin® Compact, Vent Twin Pack

Advanced Kit also includes:

- 1 x HyFin® Compact, Vent Twin Pack
- 1 x Combat Gauze® LE, Z-Fold Hemostatic Bandage

Advanced with BCD Kit also includes:

- 1 x HyFin® Compact, Vent Twin Pack
- QuikClot® Bleeding Control Dressing



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®