

# EAGLE BLS IFAK *(Individual First Aid Kit)*

Combat Casualty Response Kits® **CCRK®**

**The Basic Life Support  
Variation of the Eagle  
IFAK comes equipped  
with the Core Essentials  
for Saving Lives in  
Tactical Environments**

The Basic Life Support Variation of the Eagle IFAK is one of NAR's smallest IFAKs, complete with essential life-saving equipment for point-of-wounding and self-aid/buddy-aid care to address penetrating, blast or other traumatic injuries. It features a Quick Pull Tab, the carrying bag easily opens to reveal the removable kit with just one hand. The kit has a unique, detachable shock-cord inner panel that can be adjusted to carry different components and provides easy access to your IFAK contents by simply pulling the color-coded tab to remove the panel from the carrying bag. The compact size reduces the space requirement needed for an IFAK on your vest or belt.

A compact, versatile individual first aid kit that offers both MOLLE/PALS – style connectors to attach to the assault vest/body armor and an upper or lower mounted 2.25 in. belt loop attachment to secure the kit in the high or low position on your duty/pistol belt, all in one carrying platform.



**Available  
Colors:**



Color	Item#	NSN#
BLK	85-0177	
COY	85-0178	
MTC	85-0179	
ODG	85-0180	

*Quick Pull Tab  
easily opens the kit  
with just one hand*



**NORTH AMERICAN RESCUE®**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

PRODUCTS WITH A  
**MISSION®**

# EAGLE BLS IFAK

## Features:

- Contains the essential equipment needed for self-aid/buddy-aid to address top leading causes of potentially preventable death in a tactical environment
- Can connect to the assault vest/body armor or to a duty/pistol belt
- Quick Pull Tab easily opens the kit with just one hand
- Has a unique, detachable shock-cord inner panel that can be adjusted to your internal components
- Inner panel is easily removable for easy access to your kit

## Specifications:

- Packaged: H 8 in. x W 4 in. x D 3 in.
- Weight: 1 lb 2 oz

## Kit Contents:

- 1 x C-A-T® (Combat Application Tourniquet®) Black
- 1 x 4 in. Flat ETD™
- 1 x S-rolled Gauze (4.5 in. x 4.1 yd)
- 1 x HyFin® Vent Chest Seal Twin Pack
- 1 x pair Bear Claw® Nitrile Trauma Gloves, lg.

*The Eagle IFAK's compact size reduces the space requirement needed for an IFAK on your vest or belt.*



**NORTH AMERICAN RESCUE®**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

PRODUCTS WITH A  
**MISSION®**

