

# QUIKCLOT COMBAT GAUZE LE

REV052620



## Effective, simple and safe when trying to control life-threatening hemorrhage...

QuikClot® Combat Gauze LE™ is a soft, white, sterile, non woven 3 in. by 4 yd Z-folded gauze impregnated with kaolin, an inert mineral that does not generate heat.

QuikClot® Combat Gauze LE™ is individually wrapped in an easy rip, military grade foil pouch. Indicated for temporary external control of traumatic bleeding,

QuikClot® Combat Gauze LE™ is flexible and pliable and contours to all wounds. Effective, simple and safe to use when trying to control life-threatening hemorrhage from wounds that are not amenable to tourniquet application.

### Features:

- Z-Fold, Vacuum Packed for easy application
- Effective, Simple, and Safe – Does not generate heat
- Easy to use and remove
- Flexible and pliable
- 3 in. x 4 yd (12 ft)

### Dimensions:

- Packaged: H 5 in. x W 5.5 in. x D 0.625 in.
- Gauze Dimensions: W 3 in. x L 4 yd (12 ft)
- Weight: 0.8 oz

### QuikClot Combat Gauze LE

Item #	NSN#
30-0109	



**NORTH AMERICAN RESCUE®**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

**PRODUCTS**  
 WITH A  
**MISSION®**