

IPOK[®]

Individual Patrol Officer Kit



Compact hemorrhage control kit designed to fit into a BDU pocket...

The Individual Patrol Officer Kit (IPOK) from North American Rescue is designed to provide personnel with a compact and durable individual hemorrhage control kit to treat bleeding from penetrating and other traumatic injuries.

These kits are packaged for small cube space and designed to fit into a BDU pocket, vest pouch, or individual bag, which allows personnel to keep a compact bleeding control kit on their person, where it is needed most. The contents are vacuum sealed in a rugged and durable packaging that is easy to open.

Features:

- Compact, vacuum-sealed hemorrhage control kit designed to be carried on your person where it is needed most
- Easy-to-open streamlined kit
- Designed to fit into a BDU pocket, vest pouch, or individual bag
- Your Choice of S-Rolled Gauze or Hemostatic


Specifications:

- Packaged: H 4 in. x W 6.5 in. x D 2.75 in.
- Weight: 7.4 oz



IPOK[®] Individual Patrol Officer Kit

Hemostatic/Dressing	Item#	NSN#
S-Rolled Gauze	80-0167**	
Chito Gauze	80-0223**	
Combat Gauze LE	80-0168	



Kit Contents:

- 1 x pair Black Talon[®] Nitrile Trauma Gloves, lg.
- 1 x C-A-T[®], Black (Combat Application Tourniquet[®])
- 1 x Your Choice of Gauze (S-Rolled, ChitoGauze Pro or Combat Gauze LE)
- 1 x Flat ETD[™] 4 in. Emergency Trauma Dressing



NORTH AMERICAN RESCUE[®]
www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION[®]

IPOK® Dressing Options

REV011119



S-Rolled Gauze Option 1 for the IPOK® (Item #80-0167)

- 4.5 in. x 4.1 yd Sterile, 6-ply premium cotton gauze
- Unique and convenient packaging configuration and dispenser allows controlled application as a complete unit or in a continuous linear feed for wound packing Unused portion remains protected from contamination during application
- Used as backing gauze for hemostatic agents, in conjunction with a compression bandage to control hemorrhage



Chito Gauze PRO Gauze Option 2 for the IPOK® (Item #80-0223)

- Z-folded & packaged in 3 in. x 4 yd lengths
- Offers antimicrobial properties against 26 bacterial organisms
- Non-woven, polyester/rayon blend medical gauze is coated with chitosan to maximize its performance and control severe bleeding
- Featuring an embedded x-ray detectable strip



Combat Gauze LE Gauze Option 3 for the IPOK® (Item #80-0168)

- 3 in x 4 yd (12 ft)
- Z-Fold, Vacuum Packed for easy application
- Effective, Simple, and Safe – Does not generate heat
- Easy to use and remove
- Flexible and pliable

**Choose one
of these three
dressing options
to complete your
IPOK® Kit**



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH  **MISSION®**