

T.A.L.O.N. MANUAL NEEDLE SET

REV120117



**Caution: Federal (USA) law restricts this device to sale by or on the order of a practitioner licensed to use or order use of the device or Federal (USA) law prohibits dispensing without a prescription.*

T.A.L.O.N. Tactically Advanced Lifesaving IO Needle

Manual I/O insertion at all 7 extremity sites, including the sternum. Includes 15 gauge, 304 stainless steel needle/catheter, sharps protector, EZ-Connect® extension set, I/O wristband, quick reference card & Sternal Locator.

Seven site access includes: the sternum, bilateral proximal humerus, bilateral proximal tibia, and bilateral distal tibia. Using IO access, battlefield responders can deliver life-saving fluids and medications rapidly and effectively.

Features:

Allows for manual I/O insertion at all 7 extremity sites:

- Sternum
- Bilateral proximal humerus
- Bilateral proximal tibia
- Bilateral distal tibia




T.A.L.O.N. Manual Needle Set

Color	Size	Item #	NSN#
Brown	38.5mm	30-0132	

Specifications:

- Packaged: L 6 in. x W 3.25 in
- Weight: 1.5 oz



EZ-IO Products are restricted to Federal Purchases only. Please contact Customer Service for details.



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®