

# RAPTOR RESPONSE SHEARS (LEATHERMAN®)

## *A Must Have Multi-Tool...*

Gear up for the unexpected with the Leatherman Raptor Response Shears. These compact, portable shears cut through materials while the micro-serrations on the inner blade keep fabrics from slipping and binding. It also features a ring cutter, ruler, oxygen tank wrench, and pocket clip, so you're ready to get to work in any situation. The handle grips are contoured, finished with a durable ceramic coating, and available in two color options. Made in the USA and backed by Leatherman's 25-year warranty.

## FEATURES/SPECIFICATIONS

- Stainless Steel folding shears
- Integrated ring cutter
- 5cm ruler
- Oxygen tank wrench
- Lanyard connection fabricated into handle
- Removable pocket clip

## **DIMENSIONS**

- Closed Length: 4.3 in.
- Open Length: 7.25 in.
- Shear Length: 1.9 in.
- Width: 1.7 in | 4.3 cm
- Overall Thickness: .85 in.
- Weight: 5.5 oz

## **LEATHERMAN® RAPTOR MULTI-TOOL**

ITEM #	OPTION	NSN #
70-0328	Cement	---
70-0330	Navy	---



**Learn more** about this and other **Products with a Mission®** at **[www.NARescue.com](http://www.NARescue.com)**



**NORTH AMERICAN RESCUE®**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277