

COAG IT NP7 NASAL PACK

REV011122



The COAG IT NP7 Nasal Pack is a seven cm sterile, anatomically shaped nasal packing device for use in treating nosebleeds by medical professionals. This device was designed for patient comfort and ease of application. The smooth sides, rounded edges, and anatomical shape minimize patient discomfort during insertion and removal. While inserted the compressed sponge softly expands against the nasal cavity stopping the bleed. The blue handle removes the need for extraction strings to be taped to the patient's face and is adjustable for increased comfort.

Features:

- 7 cm single use Nasal Pack
- Anatomically shaped for maximum comfort
- Physician designed epistaxis control
- Expands to apply pressure stopping nasal bleed
- Blue handle for easy extraction

Specifications:

- Nasal Pack Length: 7 cm
- Nasal Pack Weight: 3.2 oz



Box of 10

Item #	Quantity
ZZ-1483	1
ZZ-1484	10



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®