

# TACO KYDEX TOURNIQUET HOLDER SHIRT SHIELDS

REV082020



## Protect Your Shirt...

Upgrade your Taco Kydex® Tourniquet Holders to protect your shirt from undue wear!

Constructed of aircraft-grade Kydex®, these rugged tourniquet shirt shields are designed to enhance both the Textured and Slick versions of the Racons Kydex Tourniquet Holders by protecting clothing and fabric from abrasion caused by the hook and loop of the tourniquet. The shirt shield comes with longer length screws for simple attachment to the holder using a phillips head driver (not included). The Kydex material is moisture and solvent resistant, so it doesn't add weight from rain, humidity or sweat.






## Features:

- Designed to protect clothing or fabric from abrasion when carrying windlass tourniquets on a belt, rig or bag
- Constructed of aircraft-grade Kydex® a smooth finish for easy insertion and draw (.080 thickness)
- Resistant to moisture and solvents & won't add weight from rain, humidity or sweat
- Will not damage the tourniquet or the platform it's attached to

## Specifications:

- Dimensions: 2 in. x 1 in. x 3.25 in.
- Weight: 0.06 lb

## Taco Kydex Tourniquet Holder Shirt Shield

Item #	NSN#	Color	
30-0279		Black	
30-0280		Blue	
30-0281		Coyote	
30-0282		Gray	
30-0283		OD Green	



**NORTH AMERICAN RESCUE®**

www.NARescue.com • 888.689.6277

**PRODUCTS**

WITH A

**MISSION®**