

Clear Polycarbonate  
WALL STATIONS

- Includes Individual Bleeding Control Kits packaged in red nylon bags or Vacuum-Sealed pouches.
- The 5 & 8 Pack Stations also include an easy to-open nylon carrier and one (1) Responder QuikLitter™ for moving injured casualties.

SPECIFICATIONS

**5 & 8 Pack Stations:** H 18 in. x W 15.5 in. X D 10 in.  
Empty Case Wt: 15.0 lbs (includes rack)  
Large Case Item# ZZ-0822

**Individual Stations:** H 9.25 in. x W 8.50 in. X D 4.00 in.  
Empty Case Wt: 1.68 lb  
Small Case Item# ZZ-1000

8-PACK CLEAR POLYCARBONATE WALL STATIONS

LEVEL	#KITS	BAG	WGT	ITEM #
Basic	8	Vacuum	21.8 lb	80-0471
Intermediate	8	Vacuum	22.6 lb	80-0472
Advanced BCD	8	Vacuum	23 lb	80-0526
Advanced	8	Vacuum	23 lb	80-0473
Basic	8	Nylon	23.5 lb	80-0474
Intermediate	8	Nylon	24.3 lb	80-0475
Advanced BCD	8	Nylon	27.4 lb	80-0527
Advanced	8	Nylon	27.4 lb	80-0476

5-PACK CLEAR POLYCARBONATE WALL STATIONS

LEVEL	#KITS	BAG	WGT	ITEM #
Basic	5	Vacuum	20.7 lb	80-0581
Intermediate	5	Vacuum	20.7 lb	80-0582
Advanced BCD	5	Vacuum	21.5 lb	80-0584
Advanced	5	Vacuum	21.5 lb	80-0583

IND. CLEAR POLYCARBONATE WALL STATIONS

LEVEL	#KITS	BAG	WGT	ITEM #
Basic	1	Vacuum	2.1 lb	80-0938
Intermediate	1	Vacuum	2.7 lb	80-0939
Advanced BCD	1	Vacuum	2.8 lb	80-0940
Basic	1	Nylon	2.6 lb	80-0935
Intermediate	1	Nylon	2.7 lb	80-0936
Advanced BCD	1	Nylon	2.8 lb	80-0937

Alarms (#ZZ-1075) is also available for the Clear Polycarbonate Wall Cabinet (See page 9)



8 Pack Station



Individual Kit Station



Small Individual Clear Case: ZZ-1000

Polycarbonate  
CLEAR CASE  
features

Designed to be mounted on a wall at a designated location or to be co-located with Automatic External Defibrillators (AEDs) to provide publically accessible medical supplies to control major bleeding



5 Pack Station



Large Multi-Pack Clear Case: ZZ-0822

- Durable, clear polycarbonate cabinet with rounded corners provides both physical and environmental protection for medical equipment
- Large case features wire shelf to support the contents
- External eyelet allows a secondary safety seal to easily identify if the station has been opened or tampered with
- External mounting tabs and hardware enable easy installation on any wall
- Clear design allows visibility at any angle
- Cases are water resistant. To best protect your critical life savings products cases are best mounted under cover and out of the direct elements



Semi-Recessed Metal  
WALL STATIONS

- Includes Individual Bleeding Control Kits packaged in either red nylon bags or Vacuum-Sealed pouches stored within an easy to-open nylon carrier and one (1) Responder QuikLitter™ for moving injured patients
- Outside Dimensions: H 17.5 in. x W 17.5 in. x D 1 in.  
Inside Dimensions: H 14.13 in. x W 14.13 in. x D 6.88 in.  
Empty Case Wt: 10 lb (Empty Case Item# ZZ-0921)

8-PACK SEMI-RECESSED METAL WALL STATIONS

LEVEL	#KITS	BAG	WGT	ITEM #
Basic	8	Vacuum	16.8 lb	80-0861
Intermediate	8	Vacuum	17.6 lb	80-0862
Advanced BCD	8	Vacuum	18 lb	80-0863
Advanced	8	Vacuum	18 lb	80-0864
Basic	8	Nylon	18.5 lb	80-0865
Intermediate	8	Nylon	19.3 lb	80-0866
Advanced BCD	8	Nylon	15.7 lb	80-0867
Advanced	8	Nylon	15.7 lb	80-0868

5-PACK SEMI-RECESSED METAL WALL STATIONS

LEVEL	#KITS	BAG	WGT	ITEM #
Basic	5	Vacuum	15.7 lb	80-0577
Intermediate	5	Vacuum	15.7 lb	80-0578
Advanced BCD	5	Vacuum	16.5 lb	80-0580
Advanced	5	Vacuum	16.5 lb	80-0579

All Multi-Pack Wall Stations feature:  
Highly visible "Bleeding Control Kit" letters for fast ID.  
Instructions for Use are posted on the interior back panel



8 Pack Station

Semi-Recessed  
METAL CASE  
features

- Mounts permanently inside the wall with a low profile 1 in. trim
- White powder painted corrosive resistant galvanized steel cabinet with clear plexiglass window
- Stainless steel continuous hinge & chrome plated handle designed to open with one quick pull
- Cases are water resistant. To best protect your critical life savings products cases are best mounted under cover and out of the direct elements.



Semi-Recessed Metal Case: ZZ-0921





Low Profile  
Metal Cabinet:  
ZZ-1199

Designed to be mounted on a wall with a slim 4 in. clearance and featuring easy pull-to-open design for rapid public access to major bleeding control supplies

8 Pack Nylon Station

## Low Profile Metal WALL STATIONS

- Low profile metal cabinet mounts on the wall with a slim 4 in. clearance
- Designed to safely store 8 Individual Bleeding Control Kits in public areas
- Kits are available in red nylon bags or Vacuum-Sealed pouches and stored in Two (2) easily removable Mesh Panels each equipped with 4 pockets and nylon handle
- Features highly visible "Bleeding Control Kit" lettering for fast identification and convenient Instructions for Use on the interior door panel
- Cases are water resistant. To best protect your critical life savings products cases are best mounted under cover and out of the direct elements

### SPECIFICATIONS

H 15.5 in. x W 22.5 in. X D 4 in.  
Empty Cabinet Weight: 9.365 lb

Low Pro Cabinet with ID Stickers:  
Item# ZZ-1199

### 8-PACK LOW PROFILE METAL WALL STATIONS

LEVEL	#KITS	BAG	WGT	ITEM #
Basic	8	Vacuum	17.49 lb	80-1083
Intermediate	8	Vacuum	18.19 lb	80-1084
Advanced BCD	8	Vacuum	18.69 lb	80-1085
Basic	8	Nylon	17.49 lb	80-1080
Intermediate	8	Nylon	18.19 lb	80-1081
Advanced BCD	8	Nylon	18.69 lb	80-1082

Alarms (#ZZ-1075) is also available for the Low Profile Metal Cabinet (See page 9)

### Low-Pro CABINET features

- Low profile metal cabinet with 4 in. wall clearance
- Provides convenient slimline storage for medical equipment
- White powder painted corrosive resistant galvanized steel cabinet with stainless steel continuous hinge
- Magnetic closure provides quick easy access with one pull to open - No twist and pull



## STB KITS Nylon or Vacuum-Sealed

- Small, compact and easy to use
- Kits contain the same Tourniquet and Hemostatic Dressing technology recommended by the U.S. Military & the American College of Surgeons Committee on Trauma to control bleeding
- Available in Vacuum-sealed Pouch or Red Nylon Bag



### VAC STB KIT

#### ITEM #

85-1574

### SPECIFICATIONS

H 7 in. x W 5 in. x D 2.5 in.  
Weight: 8 oz

### VACUUM-SEALED STB KIT contains:



### SPECIFICATIONS

H 7 in. x W 5 in. x D 2.5 in.  
Weight: 15 oz

### NYLON STB KIT

#### ITEM #

85-2995

### NYLON STB KIT contains:



## Accessories

### Cabinet Security ALARMS

- Alarm compatible with Clear Polycarbonate and Low Profile Metal Cabinets
- Loud 120 decibel alarm sounds when magnet separates from base more than 0.5 in.
- Easy and convenient to install
- Requires no wiring
- Mounting hardware and instructions included
- Designed for indoor use only



### SPECIFICATIONS

W 1.63 in. x H 2.41 in. x D 0.72 in.  
Weight: 6.4 oz each

### CABINET SECURITY ALARMS

ITEM#	UNIT
ZZ-1075	Pack of 2

### 3-WAY SIGN Bleeding Control



### SPECIFICATIONS

H 7 in. x W 16.5 in. x D 6.75 in.  
Weight: 0.35 oz

- Highly visible for locating lifesaving bleeding control kits immediately
- 3 way "V" shape projects from wall for easy identification from multiple directions

### 3-WAY SIGN

#### ITEM #

ZZ-1117