## NAR® RESCUE SLEEVE II™

AND ACCESSORIES

## The Ultimate in Patient Immobilization...

The NAR Rescue Sleeve II™ is the industry's ultimate immobilization system supporting both vertical and horizontal hoisting. Simple and secure, the Rescue Sleeve II™ gives confidence to operators and comfort to patients. Used extensively in the field and refined for enhanced stability and overall functionality, it is standard issue to all U.S. Navy ships.



Learn more about this and other Products with a Mission® at www.NARescue.com



## • Easy to apply

- Replaceable parts are also available: Head Block (#50-0042 | NSN: 6530-01-477-9448) *Head and Chin Straps* (#50-0001 | NSN# 6515-01-477-8921) Carrying Case (# 50-0002 | NSN# 4240-01-477-9408)
- Standard issue to all U.S. Navy ships



- Folded: L 25.5 in. x W 13 in. x D 8 in.
- Open: L 73 in. x W 18 in. x D 2 in.
- Payload: 1000 lb (453.6 kg)
- Weight: 16 lb

## FEATURES/SPECIFICATIONS

- Vertical and horizontal lifting-strap attachments
- Removable head-securing blocks
- Adjustable straps secure head and chin
- Hoist-cable storage space
- Padded head flaps
- Backboard compartment
- Padded carrying handles





NOTE: Spine board sold separately. Hoisting straps & carabiners not included.

ITEM #	COLOR	NSN
50-0011	Black	6530-01-531-6649
50-0012	Orange	6530-01-477-8525

NAR® RESCUE SLEEVE II™

SEL#: 09ME-04-SPIN SEL#: 03WA-02-BORD



