

NAR® RESCUE SLEEVE II™

AND ACCESSORIES

The Ultimate in Patient Immobilization...

The NAR Rescue Sleeve II™ is the industry's ultimate immobilization system supporting both vertical and horizontal hoisting. Simple and secure, the Rescue Sleeve II™ gives confidence to operators and comfort to patients. Used extensively in the field and refined for enhanced stability and overall functionality, it is standard issue to all U.S. Navy ships.



NOTE: Head Block and Carry Case are included and are also available separately for resupply. Head and Chin straps are sold separately.

Learn more about this and other Products with a Mission® at www.NARescue.com



FEATURES/SPECIFICATIONS

- Vertical and horizontal lifting-strap attachments
- Removable head-securing blocks
- Adjustable straps secure head and chin
- Hoist-cable storage space
- Padded head flaps
- Backboard compartment
- Padded carrying handles
- Easy to apply
- Replaceable parts are also available:
Head Block (#50-0042 | NSN: 6530-01-477-9448)
Head and Chin Straps (#50-0001 | NSN# 6515-01-477-8921)
Carrying Case (# 50-0002 | NSN# 4240-01-477-9408)
- Standard issue to all U.S. Navy ships



NOTE: Spine board sold separately. Hoisting straps & carabiners not included.



DIMENSIONS

- Folded: L 25.5 in. x W 13 in. x D 8 in.
- Open: L 73 in. x W 18 in. x D 2 in.
- Payload: 1000 lb (453.6 kg)
- Weight: 16 lb

BLK

ORG

NAR® RESCUE SLEEVE II™

ITEM #	COLOR	NSN
50-0011	Black	6530-01-531-6649
50-0012	Orange	6530-01-477-8525

SEL#: 09ME-04-SPIN

SEL#: 03WA-02-BORD



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277